AMSCO

SOFIA 2



Buy a Box of Flu Tests, Get Sofia 2 Analyzer **FREE**

Also, Buy 6 boxes for the price of 5!

Box of 25 Flu Tests

Flu Test Item Code: QUI20218





Fluorescent Immunoassay Analyzer

Sofia 2 is a small bench top analyzer that uses advanced fluorescence detection with an ultraviolet LED energy source. Sofia 2 collects hundreds of data points on the test strip contained inside of a Sofia FIA Test Cassette. Sofia 2 uses proprietary algorithms to analyze the data, interpret the signals and determine the results automatically to give you an objective result you can rely on.

Sofia 2 offers an easy-to-use interface, customizable settings, flexible workflow, barcode technology and Virena® connectivity. Results are reported automatically on the screen and can be stored in Sofia 2, optionally printed on an external printer, or sent to Virena.

Sofia 2 has two modes that can be used depending on your workflow requirements.

WALK AWAY mode may be most convenient for reading a single patient sample. In this mode, Advance Result Technology (ART) can provide results in as few as 3 minutes. The patient sample is dispensed into the Test Cassette and the Test Cassette is immediately placed in Sofia 2. Sofia 2 will automatically develop the Test Cassette for the required time period (pre-programmed for each test), scan the Test Cassette periodically, and display the test result.

READ NOW mode may be most convenient for reading several patient samples. In this mode, the patient sample is dispensed into the Test Cassette. The user then manually times the development period outside of Sofia 2. This can be done on the counter or bench top using a timer. Once the development time is complete, the Test Cassette is placed in Sofia 2 where it is immediately scanned and the test result is displayed within one (1) minute.

Influenza testing within the first 48 hours of symptoms are the most critical. Taking antivirals within the first 48 hours can help shorten the length of the flu by 30%.*

*Influenza Testing, CPT Codes and 2019 Reimbursement Information on Back of Page.







PH: 708.403.2220 FAX: 708.403.7763 925 Lambrecht Drive Frankfort, IL 60423












The ART of Flu Testing

- Accurate, objective and automated results in **as few as 3 minutes**
- Secure, remote instrument management through the Virena ecosystem
- **Excellent performance** compared to molecular and viral culture methods
- Integrated data management automatically stores test and user history
- Flexible, dual mode testing and simple, CLIA-waived procedures for a variety of laboratory environments
- User and patient ID captured with **onboard barcode scanner**

Diagnosing and treating influenza quickly, preferably within the first 48 hours, can shorten the duration of flu symptoms.



Days After Symptom Onset



©2019 Alpha Medthrift Scientific Co.