# Hippurate Disk

(English)

## INTENDED USE

Remel Hippurate Disk is a reagent-impregnated disk recommended for use in qualitative procedures to detect the hydrolysis of sodium hippurate by beta-hemolytic streptococci, *Gardnerella vaginalis*, and *Campylobacter jejuni*.

## SUMMARY AND EXPLANATION

The hippurate hydrolysis test determines the ability of group B streptococci, as well as other bacteria, to enzymatically hydrolyze sodium hippurate.<sup>1,2</sup> Hwang and Ederer employed the ninhydrin reagent in a rapid method to detect hippurate hydrolysis by the end product glycine.<sup>3</sup> Jolly later outlined minimal criteria for identification of *G. vaginalis* including hippurate hydrolysis.<sup>4</sup> In 1989, Cacho et al. further evaluated the use of the hippurate disk in the identification of *Campylobacter* spp.<sup>5</sup>

## PRINCIPLE

Hippuric acid is hydrolyzed by the enzyme hippuricase to glycine and benzoic acid. Ninhydrin evokes a five-step reaction beginning with the deamination of glycine to form hydrindantin (a reduced form of ninhydrin), carbon dioxide, and ammonia. A condensation reaction occurs with hydrindantin, ammonia, and residual ninhydrin to produce the final purple colored complex.

## REAGENTS

Hippurate Disk (25 disks/vial)  
Reactive Ingredient: Sodium Hippurate

## PRECAUTIONS

This product is For *In Vitro* Diagnostic Use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, and media after their use. Directions should be read and followed carefully.

## STORAGE

This product is ready for use and no further preparation is necessary. Store product in its original container at 2-8°C until used. Allow product to equilibrate to room temperature before use. Do not incubate prior to use.

## PRODUCT DETERIORATION

This product should not be used if (1) the color has changed from white, (2) the expiration date has passed, (3) the desiccant has changed from blue to pink, or (4) there are other signs of deterioration. Protect disks from moisture by removing from the vial only those disks necessary for testing. Promptly replace the cap and return the vial to 2-8°C.

## SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>6</sup>

## MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Plastic test tubes, (7) Forceps, (8) Ninhydrin Reagent (REF R21238) or Bactidrop™ Ninhydrin (REF R21534), (9) Demineralized water.

## PROCEDURE

1. Add 0.1 ml of sterile demineralized water to a small plastic test tube.
2. Take 1-3 colonies of a fresh isolate and emulsify in the water.
3. Using forceps drop a Hippurate Disk into the suspension.
4. Incubate aerobically for 2 hours at 35-37°C.
5. Dispense 2 drops of Ninhydrin Reagent or 3 drops Bactidrop™ Ninhydrin into each test tube, mix, and reincubate aerobically for 30 minutes at 35-37°C.
6. Observe for a blue to purple color development.

## INTERPRETATION

Positive Test - Blue to purple color development within 30 minutes  
Negative Test - Light gray color development or no color change within 30 minutes

## QUALITY CONTROL

All lot numbers of Hippurate Disk have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

### CONTROL

*Gardnerella vaginalis*  
ATCC® 14018  
*Streptococcus agalactiae*  
ATCC® 12386  
*Streptococcus pyogenes*  
ATCC® 19615

### INCUBATION

Aerobic, 2.5 h  
@ 35-37°C  
Aerobic, 2.5 h  
@ 35-37°C  
Aerobic, 2.5 h  
@ 35-37°C

### RESULTS

Positive  
Positive  
Negative

## LIMITATIONS

1. Positive and negative controls should be included with each group of tests to better correlate the color changes produced by both positive and negative tests.
2. Plastic tubes should be used, as glass may yield a false-negative result.

## BIBLIOGRAPHY

1. Greenwood, J.R. and M.J. Pickett. 1979. J. Clin. Microbiol. 9:200-204.
2. MacFaddin, J.F. 2000. Biochemical Tests for Identification of Medical Bacteria. 3<sup>rd</sup> ed. Lippincott Williams & Wilkins, Philadelphia, PA.
3. Hwang, M. and G.M. Ederer. 1975. J. Clin. Microbiol. 1:114-115.
4. Jolly, J.L.S. 1983. J. Clin. Pathol. 36:476-478.
5. Cacho, J.B., P.M. Aguirre, A. Hernanz, and A.C. Velasco. 1989. J. Clin. Microbiol. 27:359-360.
6. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.L. Landry, and M.A. Pfalter. 2007. Manual of Clinical Microbiology. 9<sup>th</sup> ed. ASM Press, Washington, D.C.

## PACKAGING

REF R21085, Hippurate Disk ..... 25 Disks/Vial


## Symbol Legend

<b>REF</b>	Catalog Number
<b>IVD</b>	In Vitro Diagnostic Medical Device
<b>LAB</b>	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
<b>LOT</b>	Batch Code (Lot Number)
	Use By (Expiration Date)
<b>EC REP</b>	Authorized European Representative

ATCC® is a registered trademark of American Type Culture Collection.  
IFU 21085, Revised December 5, 2007

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