# USA E8 CENTRIFUGES



#### **ADVANTAGES**

Accurate speeds for specific g-forces Proper separations without cell damage Quiet, vibration-free operation Auto-brake stops the rotor in less than 30 seconds Designed, Engineered, and Assembled Proudly in America

### **FEATURES**

Angled rotor, 8-tube capacity 3,500 rpm Holds 3ml to 15ml size tubes 12 volt DC maintenance-free motor Clear view port in lid for using tachometer Suction-cupped feet prevent slipping 110/220v auto-switching power adapter (CE, UL, cUL approved adapter) Lid safety shut-off switch Maximum RCF=1,534g E8 Fixed & E8 Variable: Auto-off 30-minute timer with bell

#### WEIGHT AND DIMS

E8 Touch: Touch screen, programmable

Height: 9.5" (241mm) Length: 11" (279mm) Width: 11" (279mm) Weight: 9.5 lbs (4.3kg)



The USA E8 Centrifuge is designed, engineered, and assembled at LW Scientific's ISO-certified, FDA-registered manufacturing facility in Lawrenceville, GA. Expect American quality as well as customary LW Scientific support. This model is loaded with features, such as higher speeds for faster and cleaner separations, quieter operation for peaceful work environments, lower chamber temperatures for cooler test-tube samples, shorter height to fit under any cabinet, and auto-calibrating speed function for precision results. It even has an auto-brake to stop the rotor in less than 30 seconds to save time in busy labs.



Designed, Engineered & Assembled Proudly in America

# **STANDARD E8**

### INCLUDES

Eight-place rotor Eight 15ml tube sleeves Eight 13x75mm tube sleeve inserts




		USA E8 Fixed-Speed, Angled 8-Place Centrifuge with Timer
		USA E8 Variable-Speed, Angled 8-Place Centrifuge with Timer
T	E8C-U8AT-15T3	USA E8 Digital-Speed, Touch Screen, Angled 8-Place Centrifuge

# **COMBINATION E8**

#### INCLUDES

Eight-place rotor Eight 15ml tube sleeves Eight 13x75mm tube sleeve inserts Six green tube inserts to spin 1ml-2ml Two 4-place Crit Carrier inserts (Spins 75mm capillary tubes) EZ Reader Card




F	E8C-U8AF-CR03	USA E8 Fixed-Speed, Angled 8-Place Centrifuge with Crit Carriers
V	E8C-U8AV-CR03	USA E8 Variable-Speed, Angled 8-Place Centrifuge with Crit Carriers
T	E8C-U8AT-CRT3	USA E8 Digital-Speed, Touch Screen, Angled 8-Place Centrifuge with Crit Carriers

## **PORTABLE E8**

### INCLUDES

Eight-place rotor Eight 15ml tube sleeves Eight 13x75mm tube sleeve inserts Car seat Car adapter cord



F	E8C-U8AF-150P	USA E8 Fixed-Speed, Angled 8-Place Centrifuge, Portable
V	E8C-U8AV-150P	USA E8 Variable-Speed, Angled 8-Place Centrifuge, Portable
Τ	E8C-U8AT-15TP	USA E8 Digital-Speed, Touch Screen, Angled 8-Place Centrifuge, Portable