



Magnifiers

UV Light

The ideal light for numerous applications.

Designed for a variety of uses, the Burton UV (Woods) Light is also available in a model with a magnifier and a model with a magnifier and white and UV bulbs. With handheld flexibility, large magnification and a five-year limited warranty, this UV light provides great value for your facility.

Burton's UV lights are depended on worldwide for performance and durability.

PHILIPS
burton

Features

- Special fluorescent bulbs produce UVA at approximately 365 nm wavelength
- 4-light magnifier models include optically ground and polished glass 5 diopter magnifier lens with 2.25X total power magnification

Summary of technical data

Rated life of bulb:	6000 hours
Diopter of magnifying lens:	5
Focal length:	8" (20 cm)
Type of bulb:	UV or white
Rated Life of LED Lamp:	50,000+ hours
Power:	120V
Product warranty:	5-year limited warranty
Origin:	Assembled in USA

Ordering information

Model	Description
UV501	UV Light, 2 UV Bulbs
UV502	UV Light with Magnifier, 2 UV Bulbs
UV503	UV Light with Magnifier, 2 White and 2 UV Bulbs

Accessories	Description
1007090	UV Light Wall Bracket
0001127PK	UV Replacement Bulbs, White, 4/Box
1003073PK	UV Replacement Bulbs, UV, 4/Box

UV Light product dimensions



© 2017 Koninklijke Philips N.V. All rights reserved.
Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.



Phone (800) 444-9909, Fax (800) 765-1770

LIT-0013, Rev 5

Philips Burton
11500 Melrose Avenue
Franklin Park, IL 60131
Burtonmedical.com

PHILIPS
burton