

SURGICAL



Section Index



For complete alphabetical index see pages 598 to 614 For complete numerical index see pages 616 to 639 NewProducts..... NP1-NP21


ı	Page Nos
Diagnostic Instruments	. 1-4
Anesthesia Instruments	. 5-7
Sterilizing Aids	. 8-10
Instrument Care Products	. 11
Tubing Forceps	. 12
Scalpel Blades and Handles	. 13-16
Disposable Scalpels	. 17
Operating Scissors	. 18-50
Bandage and Utility Scissors	
Cast and Utility Shears	. 53
T-C Carb-N-Sert Scissors	. 54-59
SuperCut Scissors	. 60-63
Dressing and Tissue Forceps	64-82
Splinter Forceps	. 83-87
T-C Carb-N-Sert Dressing and	
Tissue Forceps	
Hemostatic Forceps	. 90-108
Towel Clamps	. 109
Sponge Forceps	.110-111
Needle Holders	
T-C Carb-N-Sert Needle Holders	. 122-130
Surgical Needles	
Suture Sundries	. 135-143
Probes and Directors	. 144-146
Suction Tubes	. 147-149
Vein Strippers	. 150
Retractors	. 151-187
Hooks	188-194

	Page Nos
Thyroid-Goiter Instruments	195-197
Trocars and Biopsy Needles	
Gall Bladder Instruments	
Intestinal Instruments	206-214
Micro Surgery Instruments	215-224
Ophthalmic Instruments	225-268
Ear Instruments	269-288
Nasal Instruments	
Plastic Surgery Instruments	307-346
Oral Instruments	347-354
Tonsil Instruments	355-361
Laryngeal Instruments	362-364
Bronchial and Esophageal	
Instruments	
Tracheal Instruments	367-372
Atraumatic Instruments	373-387
Thoracic and Cardiovascular	
Instruments	388-408
Neurosurgical Instruments	
Orthopedic Instruments	
Cast Instruments	466-468
Rectal Instruments	469-481
Urological Instruments	482-488
Gynecological and Obstetrical	
Instruments	
Dermal Instruments	
Post Mortem Instruments	
Suggested Surgical Instrument Sets	554-597



New Products Example 1





UNIVERSAL SCISSORS

- Highest quality German stainless steel construction
- Ideal for cutting bandages, gauze, tough clothing, cast material and metal sheets
- One serrated blade and comfortable contoured plastic handles
- Available in the popular 7-1/2" design and choose from seven colored handles

Part # Description

- 5-1000 BLACK - Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, black plastic finger rings, serrated blade. Autoclavable.
- RED Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, red 5-1025 plastic finger rings, serrated blade. Autoclavable.
- ORANGE Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, 5-1026 orange plastic finger rings, serrated blade. Autoclavable.
- BLUE- Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, blue 5-1027 plastic finger rings, serrated blade. Autoclavable.
- PURPLE Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, 5-1028 purple plastic finger rings, serrated blade. Autoclavable.
- GREEN Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, 5-1029 green plastic finger rings, serrated blade. Autoclavable.
- YELLOW Bandage & Utility Universal Scissors 7-1/2" (19.1 cm), with needle destroyer, 5-1030 yellow plastic finger rings, serrated blade. Autoclavable.





MILTEX® STERILE DISPOSABLE BIOPSY PUNCHES

- Individually packaged in color-coded, transparent, sturdy blister pack to protect cutting edge
- Seamless stainless steel razor sharp cutting edge
- Ribbed handle provides maximum control and comfort
- Sterile ready to use every time
- Packed 50 Biopsy Punches to a Dispenser Box
- New 1/2 sizes

 New 1/2 sizes

 ■

Part #	Size	Part #	<u>Size</u>
33-31AA	1.0 mm	33-33	3.5 mm
33-31A	1.5 mm	33-34	4.0 mm
33-31	2.0 mm	33-35	5.0 mm
33-31B	2.5 mm	33-36	6.0 mm
33-32	3.0 mm	33-37	8.0 mm




33-38 Assortment Box 10 each of the following sizes: 2mm, 3mm, 4mm, 5mm, 6mm

See page 526 for other related products.







MILTEX® STERILE DISPOSABLE SAFETY SCALPEL

- Omplies with OSHA directive for the Needlestick Safety and Prevention Act
- Stainless steel, sterile, single-use
- User activated and controlled mechanism no springs to jam
- Audible confirmation when securely in desired position
- Metric ruler assists in determining length of incision
- Packed 10 per Box

Part #	Size
4-510	Size #10, Sterile, Single-Use
4-511	Size #11, Sterile, Single-Use
4-515	Size #15, Sterile, Single-Use






STEINMANN I.M. PINS (FOR INTRAMEDULLARY FIXATION)

- Ideal for immobilization of long bone fractures
- Manufactured from highest quality 316LVM stainless steel that is insert when implanted
- Precision ground trocar points with mirror finish
- Individually packaged and protected in clear cylinder tube

Part #	Size	
27-182	Double Trocar Smooth - 1/16" x 7"	See pages 444-445 for Kirschner
27-183	Double Trocar Smooth - 5/64" x 7"	Wire (K-Wire), Pin Cutters, Pin
27-184	Double Trocar Smooth - 3/32" x 8"	Pullers and other accessories.
27-185	Double Trocar Smooth - 7/64" x 8"	
27-186	Double Trocar Smooth - 1/8" x 9"	
27-187	Double Trocar Smooth - 9/64" x 9"	
27-188	Double Trocar Smooth - 5/32" x 10"	
27-189	Double Trocar Smooth - 3/16" x 12"	
27-190	Double Trocar Smooth - 1/4" x 12"	
27-192	Double Trocar Partial Threaded - 1/16" x 7"	
27-193	Double Trocar Partial Threaded - 5/64" x 7"	
27-194	Double Trocar Partial Threaded - 3/32" x 8"	
27-195	Double Trocar Partial Threaded - 7/64" x 8"	
27-196	Double Trocar Partial Threaded - 1/8" x 9"	
27-197	Double Trocar Partial Threaded - 9/64" x 9"	
27-198	Double Trocar Partial Threaded - 5/32" x 10"	
27-199	Double Trocar Partial Threaded - 3/16" x 12"	
27-200	Double Trocar Partial Threaded - 1/4" x 12"	
27-201	Miltex Pin and Wire Gauge - measures Steinmann Pins screws. (Includes 6" ruler)	s, K-Wire, drill bits and bone







RECTAL INSTRUMENTS


Part # Description

28-150XL McGIVNEY Hemorrhoidal Ligator, 10" shaft, includes Ligating Rings and Loading Cone - Large

28-150A McGIVNEY Hemorrhoidal Ligator Set - includes McGIVNEY Hemorrhoidal Ligator 7" shaft, McGIVNEY Grasping Forceps 7-1/2", Ligating Rings, and Loading Cone - Large

See page 475 for the complete line of McGivney patterns.





DERMAL/PODIATRY INSTRUMENT

- Double action allows for greater strength and control when cutting heavy nails
- Stainless steel construction
- Non-slip serrated grip handle
- Autoclavable

Part # Description

40-218 Double Action Nail Nipper, 6" (15.2 cm), straight jaws

See page 549 for curved pattern.



4-405

SCALPEL HANDLE

Part # Description


4-405 Self Locking Chisel Blade Handle

See pages 13-14 for a complete line of scalpel handles.









PEDERSON OPEN-SIDED VAGINAL SPECULA


Description Part #

30-63 PEDERSON (Right) Open-Sided Vag. Spec, Narrow Blades (3-1/2" x 7/8") (8.9 cm x 2.2 cm)

See page 490 for a complete line of Pederson specula.








MILTEX INSULATED INSTRUMENTS FOR ELECTROSURGICAL PROCEDURES

- Ideally suited for LEEP/LLETZ
- Highest quality German crafted instruments thermally insulated with a superior, non-conductive blue nylon coating
- Exceeds high voltage testing requirements
- Autoclavable (carefully follow Miltex Care & Handling Instructions, supplied with instrument)
- Limited 1-year warranty against defects in workmanship and material

Part #	<u>Description</u>
31-14	LEEP Dressing Forceps 8" (20.3cm) long, Serrated with Coated Tips
31-16	LEEP Dressing Forceps 10" (25.4cm) long, Serrated with Coated Tips
31-48	LEEP Tissue Forceps 8" (20.3cm) long, 1x2 Teeth with Exposed Tips
31-49	LEEP Tissue Forceps 8" (20.3cm) long, 1x2 Teeth with Coated Tips
31-50	LEEP Tissue Forceps 10" (25.4cm) long, 1x2 Teeth with Exposed Tips
31-51	LEEP Tissue Forceps 10" (25.4cm) long, 1x2 Teeth with Coated Tips








Part #	<u>Description</u>
31-949	LEEP Iris Hook 10" (25.4cm) long, Small delicate Exposed Tip with 4mm hook
31-951	LEEP Emmett Tenaculum Hook, 9" (22.9cm) long, Style 1 (half curved - small) with Exposed Tip
31-952	LEEP Emmett Tenaculum Hook, 9" (22.9cm) long, Style 2 (acute angle) with Exposed Tip
31-953	LEEP Emmett Tenaculum Hook, 9" (22.9cm) long, Style 3 (right angle) with Exposed Tip
31-954	LEEP Emmett Tenaculum Hook, 9" (22.9cm) long, Style 4 (half curved - large) with Exposed Tip
31-965	LEEP Schroeder-Braun Tenaculum Forceps 9-3/4" (24.8cm) long, straight with coated tips






Part # Description

- 31-1351 LEEP Miltex-Kogan Endocervical Speculum 10-3/4" (27.3cm) long, Narrow Tip Coated (3.2mm x 25mm), with Grip-Lock Ratchet & Gauge
- 31-1352 LEEP Miltex-Kogan Endocervical Speculum 10-3/4" (27.3cm) long, Regular Tip Coated (7.2mm x 25mm), with Grip-Lock Ratchet & Gauge
- 31-1353 LEEP Kogan Endocervical Speculum 10-3/4" (27.3cm) long, Narrow Tip Coated (3.2mm x 25mm) (Standard Ratchet Pattern)
- 31-1354 LEEP Kogan Endocervical Speculum 10-3/4" (27.3cm) long, Regular Tip Coated (7.2mm x 25mm) (Standard Ratchet Pattern)
- 31-1355 LEEP Lateral Vaginal Retractor 8-1/4" (21cm) long, with Extra Long Thumb Ratchet (22mm x 57mm blade)
- 31-70 LEEP Redundant Wall Vaginal Retractor 7-3/4" (19.5cm) long, Medium 22mm x 57mm blade
- 31-71 LEEP Redundant Wall Vaginal Retractor 7-3/4" (19.5cm) long, Large 22mm x 75mm blade








Part #	Description
31-101	LEEP Graves Speculum Small (3" x 7/8") (7.6 cm x 2.2 cm) with Smoke Evacuation Tube (2.8cm opening)
31-151	LEEP Graves Speculum Medium (4" x 1-1/4") (10.2 cm x 3.2 cm) with Smoke Evacuation Tube (2.8cm opening)
31-201	LEEP Graves Speculum Large (4-1/2" x 1-3/8") (11.4 cm x 3.5 cm) with Smoke Evacuation Tube (2.8cm opening)
31-301	LEEP Graves Speculum X-Large (6-3/4" x 1-3/8") (17.1 cm x 3.5 cm) with Smoke Evacuation Tube (2.8cm opening)
31-551	LEEP Pederson Speculum Medium (4" x $7/8$ ") (10.2 cm x 2.2 cm) with Smoke Evacuation Tube (2.8cm opening)
31-601	LEEP Pederson Speculum Large (4-3/4" x 1") (12.1 cm x 2.5 cm) with Smoke Evacuation Tube (2.8cm opening)
31-611	LEEP Pederson Speculum X-Large (6-3/4" x 1") (17.1 cm x 2.5 cm) with Smoke

Evacuation Tube (2.8cm opening)






Pan #	Description
31-62	LEEP Pederson (Right) Open-Sided Speculum Medium (4" x 7/8") (10.2 cm x 2.2 cm) with Smoke Evacuation Tube (2.8cm opening)
31-63	LEEP Pederson (Right) Open-Sided Speculum Large (4-1/2" x 1") (11.4 cm x 2.5 cm) with Smoke Evacuation Tube (2.8cm opening)
31-64	LEEP Graves Open-Sided (Right) Speculum Medium (4" x 1-1/4") (10.2 cm x 3.2 cm) with Smoke Evacuation Tube (2.8cm opening)
31-65	LEEP Graves Open-Sided (Right) Speculum Large (4-1/2" x 1-3/8") (11.4 cm x 3.5 cm) with Smoke Evacuation Tube (2.8cm opening)
31-66	LEEP Graves View-More Speculum Medium (4" x 1-1/4") (10.2 cm x 3.2 cm) with Smoke Evacuation Tube (4.0cm yoke opening)
31-67	LEEP Graves View-More Speculum Large (4-1/2" x 1-3/8") (11.4 cm x 3.5 cm) with Smoke Evacuation Tube (4.0cm yoke opening)
31-68	LEEP Graves Maxi-View Speculum Medium (4" x 1-1/4") (10.2 cm x 3.2 cm) with Smoke Evacuation Tube (7.0cm yoke opening)
31-69	LEEP Graves Maxi-View Speculum Large (4-1/2" x 1-3/8") (11.4 cm x 3.5 cm) with Smoke Evacuation Tube (7.0cm yoke opening)








Part #	<u>Description</u>
31-5000	Sterile Disposable Electrode, Ball 3mm Diameter (Box of 5)
31-5002	Sterile Disposable Electrode, Ball 5mm Diameter (Box of 5)
31-5022	Sterile Disposable Electrode, Square (10mm width x 8mm depth) (Box of 5)
31-5024	Sterile Disposable Electrode, Square (10mm width x 10mm depth) (Box of 5)
31-5054	Sterile Disposable Electrode, Loop (10mm width x 8mm depth) (Box of 5)
31-5056	Sterile Disposable Electrode, Loop (10mm width x 10mm depth) (Box of 5)
31-5070	Sterile Disposable Electrode, Loop (15mm width x 5mm depth) (Box of 5)
31-5072	Sterile Disposable Electrode, Loop (15mm width x 8mm depth) (Box of 5)
31-5074	Sterile Disposable Electrode, Loop (15mm width x 10mm depth) (Box of 5)
31-5090	Sterile Disposable Electrode, Loop (20mm width x 8mm depth) (Box of 5)
31-5092	Sterile Disposable Electrode, Loop (20mm width x 10mm depth) (Box of 5)
31-5094	Sterile Disposable Electrode, Loop (20mm width x 12mm depth) (Box of 5)
31-5096	Sterile Disposable Electrode, Loop (20mm width x 15mm depth) (Box of 5)
31-5098	Sterile Disposable Electrode, Loop (20mm width x 20mm depth).(Box of 5)



4-406



PODIATRY AND DERMATOLOGY BLADES Stainless Steel specialty surgical blades – Sterile, individually sealed in blister packs, 20 blades per box.

<u>Cat#</u> 4-601	<u>Size</u> #1	<u>Uses</u> Debridement of calloused nail groves
4-602	#2	Biopsy of small cutaneous lesions
4-603	#3	Removal of plantar warts
4-604	#15G	Debridement of heel callus
4-605	#1V	Removal of foreign bodies
4-606	#2V	Enucleation of callus with painful core
4-607	#3V	Used for a variety of work on the cuticle and corns

MEDI-KLINGEN BLADE HANDLE

4-406 MEDI - KLINGEN BLADE HANDLE with chuck, 5" (12.7cm)











19-115







19-117

30-1384

EZ-ZYME™ ALL PURPOSE ENZYME CLEANER

3/4oz. packets - box of 32. Each packet makes a gallon of solution. Proven in a Major University 3-750

Study to be the most effective enzyme cleaner against a group of competitors.

STUDENT GRADE TUNING FORKS

19-113

C-128 vibrations

19-115 19-117 C-256 vibrations C-512 vibrations

NO SCALPEL VASECTOMY INSTRUMENTS

29-805

No Scalpel Vasectomy Ring Clamp 5 1/2" (14cm), gold finger rings

29-850

No Scalpel Vasectomy Hemostat 5 1/2" (14cm), smooth sharp point, curved, gold finger rings

29-855

No Scalpel Vasectomy Instrument Set - One of each

TOWNSEND ENDOCERVICAL BIOPSY CURETTE













8-5TC

5-269TC 5-263*	KAYE Dissecting Scissors 4 1/2" (11.5cm), curved, one serrated blade, CARB-N-SERT blade KAYE Fine Dissecting Scissors 4 11/32" (11cm), serrated blades, ribbon type
5-400*	MILTEX Utility Scissors 4 1/16" (10.3cm), curved blunt tips, ribbon type
	OBJECTIVO BROWNIET - F - 4 12 (64) // 1 1 1 1


6-127 GRIFFITHS-BROWN Tissue Forceps 4 11/32" (11cm), delicate, 1.5 X 5mm Brown-type teeth
 8-3TC PAR Needle Holder 4 1/2" (11.4cm), delicate serrated jaws 16000 teeth p. sq. inch, CARR N. SERT inch.

CARB-N-SERT jaws

8-5TC PAR Needle Holder 5" (12.7cm), delicate serrated jaws 16000 teeth p. sq. inch, CARB-N-SERT jaws

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished





9-115*

REEH Stitch Scissors 3 7/8" (9.8cm), sharp tips, small hook at end of one blade,

ribbon type

9-116*

MILTEX Stitch Scissors 3 7/8" (9.8cm), curved, sharp pointed tips, ribbon type

9-117*

MILTEX Stitch Scissors 4 7/8" (12.4cm), curved, sharp pointed tips, ribbon type

18-1455

STEVENS Tenotomy Scissors 5 1/2" (14.3cm), curved

18-1456

STEVENS Tenotomy Scissors 6 1/4" (15.9cm), curved


18-SC-1652


SUPERCUT Gradle (Eye Suture) Scissors 3 3/4" (9.5cm), curved

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished


18-SC-1652











CERAM-A-GRIP® NEEDLE HOLDERS

Miltex Ceram-A-Grip® needleholders feature a patented ceramic material that is fusion bonded to the jaw to improve and increase jaw gripping power.

8-4CG	DERF 4	1.3/4"	112	1cm)

8-6CG	WEBSTER 5" (12.7cm)	, delicate, smooth jaws	, standard pattern

8-8CG	HALSEY 5	" (12.7cm),	smooth jaws
-------	----------	-------------	-------------

8-15CG	OLSEN-HEGAR 5 1/2" (14cm), needle holder/scissors combination
8-16CG	OLSEN-HEGAR 6 1/2" (16.5cm), needle holder/scissors combination









8-44CG 8-46CG 8-48CG

8-50CG 8-51CG 8-52CG

CERAM-A-GRIP® NEEDLE HOLDERS

8-40CG

Militex Ceram-A-Grip® needleholders feature a patented ceramic material that is fusion bonded to the jaw to improve and increase jaw gripping power.

·	27.01.07.11.11.21.10.11.2
8-42CG	MAYO-HEGAR 5" (12.7cm)
8-44CG	MAYO-HEGAR 6" (15.2cm)
8-46CG	MAYO-HEGAR 7" (17.8cm)
8-48CG	MAYO-HEGAR 8" (20.3cm)
8-50CG 8-51CG	CRILE-WOOD 6" (15.2cm) CRILE-WOOD 7" (17.8cm)
8-52CG	CRILE-WOOD 8" (20.3cm)

BAUMGARTNER 5 1/2" (14cm)



DISPOSABLE EAR CURETTES

Packaged 50/Box



19-320 Ring Tip 19-321 Loop Tip 19-322 Small Spoon Tip 19-323 Large Spoon Tip

FEATURES

- Four Styles of Tips
- Directional Arrow for Correct Placement
- Handle Scored for Finger Placement and Control
- Atraumatic Design
- Single Use Non Sterile
- Individually Packaged
- Packaged 50 per Dispenser Box

DISPOSABLE DERMAL CURETTES

Packaged 50/Box

33-52 2mm 33-53 3mm 33-54 4mm 33-55 5mm 33-57 7mm

FEATURES

- Five Sizes
- · Size marked on every curette
- Curette blade manufactured from beveled surgical grade stainless steel
- Sturdy handle
- Sterile, disposable
- Cost-effective
- Convenient packaging





DISPOSABLE BIOPSY PUNCHES - ASSORTED

Packaged 50/Box

10 each of the following sizes: 2mm, 3mm, 4mm, 5mm, 6mm

Color Coded by size

33-38 Assortment Box - Biopsy Punches

FEATURES

- Individually packaged in color-coded, transparent, sturdy blister pack to protect cutting edge
- Seamless stainless steel razor sharp cutting edge
- Ribbed handle provides maximum control and comfort

See page 526 for other available sizes





Miltex Edge™ Scissors Sharpener 5-5000










Figure 1



*Not intended for serrated or SuperCut Scissors

Tip-lt[™] Instrument Guards


Tip-lt™ Instrument Guards are used to protect instrument tips and edges from damage during instrument handling and sterilization. Steam and Gas Sterilizable. Radiopaque and Latex Free. 50 per pack



PRODUCT NUMBER	DIMENSION	SIZE CODE	COLOR	
3-2501	1/16" x 3/4" (1.6mm x 19mm)	1	White	
3-2502	5/64" x 3/4" (2mm x 19mm)	2	Blue	
3-2503	7/64" x 3/4" (2.8mm x 19mm)	3	Green	
3-2504	1/8" x 1" (3.2mm x 25.4mm)	1/8" x 1" (3.2mm x 25.4mm) 4		
3-2505	3/16" x 1" (4.8mm x 25.4mm)	5	Yellow	
3-2506	13/32" x 3/4" (10.3mm x 19mm)	6	Maroon	
3-2507	1/16" x 3/8" x 1" (1.6mm x 9.5mm x 25.4mm)	7	Orange	
3-2508	5/64" x 5/8" x 1" (2mm x 16mm x 25.4mm)	8	Brown	
3-2509	1/8" x 1" x 1" (3.2mm x 25.4mm x 25.4mm)	9	Black	NP 21
3-2510	Assorted, equal number of above, plus 5			







- 1-200 TAYLOR Percussion Hammer 7" (17.8 cm), Pediatric Size, solid handle
- 1-202 TAYLOR Percussion Hammer 7 ½" (19.1 cm), Standard Size, solid handle, chrome
- 1-208 TAYLOR Percussion Hammer 9" (22.9 cm), Large Head, loop handle, chrome
- 1-210 DEJERINE Percussion Hammer 9 ½" (24.1 cm), chrome
- 1-211 TROMNER Percussion Hammer 9 %*(24.8 cm), with Large and Small Rubber Tips, hollow tapered handle with sharp point, chrome
- 1-212 COMBO-5 Neuro Diagnostic Set (C-128 Tuning Fork, Buck Hammer, Pin Wheel, Brush and Needle)
- 1-213 COMBO-5 Neuro Diagnostic Set (C-256 Tuning Fork, Buck Hammer, Pin Wheel, Brush and Needle)
- 1-214 COMBO-4 Neuro Diagnostic Set (Buck Hammer, Pin Wheel, Brush and Needle)
- 1-220 BUCK Neurological Hammer 7 ¾" (19.7 cm), with Brush and Needle (screw into handle), chrome
- 1-222 RABINER Neurological Hammer 9" (22.9 cm), with Brush and Needle (screw into handle), chrome





Diagnostic Instruments






19-100

19-100





19-104









MILTEX TUNING FORKS, Aluminum Alloy

19-102


C-64 vibrations

19-102	C-128 vibrations
19-104	C-256 vibrations
19-106	C-512 vibrations
19-108	C-1024 vibrations
19-110	C-2048 vibrations
19-112	C-4096 vibrations

19-120 Set of 5 Tuning Forks, C-128, C-256,C-512, C1024 and C-2048 in


19-120A Fitted zippered case only for tuning fork set

fitted zippered case









1-226	WARTENBERG Neurological Pin Wheel 7" (17.8 cm)	
-------	--	--

1-228 COLLIN Dynamometer, Adult size, 5" x 2 1/4" (12.7 x 5.7 cm)

1-240 FRANKEL Head Band and Mirror Set. Includes Standard 3 ½" (8.9 cm) dia. mirror with 1/2" (1.3 cm) aperture and improved model Frankel Head Band No. 1-244

1-242 MILTEX Head Mirror, Standard size 3 1/2" (8.9 cm) dia. with 1/2" (1.3 cm) aperture, boilable


1-243 MILTEX Head Mirror, Small size 3" (7.6 cm) dia, with %" (1.6 cm) aperture, boilable

1-244 FRANKEL Head Band, Black fibre material, improved pattern with secure thumb screw size adjustment; includes swivel mirror connector

1-3000 EKG Caliper, with protective metal cap, graduated in mm, and ring, satin chrome finish


1-3005 EKG Caliper, with protective metal cap, graduated in mm, and pocket clip. satin chrome finish

1-3015 EKG Caliper, Economy line, with protective black plastic cap, graduated in mm. nickel plated



Diagnostic Instruments





1-3200 MILTEX Fibreglass Tape Measure, ¼" (6.4 mm) wide, graduated to 60" (152.4 cm) on one side, to 150 cm on other side. Plastic case with pushbutton tape return

1-4000 MILTEX Binocular Magnifying Loupe, lightweight, coated optics 1.8X magnification, 30 cm focal distance, 10 cm viewing area, pupillary distance adjustment, mounted on spectacle frame. With carrying case

1-4006 MILTEX Binocular Magnifying Loupe, same as No. 1-4000, but without spectacle frame and without case

1-4010 Accessory Clip-On Lenses, for No. 1-4000 and 1-4006 Loupes, increasing magnification to 2.5X, with 17 cm focal distance, 4.5 cm viewing area


1-4 FORD Stethescope with shallow bell and plastic ear tips, chrome



Anesthesia Instruments


2-104SS






Anesthesia Instruments







2-110SS



2-132 to 2-138A

2-140 to 2-147

2-100 PYNCHON Suction Tube 11"(27.9 cm), tip with large openings, with removable tubing connection, chrome

2-109 YANKAUER Pediatric Suction Tube 8 ½"(21.6 cm), with removable tip, delicate pattern

2-110 ANDREWS-PYNCHON Suction Tube 10"(25.4 cm), delicate pattern, chrome 2-110SS ANDREWS-PYNCHON Suction Tube 10"(25.4 cm), delicate pattern

2-120 DAVIS Mouth Gag 61/4" (15.9 cm), double bite, left, without blades 2-122 DAVIS Mouth Gag 61/4" (15.9 cm), double bite, right, without blades

BLADES FOR DAVIS OR McIVOR MOUTH GAGS

W/o Ether Tube	Size No.	Size	W/Ether Tube
2-130	1	34" x 1" (1.9 x 2.5 cm) Tooth Plate	
2-132	2	1" x 2 ¾" (2.5 x 7 cm) Tongue Blade	2-132A
2-134	3	1" x 3 1/4" (2.5 x 8.3 cm) Tongue Blade	2-134A
2-136	4	1" x 3 ¾" (2.5 x 9.5 cm) Tongue Blade	2-136A
2-138	5	1" x 4" (2.5 x 10.2 cm) Tongue Blade	2-138A

RING TONGUE RETRACTOR BLADES

For use with Davis or McIvor Mouth Gags. Grooved blade holds endotracheal tube firmly in place.

Left	Size No.	Size	Right
2-140	2	7/4" x 2 1/4" (2.2 x 5.7 cm)	2-144
2-141	3	11/6" x 3" (2.9 x 7.6 cm)	2-145
2-142	4	1¼" x 3 ½" (3.2 x 8.9 cm)	2-146
2-143	5	1½" x 3 ¾" (3.8 x 9.5 cm)	2-147





Anesthesia Instruments





63-20









2-162	
2-164	

63-15 63-20




Adult 51/2" (14 cm)



Sterilizing Aids






3-2 3-4	DAVIS Sterilizer Forceps 6" (15.2 cm) DAVIS Sterilizer Forceps 10" (25.4 cm)		
3-10	Three Prong Sterilizer Forceps 8" (20.3 cm)		
3-14	CHEATLE Sterilizer Forceps 11" (27.9 cm)		
3-20 3-22 3-24 3-26	Utility-Sterilizer Forceps 8" (20.3 cm), straight Utility-Sterilizer Forceps 8" (20.3 cm), curved Utility-Sterilizer Forceps 11" (27.9 cm), straight Utility-Sterilizer Forceps 11" (27.9 cm), curved		3-50
3-50	Sterilizer Spiral Rack for surgical needles		
3-60	BUNT Forceps Holder		
65-25	Grasping Forceps 6 ¼" (15.9 cm) Designed to grasp instruments from the ultrasonic and cold sterilization solution trays, containers and baskets.	65-25	3-60



Sterilizing Aids



- 3-70 MILTEX Forceps Holder 8" (20.3 cm), with secure locking device
- 3-71 MILTEX Forceps Holder 10" (25.4 cm), with secure locking device
- 3-72 MILTEX Forceps Holder 12" (30.5 cm), with secure locking device
- 3-74 MILTEX Forceps Holder 16" (40.6 cm), with secure locking device
- 3-400 Flash Autoclaving Pan 9" x 5" x 2" (22.9 x 12.7 x 5.1 cm) with drainage holes, including safety handle
- 3-401 Flash Autoclaving Pan only
- 3-402 Safety Handle only
- 3-404 Flash Autoclaving Tray 10" x 6 1/2" x 2 1/2" (25.4 x 16.5 x 6.4 cm) with perforated bottom, including safety handle
- 3-405 Flash Autoclaving Tray only

INSTRUMENT STERILIZING TRAYS


With handles. Rolled safety edges. Perforated bottoms assure complete steam penetration with no build-up of condensation.

- 3-499 Outside Dimensions 10" x 6 1/2" x 2 1/2" (25.4 x 16.5 x 6.4 cm)
- **3-500** Outside Dimensions 10" x 10 ½" x 3 ½" (25.4 x 26.7 x 8.9 cm)
- **3-501** Outside Dimensions 15" x 10 ½" x 3 ½" (38.1 x 26.7 x 8.9 cm)
- **3-502** Outside Dimensions 20" x 10 ½" x 3 ½" (50.8 x 26.7 x 8.9 cm)
- 65-28 Mesh Holder, 63/4"(17.1 cm), 21/4" (5.4 cm) inside diameter of basket, securely holds small items

such as micro clips, screws, burs etc., during cleaning, rinsing or sterilization







MICRO INSTRUMENT STERILIZING CASES

With handles and removable cover. Double rows of pins welded to bottom hold instruments in place. Strips of silicone rubber cushions on bottom and cover give added protection to instruments. Perforated cover and bottom assure complete heat penetration with no build-up of condensation. All edges rolled for safety. Easily operated lock holds cover securely.

3-550	Outside Dimensions	8"	x 8 ½"	x 1 1/4" (20.3 x 21.6 x 3.2 cm)
3-552	Outside Dimensions	8"	x 12 1/4"	x 1 1/4" (20.3 x 31.1 x 3.2 cm)
3-554	Outside Dimensions	10 16	' x 10"	x 1 1/2 (26.7 x 25.4 x 3.8 cm)

Same features as described above, except double rows of silicone rubber cushions on bottom are adjustable to lengths of instruments. Upper silicone rubber strips hold instruments securely.

3-570	Outside Dimensions 8 1/2" x 8" x 1 1/4" (21.6 x 20.3 x 3.2 cm)
3-572	Outside Dimensions 8 1/2" x 12" x 1 1/4" (21.6 x 30.5 x 3.2 cm)
3-574	Outside Dimensions 10 ½" x 10" x 1 ½" (26.7 x 25.4 x 3.8 cm)
3-576	Outside Dimensions 10 1/2" x 15" x 1 1/4" (26.7 x 38.1 x 3.2 cm)

"CTS" COLOR TRAY SYSTEM


Each tray measures 11 1/8" x 7 1/4" x 1 1/6" (28.2 x 18.4 x 4.1 cm). Lightweight anodized aluminum. Tops available in red, blue and silver. All tray bottoms are silver. Each tray supplied with 2 sets of precut foam which are steam and EO permeable. Not for use with Dry Heat or Cold Sterilizing Solutions.

3-620	Red
3-624	Blue
3-628	Silver
3-640	Pre-Cut "CTS" Foam. 6 sets of :



Instrument Care Products





- MILTEX Spray Lube, 8 fl. oz. (0.24 liter) bottles (case of 12). Pre-mixed non-silicon water soluble 3-700 formula in pump spray bottles. Lubricates all moving parts on surgical instruments.
- MILTEX Surgical Instrument Cleaner, 8 fl. oz. (0.24 liter) bottles (case of 12). Contents of one 8 oz. 3-720 bottle makes 32 gallons of solution. Neutral pH formula for use in ultrasonic or manual cleaning of surgical instruments, glassware, plastic, metal or rubber goods.
- MILTEX Surgical Instrument Cleaner, 1-gallon bottle (3.78 liter) (case of 4). Contents of 1-gallon 3-725 container makes 512 gallons of solution.
- 3-726 MILTEX Surgical Instrument Cleaner, 1/4 fl. oz. (7.4 ml) packets (dispenser box of 80). Contents of 1-box makes 80 gallons of solution.
- 3-740 MILTEX Surgical Instrument Stain Remover, 3 oz. (85 grams) (case of 12), Special powdered formula removes stains and tarnishes, restores luster to surgical instruments, autoclaves, ultrasonic cleaners. stainless steel countertops, trays and stands.
- 3-750 MILTEX EZ-Zyme All Purpose Enzyme Cleaner, 3/4 oz. (22 ml) packets (dispenser box of 32). Contents of 1-box makes 32 gallons of solution. Neutral pH multi-enzyme liquid concentrate ideal for instrument presoak, ultrasonic cleaning, laundry presoak and evacuation system cleaner.
- 3-755 MILTEX EZ-Zyme All Purpose Enzyme Cleaner, 1-gallon bottle (3.78 liter) (case of 4). Contents of 1-gallon container makes 170 gallons of solution.
- Instrument Prep Enzyme Foam, 24 oz. (0.71 liter) spray bottles (case of 12). Ready-to-use foaming 3-760 spray for pre-cleaning of soiled instruments and scopes (rigid & flexible). Ideal for applications where instrument cleaning is not immediately available.
- 3-800 MILTEX Instrument Care System Kit includes one each of the following: Spray Lube, Surgical Instrument Cleaner, Surgical Instrument Stain Remover, Instrument Care Video, Nylon Brush and Stainless Steel Brush.
- 3-810 MILTEX Instrument Care Video, 25-minute Instructional Video designed to help practitioners and staff properly care for instruments.
- 3-1000 MILTEX Instrument Cleaning Brushes 7-1/4" (18.4cm), Nylon Bristles (Set of 3). Designed to clean instruments without scratching or harming delicate surfaces. Autoclayable,
- 3-1001 MILTEX Instrument Cleaning Brushes 7-1/4" (18.4cm), Stainless Steel Bristles (Set of 3), Designed to clean instrument serrations, rasps, files and burs. Autoclavable.

Tubing Clamps






3-2000







3-2000

MILTEX Occlusion Clamp for arteriovenous cannula, 1 ¾" (4.4 cm), cross action

ST. VINCENTS Tube Occluding Clamp 6 1/2" (16.5 cm),

curved, smooth jaws

V97-590

V97-584

PRESBYTERIAN HOSPITAL Tube Occluding Clamp 7" (17.8 cm), straight, smooth jaws, heavy pattern

V97-588


English Pattern Tube Occluding Clamp 4 1/2" (11.4 cm), straight, smooth jaws

U.S. Pattern Tube Occluding Clamp 5 ¾"(14.6 cm), straight, serrated jaws, with guard

V97-585 V97-591

U.S. Pattern Tube Occluding Clamp 7 1/2" (19.1 cm), straight,


serrated jaws, with guard



V97-585 V97-591



Scalpel Handles




4-7	*KNIFE HANDLE No. 3 fitting surgical blades Nos. 10 thru 15, Extra Fine
4-7R	*Same as No. 4-7, but Economy Grade, Rustless, Nickel Plated
4-8	*KNIFE HANDLE No. 4 fitting surgical blades Nos. 20 thru 25, Extra Fine
4-8R	*Same as No. 4-8, but Economy Grade, Rustless, Nickel Plated
4-9	KNIFE HANDLE No. 7, fitting surgical blades Nos. 10 thru 15, Extra Fine
4-10	KNIFE HANDLE No. 3L, for deep surgery, fitting blades Nos. 10 thru 15, Extra Fine
4-11	KNIFE HANDLE No. 4L, for deep surgery, fitting blades Nos. 20 thru 25, Extra Fine
4-12	KNIFE HANDLE No. 3LA, angled tip, for deep surgery, fitting blades Nos. 10 thru 15, Extra Fine
4-16	KNIFE HANDLE, No. 9, fitting blades Nos, 10 thru 15, Extra Fine

^{*}Handles have mm and cm graduations

Scalpel Handles







4-17	SIEGEL Knife Handle 6" (15.2 cm) round, fitting blades Nos. 10 thru 15
4-18	Knife Handle Style No. 5, 5" (12.7 cm), green plastic w/stainless tip, fitting blades Nos. 10 thru 15
4-20	Knife Handle Style No. 6, 5% (13.3 cm), red plastic w/stainless tip, fitting blades Nos. 20 thru 25
	E-Z-OFF Reusable Plastic Knife Handles (Pat. Pend.), autoclavable, permitting quick and safe removal of blade, box of 10
4-30	Style No. 3, fitting blades Nos. 10 thru 15, grey color
4-31	Style No. 4, fitting blades Nos. 20 thru 25, blue color
4-400	Miniature Blade Handle w/chuck, 3¾" (9.5 cm), round knurled (3K type)
4-401	Miniature Blade Handle w/self-locking chuck, 3¾" (9.5 cm), round knurled (3K type)



Sterile Surgical Blades



(Fits handles #'s 3, 3L, 3LA, 5, 7, 9, SIEGEL)




MILTEX SURGICAL BLADES - Made by Kai


Individually sealed in Foil Peel Packs. Packed 100 blades to a Dispenser Box. Sterilized by Gamma Radiation 2.5M Rads.


CARBON STEEL - STERILE

	Blade Size
4-110	10
4-111	11
4-112	12
4-112B	12B
4-115	15
4-115C	15C
4-120	20
4-121	21
4-122	22
4-123	23
4-125	25

4-100








BLADE SAFE Surgical Blade Remover, will remove blades easily and safely from any style surgical handle.

Sterile Surgical Blades





MILTEX SURGICAL BLADES - Made by Kai

Individually sealed in Foil Peel Packs. Packed 100 blades to a Dispenser Box. Sterilized by Gamma Radiation 2.5M Rads.

STAINLESS STEEL - STERILE


	Blade Size
4-310	10
4-311	11
4-312	12
4-312B	12B
4-315	15
4-315C	15C
4-320	20
4-321	21
4-322	22
4-323	23
4-325	25



Sterile Disposable Scalpels

Miltex single-use Scalpels are preassembled: Stainless steel surgical blades mounted on plastic handles with fluted finger-grips. The scalpel blade is safely protected by a new type easy pull-off protector.

Each sterile scalpel is sealed in a peel-apart see-through package. 10 scalpels are packaged in a handy dispenser box for convenient use. Scalpels are supplied sterile - ready-to-use.


















MILTEX DISPOSABLE SCALPELS - Made by Kai

Packed 10 Scalpels to a Dispenser Box. Sterilized by EO.

STAINLESS STEEL - STERILE

	Blade Size
4-410	10
4-411	11
4-412	12
4-415	15
4-420	20
4-421	21
4-422	22
4-423	23






STANDARD PATTERN OPERATING SCISSORS

Straight		Curved	Straight		Curved
	Sharp-sharp points			Sharp-blunt points	
5-2	4½" (11.4 cm)	5-32	5-12	41/2" (11.4 cm)	5-42
5-4	5" (12.7 cm)	5-34	5-14	5" (12.7 cm)	5-44
5-6	51/2" (14 cm)	5-36	5-16	5½" (14 cm)	5-46
5-7	6" (15.2 cm)	5-37	5-17	6" (15.2 cm)	
5-8	6½" (16.5 cm)	5-38	5-18	61/2" (16.5 cm)	5-48
			5-19	7½" (19.1 cm)	5-49

Straight		Curved
	Blunt-blunt points	
5-22	41/2" (11.4 cm)	5-52
5-24	5" (12.7 cm)	5-54
5-26	51/2" (14 cm)	5-56
5-27	6" (15.2 cm)	******
5-28	61/2" (16.5 cm)	5-58







SERRATEX OPERATING SCISSORS

Standard Pattern Operating Scissors. One blade with fine serrations reduces tissue slippage.

Straight		Curved
5-70	5½" (14 cm) Sharp-sharp points	5-74
5-72	5½" (14 cm) Sharp-blunt points	5-76
5-73	5½" (14 cm) Blunt-blunt points	5-77
NOTE: 0		

NOTE: Serrated blades are available on all types of scissors, on special request.

LIGHTWEIGHT OPERATING SCISSORS


Delicate pattern designed for pediatric, vascular and other delicate procedures.

Straight		Curved
	Sharp-sharp points	
5-60	5" (12.7 cm)	5-63
5-80	51/2" (14 cm)	5-86
	Sharp-blunt points	
5-61	5" (12.7 cm)	5-64
5-82	51/2" (14 cm)	5-88
	Blunt-blunt points	
5-62	5" (12.7 cm)	5-65
5-84	5½" (14 cm)	5-90


DEAVER SCISSORS 51/2" (14 cm)

5-100	Straight, sharp-sharp points	5-106	Curved, sharp-sharp points
5-102	Straight, sharp-blunt points	5-108	Curved, sharp-blunt points
5-104	Straight, blunt-blunt points		






5-112 5-114	MIXTER Scissors 6¼" (15.9 cm), straight, blunt points MIXTER Scissors 6¼" (15.9 cm), curved, blunt points
5-115	MAYO Operating Scissors 51/2" (14 cm), straight, heavy pattern, sharp points
5-116	MAYO Operating Scissors 51/2" (14 cm), curved, heavy pattern, sharp points
5-120	MAYO Dissecting Scissors 51/2" (14 cm), straight, standard beveled blades
5-122	MAYO Dissecting Scissors 51/2" (14 cm), curved, standard beveled blades
5-124	MAYO Dissecting Scissors 63/4" (17.1 cm), straight, standard beveled blades
5-126	MAYO Dissecting Scissors 6¾" (17.1 cm), curved, standard beveled blades
5-128	MAYO Dissecting Scissors 9" (22.9 cm), straight, standard beveled blades
5-130	MAYO Dissecting Scissors 9" (22.9 cm), curved, standard beveled blades



5-128 5-130







5-132 5-134	SERRATEX-MAYO Dissecting Scissors 6¾" (17.1 cm), straight, one fine serrated blade SERRATEX-MAYO Dissecting Scissors 6¾" (17.1 cm), curved, one fine serrated blade
5-136 5-138	MAYO Dissecting Scissors 51/2" (14 cm), straight, rounded blades MAYO Dissecting Scissors 51/2" (14 cm), curved, rounded blades
5-140 5-142	MAYO Dissecting Scissors 6%" (17.1 cm), straight, rounded blades





5-144	MAYO-NOBLE Dissecting Scissors 61/2" (16.5 cm), straight, rounded blades
5-146	MAYO-NOBLE Dissecting Scissors 61/2" (16.5 cm), curved, rounded blades
5-148	MAYO-NOBLE Dissecting Scissors 61/2" (16.5 cm), straight, beveled blades
5-150	MAYO-NOBLE Dissecting Scissors 61/2" (16.5 cm), curved, beveled blades
5-152	SISTRUNK Scissors 51/2" (14 cm), straight
5-154	SISTRUNK Scissors 51/2" (14 cm), slightly curved




30-2575

30-2580

Operating Scissors






5-160 5-162	DOYEN Abdominal Scissors 7" (17.8 cm), straight DOYEN Abdominal Scissors 7" (17.8 cm), curved
5-166	FERGUSON Abdominal Scissors 7" (17.8 cm), angled on flat
5-172 5-174 5-176 5-178	REYNOLDS Dissecting Scissors 6" (15.2 cm), curved, tenotomy type dissecting tips REYNOLDS Dissecting Scissors 6" (15.2 cm), curved, tenotomy type dissecting tips, one serrated blade REYNOLDS Dissecting Scissors 7" (17.8 cm), curved, tenotomy type dissecting tips REYNOLDS Dissecting Scissors 7" (17.8 cm), curved, tenotomy type dissecting tips, one serrated blade
30-2190 30-2195	BRAUN Episiotomy Scissors 5½" (14 cm), angled to side, guarded lower blade BRAUN Episiotomy Scissors 8½" (21.6 cm), angled to side, guarded lower blade

Umbilical Scissors 41/8" (10.5 cm), American Pattern

BUSCH Scissors 61/2" (16.5 cm), curved on side

30-2580






METZENBAUM SCISSORS – STANDARD PATTERN		METZENBAU	M SCISSORS - DELICAT	E PATTERN		
Straight		Length	Curved	Straight	Length	Curved
5-179	5½"	(14 cm) (Lahey)	5-180	5-283	5" (12.7 cm)	5-284
5-181	7"	(17.8 cm)	5-182	_	6" (15.2 cm)	5-286
	8"	(20.3 cm)	5-184		7" (17.8 cm)	5-288
5-186	9"	(22.9 cm) (Nelson)	5-188	5-183	8" (20.3 cm)	5-185
5-190	11"	(27.9 cm) (Nelson)	5-192	5-187	9" (22.9 cm)	5-189
5-194	141/2	(36.8 cm) (Nelson)	5-196		11" (27.9 cm)	5-193

5-200 HARRINGTON Scissors 11½" (29.2 cm), straight 5-202 HARRINGTON Scissors 11½" (29.2 cm), curved








5-220	SIMS Scissors 8" (20.3 cm), straight, sharp-sharp points
5-222	SIMS Scissors 8" (20.3 cm), straight, sharp-blunt points
5-224	SIMS Scissors 8" (20.3 cm), straight, blunt-blunt points
5-226	SIMS Scissors 8" (20.3 cm), curved, sharp-sharp points
5-228	SIMS Scissors 8" (20.3 cm), curved, sharp-blunt points
5-230	SIMS Scissors 8" (20.3 cm), curved, blunt-blunt points
5-243	JORGENSON Scissors 9" (22.9 cm), strongly curved blades
5-244	THOREK Scissors 71/4" (18.4 cm), full curve
5-245	THOREK Scissors 10" (25.4 cm), full curve
27-1000	MARTIN Cartilage Scissors 8" (20.3 cm), serrated non-slip cutting edges
30-2605	DUBOIS Scissors 101/2" (26.7 cm), curved, extra heavy pattern, blunt points
	and the second s







5-246	POTTS-SMITH Scissors 7½" (19.1 cm), angled on side 25°, 20 mm blades	
5-247	POTTS-SMITH Scissors 7½" (19.1 cm), angled on side 45°, 20 mm blades	
5-248	POTTS-SMITH Scissors 7½" (19.1 cm), angled on side 60°, 20 mm blades	\mathcal{O}
5-249A	POTTS-SMITH Scissors 7" (17.8 cm), delicate, angled on side 25°,13 mm blades	25-1240B 25-1240C
5-249B	POTTS-SMITH Scissors 7" (17.8 cm), delicate, angled on side 45°,13 mm blades	
5-249C	POTTS-SMITH Scissors 7" (17.8 cm), delicate, angled on side 60°,13 mm blades	
25-1240B	DeBAKEY-POTTS Scissors 6¾" (17.1 cm), dull tips, angled on side 45°	
25-1240C	DeBAKEY-POTTS Scissors 6%" (17.1 cm), dull tips, angled on side 60°	









5D-329







25-1280	SATINSKY Scissors 10" (25.4 cm), S-curved
28-240	MILLER Rectal Scissors 6%" (17.1 cm), bayonet shape, delicate serrated blades with blunt points
5-250 5 D- 251	KELLY Fistula Scissors 6¼" (15.9 cm), sharp points, straight KELLY Fistula Scissors 6¼" (15.9 cm), sharp points, straight, with one serrated blade
5-252 5D-253	KELLY Fistula Scissors 61/4" (15.9 cm), sharp points, curved KELLY Fistula Scissors 61/4" (15.9 cm), sharp points, curved,
5-254	with one serrated blade KELLY Fistula Scissors 61/4" (15.9 cm), sharp points, angled on side
5D-255	KELLY Fistula Scissors 6¼" (15.9 cm), sharp points, angular, with one serrated blade
5D-328	LOCKLIN Scissors 61/4" (15.9 cm), angled on side 25°, with one serrated blade


LOCKLIN Scissors 6% (15.9 cm), curved shanks, angled on side 25° , with one serrated blade














5-256	BOETTCHER Dissecting Scissors 71/4" (18.4 cm) double edged, curved
22-962	LILLY Dissecting Scissors 7½" (19.1 cm), strong curve
5-260	SWEET Scissors 8½" (21.6 cm), angled on flat
5-262	STRULLY Dissecting Scissors 8" (20.3 cm), slightly curved blades with probe tips
5-264	DEAN Dissecting Scissors 6%" (17.1 cm), blades angled on flat, one serrated blade
5-265	KAHN Dissecting Scissors 51/2" (14 cm), straight
5-266	KAHN Dissecting Scissors 51/2" (14 cm), curved
5-267	KAHN Dissecting Scissors 7" (17.8 cm), straight
5-268	KAHN Dissecting Scissors 7" (17.8 cm), curved



















PLASTIC SURGERY SCISSORS 4¾" (12.1cm)

		,
Straight		Curved
5-270	Sharp-sharp points	5-276
5-271	Sharp-sharp points, serrated blade	5-277
5-272	Sharp-blunt points	5-278
5-273	Sharp-blunt points, serrated blade	5-279
5-274	Blunt-blunt points	5-280
5-275	Blunt-blunt points, serrated blade	5-281
	Angled on side	

5-282 Sharp-sharp points Sharp-sharp points, serrated blade 5-282A








flat tip, curved	secting Scissors 7" (17.8 cm),	6	
5-302 Iris Scissors 4 5D-303 Iris Scissors 4 5-304 Iris Scissors 4 5-306 Iris Scissors 4	" (10.2 cm), straight " (10.2 cm), curved ½" (11.4 cm), curved sideways ½" (11.4 cm), straight ½" (11.4 cm), curved ½" (11.4 cm), angled on side	5-306	5-308







18-1392	BONN Miniature Iris Scissors 3½" (8.9 cm), with 15 mm blades, straight
18-1394	BONN Miniature Iris Scissors 3½" (8.9 cm), with 15 mm blades, curved

18-1396 Iris Scissors 31/2" (8.9 cm), with 20 mm blades, straight, delicate

18-1398 Iris Scissors 31/2" (8.9 cm), with 20 mm blades, curved, delicate

18-1404 MILTEX Iris Scissors 4" (10.2 cm), straight, heavy pattern

18-1405 MILTEX Iris Scissors 4" (10.2 cm), curved, heavy pattern


18-1415* Iris Scissors 4" (10.2 cm), straight, standard pattern, ribbon type

18-1416* Iris Scissors 4" (10.2 cm), curved, standard pattern, ribbon type

Iris Scissors 3¾" (9.5 cm), straight, miniature pattern 18-1417* with 14 mm blades, ribbon type


18-1418* Iris Scissors 33/4" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon type

18-1419* Utility Scissors 3¾" (9.5 cm), curved blunt blades, 15 mm long, ribbon type



^{*}Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished





18-1426 18-1428







18-1516

18-1420	KNAPP Iris Scissors 4" (10.2 cm), straight, sharp-blunt points
18-1422	KNAPP Iris Scissors 4" (10.2 cm), straight, blunt-blunt points
18-1426	KNAPP Iris Scissors 4" (10.2 cm), curved, sharp-blunt points
18-1428	KNAPP Iris Scissors 4" (10.2 cm), curved, blunt-blunt points
18-1440	WILMER Conjunctival and Utility Scissors 4" (10.2 cm), angled on flat
18-1510	NOYES Iris Scissors 41/2" (11.4 cm), straight, sharp-sharp points
18-1512	NOYES Iris Scissors 41/2" (11.4 cm), straight, sharp-blunt points
18-1514	NOYES Iris Scissors 4½" (11.4 cm), straight, blunt-blunt points
18-1516	NOYES Iris Scissors 41/2" (11.4 cm), curved, sharp-sharp points
18-1520 18-1521	BARRAQUER Iris Scissors 2% " (5.7cm), angled 7 mm blades, sharp points BARRAQUER Iris Scissors 2% " (5.7cm), angled 7 mm blades, blunt points



18-1520 18-1521













18-1524 18-1526

18-1550 18-1552

18-1554


18-1619

18-1524 18-1526	DeWECKER Iris Scissors 4¼" (10.8 cm), angled 7 mm blades, sharp-sharp points DeWECKER Iris Scissors 4¼" (10.8 cm), angled 7 mm blades, sharp-blunt points
18-1550	McCLURE Iris Scissors 3%" (9.5 cm), angled over edge, 4 mm blades with very delicate sharp points
18-1552	McCLURE Iris Scissors 3¾" (9.5 cm), angled on flat, 4 mm blades with very delicate sharp points
18-1554	McCLURE Iris Scissors 4" (10.2 cm), angled on flat, regular 10 mm blades with sharp points
18-1618 18-1619	Micro Iris Scissors 4" (10.2 cm), straight, sharp points Micro Iris Scissors 4" (10.2 cm), curved, sharp points
18-1625 18-1626	McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), straight, sharp points McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), curved, sharp points



18-1625 18-1626





18-1430	Eye Scissors with probe points 4" (10.2 cm), straight
18-1432	Eye Scissors with probe points 4" (10.2 cm), curved
5-310	Eye Scissors with probe points 41/2" (11.4 cm), angled on side
5-312	Strabismus Scissors 4" (10.2 cm), straight
5-314	Strabismus Scissors 4" (10.2 cm), curved
18-1453*	Strabismus Scissors 41/4" (10.8 cm), curved, ribbon type
18-1454*	Strabismus Scissors 41/4" (10.8 cm), curved, ribbon type



*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished
























21-605



21-558 21D-559	QUIMBY Scissors 5" (12.7 cm), blades 15 mm, strongly curved and pointed QUIMBY Scissors 5" (12.7 cm), blades 15 mm, strongly curved and pointed, one serrated blade
21-560	TAYLOR Dural Scissors 51/2" (14 cm), with probe tip
21-600 21-602	FOMON Saber-Back 5" (12.7 cm), slightly curved blades with semi-sharp outer edges used for dissection FOMON Lower Lateral Scissors 5" (12.7 cm), strong curve, dissecting blades
21-604	FOMON Dorsal Scissors 5½" (14 cm), light pattern, angular shanks
21-605	MILTEX Saber-Back Rhytidectomy Scissors 6½" (16.5 cm), curved with semi-sharp outer edges used for dissection







21-606 CASTANARES Facelift Scissors 61/4" (15.9 cm), straight, short strong blades 21-607 GORNEY Facelift Scissors 7¾" (19.7 cm), straight, one serrated blade GORNEY Facelift Scissors 73/4" (19.7 cm), slightly curved, one serrated blade 21-608 COTTLE Dorsal Scissors 61/2" (16.5 cm), heavy pattern with rounded blades, 21-610 angular shanks 21-612 COTTLE Bulldog Nasal Scissors 41/2" (11.4 cm), slightly curved blunt blades











21-621 GORNEY Septal Scissors 7¾" (19.7 cm), double action, angled handle, serrated blades with 1¾" (4.4 cm) cut

21-622 CAPLAN Nasal Bone Scissors 8" (20.3 cm), double action, serrated blades

21-630 JOSEPH Nasal Scissors 5¾" (14.6 cm), straight, sharp points
21-632 JOSEPH Nasal Scissors 5¾" (14.6 cm), curved, sharp points
33-240 BIRO Dermal Naevus Scissors 4" (10.2 cm), angled on flat, with sharp outside edges


33-242 MILTEX Undermining Scissors 4½" (11.4 cm), curved blades,

with flat pointed tip











9-96	HEATH Suture Scissors 61/4" (15.9 cm), serrated blade
9-98	NEW'S Suture Scissors 5½" (14 cm), angled on flat
9-100	SPENCER Stitch Scissors 31/2" (8.9 cm), delicate
9-101	SHORTBENT Stiftch Scissors $3\%"$ (8.9 cm), curved, delicate
9-102	LITTAUER JR. Stitch Scissors 41/2" (11.4 cm), light pattern
9-103	NORTHBENT Stitch Scissors 4¾" (12.1 cm), curved, light pattern
9-104	LITTAUER Stitch Scissors 5½" (14 cm), standard pattern







9-106	LONG Oral Surgery Stron Scissors 6" (15.2 cm)	
9-108	MILTEX Stitch Scissors 41/2" (11.4 cm), angled, delicate	
9-110	O'BRIEN Stitch Scissors 3¾" (9.5 cm), angled, sharp points	
	WIRE CUTTING SCISSORS 43/4" (12.1 cm)	

9-120 Straight, one serrated blade

9D-121 Straight

9-122 Curved, one serrated blade

9D-123 Curved

9-124 Angled to side, one serrated blade
















WIRE CUTTING SCISSORS 41/4" (10.8 cm)

9D-132	Straight, one serrated blade
--------	------------------------------

9D-133 Straight

9D-134 Curved, one serrated blade

9D-135 Curved

18-1486 WESTCOTT Stitch Scissors 41/2" (11.4 cm), curved, sharp points

Eye Suture (GRADLE) Scissors 3¾" (9.5 cm), slightly curved, sharp points 18-1652 Eye Suture Scissors 4" (10.2 cm), slightly curved, sharp points, ribbon type 18-1653*

21-536 LITTLER Suture Carrying Scissors 45/8" (11.8 cm), with suture hole in blades, curved

21-536

^{*} Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished







17-500

PENCIL-TYPE MICRO SCISSORS

7¼" (18.4 cm) overall length; extremely delicate working ends

17-115	Scissors, 4 mm long blades, straight shaft
17-120	Scissors, 8 mm long blades, straight shaft
17-215	Scissors, 4 mm long blades, angled shaft

17-500 MILTEX Micro Dissecting Scissors 5%" (14.6 cm), double action, sharp points, spring handle

17-115 17-120 17-215







17-2060 17-2040



17-2150 17-2100



17-2110



18-1559

MILTEX MICRO SURGERY SCISSORS, Sharp points

17-2060	6" (15.2 cm) straight, 8 mm blades, flat handles	
17-2040	6" (15.2 cm) curved, 8 mm blades, flat handles	
17-2150	6" (15.2 cm) straight, 6 mm blades, round handles	
17-2100	6" (15.2 cm) curved, 6 mm blades, round handles	
17-2160	71/8" (18.1 cm) straight, 6 mm blades, round handles	
17-2110	71/6" (18.1 cm) curved, 6 mm blades, round handles	
18-1558	CASTROVIEJO Corneal Scissors 4" (10.2 cm), right, sub-miniature for micro surgery, 6 mm blades	
18-1559	CASTROVIEJO Corneal Scissors 4" (10.2 cm),	











18-1618	Micro Iris Scissors 4" (10.2 cm), straight, sharp points
18-1619	Micro Iris Scissors 4" (10.2 cm), curved, sharp points
18-1625 18-1626	McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), straight, sharp points McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), curved, sharp points
18-1630	Micro-VANNAS Scissors 31/4" (8.3 cm), straight, ultrafine blades
18-1631	Micro-VANNAS Scissors 31/4" (8.3 cm), curved, ultrafine blades
18-1633	Micro Scissors 3¼" (8.3 cm), straight, superfine blades
18-1634	Micro Scissors 3¼" (8.3 cm), curved, superfine blades





18-1460 to 18-1476

18-1478*



18-1460 to 18-1466



18-1470 to 18-1476

STEVENS TENOTOMY SCISSORS

18-1460	Straight sharp points 41/8" (10.5 cm), short blades
18-1462	Straight blunt points 41/6" (10.5 cm), short blades
18-1464	Curved sharp points 41/6" (10.5 cm), short blades
18-1466	Curved blunt points 41/8" (10.5 cm), short blades
18-1470	Straight sharp points 41/2" (11.4 cm), long blades
18-1472	Straight blunt points 41/2" (11.4 cm), long blades
18-1474	Curved sharp points 41/2" (11.4 cm), long blades
18-1476	Curved blunt points 41/2" (11.4 cm), long blades
18-1477*	STEVENS Tenotomy Scissors 3¾" (9.5 cm), straight, blunt points, ribbon type

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished

STEVENS Tenotomy Scissors 3¾" (9.5 cm), curved, blunt points, ribbon type













18-1480 WESTCOTT Utility Scissors 4½" (11.4 cm), curved, blunt tips
 18-1482 SCHEIE-WESTCOTT Corneal Section Scissors 5" (12.7 cm), right, curved, blunt tips
 18-1484 SCHEIE-WESTCOTT Corneal Section Scissors 5"(12.7 cm), left, curved, blunt tips
 18-1487 WESTCOTT Tenotomy Scissors 5½" (13.3 cm), right, wide handles
 18-1488 WESTCOTT Tenotomy Scissors 5½" (13.3 cm), left, wide handles

Enucleation Scissors 5"(12.7 cm), light curve

Enucleation Scissors 5"(12.7 cm), full curve

Enucleation Scissors 5"(12.7 cm), medium curve

18-1494 to 18-1498

18-1494

18-1496

18-1498













18-1570

18-1572





18-1540

18-1540

NOYES Iridectomy Scissors 51/2" (14 cm), straight, 7 mm blades, sharp-sharp points

18-1560 CASTROVIEJO Corneal Section Scissors 41/2" (11.4 cm), right,

inner blade 1 mm longer, blunt tips

18-1562 CASTROVIEJO Corneal Section Scissors 41/2" (11.4 cm), left, inner blade 1 mm longer, blunt tips


18-1570 CASTROVIEJO Keratoplasty Scissors 3¾" (9.5 cm), angled blades 11 mm, sharp points

CASTROVIEJO Corneal Scissors 3¾" (9.5 cm), angled blades 18-1572 11 mm, blunt points

CASTROVIEJO Corneal Scissors 33/4" (9.5 cm), curved, sharp points 18-1576 CASTROVIEJO Corneal Scissors 3¾" (9.5 cm), curved, blunt points

18-1578

McGUIRE Corneal Scissors 41/8" (10.5 cm), right 18-1594 18-1596 McGUIRE Corneal Scissors 41/8" (10.5 cm), left











18-1600 to 18-1604

very delicate

18-1620 18-1622 18-1623 18-1624


18-1600 18-1602 18-1604	AEBLI Corneal Scissors 4" (10.2 cm), straight AEBLI Corneal Scissors 4" (10.2 cm), angled on flat, right AEBLI Corneal Scissors 4" (10.2 cm), angled on flat, left	
18-1620	VANNAS Capsulotomy Scissors 31/4" (8.3cm), straight, very delicate sharp points	
18-1622	VANNAS Capsulotomy Scissors 31/4" (8.3 cm), curved, very delicate sharp points	
18-1623	VANNAS Capsulotomy Scissors 31/4" (8.3 cm), angled to side, very delicate sharp points	19-2060
18-1624	VANNAS Capsulotomy Scissors 31/4" (8.3 cm), angled on flat, very delicate sharp points	

WULLSTEIN Ear Scissors 3" (7.6 cm) shaft, straight,

19-2060







BELLUCCI Micro Ear Scissors 2¾" (7 cm) shaft, 4 mm blades, ebony finish 19-2150-B 20-160 KNIGHT Nasal Scissors 6¾" (17.1 cm), angled on side, standard pattern 20-164 HEYMAN Nasal Scissors 6¾" (17.1 cm), angled on side, rounded blades 20-168 SEILER Turbinate Scissors 61/4" (15.9 cm), angled on side

Operating Scissors













34-80	Dissecting Scissors 4½" (11.4 cm), straight, two sharp points economy grade
34-92	RICHTER Dissecting Scissors 5" (12.7 cm), angular blades, one sharp and one blunt point
34-98	Enterotomy Scissors 8" (20.3 cm), with hook blade
40-260	Ingrown Naii Splitting Scissors 6" (15.2 cm), straight one serrated blade.
40-425	Nail Scissors 31/2" (8.9 cm), curved blades, chrome
40-435	Cuticle Scissors 3½" (8.9 cm), curved blades, standard pattern, chrome
40-445	Cuticle Scissors 3½" (8.9 cm), curved blades, extra delicate, chrome









Bandage and Utility Scissors





5-510	LISTER Bandage Scissors 31/2" (8.9 cm), extra fine
5-512	LISTER Bandage Scissors 41/2" (11.4 cm), extra fine
5-514	LISTER Bandage Scissors 51/2" (14 cm), extra fine
5-516	LISTER Bandage Scissors 71/4" (18.4 cm), extra fine
5-550	LISTER Bandage Scissors 8" (20.3 cm), with one large finger ring
5-560 5-561	KNOWLES Bandage Scissors 5½" (14 cm), straight KNOWLES Bandage Scissors 5½" (14 cm), shanks angled on side, serrated blade
5-547	HERCULES Heavy Duty Bandage and Plaster Shears $7\%"$ (19.1 cm), serrated blade
5-564	U.S.A. Gauze Scissors 8" (20.3 cm), chrome





Bandage and Utility Scissors













5-568 ESMARCH Heavy Duty Bandage and Cast Shears 8" (20.3 cm)

5-600 MILTEX First Aid and Utility Shears 8½" (21.6 cm), with wire cutter and serrated blade

Nurse's Scissors 5½" (14 cm). Fluoride coated blades help prevent adhesive from sticking to the blade surfaces. Scissors have plastic finger rings and safety guard. Supplied in a plastic carrying pouch. Autoclavable

5-700	Pink
5-702	Blue
5-800	Bandage and Utility Scissors 8½" (21.6 cm), with fluoride coated blades to help prevent adhesive or fiberglass casting material from sticking to the blade surfaces. Neon green plastic finger rings, serrated blade, Autoclauable


5-998 Bandage and Utility Scissors 6½" (16.5 cm), with needle destroyer, black plastic finger rings, serrated blade. Autoclavable

5-1000 Bandage and Utility Scissors 7½" (19.1 cm), with needle destroyer, black plastic finger rings, serrated blade. Autoclavable



Cast and Utility Shears








27-3059 STILLE Plaster Shears 8" (20.3 cm), Pediatric Size, serrated blades STILLE Plaster Shears 11" (27.9 cm), Medium Size, serrated blades 27-3060 27-3062 STILLE Plaster Shears 15" (38.2 cm), Standard Size, serrated blades

27-3070 BRUNS Plaster Shears 91/4" (23.5 cm), serrated blade

Carb-N-Sert Scissors



OUR TOP QUALITY LINE OF SURGICAL SCISSORS



Manufactured to traditionally high Miltex standards, these tungsten carbide blades are bonded to the finest quality stainless steel surgical scissors, then precision ground and honed on diamond wheels.

- Precision screws, fully tempered, will keep these scissors in perfect alignment.
- Scissors shanks are tempered to a hardness that permits proper flexibility.
- Finger rings are well rounded and polished for better "feel" and prevention of finger fatigue, and each instrument is satin finished to reduce glare.

Gold plated finger rings identify these tungsten carbide instruments, and the Millex catalog number has been etched on each scissors to facilitate re-ordering.


Carb-N-Sert Scissors will rarely require sharpening. While the initial cost of Carb-N-Sert Scissors is higher than standard stainless scissors, Carb-N-Sert Scissors represent a considerable savings to the user in the long run because of superior performance and longer service life.



Carb-N-Sert Scissors



CARB-N-SERT SCISSORS with tungsten carbide insert blades will stay sharp for years




5-6TC	Operating Scissors 51/2" (14 cm), straight, sharp-sharp points
5-16TC	Operating Scissors 51/2" (14 cm), straight, sharp-blunt points
5-26TC	Operating Scissors 51/2" (14 cm), straight, blunt-blunt points
5-36TC	Operating Scissors 51/2" (14 cm), curved, sharp-sharp points
5-46TC	Operating Scissors 51/2" (14 cm), curved, sharp-blunt points
5-56TC	Operating Scissors 51/2" (14 cm), curved, blunt-blunt points





CARB-N-SERT SCISSORS with tungsten carbide insert blades will stay sharp for years



5-120TC	MAYO Dissecting Scissors 51/2" (14 cm), straight, standard beveled blades
5-122TC	MAYO Dissecting Scissors 5½" (14 cm), curved, standard beveled blades
5-124TC	MAYO Dissecting Scissors 6¾" (17.1 cm), straight, standard beveled blades
5-126TC	MAYO Dissecting Scissors 6¾" (17.1 cm), curved, standard beveled blades
5-128TC	MAYO Dissecting Scissors 9" (22.9 cm), straight, standard beveled blades
5-130TC	MAYO Dissecting Scissors 9" (22.9 cm), curved, standard beveled blades
5-136TC	MAYO Dissecting Scissors 5½" (14 cm), straight, rounded blades
5-138TC	MAYO Dissecting Scissors 5½" (14 cm), curved, rounded blades
5-140TC	MAYO Dissecting Scissors 6¾" (17.1 cm), straight, rounded blades
5-142TC	MAYO Dissecting Scissors 6¾" (17.1 cm), curved, rounded blades



Carb-N-Sert Scissors



CARB-N-SERT SCISSORS with tungsten carbide insert blades will stay sharp for years



5-182TC



5-188TC 5-193A-TC 5-189A-TC



5-180A-TC



5-180B-TC




5-182A-TC

METZENBAUM SCISSORS STANDARD PATTERN — BLUNT POINTS

Straight	Length	Curved
5-179TC 5-181TC	5½" (14 cm) 7" (17.8 cm) 8" (20.3 cm)	5-180TC 5-182TC 5-184TC
5-186TC —	9" (22.9 cm) 10" (25.4 cm)	5-188TC 5-189A-TC
_	11" (27.9 cm) 12" (30.5 cm)	5-192TC 5-193A-TC


METZENBAUM SUISSONS - DELIGATE PATTERN				
Straight	Length	Curved		
5-179A-TC 5-179B-TC 5-181A-TC 5-181B-TC 5-187TC	5½" (14 cm) blunt points 5½" (14 cm) sharp points 7" (17.8 cm) blunt points 7" (17.8 cm) sharp points 9" (22.9 cm) blunt points	5-180A-TC 5-180B-TC 5-182A-TC 5-182B-TC 5-189TC		







CARB-N-SERT SCISSORS with tungsten carbide insert blades will stay sharp for years



5-304TC Iris Scissors 4½" (11.4 cm), straight **5-306TC** Iris Scissors 4½" (11.4 cm), curved


5-312TC Strabismus Scissors 4½" (11.4 cm), straight 5-314TC Strabismus Scissors 4½" (11.4 cm), curved



Carb-N-Sert Scissors



CARB-N-SERT SCISSORS with tungsten carbide insert blades will stay sharp for years



GOLDMAN-FOX Wound Debridement Scissors 5" (12.7 cm), curved fine tips, 5-320TC

one serrated blade

Wire Cutting Scissors 43/4" (12.1 cm), angled to side, one serrated blade 9-124TC

MILTEX Blepharoplasty Scissors 4½" (11.4 cm), curved flat tips, one serrated blade 21-537TC

21-604TC FOMON Dorsal Scissors, 5½" (14 cm), angled shanks

C-N-S Pin and Wire Cutter 7" (17.8 cm), double action, cuts up to 27-159TC

2 mm (5/4") soft wire or 1.5 mm (1/16") hard wire











5-SC-126

Manufactured with the highest quality stainless steel and German craftsmanship, Miltex SuperCut Scissors are ideal for various surgical and dental procedures or applications. SuperCut Scissors have a specially designed razor-sharp upper blade edge that cuts effortlessly through tissue. The lower blade has micro-fine serrations to prevent tissue slippage. For easy identification, Miltex SuperCut Scissors have two black finger rings.


Miltex also offers SuperCut/Carb-N-Serte Scissors, which have precision SuperCut performance with one lower serrated blade and one upper razor-sharp blade, as well as the durability of Tungsten Carbide cutting edges. These scissors are easily identified with one short and one long gold-plated finger ring.



5-SC-16 5-SC-46	Operating Scissors 5½" (14 cm), straight, sharp/blunt Operating Scissors 5½" (14 cm), curved, sharp/blunt
0 00 40	Operating deleters 572 (11 am), early of charpes and
5-SC-120	MAYO Scissors 5½" (14 cm), straight
5-SC-122	MAYO Scissors 51/2" (14 cm), curved
5-SC-124	MAYO Scissors 6¾" (17.1 cm), straight
5-SC-126	MAYO Scissors 6¾" (17.1 cm), curved







5-SC-128 5-SC-130	MAYO Scissors 9" (22.9 cm), straight MAYO Scissors 9" (22.9 cm), curved
5-SC-180	METZENBAUM Scissors 51/2" (14 cm), curved
5-SC-182	METZENBAUM Scissors 7" (17.8 cm), curved
5-SC-18B	METZENBAUM Scissors 9" (22.9 cm), curved






5-SC-283 5-SC-284	METZENBAUM Scissors 4½" (11.4 cm), straight, delicate METZENBAUM Scissors 4½" (11.4 cm), curved, delicate
5-SC-304 5-SC-306	Iris Scissors 4½" (11.4 cm), straight Iris Scissors 4½" (11.4 cm), curved
5-SC-320	GOLDMAN-FOX Scissors 5" (12.7 cm), curved







18-SC-1472 18-SC-1476 STEVENS Tenotomy Scissors 41/2" (11.4 cm), straight STEVENS Tenotomy Scissors 41/2" (11.4 cm), curved

5-SC-182TC


METZENBAUM Scissors 7" (17.8 cm), curved, SuperCut Carb-N-Sert Blades

5-SC-304TC 5-SC-306TC

Iris Scissors 41/2" (11.4 cm), straight, SuperCut Carb-N-Sert Blades Iris Scissors 41/2" (11.4 cm), curved, SuperCut Carb-N-Sert Blades



18-781



DRESSING FORCEPS, SERRATED TIPS


Standard Pattern Serrated Handles	Standard Pattern Lightweight with Fenestrated Handles	Length	Delicate Pattern Fluted Handles	Delicate Pattern Lightweight with Fenestrated Handles	
6-2	-	4" (10.2 cm)		-	
6-4		41/2" (11.4 cm)	6-24		/ / [2]
6-6		5" (12.7 cm)	6-26		
6-8	6-8XL	51/2" (14 cm)	6-28	6-28XL	1/2
6-10	_	6" (15.2 cm)	6-30		WI /
6-12	_	7" (17.8 cm)	_		M
6-14		8" (20.3 cm)	-		8
6-16		10" (25.4 cm)		_	
6-18	_	12" (30.5 cm)	_	_	\(\frac{\psi}{y}\)
					V

EYE DRESSING FORCEPS 4" (10.2cm), SERRATED TIPS

Standard Pattern 0.8 mm wide tips		Extra Delicate Pattern 0.5 mm wide tips	
18-780	Straight	18-779	Straight
18-782	Half curved	18-781	Half curved
10 704	Eull augusd	10 702	Cull augrad








TISSUE FORCEPS, 1 X 2 TEETH

Standard Pattern Serrated Handles	Standard Pattern Lightweight with Fenestrated Handles	Length	Delicate Pattern Fluted Handles
6-40	_	41/2" (11.4 cm)	6-90
6-42		5" (12.7 cm)	6-92
6-44	6-44XL	51/2" (14 cm)	6-94
6-46	_	6" (15.2 cm)	6-96
6-47	_	7" (17.8 cm)	_
6-48		8" (20.3 cm)	******
6-50	******	10" (25.4 cm)	
6-52	_	12" (30.5 cm)	
6-60	Tissue Forceps	4½" (11.4 cm)	2 x 3 teeth
6-62	Tissue Forceps	5" (12.7 cm)	2 x 3 teeth
6-64	Tissue Forceps	51/2" (14 cm)	2 x 3 teeth
6-66	Tissue Forceps	6" (15.2 cm)	2 x 3 teeth
6-72	Tissue Forceps	5" (12.7 cm)	3 x 4 teeth
6-74	Tissue Forceps	5½" (14 cm)	3 x 4 teeth
6-76	Tissue Forceps	6" (15.2 cm)	3 x 4 teeth
6-82	Tissue Forceps	5" (12.7 cm)	4 x 5 teeth
6-84	Tissue Forceps	51/2" (14 cm)	4 x 5 teeth
6-86	Tissue Forceps	6" (15.2 cm)	4 x 5 teeth






IRIS TISSUE FORCEPS 4" (10.2 cm), 1 x 2 TEETH

	ard Pattern n wide tips		licate Pattern n wide tips
18-786	Straight	18-785	Straight
18-787	Straight (Non Magnetic)	18-789	Half curved
18-788	Half curved	18-791	Full curved
18-790	Full curved		








6-104 SEMKEN Dressing Forceps 5" (12.7 cm), serrated	
6-106 SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth 6-107 SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth, curved tip 6-108 SEMKEN Tissue Forceps 5" (12.7 cm), 2 x 3 teeth 6-109 SEMKEN Dressing Forceps 6" (15.2 cm), serrated 6-109A SEMKEN Tissue Forceps 6" (15.2 cm), 1 x 2 teeth	os





Lightweight Forceps with Fenestrated Handles provide the user with a firm, non-slip grip, reducing hand fatigue.


6-8XL 6-28XL	Dressing Forceps 5½" (14 cm), standard pattern Dressing Forceps 5½" (14 cm), delicate pattern
6-44XL	Tissue Forceps 51/2" (14 cm), 1 x 2 teeth, standard pattern
6-118XL 6-120XL 6-123XL	ADSON Dressing Forceps 4¾" (12.1 cm), serrated ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth ADSON Tissue Forceps 4¾" (12.1 cm), w/platform, 1 x 2 teeth
6-124XL	BROWN-ADSON Tissue Forceps 4%" (12.1 cm), 7 x 7 side grasping teeth
24-570XL	DeBAKEY Tissue Forceps 6" (15.2 cm), atraumatic teeth, 1.5 mm wide tips



6-110

Dressing and Tissue Forceps






HUDSON (EWALD) Dressing Forceps 4¾" (12.1 cm), serrated 6-112 HUDSON (EWALD) Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth 6-114 6-118 ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate ADSON Dressing Forceps 4%" (12.1 cm), serrated, delicate, angled 6-118A ADSON Dressing Forceps 4¾" (12.1 cm), serrated, lightweight, fenestrated handles 6-118XL

17-2510 Micro ADSON Dressing Forceps 43/4" (12.1 cm), serrated tips, 0.5 mm wide

STILLE Pattern Tissue Forceps 5" (12.7 cm), 4 x 5 teeth






6-120	ADSON Tissue Forceps 4%" (12.1 cm), 1 x 2 teeth, delicate			
6-120A	ADSON Tissue Forceps 49/4" (12.1 cm), 1 x 2 teeth, delicate, angled			
6-120XL	ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, lightweight, fenestrated handles			
6-121	ADSON Tissue Forceps 4% (12.1 cm), 1 x 2 teeth and cross serrated tips, delicate		no er	
17-2500	Micro ADSON Tissue Forceps 4% " (12.1 cm), tips 0.5 mm wide, 1 x 2 teeth			
6-122	ADSON Tissue Forceps 4¾" (12.1 cm), 2 x 3 teeth, delicate		(900)	
6-123	ADSON Tissue and Suture Forceps 4%4" (12.1 cm), 1 x 2 teeth with tying platform			
6-123XL	ADSON Tissue Forceps 4%." (12.1 cm), w/platform, 1 x 2 teeth, lightweight, fenestrated handles		Ō	sis.
6-124	BROWN-ADSON Tissue Forceps 4%" (12.1 cm), 7 x 7 side grasping teeth			
6-124XL	BROWN-ADSON Tissue Forceps 4%" (12.1 cm), 7 x 7 side grasping teeth, lightweight, fenestrated handles			
6-125	DeBAKEY-ADSON Tissue Forceps 4¾" (12.1 cm), atraumatic tips 1.5 mm wide	6-124	6-124XL	6-125
	•			








6-159



6 - 128





6-134

6-136





6-126 BROWN-ADSON Tissue Forceps 43/4" (12.1 cm), 7 x 7 side grasping teeth, angular HEANEY Tissue Forceps 51/2" (14 cm), 2 x 3 teeth and cross serrated tips 6-128 STONE Tissue Forceps 6" (15.2 cm), 4 x 5 teeth 6-132 6-134 MARTIN Tissue Forceps 6" (15.2 cm), 7 x 8 teeth

6-136 MARTIN Tissue Forceps 71/2" (19.1 cm), 7 x 8 teeth 6-138 BROWN Tissue Forceps 6" (15.2 cm), 8 x 8 side grasping teeth 6-140 BROWN Tissue Forceps 8" (20.3 cm), 8 x 8 side grasping teeth BROWN Tissue Forceps 10" (25.4 cm), 8 x 8 side grasping teeth 6-141 RUSSIAN Tissue Forceps 6" (15.2 cm) 6-142

RUSSIAN Tissue Forceps 8" (20.3 cm) 6-144 6-145 RUSSIAN Tissue Forceps 10" (25.4 cm) 6-148 BONNEY Tissue Forceps 7" (17.8 cm), 1 x 2 teeth and serrated tips

6-150 BONNEY Tissue Forceps 7" (17.8 cm), 2 x 3 teeth and serrated tips POTTS-SMITH Dressing Forceps 7" (17.8 cm), serrated 6-154 POTTS-SMITH Dressing Forceps 81/4" (21 cm), serrated 6-156 6-158 POTTS-SMITH Dressing Forceps 91/2" (24.1 cm), serrated


POTTS-SMITH Dressing Forceps 12" (30.5 cm), serrated

6-148

6-150








6-160 6-162 6-164 6-164A	POTTS-SMITH Tissue Forceps 7" (17.8 cm), 1 x 2 teeth and serrated tips POTTS-SMITH Tissue Forceps $8\frac{1}{4}$ " (21 cm), 1 x 2 teeth and serrated tips POTTS-SMITH Tissue Forceps $9\frac{1}{2}$ " (24.1 cm), 1 x 2 teeth and serrated tips POTTS-SMITH Tissue Forceps 12" (30.5 cm), 1 x 2 teeth and serrated tips	T Control of the second	a di
6-166	CUSHING Dressing Forceps 7" (17.8 cm), serrated, Gutch handle, scraper end		
6-168	CUSHING Tissue Forceps 7" (17.8 cm), 1 x 2 teeth, Gutch handle, scraper end		
6-180 6-182 6-184 6-186	GERALD Dressing Forceps 7" (17.8 cm), straight, serrated GERALD Tissue Forceps 7" (17.8 cm), straight, 1 x 2 teeth GERALD Dressing Forceps 7" (17.8 cm), curved, serrated GERALD Tissue Forceps 7" (17.8 cm), curved, 1 x 2 teeth		
6-188 6-189	CUSHING Dressing Forceps 7" (17.8 cm), serrated, scraper end CUSHING Tissue Forceps 7" (17.8 cm), 1 x 2 teeth, scraper end	6-188	6-189








6-190 6-192	CUSHING Dressing Forceps 7% " (18.4 cm), bayonet, serrated, scraper end CUSHING Tissue Forceps 7% " (18.4 cm), bayonet, 1 x 2 teeth, scraper end
6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps 61/2" (16.5 cm), serrated
6-198	ADSON Bayonet Dressing Forceps 81/4" (21 cm), serrated
6-200 6-202 6-204	KELLY Tissue Forceps 9" (22.9 cm), 1 x 2 teeth, slotted Gutch handle KELLY Tissue Forceps 9" (22.9 cm), 2 x 3 teeth, slotted Gutch handle KELLY Tissue Forceps 9" (22.9 cm), 3 x 4 teeth, slotted Gutch handle
6-210	MAYO RUSSIAN Tissue Forceps 9" (22.9 cm), fenestrated jaws, slotted grooved handles
6-212	DAVIS Tissue Forceps 10" (25.4 cm)
6-214	SINGLEY Tissue Forceps 9" (22.9 cm), fenestrated serrated jaws
16-80	WANGENSTEEN Tissue Forceps 9" (22.9 cm), rounded jaws with fine cross serrations





18-795A





RHOTON Micro Suture Forceps 7" (17.8 cm), straight, with 0.7 mm RHOTON Micro Forceps 7" (17.8 cm), straight, 0.3 mm tips 18-620 KIRBY Iris Tissue Forceps 4" (10.2 cm), straight, 1 x 2 teeth KIRBY Iris Tissue Forceps 4" (10.2 cm), curved, 1 x 2 teeth 18-624 18-632 KIRBY Tissue Forceps 4" (10.2 cm), angled tips with 1 x 2 teeth 18-793 O'CONNOR Iris Tissue Forceps 3¾" (9.5 cm), angular jaws 0.5 x 9 mm with 1 x 2 teeth 18-794 BRACKEN Iris Dressing Forceps 4" (10.2 cm), with fine cross-serrated tips 18-795 FOERSTER Eye Dressing Forceps 3¾" (9.5 cm), octagonal grip, straight, serrated tips 0.5 mm wide FOERSTER Eye Dressing Forceps 3¾" (9.5 cm), octagonal grip, full curved, 18-795A serrated tips 0.5 mm wide FOERSTER Iris Tissue Forceps 3¾" (9.5 cm), octagonal grip, straight, 18-796 1 x 2 teeth, 0.5 mm wide 18-798 FOERSTER Iris Tissue Forceps 3¾" (9.5cm), octagonal grip, full curved, 1 x 2 teeth, 0.5 mm wide




18-904

18-835

Dressing and Tissue Forceps






10-004	1 mm wide	n sn	11	u
18-814 18-816	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,			
18-818 18-820				
18-825	GILL (CHANDLER) Iris Dressing Forceps 3" (7.6 cm), fine cross serrated tips, 0.6 mm wide	Mann		
18-834	GILL Iris Tissue Forceps 21/6" (7.3 cm), tips 0.6 mm wide, with 1 x 2 teeth	-	1	
18-827	TROUTMAN-BARRAQUER Colibri Type Tissue Forceps 3" (7.6 cm), with 1 x 2 teeth, 0.4 mm, set at 45° angle, with tying platform	18-830	₩ 18-835	
18-828	TROUTMAN-BARRAQUER Colibri Type Tissue Forceps 3" (7.6 cm), with 1 x 2 teeth, 0.12 mm, set at 45° angle, with tying platform	10-030	10-000	
18-829	BOTVIN Iris Dressing Forceps 21/6" (7.3 cm), fine cross serrated tips, 0.5 mi	m wide		
18-830	GILL HESS Iris Forceps 2%" (7 cm), with 1 x 2 teeth forward 0.5 mm			

BONN Micro Tissue Forceps 23/4" (7 cm), 1 x 2 teeth 0.12 mm, with tying platform

STEVENS Iris Tissue Forcens 4" (10.2 cm), straight, 1 v 2 teeth






18-836	BONN Micro Tissue Forceps 3¾" (9.5 cm), 1 x 2 teeth 0.12 mm, with tying platform
18-837	McPHERSON Micro Tissue Forceps $3\frac{1}{2}$ " (8.9 cm), straight, 1 x 2 teeth 0.2 mm, with tying platform
18-838	McPHERSON Micro Tissue Forceps 3½" (8.9 cm), angled, 1 x 2 teeth 0.2 mm, with tying platform
18-866	BISHOP-HARMON Dressing Forceps 3%" (8.6 cm), fine cross serrated tips 0.5 mm wide
18-867	BISHOP-HARMON Micro Tissue Forceps 3%" (8.6 cm), delicate 0.3 mm 1 x 2 teeth
18-868	BISHOP-HARMON Tissue Forceps 3%" (8.6 cm), standard 0.7 mm 1 x 2 teeth
18-840	Strabismus Tissue Forceps 4" (10.2 cm), 1 mm wide tips with 1 x 2 teeth
18-850	GRAEFE Tissue Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth
18-852	GRAEFE Tissue Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth, Non-Magnetic
18-854	GRAEFE Tissue Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth, with catch
















18-880 18-881	PIERSE Colibri Type Tissue Forceps 31/4" (8.3 cm), tips 0.12 mm PIERSE Colibri Type Tissue Forceps 31/4" (8.3 cm), tips 0.3 mm
18-883	CASTROVIEJO Colibri Type Tissue Forceps 4" (10.2 cm), tips 0.12 mm, with tying platform
18-896 18-897	GREEN Tissue Forceps 4" (10.2 cm), T-shape, 12 mm wide 19 x 20 teeth GREEN Tissue Forceps 4" (10.2 cm), T-shape, 5 mm wide 8 x 9 teeth
18-912 18-914	LESTER Tissue Forceps 3%" (9.5 cm), delicate 1 x 2 teeth, 0.6 mm wide LESTER Tissue Forceps 3%" (9.5 cm), delicate 2 x 3 teeth, 1 mm wide





18-916	HARMON Tissue Forceps 3" (7.6 cm), 4 x 5 teeth, delicate 2.5 mm wide tips
18-946 18-947	HARMS Suturing Forceps 4% " (10.5 cm), straight, with 11 mm wide handles, tying platform, 0.6 mm wide at til HARMS Suturing Forceps 4% " (10.5 cm), curved, with 11 mm wide handles, tying platform, 0.6 mm wide at til 10 mm wide handles, tying platform, 0.6 mm wide at til 10 mm wide handles, tying platform, 0.6 mm wide handles, 0.6 mm wide h
18-948 18-949	McPHERSON Micro Iris Suturing Forceps 3% " (8.9 cm), straight, tying platform 0.3×5 mm, no teeth McPHERSON Micro Iris Suturing Forceps 3% " (8.9 cm), angled, tying platform 0.3×5 mm, no teeth
18-956 18-958	NUGENT Utility Forceps 4¼" (10.8 cm), smooth tips 1.2 mm wide NUGENT Utility Forceps 4¼" (10.8 cm), cross serrated tips 1.2 mm wide
18-962 18-964 18-966	McCULLOUGH Utility Forceps 4" (10.2 cm), smooth tips 1.5mm wide McCULLOUGH Utility Forceps 4" (10.2 cm), cross serrated tips 1.5 mm wide McCULLOUGH Suturing Forceps 4" (10.2 cm), cross serrated tips 1.5 mm wide with 1 x 2 teeth











18-978







WILLS HOSPITAL Utility Forceps 4" (10.2 cm), with cross serrated tying 18-970 platform 1.2 mm wide

BONACCOLTO Utility Forceps 4" (10.2 cm), longitudinal serrations with 18-971 cross serrations on tip, 1.2 mm wide with tapered tips


18-978 SCHAAF Foreign Body Forceps 3¾" (9.5 cm), grooved tips

18-1100 LITTAUER Cilia Forceps 31/2" (8.9 cm), 4 mm wide jaws with fine horizontal serrations

BEER Cilia Forceps 31/2" (8.9 cm), 4 mm wide smooth jaws 18-1102

ZIEGLER Cilia Forceps 31/2" (8.9 cm), 2 mm wide smooth jaws 18-1104






18-1106 DOUGLAS Cilia Forceps 31/2" (8.9 cm), thin tips 2.5 mm wide, with fine horizontal serrations 18-1107 Swiss Cilia and Suture Forceps 33/4" (9.5 cm), precision fitted diagonal jaws 3 mm wide 18-1108 BERGH Cilia Forceps 41/4" (10.8 cm), 10 mm wide jaws with horizontal serrations BERGH Cilia Forceps 31/2" (8.9 cm), 5 mm wide jaws with horizontal serrations 18-1110 18-1111 BARRAQUER Cilia and Suture Forceps 41/2" (11.4 cm), with 5 mm smooth platform 18-1112 BENNETT Cilia Forceps 3" (7.6 cm), cup shaped jaws 3 mm diameter BEAUPRE Cilia Forceps 41/2" (11.4 cm), smooth jaws 12 mm long, 45° angle 18-1114 GRADLE Cilia Forceps 4" (10.2 cm), precision ground jaws 1.4 mm wide 18-1116










19-350 19-356	WILDE Dressing Forceps 5" (12.7 cm), angular with serrated tips WILDE Dressing Forceps 5¼" (13.3 cm), cross action, angular with serrated tips	
19-360	BLAKE Dressing Forceps 41/2" (11.4 cm), very delicate, angular with serrated tips	
19-370	LUCAE Dressing Forceps 51/2" (14 cm), bayonet shape with serrated tips	
19-371	JANSEN Dressing Forceps 6% " (15.9 cm), bayonet shape, extra delicate pattern with 2 mm wide serrated tips	
20-120	WILDE Dressing Forceps 6½" (16.5 cm), angular	

20-120





DeBAKEY ATRAUMATIC TISSUE FORCEPS, 2.5 cm wide tips

24-550	6" (15.2 cm)	24-554	91/2" (24.1 cm)
24-552	73/4" (19.7 cm)	24-556	12" (30.5 cm)

DeBAKEY ATRAUMATIC TISSUE FORCEPS, 1.5 mm wide jaw tips

	24-570 24-570XL 24-572 24-574	6" 8"	(15.2 cm) lightweight with fenestrated handles (20.3 cm)
500	OCCUEVA: CT		01/450 17 0 21

24-582	COOLEY Atraumatic Tissue Forceps 8 (10.3 cm), tips 2 mm wide
26-956 26-958	FERRIS SMITH Tissue Forceps 7" (17.8 cm), serrated platform tips with 1x 2 teeth FERRIS SMITH Tissue Forceps 7" (17.8 cm), serrated platform tips with 2 x 3 teeth
33-500	MILTEX Depilatory Forceps 51/2" (13.3 cm), flat square tips with deep serrations




Splinter Forceps















6-300 Plain Splinter Forceps 31/2" (8.9 cm) 6-302 Plain Splinter Forceps 4" (10.2 cm) Plain Splinter Forceps 41/2" (11.4 cm) 6-304 CARMALT Splinter Forceps 41/4" (10.8 cm), straight 6-308 CARMALT Splinter Forceps 41/4" (10.8 cm), curved 6-310 6-312 WALTER Splinter Forceps 41/4" (10.8 cm), straight 6-314 WALTER Splinter Forceps 41/4" (10.8 cm), curved 6-316 FEILCHENFELD Splinter Forceps 3" (7.6 cm) 6-318 FEILCHENFELD Splinter Forceps 41/2" (11.4 cm)

Splinter Forceps










V96-329	Splinter Forceps 3" (7.6 cm), with magnifying glass, in carrying case
6-330	Physicians' Splinter Forceps 5" (12.7 cm)
6-336 6-337	VIRTUS Splinter Forceps 6" (15.2 cm), straight VIRTUS Splinter Forceps 6" (15.2 cm), 45° angle
6-338 6-340 6-342	WILLIAMS Splinter Forceps 5½" (14 cm), straight WILLIAMS Splinter Forceps 5½" (14 cm), 45° angle WILLIAMS Splinter Forceps 5½" (14 cm), right angle






Splinter Forceps









PEET Splinter Forceps 43/4" (12.1 cm), double curve 6-350

6-400 MILTEX Splinter Liberator 5" (12.7 cm), trocar point in reversible handle

V96-410

Professional Splinter Removal Kit. The Liberator is used to approach the embedded fragment by elevating the epidermis. After the fragment is exposed, the Walter Forceps is used for quick and easy removal of splinter.

Display carton with 12 ea. V96-410 Professional Splinter V96-415

Removal Kits (not shown)

Splinter Forceps Splinter







17-302









17-301

SWISS JEWELER STYLE FORCEPS Non-Magnetic Stainless Steel

17-301	Style 1, fine, 4¾" (12.1 cm)
17-302	Style 2, strong, fine, 4¾" (12.1 cm)
17-303	Style 3, narrow, fine, 4¾" (12.1 cm)
17-303C	Style 3C, narrow, fine, 4%" (11.2 cm
17-303X	Style 3, micro-fine, 4¾" (12.1 cm)
17-304	Style 4, fine, 4%" (11.2 cm)
17-305	Style 5, super fine, 4%" (11.2 cm)





Splinter Forceps







17-306



17-307





SWISS JEWELER STYLE FORCEPS Non-Magnetic Stainless Steel


Style 5, micro-fine, 4%" (11.2 cm) 17-305X 17-306 Style 6, pickup, 41/2" (11.4 cm) 17-307 Style 7, fine curved, 41/2" (11.4 cm) 17-307X Style 7, micro-fine, curved, 41/2" (11.4 cm)

THORPE Splinter Forceps 4" (10.2 cm), very sharp points, fine serrated tips $\,$ 18-980

Carb-N-Sert Dressing and Tissue Forceps



CARB-N-SERT DRESSING & TISSUE FORCEPS with tungsten carbide insert tips permit atraumatic and firm grasp of tissue. Tips have 4000 teeth per square inch



6-30TC Dressing Forceps 6" (15.2 cm), cross serrated tips

6-118TC ADSON Dressing Forceps 4¾" (12.1 cm), cross serrated tips 6-123TC ADSON Tissue and Suture Forceps 4¾" (12.1 cm), 1 x 2 teeth


and cross serrated tying platform

BROWN-ADSON Tissue and Suture Forceps 4%" (12.1 cm), 6-124TC

7 x 7 side grasping teeth and cross serrated tying platform




CARB-N-SERT DRESSING & TISSUE FORCEPS with tungsten carbide insert tips permit atraumatic and firm grasp of tissue. Tips have 4000 teeth per square inch



6-154TC POTTS-SMITH Dressing Forceps 7" (17.8 cm), cross serrated tips
6-156TC POTTS-SMITH Dressing Forceps 81/4" (21 cm), cross serrated tips
6-158TC POTTS-SMITH Dressing Forceps 93/4" (24.8 cm), cross serrated tips
6-165TC EVANS Pelvic Tissue Forceps 10" (25.4 cm), cross serrated tips
6-166TC CUSHING Dressing Forceps 7" (17.8 cm), Gutch handle, cross serrated tips
6-198TC ADSON Bayonet Dressing Forceps 81/2" (21.6 cm), cross serrated tips






7-2 7-4	HALSTED Mosquito Forceps 5" (12.7 cm), straight HALSTED Mosquito Forceps 5" (12.7 cm), curved	
18-1934 18-1936	Non-Magnetic HALSTED Mosquito Forceps 5" (12.7 cm) Non-Magnetic HALSTED Mosquito Forceps 5" (12.7 cm)	
7-8 7-10	HALSTED Mosquito Forceps 5" (12.7 cm), straight, extra delicate HALSTED Mosquito Forceps 5" (12.7 cm), curved, extra delicate	/ · \
17-2600	JACOBSON Micro Mosquito Forceps 5" (12.7 cm), straight, extremely delicate	
17-2602	JACOBSON Micro Mosquito Forceps 5" (12.7 cm), curved, extremely delicate	
7-11 7-12	PACKER Mosquito Forceps 5" (12.7 cm), straight, flat jaws PACKER Mosquito Forceps 5" (12.7 cm), curved, flat jaws	7-11 7-12









7-14 7-16 7-20	HALSTED Mosquito Forceps 5" (12.7 cm), straight, 1 x 2 teeth HALSTED Mosquito Forceps 5" (12.7 cm), curved, 1 x 2 teeth HALSTED Mosquito Forceps 5" (12.7 cm), jaws 45° angled on flat, extra delicate
7-24 7-26	HARTMAN Mosquito Forceps 3% " (8.9 cm), straight HARTMAN Mosquito Forceps 3% " (8.9 cm), curved
7-24A	HARTMAN LEE Mosquito Forceps $3\% \mbox{''}$ (8.9 cm), straight jaws, shanks angled on flat
7-30 7-32	PROVIDENCE HOSPITAL Forceps 5½" (14 cm), straight PROVIDENCE HOSPITAL Forceps 5½" (14 cm), curved

































7-102	DANDY Forceps 51/2" (14 cm), curved sideways
7-106	PRATT T-Shape Forceps 6" (15.2 cm)
7-110 7-112	Baby PEAN Forceps 51½" (14 cm), straight, extra delicate Baby PEAN Forceps 5½" (14 cm), curved, extra delicate


ROCHESTER PEAN FORCEPS

Straight	Length	Curved
7-118	5½" (14 cm)	7-136
7-120	61/4" (15.9 cm)	7-138
7-122	71/4" (18.4 cm)	7-140
7-124	8" (20.3 cm)	7-142
7-126	9" (22.9 cm)	7-144
7-128	101/4" (26 cm)	7-146
7-130	12" (30.5 cm)	7-148









ROCHESTER OCHSNER FORCEPS, 1 x 2 TEETH

Straight	Length	Curved
7-150	6¼" (15.9 cm)	7-160
7-152	7¼" (18.4 cm)	7-162
7-154	8" (20.3 cm)	7-164
7-156	9" (22.9 cm)	7-166
7-157	10" (25.4 cm)	7-167
7-170	ROCHESTER CARMALT Forceps 61/4"	(15.9 cm), straight
7-172	ROCHESTER CARMALT Forceps 61/4"	(15.9 cm), curved
7-180	ROCHESTER CARMALT Forceps 8"	(20.3 cm), straight
7-182	ROCHESTER CARMALT Forceps 8"	(20.3 cm), straight
7-190	MICULICZ Hemostatic Forceps 8" (20.3 1 x 2 teeth	cm), slightly curved,








30-1700

30-1718






HEANEY Forceps 81/4" (21 cm), heavy pattern, single tooth, curved HEANEY Forceps 81/4" (21 cm), heavy pattern, 30-1710 double tooth, curved 30-1715 HEANEY Forceps 93/4" (24.8 cm), extra heavy pattern, double tooth, curved HEANEY-BALLENTINE Forceps 81/4" (21 cm), longitudinal serrations, single tooth, curved 30-1716

> BALLENTINE Forceps 81/2" (21.6 cm), longitudinal serrations, single tooth, straight

GARLAND Forceps 81/4" (21 cm), longitudinal serrations, 30-1720 single tooth, angled on flat














30-1730 LONG Forceps 71/2" (19.1 cm), longitudinal serrations, 1 x 2 teeth, curved

PHANEUF Forceps 8" (20.3 cm), horizontal serrated jaws, 30-1735 1 x 2 teeth, straight

30-1740


PHANEUF Forceps 8" (20.3 cm), horizontal serrated jaws, 1 x 2 teeth, angled on flat


PEAN Forceps 9" (22.9 cm), curved, with longitudinal serrations 30-1805 30-1815 PEAN Forceps 101/4" (26 cm), curved, with longitudinal serrations












29-280 WERTHEIM Pedicle Clamp 9¾" (24.8 cm), curved jaws 2" (5.1 cm) long


29-282 WERTHEIM-CULLEN Pedicle Clamp 8½" (21.6 cm), right angle jaws 2" (5.1 cm) long

7-200 MIXTER Forceps 6¼" (15.9 cm), fully curved 7-202 MIXTER Forceps 7¼" (18.4 cm), fully curved 7-206 MIXTER Forceps 9" (22.9 cm), fully curved











25-812	MIXTER Forceps 9" (22.9 cm), right angle jaws, longitudinal serrations
25-814	MIXTER Forceps 11" (27.9 cm), right angle jaws, longitudinal serrations
25-821	MILTEX Forceps 11" (27.9 cm), a delicate light-weight instrument with fine serrated 45 mm long jaws
25-834 25-835	BRIDGE Forceps 11" (27.9 cm), straight, very delicate BRIDGE Forceps 11" (27.9 cm), curved, very delicate
25-836 25-837 25-838	KANTROWITZ Forceps 7½" (19.1 cm), delicate right angle jaws KANTROWITZ Forceps 9½" (24.1 cm), delicate right angle jaws KANTROWITZ Forceps 11" (27.9 cm), delicate right angle jaws
25-840	FINOCHIETTO Forceps 9" (22.9 cm), slightly curved jaws







7-210	Baby MIXTER Forceps 51/4" (13.3 cm), fully curved, extra delicate	
7-210A	MIXTER Forceps, Lee modif. 51/4" (13.3 cm), curved jaws and angled shanks, delicate	
14-46	LOWER Forceps 7" (17.8 cm)	
14-48 14-49	LAHEY Forceps 7½" (19.1cm) LAHEY Forceps 9" (22.9 cm)	
14-50	MOYNIHAN Forceps 81/2" (21.6 cm)	














14-51 GREEN Forceps 81/2" (21.6 cm), delicate 14-52 JOHNS HOPKINS Forceps 8" (20.3 cm)

14-56 GRAY Forceps 9" (22.9 cm), serrated, set of two

14-62 SHALLCROSS Forceps 7" (17.8 cm)













22-660 SCHNIDT Hemostatic Forceps 7½" (19.1 cm), slight curve, open ring
22-662 SCHNIDT Hemostatic Forceps 7½" (19.1 cm), slight curve, closed rings


22-666 SAWTELL Hemostatic Forceps 7½" (19.1 cm). full curve, open ring
22-668 SAWTELL Hemostatic Forceps 7½" (19.1 cm), full curve, closed rings

7-220 PENNINGTON Hemostatic Forceps 6" (15.2 cm)









14-66	LOVELACE Gall Bladder F 1"(2.5 cm) wide jaws	orceps 71/41	' (18.4 cm	ı), angular,
= 000	DAMADDIDOF F	(4E0 \	1. 2.10	L 9 P 1

7-230 BAINBRIDGE Forceps 6' (15.2 cm), straight, longitudinal serrations with cross tip serrations

7-232 BAINBRIDGE Forceps 6" (15.2 cm), curved, longitudinal serrations with cross tip serrations

7-234 BAINBRIDGE Forceps 71/4" (18.4 cm), straight, longitudinal serrations with cross tip serrations

7-236 BAINBRIDGE Forceps 71/4" (18.4cm), curved, longitudinal serrations with cross tip serrations

7-240 ADSON Hemostatic Forceps 71/4" (18.4 cm), straight 7-242 ADSON Hemostatic Forceps 71/4" (18.4 cm), curved


17-2603 JACOBSON Hemostatic Forceps 7" (17.8 cm), curved, extremely delicate

7-250 MOSQUITO Hemostatic Forceps 81/4" (21cm), straight, delicate 7-252 MOSQUITO Hemostatic Forceps 81/4" (21cm), curved, delicate 7-254 MOSQUITO Hemostatic Forceps 81/4" (21cm), right angle, delicate













16-99 DENNIS Clamp 9" (22.9 cm), longitudinal serrations
 16-100 ALLEN Clamp 8" (20.3 cm), longitudinal serrations with 1 x 2 teeth
 16-102 Baby ALLEN Clamp 6½" (16.5 cm), delicate, longitudinal serrations with 1 x 2 teeth
 16-108 Baby DOYEN Forceps 6½" (16.5 cm), straight, flexible blades

Baby DOYEN Forcepts 6½ (16.5 cm), straight, hexible blade with diagonal serrations

16-110 Baby DOYEN Forceps 6½" (16.5 cm), curved, flexible blades with diagonal serrations


16-120KOCHER Forceps 10½" (26.7 cm), straight, longitudinal serrations16-122KOCHER Forceps 10½" (26.7 cm), curved, longitudinal serrations






Hemostatic Forceps _____









- 16-130 MAYO-ROBSON Forceps 10" (25.4 cm), straight, longitudinal serrations
- MAYO-ROBSON Forceps 10" (25.4 cm), curved, longitudinal 16-132 serrations
- SCUDDER Forceps 91/2" (24.1 cm), straight, smooth solid blades 16-140
- DOYEN Forceps 9" (22.9 cm), straight, flexible blades 16-160 with longitudinal serrations
- 16-162 DOYEN Forceps 9" (22.9 cm), curved, flexible blades with longitudinal serrations
- 7-255 GEMINI-MIXTER Forceps 51/2" (14 cm), delicate, full curved jaws
- 7-256 GEMINI-MIXTER Forceps 7" (17.8 cm), delicate, full curved jaws
- 7-257 GEMINI-MIXTER Forceps 8" (20.3 cm), delicate, full curved jaws
- 7-258 GEMINI-MIXTER Forceps 9" (22.9 cm), delicate, full curved jaws
- 7-259 GEMINI-MIXTER Forceps 11" (27.9 cm), delicate, full curved jaws











25-780 RUMEL Forceps 9" (22.9 cm), slight curve, 2%" (6.1 cm), long jaws
25-782 RUMEL Forceps 9" (22.9 cm), medium curve, 2%" (6.1 cm), long jaws
25-784 RUMEL Forceps 9" (22.9 cm), full curve, 2%" (6.1 cm), long jaws
25-786 RUMEL Forceps 9" (22.9 cm), angular, 2%" (6.1 cm), long jaws
25-800 SAROT Artery Forceps 91/2" (24.1 cm), with 21/2" (6.4 cm) long jaws











7-260 7-262 7-264 7-266	FERGUSON Angiotribe Forceps 6½" (16.5 cm), straight FERGUSON Angiotribe Forceps 6½" (16.5 cm), curved FERGUSON Angiotribe Forceps 7½" (19.1 cm), straight FERGUSON Angiotribe Forceps 7½" (19.1 cm), curved
30-1025	JACOBS Forceps 8½" (21.6 cm), straight, serrated jaws, 2 x 2 teeth
30-1030	JACOBS Forceps 8½" (21.6 cm), curved sideways, serrated jaws, 2 x 2 teeth
30-1050-6	MUSEUX Forceps 9½" (24.1 cm), straight, 2 x 2 teeth, 6 mm wide jaws
30-1050-8	MUSEUX Forceps 91/2" (24.1 cm), straight, 2 x 2 teeth, 8 mm wide jaws
30-1050-10	MUSEUX Forceps 9½" (24.1 cm), straight, 2 x 2 teeth, 10 mm wide jaws



Hemostatic Forceps Management











7-280 to 7-298

7-280	JOHNS HOPKINS Bulldog Clamp 11/2" (3.8 cm), straight
7-282	JOHNS HOPKINS Bulldog Clamp 2" (5 cm), straight
7-284	JOHNS HOPKINS Bulldog Clamp 21/2" (6.3 cm), straight
7-286	JOHNS HOPKINS Buildog Clamp 3" (7.6 cm), straight
7-288	JOHNS HOPKINS Bulldog Clamp 31/2" (8.9 cm), straight
7-290	JOHNS HOPKINS Bulldog Clamp 11/2" (3.8 cm), curved
7-292	JOHNS HOPKINS Buildog Clamp 2" (5 cm), curved
7-294	JOHNS HOPKINS Bulldog Clamp 21/2" (6.3 cm), curved
7-296	JOHNS HOPKINS Buildog Clamp 3" (7.6 cm), curved
7-298	JOHNS HOPKINS Bulldog Clamp 31/2" (8.9 cm), curved
7-300	DIEFFENBACH Serrefine 21/4" (5.7 cm), straight
7-302	DIEFFENBACH Serrefine 21/4" (5.7 cm), curved
7-304	DIEFFENBACH Serrefine 1½" (3.8 cm), straight
7-305	DIEFFENBACH Serrefine 11/2" (3.8 cm), curved
7-310	DIETHRICH Micro Bulldog Clamp 1%" (4.7 cm), straight jaws 8 mm x 1.2 mm
7-315	DIETHRICH Micro Bulldog Clamp 1¾" (4.4 cm), angled jaws 8 mm x 1.2 mm


For Atraumatic Bulldog clamps see pages 374 - 376



Towel Clamps











7-500

7-501


7-504

7-506 7-510

7-501

JONES Towel Clamp 31/2" (8.9 cm), standard size JONES Towel Clamp 21/4" (5.7 cm), lightweight BACKHAUS Towel Clamp 31/2" (8.9 cm) BACKHAUS Towel Clamp 51/4" (13.3 cm) ROEDER Towel Clamp 51/4" (13.3 cm), with ball stops

7-514 LORNA Non-Perforating Towel Clamp 31/2" (8.9 cm) 7-515 LORNA Non-Perforating Towel Clamp 51/4" (13.3 cm)




7-515

Sponge Forceps





30-2505



7-610 to 7-616

FOERSTER SPONGE FORCEPS

7-600	91/2" (24.1 cm), straight, serrated jaws
7-602	91/2" (24.1 cm), straight, smooth jaws
7-604	7" (17.8 cm), straight, serrated jaws
7-606	7" (17.8 cm), straight, smooth jaws
7-610	91/2" (24.1 cm), curved, serrated jaws
7-612	91/2" (24.1 cm), curved, smooth jaws
7-614	7" (17.8 cm), curved, serrated jaws
7-616	7" (17.8 cm), curved, smooth jaws


KELLY Sponge Forceps 12½" (31.8 cm), serrated ring jaws, slightly curved






Sponge Forceps







7-620


7-622

7-624

BOZEMAN Sponge Forceps 101/2" (26.7 cm), straight 7-620


7-622 BOZEMAN Sponge Forceps 101/2" (26.7 cm), double curved 7-624 BOZEMAN Sponge Forceps $101\!\!/\!\!2''$ (26.7 cm), double curved, with one large finger ring


FLETCHER Sponge Forceps 91/2" (24.1 cm), straight 7-630 FLETCHER Sponge Forceps 91/2" (24.1 cm), curved 7-632




7-632




























8-22 BROWN Needle Holder 5¼" (13.3 cm), convex jaws, fine serrations
8-23 BROWN Needle Holder 6¾" (17.1 cm), convex jaws, fine serrations
8-26 BOYNTON Needle Holder 5" (12.7 cm), delicate
8-28 MATHIEU Needle Holder 5½" (14 cm), with spring handles
8-30 MATHIEU Needle Holder 6¼" (15.9 cm), with spring handles
8-32 MATHIEU Needle Holder 7½" (19.1 cm), with spring handles
8-34 BAUMGARTNER Needle Holder 5½" (14 cm)






8-59

Needle Holders









8-42	MAYO-HEGAR Needle Holder 5" (12.7 cm)
8-44	MAYO-HEGAR Needle Holder 6" (15.2 cm)
8-46	MAYO-HEGAR Needle Holder 7" (17.8 cm)
8-48	MAYO-HEGAR Needle Holder 8" (20.3 cm)
8-130	MAYO-HEGAR (MASSON) Needle Holder 101/2" (26.7 cm)
8-50	CRILE-WOOD Needle Holder 6" (15.2 cm)
8-51	CRILE-WOOD Needle Holder 7" (17.8 cm)
8-52	CRILE-WOOD Needle Holder 8" (20.3 cm)
8-53	CRILE-WOOD Needle Holder 9" (22.9 cm)
	GILLIES-SHEEHAN Needle Holders combined with Suture Scissors, one fenestrated jaw and one angled finger ring
8-58	61/2" (16.5 cm), straight

6½" (16.5 cm), curved








8-70	ADSON Needle Holder 7" (17.8 cm), fenestrated jaws
8-74	METZENBAUM Needle Holder 744° (18.4 cm), straight, fenestrated jaws
8-76	METZENBAUM Needle Holder 71/4" (18.4 cm), curved, fenestrated jaws
8-80	NEW ORLEANS Needle Holder 7" (17.8 cm)
8-84	M.G.H. Needle Holder 7" (17.8 cm)
8-94	CRILE Needle Holder 6" (15.2 cm)








8-96 SAROT Needle Holder 71/8" (18.1 cm), light model

8-114 SAROT Needle Holder 101/2" (26.7 cm)

8-100 HEANEY Needle Holder 81/2" (21.6 cm), curved jaws

STRATTE Needle Holder 9" (22.9 cm), double curved jaws and bent shanks $\,$ 8-102








8-122 FINOCHIETTO Needle Holder 10½" (26.7 cm), angled jaws
8-134 WANGENSTEEN Needle Holder 10¾" (27.3 cm), long jaws
8-136 WANGENSTEEN Needle Holder 10½" (26.7 cm), short jaws















MILTEX Micro Surgery Needle Holders with round handles, 0.6 mm tips

17-1000	5¼" (13.3 cm), curved jaws
17-980	71/8" (18.1 cm), straight jaws
17-1010	71/s" (18.1 cm), curved jaws
17-1020	51/4" (13.3 cm), curved jaws, with lock
17-1024	71/s" (18.1 cm), straight jaws, with lock
17-1025	71/a" (18.1 cm), curved jaws, with lock





18-1810 GREEN Needle Holder and Suturing Forceps 41/4" (10.8 cm), angled jaws, double end 18-1820 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, without lock 18-1824 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, without lock 18-1828 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock 18-1832 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock 18-1831 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, extra delicate jaws, with lock 18-1833 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, extra delicate jaws, with lock	18-1800	KALT Needle Holder 5¼" (13.3 cm), standard model, delicate jaws
18-1824 CASTROVIEJO Needle Holder 5½" (14 cm), curved, without lock 18-1828 CASTROVIEJO Needle Holder 5½" (14 cm), straight, with lock 18-1832 CASTROVIEJO Needle Holder 5½" (14 cm), curved, with lock 18-1831 CASTROVIEJO Needle Holder 5½" (14 cm), straight, extra delicate jaws, with lock	18-1810	GREEN Needle Holder and Suturing Forceps 41/4" (10.8 cm), angled jaws, double end
18-1828 CASTROVIEJO Needle Holder 5½" (14 cm), straight, with lock 18-1832 CASTROVIEJO Needle Holder 5½" (14 cm), curved, with lock 18-1831 CASTROVIEJO Needle Holder 5½" (14 cm), straight, extra delicate jaws, with lock	18-1820	CASTROVIEJO Needle Holder 51/2" (14 cm), straight, without lock
18-1832 CASTROVIEJO Needle Holder 5½" (14 cm), curved, with lock 18-1831 CASTROVIEJO Needle Holder 5½" (14 cm), straight, extra delicate jaws, with lock	18-1824	CASTROVIEJO Needle Holder 51/2" (14 cm), curved, without lock
18-1831 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, extra delicate jaws, with lock	18-1828	CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock
	18-1832	CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock
18-1833 CASTROVIEJO Needle Holder 5½" (14 cm), curved, extra delicate jaws, with lock	18-1831	CASTROVIEJO Needle Holder 51/2" (14 cm), straight, extra delicate jaws, with lock
	18-1833	CASTROVIEJO Needle Holder 51/2" (14 cm), curved, extra delicate jaws, with lock
















18-1834 PATON Needle Holder 41/2" (11.4 cm), straight, without lock 18-1837 McPHERSON Needle Holder for microsurgery 4" (10.2 cm), curved, smooth tapered jaws, with lock 18-1838 BARRAQUER Needle Holder 5" (12.7 cm), curved, without lock, 10 mm wide hollow round handles BARRAQUER Needle Holder 5" (12.7cm), curved, with lock, 10 mm wide hollow round handles 18-1839

18-1840 BARRAQUER Needle Holder 51/8" (13 cm), curved, delicate tapered jaws, 8 mm wide round handles. without lock

STEVENS Needle Holder 61/2" (16.5 cm), with ratchet lock 18-1858



OUR TOP QUALITY LINE OF NEEDLE HOLDERS



Carb-N-Sert jaws are made of tungsten carbide, the hardest steel made for surgical instruments. They are then bonded to needle holders made of the finest quality stainless steel.

- Teeth are precision milled into the tungsten carbide, forming sharp pyramid-shapes to provide a positive non-slip grip on needles and sutures.
- Carb-N-Sert jaws are custom styled with teeth ranging from 2600 p.sq. inch to 16000 p.sq. inch, according to the size of the instrument and the sizes of needles and sutures used. Some needle holder patterns are also available with smooth Carb-N-Sert jaws, designed to hold the smallest sutures in micro-surgical and plastic surgery procedures.
- The edges around the outside of the jaws are chamfered to prevent

accidental cutting of suture material when held in the jaws.

- The edges around the box lock have distinct non-catch bevels, preventing fraying of suture material.
- Finger Rings are well rounded and polished for better "feel" and prevention of finger fatigue.
- Each instrument is satin finished to reduce glare.


Gold plated finger rings identify these tungsten carbide instruments, and the Miltex catalog number has been etched on each needle holder to facilitate re-ordering.

Carb-N-Sert Needle Holders are specifically designed for precise suturing, thereby facilitating the surgeon's task.





CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures



8-4TC DERF Needle Holder 4¾" (12.1 cm), serrated jaws 4000 teeth p. sq. inch

8-6TC WEBSTER Needle Holder 51/4" (13.3 cm), smooth jaws, standard pattern, delicate


8-7TC WEBSTER Needle Holder 4¾" (12.1 cm), smooth jaws, extra delicate WEBSTER (JOHNSON) Needle Holder 4¾" (12.1cm), serrated jaws 8-7A-TC

4000 teeth p. sq. inch, extra delicate

HALSEY Needle Holder 5" (12.7 cm), smooth jaws, standard pattern 8-8TC 8-8A-TC HALSEY Needle Holder 5" (12.7 cm), serrated jaws 4000 teeth p. sq. inch



CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures





8-15TC 8-16TC 8-17TC



8-28TC 8-30TC

OLSEN-HEGAR NEEDLE HOLDERS combined with Suture Scissors

8-14TC 4¾" (12.1 cm) extra delicate, serrated jaws

4000 teeth p. sq. inch

8-14A-TC 4¾" (12.1 cm), extra delicate, smooth jaws 8-15TC 5½" (14 cm) serrated laws 2600 teeth p. sq. inch 61/2" (16, 5 cm) serrated jaws 2600 teeth p. sq. inch 8-16TC 71/2" (19.1 cm) serrated jaws 2600 teeth p. sq. inch 8-17TC

8-28TC MATHIEU Needle Holder 51/2" (14 cm), serrated jaws


2600 teeth p. sq. inch

8-30TC MATHIEU Needle Holder 61/4" (15.9 cm), serrated jaws

2600 teeth p. sq. inch

8-40TC BAUMGARTNER Needle Holder 51/2" (14 cm), serrated


jaws 4000 teeth p. sq. inch









CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures








8-42TC MAYO-HEGAR Needle Holder 5" (12.7 cm), serrated jaws 2600 teeth p. sq. inch MAYO-HEGAR Needle Holder 6" (15.2 cm), serrated jaws 2600 teeth p. sq. inch 8-44TC MAYO-HEGAR Needle Holder 7" (17.8 cm), serrated jaws 2600 teeth p. sq. inch 8-46TC MAYO-HEGAR Needle Holder 8" (20.3 cm), serrated jaws 2600 teeth p. sq. inch 8-48TC MAYO-HEGAR (MASSON) Needle Holder 101/2" (26.7 cm), serrated jaws 2600 8-130TC teeth p. sq. inch 8-49TC MAYO-HEGAR Needle Holder 12" (30.5 cm), serrated jaws 2600 teeth p. sq. inch 8-50TC CRILE-WOOD Needle Holder 6" (15.2 cm), serrated jaws 4000 teeth p. sq. inch 8-50A-TC CRILE-WOOD Needle Holder 6" (15.2 cm), smooth jaws



CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures



8-51TC 8-52TC 8-53TC 8-55TC 8-57TC	CRILE-WOOD (DeBAKEY) Needle Holder 7" (17.8 cm), serrated jaws 4000 teeth p. sq. inch CRILE-WOOD (DeBAKEY) Needle Holder 8" (20.3 cm), serrated jaws 2600 teeth p. sq. inch CRILE-WOOD (DeBAKEY) Needle Holder 9" (22.9 cm), serrated jaws 2600 teeth p. sq. inch CRILE-WOOD (DeBAKEY) Needle Holder 10½" (26.7 cm), serrated jaws 2600 teeth p. sq. inch CRILE-WOOD (DeBAKEY) Needle Holder 12" (30.5 cm), serrated jaws 2600 teeth p. sq. inch
8-62TC 8-63TC	COOLEY Needle Holder 6" (15.2 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch COOLEY Needle Holder 6%" (17.1 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch
8-65TC 8-66TC 8-67TC 8-69TC	RYDER-Type Intracardiac Needle Holder 5" (12.7 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch RYDER-Type Intracardiac Needle Holder 6" (15.2 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch RYDER-Type Intracardiac Needle Holder 7" (17.8 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch RYDER-Type Intracardiac Needle Holder 9" (22.9 cm), extra delicate, serrated jaws 4000 teeth p. sq. inch
8-81TC	NEW ORLEANS Needle Holder 7" (17.8cm), serrated jaws 2600 teeth p. sq. inch





CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures



8-96TC





8-83TC

8-82TC

JULIAN PATTERN Needle Holder 81/4" (21cm), straight fine tips. serrated jaws 2600 teeth p. sq. inch

8-83TC JULIAN PATTERN Needle Holder 81/4" (21cm), curved fine tips.

8-114TC

serrated jaws 2600 teeth p. sq. inch

8-96TC SAROT Needle Holder 71/4" (18.4cm), serrated jaws

2600 teeth p. sq. inch

8-114TC SAROT Needle Holder 101/2" (26.7cm), serrated jaws

2600 teeth p. sa. inch

8-100TC HEANEY Needle Holder 81/2" (21.6 cm), curved, serrated jaws

2600 teeth p. sq. inch

STRATTE Needle Holder 9" (22.9cm), serrated jaws 8-102TC


2600 teeth p. sq. inch

22-500TC Wire Twister 6" (15.2 cm) 22-502TC Wire Twister 8" (20.3 cm)





CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures





8-101TC	RYDER Needle Holder 5" (12.7 cm), serrated jaws 4000 teeth p. sq. inch
8-103TC	RYDER Needle Holder 6" (15.2 cm), serrated jaws 4000 teeth p. sq. inch
8-104TC	RYDER Needle Holder 7" (17.8 cm), serrated jaws 4000 teeth p. sq. inch
8-105TC	RYDER Needle Holder 8" (20.3 cm), serrated jaws 4000 teeth p. sq. inch
8-106TC	RYDER Needle Holder 9" (22.9 cm), serrated jaws 4000 teeth p. sq. inch
8-107TC	RYDER Needle Holder 101/4" (26 cm), serrated jaws 4000 teeth p. sq. inch
8-134TC	WANGENSTEEN Needle Holder 10%" (27.3 cm), serrated jaws 2600 teeth p. sq. inch
18-1800TC	KALT Needle Holder 51/2" (14 cm), serrated jaws 16000 teeth p. sq. inch

18-1816TC MILTEX-CASTROVIEJO Needle Holder 41/6" (10.5 cm), round handles,

with lock, delicate smooth jaws

BIHRLE Needle Holder/Clamp 11" (27.9 cm)

18-1816TC

29-310TC





CARB-N-SERT NEEDLE HOLDERS

with tungsten carbide inserts insure firm holding of needles and sutures







18-1828A-TC

18-1832A-TC

17-1090TC to 17-1092TC

CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock, smooth jaws 18-1828TC

CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock, serrated jaws 16000 teeth p. sg. inch 18-1828A-TC

CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock, smooth jaws 18-1832TC

CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock, serrated jaws 16000 teeth p. sq. inch 18-1832A-TC

17-1090TC JACOBSON Microvascular Needle Holder 7" (17.8 cm), straight, with lock, serrated jaws

16000 teeth p. sq. inch

17-1091TC JACOBSON Microvascular Needle Holder 7" (17.8 cm), curved, with lock, serrated jaws

16000 teeth p. sq. inch


JACOBSON Microvascular Needle Holder 81/2" (21.6 cm), straight, with lock, serrated jaws 17-1092TC

16000 teeth p. sq. inch



CARB-N-SERT NEEDLE HOLDERS with tungsten carbide inserts insure firm holding of needles and sutures

With QUICK-RELEASE RATCHET (suitable for Left Handers)



8-706TC

8-753TC






8-744TC	MAYO-HEGAR Quick-Release Needle Holder 6" (15.2 cm), serrated jaws 2600 teeth p. sq. inch
8-746TC	MAYO-HEGAR Quick-Release Needle Holder 7" (17.8 cm), serrated jaws 2600 teeth p. sq. inch
8-748TC	MAYO-HEGAR Quick-Release Needle Holder 8" (20.3 cm), serrated jaws 2600 teeth p. sq. inch
8-750TC	CRILE-WOOD Quick-Release Needle Holder 6" (15.2 cm), serrated jaws 4000 teeth p. sq. inch
8-751TC	CRILE-WOOD Quick-Release Needle Holder 7" (17.8 cm), serrated jaws 4000 teeth p. sq. inch
8-752TC	CRILE-WOOD Quick-Release Needle Holder 8" (20.3 cm), serrated jaws 4000 teeth p. sq. inch

CRILE-WOOD Quick-Release Needle Holder 9" (22.9 cm), serrated jaws 4000 teeth p. sq. inch


WEBSTER Quick-Release Needle Holder 5" (12.7 cm), smooth jaws





All needles shown actual size






3/8 Circle Cutting Edge

> MS-140S Same as above with Spring Eye



MS-153 MAYO INTESTINAL 3/8 Circle Taper Point


All needles shown actual size.



FERGUSON 1/2 Circle Taper Point



FISTULA 1/2 Circle Cutting Edge



MS-140A REGULAR SURGEON'S 3/8 Circle Taper Point




MS-179 MURPHY'S INTESTINAL 1/2 Circle Taper Point







MS-431 KEITH ABDOMINAL Straight Triangular Point



MARTIN'S UTERINE 1/2 Circle Cutting Edge

Eye Needles



MS-172 EYE NEEDLE 3/8 Circle Cutting Edge



MS-175 EYE NEEDLE 1/2 Circle Cutting Edge



MS-177 JAMESON STRABISMUS Cutting Edge



MS-176 EYE NEEDLE 1/2 Circle Taper Point

Dental Needles



MS-155 LANE CLEFT PALATE 3/8 Circle Cutting Edge




MS-156 LANE CLEFT PALATE 1/2 Circle Cutting Edge



MSRC-14 MSRC-18 DENTAL REVERSE CUTTING EDGE 1/2 Circle






All needles have Regular eyes, unless otherwise noted.








9-1 to 9-10







MICHEL WOUND CLIPS (100 to an envelope, 1000 to a box), stainless

9-1	size 7.5 mm	9-6	size 16 mm
9-2	size 11 mm	9-8	size 18 mm
9-4	size 14 mm	9-10	size 22 mm

9-32 McKENZIE Brain Clip Holding Rack, for 20 clips

3-50 Sterilizing Spiral Rack for Surgical Needles

BABCOCK SUTURE WIRE, Stainless Steel Type AISI 316-L

B&S Gauge	Diameter	1 oz. (28 g) Spools	4 oz. (113 g) Spools
18	.040" (1 mm)	9-40-18	9-44-18
20	.032" (.8 mm)	9-40-20	9-44-20
22	.025" (.6 mm)	9-40-22	9-44-22
24	.020" (.5 mm)	9-40-24	9-44-24
25	.018" (.45 mm)	9-40-25	9-44-25
26	.016" (.4 mm)	9-40-26	9-44-26
28	.012" (.3 mm)	9-40-28	9-44-28
30	.010" (.25 mm)	9-40-30	9-44-30
32	.008" (.2 mm)	9-40-32	9-44-32
34	.006" (.15 mm)	9-40-34	9-44-34
36	.005" (.13 mm)	9-40-36	














9-50	MICHEL Clip Applying Forceps 4%" (12.1 cm)
9-52	HEGENBARTH Clip Applying Forceps 4¾" (12.1 cm), cross action
9-54	WACHENFELD Clip Applying-Removing Forceps $4\%^{\circ}$ (12.1 cm), cross action
9-60	MICHEL Clip Applying-Removing Forceps 5" (12.7 cm), double ended
9-62	MILTEX Skin Staple Removing Forceps 51/2" (14 cm)






9-74	SMITHWICK Clip Applying Forceps 9" (22.9 cm), curved jaws, straight handles
9-80	McKENZIE Clip Applying Forceps 6" (15.2 cm), curved jaws, straight handles
9-90	SCOVILLE Clip Applying Forceps 6" (15.2 cm) shaft
9-94	RICHTER HEATH Clip Removing Forceps 51/2" (14 cm)








18-1107







17-303 Suture Forceps (style 3), narrow, fine 4¾" (12.1 cm)

18-1107 Swiss Cilia and Suture Forceps 3¾" (9.5 cm), precision fitted diagonal jaws 3 mm wide


9-96 HEATH Suture Scissors 6¾" (15.4 cm), serrated blade

9-98 NEW'S Suture Scissors 5½" (14 cm), angled on flat

9-100 SPENCER Stitch Scissors 3½" (8.9 cm), delicate

9-101 SHORTBENT Stitch Scissors 3½" (8.9 cm), curved, delicate

9-102 LITTAUER JR. Stitch Scissors 4½" (11.4 cm), light pattern
















9-103 NORTHBENT Stitch Scissors 4¾" (12.1cm), curved, light pattern LITTAUER Stitch Scissors 51/2" (14 cm), standard pattern 9-104 9-106 LONG Oral Surgery Stitch Scissors 6" (15.2 cm) MILTEX Stitch Scissors 41/2" (11.4 cm), angled, delicate 9-108 9-110 O'BRIEN Stitch Scissors 3¾" (9.5 cm), angled, sharp points













18-1486 WESTCOTT Stitch Scissors 41/2" (11.4 cm), curved, sharp points

Eye Suture (GRADLE) Scissors 3¾" (9.5 cm), slightly curved, sharp points 18-1652 18-1653* Eye Suture Scissors 4" (10.2 cm), slightly curved, sharp points, ribbon type

21-536 LITTLER Suture Carrying Scissors 45/8" (11.8 cm), with suture hole


in blades, curved


For a complete selection of Carb-N-Sert Eye Scissors with Tungsten Carbide blades and Super Cut Scissors see pages 58, 62 and 63.

*Ribbon type scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.














WIRE CUTTING SCISSORS 4%" (12.1 cm)

9-120 Straight, one serrated blade

9D-121 Straight

9-122 Curved, one serrated blade

9D-123 Curved

Angled to side, one serrated blade 9-124 9-124TC

Angled to side, one serrated blade, Carb-N-Sert inserts


WIRE CUTTING SCISSORS 41/4" (10.8 cm)

Straight, one serrated blade 9D-132

9D-133 Straight

9D-134 Curved, one serrated blade

9D-135 Curved












9-128 Double Action Parallel Pliers and Wire Cutter 4½" (11.4 cm), nickel plated

22-500TC Wire Twister 6" (15.2 cm), CARB-N-SERT jaws 22-502TC Wire Twister 8" (20.3 cm), CARB-N-SERT jaws

27-176TC MILTEX Pin Puller 5½" (14 cm), CARB-N-SERT cross serrated top jaw measures 1 cm in length. Grooved cross serrated opposing jaw (groove measures 1 cm in length).

9-140 DESCHAMPS Ligature Carrier 8" (20.3 cm), blunt, left
 9-142 DESCHAMPS Ligature Carrier 8" (20.3 cm), blunt, right
 9-146 DESCHAMPS Ligature Carrier 8" (20.3 cm), sharp, left
 9-148 DESCHAMPS Ligature Carrier 8" (20.3 cm), sharp, right




For additional Orthopedic Pin & Wire instruments see pages 443,444









9-156



12-90

12-92



21-900









DUPUYS-WEISS Ligature Carrier 81/2" (21.6 cm), right 9-154 9-156 DUPUYS-WEISS Ligature Carrier 81/2" (21.6 cm), left 12-90 LAHEY Ligature Carrier 71/2" (19.1 cm), curved left LAHEY Ligature Carrier 71/2" (19.1 cm), curved right 12-92 21-900 OBWEGESER Mandibular Awl for wire suturing 51/2" (14 cm), slightly curved 21-901 OBWEGESER Maxillary Awl for wire suturing 7" (17.8 cm), strong curve 21-902 OBWEGESER Zygomatic Arch Awl for wire suturing 9" (22.9 cm), slightly curved RANEY Scalp Clip Applying Forceps 26-590

61/4" (15.9 cm)








Probes and Directors





PROBES, DOUBLE ENDED, Malleable

PROBES WITH EYE, Malleable

Nickel-Silver		Sterling	Stainless		Sterling
10-2-NS	4½" (11.4 cm)	10-2-ST	10-22-SS	4½" (11.4 cm)	10-22-ST
10-4-NS	5" (12.7 cm)	10-4-ST	10-24-SS	5" (12.7 cm)	10-24-ST
10-6-NS	5½" (14 cm)	10-6-ST	10-26-SS	5½" (14 cm)	10-26-ST
10-8-NS	6" (15.2 cm)	10-8-ST	10-28-SS	6" (15.2 cm)	10-28-ST
10-10-NS	7" (17.8 cm)	10-10-ST	10-30-SS	7" (17.8 cm)	10-30-ST
10-12-NS	8" (20.3 cm)	10-12-ST	10-32-SS	8" (20.3 cm)	10-32-ST
10-14-NS	10" (25.4 cm)	10-14-ST	10-34-SS	10" (25.4 cm)	10-34-ST

10-38

Probe with Myrtle Leaf 5" (12.7 cm)

GROOVED DIRECTORS

With Tongue Tie		With Probe Tip and Tongue Tie
10-69 10-70	4½" (11.4 cm) 5" (12.7 cm)	10-79 10-80
10-72	5½" (14 cm)	10-82
10-74	6" (15.2 cm)	10-84
10-76	8" (20.3 cm)	_

28,100


LARRY Probe and Grooved Director 7" (17.8 cm)














FRENCH PATTERN LACRIMAL PROBES 5" (12.7 cm)

Bayonet shape


18-710	size 00	18-716	size 2
18-712	size 0	18-718	size 3
18-714	size 1	18-720	size 4


WILLIAMS LACRIMAL PROBES 5" (12.7 cm), sterling


18-722	size 0000-000	18-728	size 3-4
18-724	size 00-0	18-730	size 5-6
18-726	size 1-2	18-732	size 7-8


BOWMAN LACRIMAL PROBES 5" (12.7 cm), sterling

18-740	size 0000-000	18-746	size 3-4
18-742	size 00-0	18-748	size 5-6
18-744	size 1-2	18-750	size 7-8











19-210





28-106	BUIE Fistula Probe, with 51/2" (14 cm) sterling shaft, 63/4" (17.1 cm) overall
28-108	PRATT Rectal Probe, with 7¾" (19.7 cm) sterling shaft, 11" (27.9 cm) overall
28-109	BARR Double Ended Probe 7" (17.8 cm), rigid
14-2	OCHSNER (FENGER) Gall Duct Probe, cannulated flexible spiral 14" (35.6 cm) with stylet
14-4-15 14 - 4-18	MAYO Common Duct Probe 10" (25.4 cm), malleable shaft, 15 French (5 mm) MAYO Common Duct Probe 10" (25.4 cm), malleable shaft, 18 French (5.9 mm)
14-12	MIXTER Irrigating Dilaprobe 12" (30.5 cm), malleable, with Luer Lock, silver plated



Suction Tubes







2-110

2-110SS





2-106 2-106SS



2-100 PYNCHON Suction Tube 11" (27.9 cm), tip with large openings, with removable tubing connection, chrome

2-109 YANKAUER Pediatric Suction Tube 8½" (21.6 cm), with removable tip, delicate pattern

2-110 ANDREWS-PYNCHON Suction Tube 10" (25.4 cm), delicate pattern, chrome

2-110SS ANDREWS-PYNCHON Suction Tube 10" (25.4 cm), delicate pattern

2-104 YANKAUER Suction Tube 11% (29.8 cm), with removable tip and tubing connection, chrome

2-104SS YANKAUER Suction Tube 11%" (29.8 cm), with removable tip and tubing connection

2-106 Tip for Yankauer Suction Tube, chrome

2-106SS Tip for Yankauer Suction Tube

10-220 ANTHONY Suction Tube 10" (25.4 cm), lumen 2 mm, chrome



Suction Tubes Succession













20-466

COUPLAND Suction Tube 10" (25.4 cm), with 2 interchangeable

tips 1.5 and 2 mm, chrome

10-310 10-312 POOLE Suction Tube 23 Fr. (7.6 mm), curved POOLE Suction Tube 30 Fr. (9.9 mm), straight

HOUSE SUCTION-IRRIGATION TUBES, 9 cm working length

19-562 1.8 mm Suction Tube and 1.8 mm Irrigating Tube
19-563 2.5 mm Suction Tube and 2 mm Irrigating Tube
19-564 3 mm Suction Tube and 2 mm Irrigating Tube
19-565 4.5 mm Suction Tube and 2 mm Irrigating Tube

FRAZIER FERGUSON SUCTION TUBES, with finger cut-off, angled

 19-570
 6 French (2 mm)
 19-575
 10 French (3.3 mm)

 19-572
 7 French (2.3 mm)
 19-576
 11 French (3.6 mm)

 19-573
 8 French (2.6 mm)
 19-577
 12 French (4 mm)

 19-574
 9 French (2 mm)
 19-577
 12 French (4 mm)


19-574 9 French (3 mm)

20-460 FRAZIER FERGUSON Suction Tube 9 French (3 mm), straight, with finger cut-off

20-460



Suction Tubes



19-580 19-582 19-584	BARON Suction Tube 3 French (1 mm), with finger cut-off, 7.5 cm tube BARON Suction Tube 5 French (1.7 mm), with finger cut-off, 7.5 cm tube BARON Suction Tube 7 French (2.3 mm), with finger cut-off, 7.5 cm tube
19-1820	Ear Suction Tube with 2 interchangeable tips 20 ga. and 22 ga.
19-589	ANTHONY Suction Tube 61/2" (16.5 cm), with 3 mm lumen
26-780 26-782	ADSON Suction Tube 6" (15.2 cm), with finger valve and obturator, 11 Fr. (3.6 mm) ADSON Suction Tube 8" (20.3 cm), with finger valve and obturator, 15 Fr. (5 mm)
28-370	BUIE Suction Tube 16" (40.6 cm), with finger valve, 15 Fr. (5 mm)





10-290 10-291

NABATOFF Vein Stripper Set, including 2 probe tips, 9 mm, 12 mm and 15 mm olives, and handle with screws on either end of 34" (86.4 cm) cable. Complete in sterilizer case. 10-290

NABATOFF Vein Stripper Set, without sterilizer case. 10-291




Retractors











11-2 JANSEN Retractor 41/4" (10.8 cm), spread 3" (7.6 cm), 3 x 3 blunt prongs 10 mm deep, standard model JANSEN Retractor 41/2" (11.4 cm), spread 3" (7.6 cm), 4 x 4 blunt prongs 11-4 20 mm deep 11-6 ALLPORT Retractor 4" (10.2 cm), spread 21/2" (6.4 cm), 4 x 4 sharp prongs

11-8 JANSEN-WAGNER Retractor 5" (12.7 cm), 5 x 5 sharp prongs, self-adjusting

11-10 DAVIS Retractor 41/2" (11.4 cm), 3 x 3 blunt prongs



11-10

Retractors Retractors











11-12 11-13	ALM Retractor 2¾" (7 cm), spread 2¼" (5.7 cm), prongs ½" (.32 cm) deep ALM Retractor 3%" (9.8 cm), spread 3" (7.6 cm), prongs ½" (1.3 cm) deep
11-14	MILTEX Cross Action Retractor 4" (10.2 cm), sharp prongs
11-16	MILTEX Cross Action Retractor 4" (10.2 cm), blunt prongs
11-17	FARR Wire Retractor 5½" (14 cm), prongs ¾" (1.9 cm)
11-18	FARR Wire Retractor 5½" (14 cm), prongs 1½" (2.8 cm)
11-19	Spring Wire Retractor 31/4" (8.3 cm), with blades 10 mm wide x 13 mm deep



Retractors











18-80 AGRICOLA Retractor 11/2" (3.8 cm), 3 x 3 pointed prongs 4 mm deep, non-magnetic 18-84 MELLER Retractor 27%" (7.3 cm), 3 x 3 prongs 3 mm deep, non-magnetic 18-88 GOLDSTEIN Retractor 11/6" (2.8 cm), 3 x 3 pointed prongs 5 mm deep, non-magnetic 18-96 STEVENSON Retractor 31/4" (8.3 cm), 3 x 3 prongs 3 mm deep, non-magnetic

Retractors East










11-23





11-20 LAHEY Retractor 73/4" (19.7 cm) 11-22 GREEN Retractor 81/2" (21.6 cm) 11-23 LITTLE Retractor 73/4" (19.7 cm) 18-110 DESMARRES Retractor 51/2" (14 cm), size 1, 11 mm 18-112 DESMARRES Retractor 51/2" (14 cm), size 2, 13 mm 18-114 DESMARRES Retractor 51/2" (14 cm), size 3, 15 mm 18-116 DESMARRES Retractor 51/2" (14 cm), size 4, 17 mm 11D-27 Flexible Neck Rake Retractor 5" (12.7 cm), 3 sharp prongs 11-30 Flexible Neck Rake Retractor 6" (15.2 cm), 1 sharp prong 11-32 Flexible Neck Rake Retractor 6" (15.2 cm), 2 sharp prongs Flexible Neck Rake Retractor 6" (15.2 cm), 3 sharp prongs 11-34 11-36 Flexible Neck Rake Retractor 6" (15.2 cm), 1 blunt prong 11-38 Flexible Neck Rake Retractor 6" (15.2 cm), 2 blunt prongs 11-40 Flexible Neck Rake Retractor 6" (15.2 cm), 3 blunt prongs



11-30 to 11-40




20-838

20-840

Retractors Emili





11-50 11-52	Rigid Neck Rake Retractor 6" (15.2 cm), 1 sharp prong Rigid Neck Rake Retractor 6" (15.2 cm), 2 sharp prongs	
11-54 11-56 11-58 11-60	Rigid Neck Rake Retractor 6" (15.2 cm), 3 sharp prongs Rigid Neck Rake Retractor 6" (15.2 cm), 1 blunt prong Rigid Neck Rake Retractor 6" (15.2 cm), 2 blunt prongs Rigid Neck Rake Retractor 6" (15.2 cm), 3 blunt prongs	
18-100 18-102	KNAPP (ROLLET) Retractor 51/4" (13.3 cm), with 4 sharp prongs KNAPP Retractor 51/4" (13.3 cm), with 4 blunt prongs	
18-439	McGANNON Retractor 5" (12.7 cm), fenestrated blade 2.5 mm wide	
20-835	LOTHROP Retractor 8" (20.3 cm) shaft, grooved blade 40 x 18 mm with 5 mm lip	
20-836	LOVE Retractor, 18 mm wide blade	

LOVE Retractor, 20 mm wide blade


LOVE Retractor, 22 mm wide blade

20-836 to 20-840















21-107 MILLARD Thimble Retractor, %" (1.6 cm) working length, sharp single hook 21-108 COTTLE Thumb Retractor, 11/2" (3.8 cm) working length, sharp double hook 8 mm wide 21-148 FOMON Retractor 61/2" (16.5 cm), two ball end prongs 11 mm wide 21-150 COTTLE Retractor 51/2" (14 cm), four blunt prongs, 2 center prongs recessed CONVERSE Alar Retractor 4" (10.2 cm), double end, with blunt double hooks 10 mm and 14 mm wide 21-163

COTTLE Knife Guide and Retractor 8" (20.3 cm), double end, ball tip prongs 10 mm wide, and 5 mm wide slotted hook

21-165


21-165



Retractors








21-176	NEIVERT Knife Guide and Retractor 6% (16.5 cm), double end, open blade 14 mm wide
21-182	CONVERSE Blade Retractor 3¾" (9.5 cm), size 1
21-183	CONVERSE Blade Retractor 3¾" (9.5 cm), size 2
21-184	CONVERSE Blade Retractor 3¾" (9.5 cm), size 3
21-185	CONVERSE Blade Retractor 3¾" (9.5 cm), size 4
21-188	AUFRICHT Retractor-Speculum 61/2" (16.5 cm), solid blade 45 mm long
21-190	AUFRICHT Retractor 7" (17.8 cm), fenestrated blade 45 mm long



11-79




VOLKMAN FINGER RETRACTORS 41/2" (11.4 cm)

Sharp		Blunt
11-42	2 prongs	11-46
11-43	3 prongs	11-47
11-44	4 prongs	11-48
11-45	6 prongs	11-49

- 11-73 RAGNELL Retractor 6" (15.2 cm), double end 3 x 8 mm and 5 x 15 mm blades, delicate
 11-74 SENN Retractor 6%" (16.2 cm), sharp, double end
 11-76 SENN Retractor 6%" (16.2 cm), blunt, double end
 11-78 CRILE Retractor 4½" (11.4 cm), double end 3 x 11 mm and 6 x 20 mm blades
- 11-79 LUER S-Shape Retractor $5\frac{1}{2}$ " (14 cm), double end 10 x 5 mm and 14 x 9 mm








21-180 LOVE Nerve Retractor 81/2" (21.6 cm), straight, 7 mm blade 26-1352 LOVE Nerve Root Retractor 81/2" (21.6 cm), angled 45°, 7 mm blade 26-1354 LOVE Nerve Root Retractor 81/2" (21.6 cm), angled 90°, 7 mm blade 26-1360 SCOVILLE Nerve Root Retractor 9" (22.9 cm), straight, 8 mm wide blade SCOVILLE Nerve Root Retractor 9" (22.9 cm), angular, 8 mm wide blade 26-1362 11-80 MATHIEU Retractor 61/2" (16.5 cm), sharp, double end MATHIEU Retractor 61/2" (16.5 cm), blunt, double end 11-81









11-82	MEYERDING Skin Hook and Retractor 61/4" (15.9 cm), double end
11-88	MEYERDING Finger Retractor 7" (17.8 cm), style A. Blade ¼" (.64 cm) x 1½ 6" (1.7 cm)
11-90	MEYERDING Finger Retractor 7" (17.8 cm), style B. Blade 3/16" (.48 cm) x 5/1" (1.6 cm)
11-92	MEYERDING Finger Retractor 7" (17.8 cm), style C. Blade 1/8" (.32 cm) x 1/4" (.64 cm)
11-94	MEYERDING Finger Retractor 7" (17.8 cm), style D. Blade 32" (.24 cm) x 36" (1 cm)
11-96	MEYERDING Finger Retractor 7" (17.8 cm), style E. Blade $\frac{1}{4}$ " (.64 cm) x $\frac{1}{4}$ " (2.2 cm)
23-1047	LUKENS Retractor 6¾" (17.1 cm), 10 mm wide curved and 7 mm wide angular blades
23-1049	JACKSON Retractor 7" (17.8 cm), double end, 2 blunt prongs 5%" (2.4 cm) wide, and solid blade %" (1 cm) wide







25-274





11-102





SAUERBRUCH Retractor 9" (22.9 cm), blade 1/8" (2.2 cm) x 21/4" (5.7 cm) 25-272 SAUERBRUCH Retractor 9" (22.9 cm), blade 3/4" (1.9 cm) x 27/6" (7.3 cm) 25-274 11-100 CUSHING Decompression Retractor 81/2" (21.6 cm) CUSHING Vein Retractor 9" (22.9 cm), blade 11 x 13 mm 11-102 11-110 FARABEUF Retractors 6" (15.2 cm), double end, set of two Baby FARABEUF Retractors 43/4" (12.1 cm), double end, set of two 11-112 PARKER Retractors 5" (12.7 cm), double end, set of two 11-114 11-116 PARKER Retractors 71/4" (18.4 cm), double end, set of two


Retractors Extractors











22-200 MILTEX Lip Retractor 6" (15.2 cm), double end 25 mm and 36 mm wide
22D-202 COLUMBIA Lip Retractor 5½" (14 cm), double end 37 mm and 50 mm wide
22D-209 UNIVERSITY OF MINNESOTA Cheek Retractor 5½" (14 cm), double end
14 mm and 20 mm wide
22-210 Cheek Retractor 7¾" (19.7 cm), with swivel blade 35 mm wide
22-854 HURD Tonsil Dissector and Pillar Retractor 9" (22,9 cm), standard pattern,
9 x 25 mm blade




22-854



Retractors Esse





11-144

11-146











22-856 MORRISON-HURD Tonsil Dissector and Pillar Retractor 9" (22.9 cm), flat handle, 8 x 25 mm serrated blade 22-890 COLVER Tonsil Examining Retractor 8" (20.3 cm), blunt tip 11-118 MAYO-COLLINS Retractor 6" (15.2 cm), double end, set of two 11-120 GOELET Retractor 71/2" (19.1 cm), double end 11-122 U.S. ARMY Retractor 81/2" (21.6 cm), double end, set of two 11-130 Ribbon Retractor 3/4" (1.9 cm) x 13" (33 cm), malleable 11-132 Ribbon Retractor 1" (2.5 cm) x 13" (33 cm), malleable 11-134 Ribbon Retractor 11/4" (3.2 cm) x 13" (33 cm), malleable 11-136 Ribbon Retractor 11/2" (3.8 cm) x 13" (33 cm), malleable 11-140 Ribbon Retractor 2" (5.1 cm) x 13" (33 cm), malleable

> Pediatric Ribbon Retractor %" (1.6 cm) x 8" (20.3 cm), malleable Pediatric Ribbon Retractor 1" (2.5 cm) x 8" (20.3 cm), malleable



11-130 to 11-146






11-150	PARKER-MOTT Retractors 6" (15.2 cm), double end, set of two
11-152 11-153 11-154	ROUX Retractors 6%" (16.8 cm), double end, set of three ROUX Retractor 5½" (13.3 cm), small %" (1.9 cm) x 1" (2.5 cm) ROUX Retractor 5¾" (14.6 cm), medium %" (2.2 cm) x 1%" (3.5 cm)
11-155	ROUX Retractor 6%" (16.8 cm), large 1" (2.5 cm) x 1%" (4.1 cm)
11-160 11-162 11-164 11-166	LANGENBECK Retractor 8" (20.3 cm), blade %" (1 cm) x 1% " (3.2 cm) LANGENBECK Retractor 8" (20.3 cm), blade %" (1 cm) x 1% " (3.8 cm) LANGENBECK Retractor 8" (20.3 cm), blade $\%$ " (1.3 cm) x 1% " (3.8 cm) LANGENBECK Retractor 8" (20.3 cm), blade $\%$ " (1.6 cm) x 1% " (3.8 cm)
11-180 11-182 11-184	HEANEY Retractor 10½" (26.7 cm), blade 1" (2.5 cm) x 1¾" (4.4 cm) HEANEY Retractor 10½" (26.7 cm), blade 1½" (2.8 cm) x 2½" (7.3 cm) HEANEY Retractor 10½" (26.7 cm), blade 1½" (2.8 cm) x 3½" (8.3 cm)








11-185





HEANEY-SIMON Retractor 111/2 " (29.2 cm), blade 1" (2.5 cm) x 41/2" (11.4 cm)


30-320 to 30-330

30-300 30-305 30-310	JACKSON Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 3" (7.6 cm) JACKSON Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 3½" (8.9 cm) JACKSON Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 4" (10.2 cm)
30-320 30-325 30-330	EASTMAN Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 3" (7.6 cm) EASTMAN Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 3½" (8.9 cm) EASTMAN Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 4" (10.2 cm)
30-340	DE LEE Universal Retractor 91/2" (24.1 cm), Blade 23/4" (7 cm) x 2" (5.1 cm)



30-340





30-395 30-406 30-415	DOYEN Retractor 9" (22.9 cm), Blade 1%" (4.4 cm) x 21/8" (6.4 cm) DOYEN Retractor 9" (22.9 cm), Blade 21/8" (5.4 cm) x 3" (7.6 cm) DOYEN Retractor 9" (22.9 cm), Blade 21/4" (5.7 cm) x 31/8" (8.9 cm)
11-192 11-194 11-196	JUNG Mini-Chole Retractor 10½" (26.7 cm), 3 cm width blade JUNG Mini-Chole Retractor 9" (22.9 cm), 4 cm width blade JUNG Mini-Chole Retractors, set of 2 each of 3 cm width blade and 1 each of 4 cm width blade
11-200 11-202 11-204	DEAVER Retractor 1" (2.5 cm) x 9" (22.9 cm) DEAVER Retractor 1" (2.5 cm) x 12" (30.5 cm) DEAVER Retractor 1" (2.5 cm) x 13" (33 cm)
















11-206 11-208 11-210	DEAVER Retractor 11½" (3.8 cm) x 12" (30.5 cm) DEAVER Retractor 2" (5.1 cm) x 12" (30.5 cm) DEAVER Retractor 3" (7.6 cm) x 12" (30.5 cm)
11-218 11-219 11-220	Pediatric DEAVER Retractor %" (1.6 cm) x 8½" (21.6 cm) Pediatric DEAVER Retractor ¾" (1.9 cm) x 7" (17.8 cm) Pediatric DEAVER Retractor %" (2.2 cm) x 8" (20.3 cm)
11-221 11-222 11-223	DEAVER Retractor with Hollow Grip Handle 1" (2.5 cm) x 12" (30.5 cm) DEAVER Retractor with Hollow Grip Handle 1" (2.5 cm) x 13" (33 cm) DEAVER Retractor with Hollow Grip Handle 1½" (3.8 cm) x 12" (30.5 cm)












11-224 DEAVER Retractor with Hollow Grip Handle, 2" (5.1 cm) x 12" (30.5 cm) 11-224A DEAVER Retractor with Hollow Grip Handle, 3" (7.6 cm) x 12" (30.5 cm) 11-225 DEAVER Retractor with Hollow Grip Handle, 4" (10.2 cm) x 12" (30.5 cm) 11-240 MEYERDING Retractor 91/2" (24.1 cm), large blade, 31/2" (8.9 cm) x 2" (5.1 cm) 11-242 MEYERDING Retractor 9½" (24.1 cm), medium blade, 3" (7.6 cm) x 1" (2.5 cm) 11-244 MEYERDING Retractor 9½" (24.1 cm), small blade, 2" (5.1 cm) x %" (1.6 cm) 11-246 RICHARDSON Retractor 91/2" (24.1 cm), Loop Handle, 1" (2.5 cm) x 3/4" (1.9 cm) RICHARDSON Retractor 91/2" (24.1 cm), Loop Handle, 11/4" (3.2 cm) x 1" (2.5 cm) 11-247 11-248 RICHARDSON Retractor 91/2" (24.1 cm), Loop Handle, 11/2" (3.8 cm) x 11/2" (3.8 cm) 11-249 RICHARDSON Retractor 91/2" (24.1 cm), Loop Handle, 2" (5.1 cm) x 3/4" (1.9 cm) 11-250 RICHARDSON Retractor 9½" (24.1 cm), Hollow Grip Handle, 1" (2.5 cm) x ¾" (1.9 cm) 11-252 RICHARDSON Retractor 9½" (24.1 cm), Hollow Grip Handle, 1¼" (3.2 cm) x 1" (2.5 cm) RICHARDSON Retractor 91/2" (24.1 cm), Hollow Grip Handle, 11/2" (3.8 cm) x 1" (2.5 cm) 11-254 11-256 RICHARDSON Retractor 91/2" (24.1 cm), Hollow Grip Handle, 2" (5.1 cm) x 3/4" (1.9 cm)













11-260 to 11-260B

11-260

11-266 to 11-269


11-270 to 11-274

11-276

11-260A	RICHARDSON-EASTMAN Retractor, double end, 11" (27.9 cm), large size
11-260B	RICHARDSON-EASTMAN Retractor, double end, 10" (25.4 cm), small size
11-266	KELLY Retractor 9¾" (24.8 cm), Loop Handle, 2" (5.1 cm) x 1½" (3.8 cm)
11-267	KELLY Retractor 9¾" (24.8 cm), Loop Handle, 2½" (6.4 cm) x 2" (5.1 cm)
11-268	KELLY Retractor 9¾" (24.8 cm), Loop Handle, 3" (7.6 cm) x 2½" (6.4 cm)
11-269	KELLY Retractor 9¾" (24.8 cm), Loop Handle, 3½" (8.9 cm) x 2¾" (7 cm)
11-270	KELLY Retractor 91/2" (24.1 cm), Hollow Grip Handle, 2" (5.1 cm) x 11/2" (3.8 cm)
11-272	KELLY Retractor 10" (25.4 cm), Hollow Grip Handle, 21/2" (6.4 cm) x 2" (5.1 cm)
11-274	KELLY Retractor 10" (25.4 cm), Hollow Grip Handle, 3" (7.6 cm) x 2½" (6.4 cm)
11-276	KELLY Retractor 10" (25.4 cm), Hollow Grip Handle, 31/2" (8.9 cm) x 23/4" (7 cm)

RICHARDSON-EASTMAN Retractor, double end, set of two





11-290 11-292	MURPHY Gall Bladder Retractor 10½" (26.7 cm), blade 1½" (3.2 cm) wide MURPHY Gall Bladder Retractor 10½" (26.7 cm), blade 2" (5.1 cm) wide
11-300	FOSS Gall Bladder Retractor 10" (25.4 cm), blade 1%" (4.4 cm) wide
11-314 11-316 11-318	HARRINGTON Retractor 12" (30.5 cm), blade 2½" (6.4 cm) wide HARRINGTON Retractor 12" (30.5 cm), blade 1½" (3.8 cm) wide HARRINGTON Retractor 9½" (24.1 cm), blade 1" (2.5 cm) wide
11-342	FRITSCH Abdominal Retractor 91/2" (24.1 cm), Hollow Grip Handle, 3" (7.6 cm) blade




Retractors ____





25-302





 11-346
 MAYO Abdominal Retractor 10" (25.4 cm), blade 2¾" (7 cm) wide


 25-302
 ALLISON Lung Retractor 12¾" (32.4 cm), standard size 2" (5.1 cm) wide blade

 25-306
 DAVIDSON Scapula Retractor, blade 3½" (8.9 cm) wide x 3½" (8.3 cm) deep

 26-91
 TAYLOR Spinal Retractor 7½" (18.4 cm), blade 1½" (3.2 cm) x 3" (7.6 cm)

 26-92
 TAYLOR Spinal Retractor 8" (20.3 cm), blade 1½" (3.2 cm) x 4" (10.2 cm)






28-70 28-72 28-74	SAWYER Retractor 11" (27.9 cm), blade 7%" (2.2 cm) x 21%" (6.4 cm) SAWYER Retractor 11" (27.9 cm), blade 11%" (2.8 cm) x 3" (7.6 cm) SAWYER Retractor 11" (27.9 cm), blade 11%" (3.8 cm) x 31%" (8.9 cm)	
11-349	OLLIER Retractor 8½" (21.6 cm), blunt, 4 prongs 2½" (5.7 cm) x 1¾" (4.4 cm)	
11-350 11-352	ISRAEL Retractor 9" (22.9 cm), blunt, 4 prongs 11 %" (3.8 cm) x 15 %" (4.1 cm) ISRAEL Retractor 9" (22.9 cm), blunt, 6 prongs 21 %" (5.7 cm) x 15 %" (4.4 cm)	
11-360 11-362	MAYO-ADAMS Retractor 6¾" (17.1 cm), with grip lock MAYO-ADAMS Retractor 6¾" (17.1 cm), with center blade and grip lock	







26-1973 BECKMAN-EATON Laminectomy Retractor 12¾" (32.4 cm), hinged, sharp 7 x 7 teeth 1¾" (4.4 cm) wide x 2" (5.1 cm) deep


26-2000 INGE Lamina Spreader 6" (15.2 cm), jaws open to 1¾" (2.8 cm), without teeth

1NGE Lamina Spreader 6" (15.2 cm), jaws open to 1¾" (2.8 cm), without teeth

26-2005 INGE Lamina Spreader 91/2" (24.1 cm), jaws open to 11/4" (4.4 cm), smooth handles


11-365 Pediatric GELPI Retractor 5%" (14 cm), with grip lock
11-366 GELPI Retractor 6%" (17.1 cm), with grip lock
11-368 GELPI Retractor 6%" (17.1 cm), with ball stops and grip lock


11-374 RIGBY Appendectomy Retractor 6¾" (17.1 cm), with grip lock



Retractors Example 2









11-390 MASSON-JUDD Bladder Retractor with 2 swivel blades 1" (2.5 cm) x 2¾" (7 cm)
11-392 MASSON-JUDD Bladder Retractor with 2 swivel blades and center blade
1" (2.5 cm) x 2¾" (7 cm)

BARR Anal Retractor 8" (20.3 cm), self-retaining, blades 2¾" (7 cm) x ¾" (2.2 cm)

COOK Rectal Retractor 8" (20.3 cm), triple bladed, self-retaining, 4" (10.2 cm) deep blades, lower blade $\frac{1}{2}$ " (1.3 cm) wide

28-8

28-6


28-8




11-404

Retractors ____











SMITH (BUIE) Anal Retractor, self-retaining, swivel blades 21/2" (6.4 cm) x 7/6" (2.2 cm) 28-64 open to 4" (10.2 cm) SMITH (BUIE) Anal Retractor, self-retaining, extra long swivel blades 4" (10.2 cm) x 7/6" (2.2 cm) 28-65 open to 4" (10.2 cm) HILL-FERGUSON Retractor 81/2" (21.6 cm), small blade 2.2 x 6.4 cm 28-68 HILL-FERGUSON Retractor 8¾" (21 cm), large blade 3 x 7.6 cm 28-69 11-400 GOSSET Retractor 4" (10.2 cm) spread, 1" (2.5 cm) x 11/2" (3.8 cm) blades GOSSET Retractor 51/2" (14 cm) spread, 1" (2.5 cm) x 21/2" (6.4 cm) blades 11-402

Pediatric GOSSET Retractor 31/2" (8.9 cm) spread, 3/4" (1.9 cm) x 1" (2.5 cm) blades, light pattern







11-410 O'SULLIVAN-O'CONNOR Abdominal Retractor, with 3 interchangeable blades, two small size 1¾" (4.4 cm) x 3¼" (8.3 cm), one large size 2¾" (7 cm) x 4¼" (10.8 cm)


30-370 O'SULLIVAN-O'CONNOR Vaginal Retractor, with 3 interchangeable blades, ½" (1.3 cm) x 1½" (4.7 cm), 1½" (3.8 cm) x 3½" (8.9 cm) and 1¾" (4.4 cm) x 3½" (8.9 cm)







11-414



COLLIN Abdominal Retractor 9" (22.9 cm), 11-414


with 2 swivel blades

BALFOUR ABDOMINAL RETRACTORS with fenestrated side blades 2½" (6.4 cm) deep and solid center blade 3" (7.6 cm) x 1¾" (4.4 cm) deep

7" (17.8 cm) spread (standard size) 11-420


10" (25.4 cm) spread 11-421

BALFOUR Abdominal Retractor, with solid side 11-424 blades 2%" (7 cm) deep and solid center blade 3" (7.6 cm) x 1%" (4.4 cm) deep, 7" (17.8 cm) spread



Retractors Exercise







Retractors ____










11-440	BALFOUR Center Blade, standard size 3" (7.6 cm) x 1¾" (4.4 cm) deep
11-442	BALFOUR Center Blade, extra deep model, 2¾" (7 cm) x 4" (10.2 cm) deep
11-444	BALFOUR Center Blade, 2%" (6.7 cm) x 1¾" (4.4 cm) deep
11-446	BALFOUR Center Blade, 31/4" (8.3 cm) x 23/4" (7 cm) deep
11-448	BALFOUR Center Blade, 2" (5.1 cm) x 5" (12.7 cm) deep, malleable

Retractors Reserved





11-470 FRANZ Abdominal Retractor, with 8 interchangeable blades, folding for easy storage.

Blades No. 1 and 2 - 2" (5.1 cm) x 134" (4.4 cm) deep

Blades No. 3 and 4 $-\,2\%$ (6.4 cm) x 1½" (3.2 cm) deep

Blades No. 5 and 6 - 3" (7.6 cm) x 2" (5.1 cm) deep

Blades No. 7 and 8 - 31/2" (8.9 cm) x 11/4" (3.2 cm) deep

















11-302

27-230 27-232

11-504 to 11-516

11-530 to 11-544

11-500	VOLKMAN Hook 8½" (21.6 cm), sharp, hollow handle
11-502	VOLKMAN Hook 81/2" (21.6 cm), sharp, ring handle
27-230	VOLKMAN Hook Retractor 9" (22.9 cm), sharp point 15 mm deep
27-232	VOLKMAN Hook Retractor 9" (22.9 cm), blunt tip 15 mm deep

VOLKMAN RETRACTORS 81/2" (21.6 cm), Hollow Handle

Sharp		Blunt
11-504	2 prongs	11-512
11-506	3 prongs	11 - 514
11-508	4 prongs	11-51 6
VOLKMAN RET	TRACTORS 8" (20.3 cr	n), Ring Handle
Sharp		Blunt
11-530	2 prongs	11-538
11-532	3 prongs	11-540
11-534	4 prongs	11-542
11-536	6 prongs	11-544











11-580 11-582


VOLKMAN RETRACTORS 8%" (22.2 cm), USA Pattern Ring Handle

Sharp		Blunt
11-546	2 prongs	11-551
11-547	3 prongs	11-552
11-548	4 prongs	11-553
11-549	6 prongs	11-554

MURPHY RETRACTORS 7½" (19.1 cm)

		Blunt
2 prongs		11-568
3 prongs		11-570
4 prongs		11-572
6 prongs		11-574
	3 prongs 4 prongs	3 prongs 4 prongs

11-580	BENNETT Tibia Retractor 91/2" (24.1 cm), small blade 13/4" (4.4 cm) wide
11-582	BENNETT Tibia Retractor 91/2" (24.1 cm), large blade 21/2" (6.4 cm) wide
11-590	HIBBS Betractor 91/4" (23.5 cm), blade 1" (2.5 cm) x 3" (7.6 cm) deep




11-590















27-213 MINI-HOHMANN Retractor 61/4" (15.9 cm), blade 6 mm wide 27-214 MINI-HOHMANN Retractor 61/4" (15.9 cm), blade 8 mm wide 27-215 HOHMANN Retractor 91/4" (23.5 cm), blade 10 mm wide 27-216 HOHMANN Retractor 91/4" (23.5 cm), blade 18 mm wide 27-217 HOHMANN Retractor 91/4" (23.5 cm), blade 25 mm wide 27-218 HOHMANN Retractor 91/4" (23.5 cm), blade 43 mm wide 27-219 HOHMANN Retractor 93/4" (24.8 cm), blade 70 mm wide 27-223 BLOUNT Knee Retractor 7" (17.8 cm), blade 1¾" (4.4 cm) deep. %" (1.6 cm) wide tapering to 3/16" (.48 cm) 27-224 SMILLIE Meniscus Hook Retractor 53/4" (14.6 cm), blade 34" (1.9 cm) x 1%" (3.5 cm) HAYS Hand and Finger Retractor 6" (15.2 cm), angled blade 27-225 11/2" (3.8 cm) long, 3/6" (1 cm) wide tapering to 1/6" (.32 cm)




Retractors Retractors









11-600 to 11-602

11-608 to 11-610 11-614 to 11-622

11-611


WEITLANER RETRACTORS

Sharp	Size	Blunt
11-600	4" (10.2 cm) 2 x 3 teeth	11-602
11-608	51/2" (14 cm) 3 x 4 teeth	11-610
	5½" (14 cm) solid blades	11-611
11-614	61/2" (16.5 cm) 3 x 4 teeth	11-616
11-617	8" (20.3 cm) 3 x 4 teeth	11-618
11-620	9½" (24.1 cm) 3 x 4 teeth	11-622

BECKMAN-WEITLANER RETRACTORS with Hinged Blades

	with ranged blades	
Sharp	Size	Blunt
11-630-SH	5½" (14 cm) 3 x 4 teeth	11-630-BL
11-634-SH	6½" (16.5 cm) 3 x 4 teeth	11-634-BL
11-642-SH	9" (22.9 cm) 3 x 4 teeth	11-642-BL


40-1235 SCHINK Retractor 41/2" (11.4 cm)


















MEYERDING SELF-RETAINING RETRACTORS

11-650	1" (2.5 cm) x 2%" (6.1 cm) blades
11-652	1" (2.5 cm) x 21/6" (5.4 cm) blades
11-654	1" (2.5 cm) x 13/4" (4.4 cm) bladeş
11-656	3/4" (1.9 cm) x 23/6" (6.1 cm) blades
40.055	DOLLHIJAROUT D I. D


19-655	SCHUHKNECHT Postauricular Retractor 5" (12.7 cm),
	3 x 3 teeth, 3" (7.6 cm) opening


19-660	WULLSTEIN Endaural Retractor 5" (12.7 cm), 3 x 3 sharp teeth, ebony non-glare finish
19-668	PERKINS Endaural Retractor 5¼" (13.3 cm), left serrated solid blade, right side 3 teeth

19-669	PERKINS Endaural Retractor 51/4" (13.3 cm),
	right serrated solid blade, left side 3 teeth










11-660 FRAZIER Retractor, 3 x 4 blunt teeth, 31/4" (8.3 cm) spread

11-672 BADGLEY Retractor, 4 x 4 blunt teeth, 2" (5.1 cm) x 1¾" (4.4 cm), 6½" (16.5 cm) spread

11-676 TUFFIER Retractor, 2" (5.1 cm) x 1¾" (4.4 cm), grooved blades, 6½" (16.5 cm) spread






Retractors Extended



30-1355








BECKMAN Retractor 12½" (31.8 cm), hinged, sharp 4 x 4 teeth %" (1.9 cm) x 1" (2.5 cm) 11-684

BECKMAN Retractor 12½" (31.8 cm), hinged, sharp 4 x 4 teeth 1¾" (4.4 cm) x 1¼" (3.2 cm) 11-686


11-700 ADSON Retractor 7½" (19.1 cm), sharp, 4 x 4 teeth ½" (1.9 cm) x ¾" (1.9 cm), straight arms

ADSON Retractor $7\frac{1}{2}$ " (19.1 cm), sharp, 4 x 4 teeth $\frac{3}{4}$ " (1.9 cm) x $\frac{3}{4}$ " (1.9 cm), angled arms 11-702

> Lateral Vaginal Retractor 814" (21 cm), with extra long thumb ratchet retracts up to a full 214" (5.7 cm) for an extra wide viewing area







11-82	MEYERDING Skin Hook and Retractor 61/4" (15.9 cm), double end
11-500	VOLKMAN Hook 8½" (21.6 cm), sharp, hollow handle
11-502	VOLKMAN Hook 8½" (21.6 cm), sharp, ring handle
27-230	VOLKMAN Hook Retractor 9" (22.9 cm), sharp point 15 mm deep
27-232	VOLKMAN Hook Retractor 9" (22.9 cm), blunt tip 15 mm deep
18-412	TYRELL Hook 5" (12.7 cm), sharp
18-414	TYRELL Hook 5" (12.7 cm), blunt
18-416	GUTHRIE Fixation Hook 5" (12.7 cm), sharp double pronos 1.5 mm wide








WIENER Hook 5" (12.7 cm), sharp 18-418 18-420 WIENER Hook 5" (12.7 cm), blunt 18-424 O'CONNOR Hook 43/4" (12.1 cm), sharp 18-428 NUGENT Hook 43/4" (12.1 cm), right angle, sharp 18-440 STEVENS Hook 5" (12.7 cm), right angle 5 mm deep 18-442 GREEN Hook 5" (12.7 cm), 9 mm deep 18-446 JAMESON Hook 5" (12.7 cm), 10 mm deep 18-450 VON GRAEFE Hook 51/2" (14 cm), small size 8 mm deep 18-452 VON GRAEFE Hook 51/2" (14 cm), medium size 10 mm deep VON GRAEFE Hook 51/2" (14 cm), large size 11 mm deep 18-454





19-220 19-222	DAY Hook 6½" (16.5 cm), small size DAY Hook 6½" (16.5 cm), large size	
19-226	ALLPORT Hook 61/2" (16.5 cm), with ball tip	
19-228	GROSS Hook and Spoon 5" (12.7 cm)	S
21-85 21-86 21-87	KLEINERT-KUTZ Hook 5" (12.7 cm), small size 3 mm dia. KLEINERT-KUTZ Hook 5" (12.7 cm), medium size 5 mm dia. KLEINERT-KUTZ Hook 5" (12.7 cm), large size 7 mm dia.	वार के दुरावार है।
21-88 21-89	GILLIES (CONVERSE) Skin Hook 7" (17.8 cm), large size GILLIES (CONVERSE) Skin Hook 7" (17.8 cm), small size	
21-90 21-92 21-94	COTTLE Skin Hook 5½" (14 cm), small, shallow curve COTTLE Skin Hook 5½" (14 cm), small, deep curve COTTLE Skin Hook 5½" (14 cm), large, deep curve	26-1084 26-1086
26-1084 26-1086	FRAZIER Dura Hook 5" (12.7 cm), sharp point FRAZIER Dura Hook 5" (12.7 cm), blunt point	



Hooks 🔛



21-158

21-160












21-102 FREER Skin Hook 6" (15.2 cm), sharp double prongs 2.5 mm wide 21-103 FREER Skin Hook 6" (15.2 cm), blunt double prongs 2.5 mm wide 21-107 MILLARD Thimble Hook Retractor, %" (1.6 cm) working length, sharp single hook COTTLE Thumb Hook Retractor, 11/2" (3.8 cm) working length, sharp 21-108 double hook 8 mm wide COTTLE Double Hook 51/2" (14 cm), two sharp prongs 12 mm wide 21-144 21-152 MILTEX Skin Hook 43/4" (12.1 cm), sharp prong 2 mm dia. 21-153 JOSEPH Single Hook 61/4" (15.9 cm), one sharp prong JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 2 mm 21-154 21-156 JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 5 mm

JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 7 mm JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 10 mm



21-154 to 21-160





21-430	MILTEX Tissue Twist Hook 5½" (14 cm)
23-1048	HUPP Hook 61/2" (16.5 cm), sharp point
23-1050	JACKSON Tenaculum Hook 5%" (14.6 cm), sharp point
23-1051	NEWS Hook 6" (15.2 cm), sharp point
26-1080	CUSHING Dura Hook 61/2" (16.5 cm), blunt
26-1090 26-1092	ADSON Dura & Dissecting Hook 8" (20.3 cm), sharp blade 5 mm at right angle ADSON Dura Hook 8" (20.3 cm), blunt 5 mm hook at right angle



Hooks Example












26-1150 DANDY Hook 9" (22.9 cm), straight GRAHAM Hook 7" (17.8 cm), blunt, 6 mm deep 26-1156 SMITHWICK Hook 7" (17.8 cm), blunt, 10 mm deep 26-1170 SMITHWICK Hook and Dissector 12" (30.5 cm) 26-1174 .28-120 PRATT Crypt Hook 81/2" (21.6 cm), narrow style 28-122 PRATT Crypt Hook 81/2" (21.6 cm), short style PRATT Crypt Hook 81/2" (21.6 cm), long style 28-124 STEWART Crypt Hook 81/2" (21.6 cm), short style 28-128 28-130 STEWART Crypt Hook 81/2" (21.6 cm), long style



28-131	ROSSER Crypt Hook 6¾" (17.1 cm)			
28-140	BARR Crypt Hook 10" (25.4 cm), short and long hook ends			
30-940 30-950	MILTEX I.U.D. Extractor Hook MILTEX Tenaculum Hook 8" (2 handle graduated in cm	, ,	delicate pointed hook,	
	EMMETT TENACULUM HOO	KS 9" (22.9 cm), Sharp Points	
30-951	Style 1 Half-Curved, small	30-954	Style 4 Half-Curved, large	
30-952 30-953	Style 2 Acute Angle Style 3 Right Angle	30-955	Style 5 Double Angle	
33-121	SELTZER Dermal Biopsy Hoo	k 4½" (11.4 cm	ı), sharp point	(Contract)
34-250	Dissecting Tenaculum Hook 6	½" (16.5 cm), s	harp point	34-250



Thyroid-Goiter Instruments





11-20









11-81

11-20	LAHEY Thyroid Retractor 7%" (19.7 cm)
44.00	ODEEN The 14 Debender 01/11/04 0

GREEN Thyroid Retractor 81/2" (21.6 cm) 11-22

SENN Retractor 6%" (16.2 cm), sharp, double end 11-74 11-76 SENN Retractor 6%" (16.2 cm), blunt, double end


CRILE Thyroid Retractor 4%" (11.4 cm), double end 3 x 11 mm and 6 x 20 mm blades 11-78

11-80 MATHIEU Retractor 61/2" (16.5 cm), sharp, double end 11-81 MATHIEU Retractor 61/2" (16.5 cm), blunt, double end

11-642-SH BECKMAN-WEITLANER Retractor 9" (22.9 cm),

hinged blades, sharp 3 x 4 teeth 11-642-BL BECKMAN-WEITLANER Retractor 9" (22.9 cm),

hinged blades, blunt 3 x 4 teeth



11-642-BL

Thyroid-Goiter Instruments













11-684	BECKMAN Thyroid Retractor 121/2" (31.8 cm), hinged, sharp 4 x 4 teeth 3/4" (1.9 cm) x 1" (2.5 cm)
12-80	LAHEY Traction Forceps 61/4" (15.9 cm), 3 x 3 teeth
12-90 12-92	LAHEY Ligature Carrier 71/2" (19.1 cm), curved left LAHEY Ligature Carrier 71/2" (19.1 cm), curved right
5-179 5-180	LAHEY Goiter Scissors 5½" (14 cm), straight LAHEY Goiter Scissors 5½" (14 cm), curved
7-30 7-32	PROVIDENCE HOSPITAL Hemostatic Forceps 5½" (14cm), straight PROVIDENCE HOSPITAL Hemostatic Forceps 5½" (14cm), curved





Thyroid-Goiter Instruments






7,60

7-236

12-170





7-62	LAHEY Hemostatic Forceps 5½" (14 cm), curved
7-98	LOVELACE Hemostatic Forceps 6% " (15.9 cm), straight, cross serrations, 1 x 2 teeth
7-230	BAINBRIDGE Hemostatic Forceps 6" (15.2 cm), straight, longitudinal serrations with cross tip serrations
7-232	BAINBRIDGE Hemostatic Forceps 6" (15.2 cm), curved, longitudinal serrations with cross tip serrations
7-234	BAINBRIDGE Hemostatic Forceps 7¼" (18.4 cm), straight, longitudinal serrations with cross tip serrations

and one serrated jaw, with four pins

BAINBRIDGE Hemostatic Forceps 7% " (18.4 cm), curved, longitudinal serrations with cross tip serrations

MASTIN Muscle Clamp 7" (17.8 cm), shanks curved to side, one smooth


LAHEV Hemostatic Forcers 51/6" (14 cm), straight

12-170

197

Trocars and Biopsy Needles





13-2 NESTED TROCAR SET of 4, sizes 6 French (2 mm), 9 French (3 mm), 13 French (4.3 mm) and 17 French (5.7 mm), in flat chrome plated handle. Trocar points and cannula are stainless.

PATTERSON TROCARS 81/2" (21.6 cm). Trocar points, Cannula and handle are stainless

13-8	size 1	6 mm dia.
13-9	size 2	7.5 mm dia.
13-10	sixe 3	10.5 mm dia.
13-11	size 4	12 mm dia.

13-20 DUKE Trocar and Cannula 4½" (11.4 cm), with perforated tube, 5½" (13 cm) long, size 17 French (5.7 mm), all stainless

OCHSNER TROCARS, Original Model 6" (15.2 cm), with tubing connection

	0 (101 01.11/1 111111		
3-40	12 French (4 mm)	13-46	20 French (6.7 mm)
3-42	14 French (4.7 mm)	13-48	22 French (7.3 mm)
3-44	17 French (5.7 mm)	13-50	24 French (8 mm)

OCHSNER TROCARS, Improved Philadelphia Pattern,

		61	15.9 "4"	cm) for in	troduc	tion o	f catheters
1	3-60	12	French	(4 mm)	for 8	Fr.	(2.7 mm) catheters
1	3-62	14	French	(4.7 mm)	for 10	Fr.	(3.3 mm) catheters
1	3-64	16	French	(5.3 mm)	for 12	Fr.	(4 mm) catheters
1	3-66	18	French	(6 mm)	for 14	Fr.	(4.7 mm) catheters
- 13	3-68	20	French	(6.7 mm)	for 16	Fr.	(5.3 mm) catheters
							(6 mm) catheters
							(6.7 mm) catheters
1	3-74	28	French	(9.3 mm)	for 22	Fr.	(7.3 mm) catheters


MENGHINI BIOPSY NEEDLES,

Improved Model, with built-in stop plug that never requires replacement

For Liver Biopsy Child Size – 40 mm long			For Liver Biopsy Adult Size – 70 mm long		Fo	For Kidney Biopsy 120 mm long		
13-120 13-122 13-124 13-126	1 1.2 1.4	side diameter mm – 19 gauge mm – 18 gauge mm – 17 gauge mm – 16 gauge	13-128 13-130 13-132 13-134	1 1.2 1.4	nside diameter mm – 19 gauge mm – 18 gauge mm – 17 gauge mm – 16 gauge	13-150 13-152 13-154 13-156	1 1.2 1.4	nside diameter mm – 19 gauge mm – 18 gauge mm – 17 gauge mm – 16 gauge







14-2 OCHSNER (FENGER) Gall Duct Probe, cannulated flexible spiral 14" (35 cm) with stylet

MAYO Common Duct Probe 10" (25.4 cm), malleable shaft, 15 French (5 mm) 14-4-15 14-4-18 MAYO Common Duct Probe 10" (25.4 cm), malleable shaft, 18 French (6 mm)


BAKES Common Duct Dilator Set; fitted plastic case containing 9 dilators 14-8 sizes 3 to 11 mm, malleable shafts

BAKES COMMON DUCT DILATORS 8¾" (22.2 cm), single sizes, malleable shafts

14-10-3	3 mm dia.	14-10-8	8 mm dia.
14-10-4	4 mm dia.	14-10-9	9 mm dia.
14-10-5	5 mm dia.	14-10-10	10 mm dia.
14-10-6	6 mm dia.	14-10-11	11 mm dia.
14-10-7	7 mm dia		

14-12 MIXTER Irrigating Dilaprobe 12" (30.5 cm), malleable, with Luer Lock, silver plated






14-22	MAYO Gall Stone Scoop 11" (27.9 cm), double end
14-24-17 14-24-21 14-24-24 14-24-27	DESJARDIN Gall Stone Scoop 11" (27.9 cm), malleable shaft, 17 French (5.7 mm) DESJARDIN Gall Stone Scoop 11" (27.9 cm), malleable shaft, 21 French (7 mm) DESJARDIN Gall Stone Scoop 11" (27.9 cm), malleable shaft, 24 French (8 mm) DESJARDIN Gall Stone Scoop 11" (27.9 cm), malleable shaft, 27 French (9 mm)
14-26-L	FERGUSON Gall Stone Scoop 9½" (24.1 cm), double end, large size scoops ½" (1.3 cm) x 1½" (3.2 cm) and $\frac{1}{2}$ " (1.1 cm) x 1" (2.5 cm)
14-26-M	FERGUSON Gall Stone Scoop 9½" (24.1 cm), double end, medium size scoops ½" (1.3 cm) x 1½" (2.8 cm) and $\%$ " (.95 cm) x 1" (2.5 cm)
14-26-S	FERGUSON Gall Stone Scoop 9½" (24.1 cm), double end, small size scoops %" (.95 cm) x 1" (2.5 cm) and %6" (.79 cm) x ¾" (1.9 cm)
14-27 14-28 14-29	MAYO Common Duct Scoop 10½" (26.7 cm), malleable shaft, small size, silver plated MAYO Common Duct Scoop 10½" (26.7 cm), malleable shaft, medium size, silver plated MAYO Common Duct Scoop 10½" (26.7 cm), malleable shaft, large size, silver plated










7-200 7-202	MIXTER Gall Duct Forceps 6¼" (15.9 cm) MIXTER Gall Duct Forceps 7¼" (18.4 cm)
14-46	LOWER Gall Duct Forceps 7" (17.8 cm)
14-48 14-49	LAHEY Gall Duct Forceps 7½" (19.1 cm) LAHEY Gall Duct Forceps 9" (22.9 cm)
14-50	MOYNIHAN Gall Duct Forceps 81/2" (21.6 cm)
14-51	GREEN Cystic Duct Forceps 8½" (21.6 cm), delicate













14-52 JOHNS HOPKINS Gall Bladder Forceps 8" (20.3 cm)
 7-206 MIXTER Gall Bladder Forceps 9" (22.9 cm)
 14-56 GRAY Cystic Duct Forceps 9" (22.9 cm), serrated, set of two








7-255 GEMINI-MIXTER Forceps 51/2" (14 cm), delicate 7-256 GEMINI-MIXTER Forceps 7" (17.8 cm), delicate 7-257 GEMINI-MIXTER Forceps 8" (20.3 cm), delicate 7-258 GEMINI-MIXTER Forceps 9" (22.9 cm), delicate 7-259 GEMINI-MIXTER Forceps 11" (27.9 cm), delicate 14-62 SHALLCROSS Cystic Duct Forceps 7" (17.8 cm) 14-66 LOVELACE Gall Bladder Forceps 71/4" (18.4 cm), angular, 1" (2.5 cm) wide jaws





2-50	COLLIN (JUDD-DeMARTEL) Gall Bladder Forceps 6" (15.2 cm)
14-69	MALIK Cystic Duct Catheter Clamp for Cholangiography 8% (21.6 cm), jaws are designed to hold different sizes of catheters.
11-192	JUNG Mini-Chole Retractor 101/2" (26. cm), 3 cm width blade
11-194	JUNG Mini-Chole Retractor 9" (22.9 cm), 4 cm width blade
11-196	JUNG Mini-Chole Retractors, set of 2 each of 3 cm width blade and 1 each of 4 cm width blade

For a full line of Retractors see pages 151 - 187














BLAKE Gall Stone Forceps 81/2" (21.6 cm), straight 14-84 14-86 BLAKE Gall Stone Forceps 81/2" (21.6 cm), curved DESJARDIN (ROCHESTER) Gall Stone Forceps 9" (22.9 cm) 14-88 MIXTER Gall Stone Forceps 9" (22.9 cm), length from tip to angle 3" (7.6 cm) 14-90 THOREK Gall Bladder Scissors 71/4" (18.4 cm), full curve 5-244 THOREK Gall Bladder Scissors 10" (25.4 cm), full curve 5-245











16-2	LOCKWOOD-ALLIS Tissue Forceps 61/4" (15.9 cm)
16-3 16-4	Baby ALLIS Tissue Forceps 5" (12.7 cm), 4 x 5 teeth, extra delicate jaws 2.5 n Baby ALLIS Tissue Forceps 5% " (14 cm), 4 x 5 teeth, delicate jaws 4 mm wide
16-6	ALLIS Tissue Forceps 6" (15.2 cm), 4 x 5 teeth

 16-6
 ALLIS Tissue Forceps
 6" (15.2 cm), 4 x 5 teeth

 16-8
 ALLIS Tissue Forceps
 6" (15.2 cm), 5 x 6 teeth

 16-10
 ALLIS Tissue Forceps
 7¼" (18.4 cm), 5 x 6 teeth


 16-12
 ALLIS Tissue Forceps
 10" (25.4 cm), 5 x 6 teeth

 16-20
 ALLIS Tissue Forceps
 6" (15.2 cm), 4 x 5 teeth, angular jaws

 16-22
 ALLIS Tissue Forceps
 6" (15.2 cm), 5 x 6 teeth, angular jaws

16-24 ALLIS Tissue Forceps 7½" (18.4 cm), 5 x 6 teeth, angular jaws
 16-25 ALLIS-NC Tissue Forceps 6" (15.2 cm), with double row of non-traumatic teeth
 16-26 ALLIS-NC Tissue Forceps 7½" (18.4 cm), with double row of non-traumatic teeth

16-30 ALLIS-ADAIR Tissue Forceps 61/2" (16.5 cm), 9 x 10 teeth

















16-31	BOYS-ALLIS Tissue Forceps 6" (15.2 cm), 5 x 6 teeth, heavy pattern
16-33 16-34	JUDD-ALLIS Tissue Forceps 6" (15.2 cm), 3 x 4 teeth, delicate JUDD-ALLIS Tissue Forceps $71/2$ " (19.1 cm), 3 x 4 teeth, delicate
16-38	THOMS-ALLIS Tissue Forceps 8" (20.3 cm), 6 x 7 teeth
16-40	LOCKWOOD Tissue Forceps 73/4" (19.7 cm)
16-42	Baby BABCOCK Intestinal Forceps 51/2" (14 cm), delicate jaws 6 mm wide
16-44	BABCOCK Intestinal Forceps 61/4" (15.9 cm), jaws 9 mm wide
16-46	BABCOCK Intestinal Forceps 81/4" (21 cm), jaws 10 mm wide
16-48	BABCOCK Intestinal Forceps 91/2" (24.1 cm), jaws 15 mm wide













7-106 PRATT T-Clamp 6" (15.2 cm)

16-51 ADAIR Tissue Holding Tenaculum Forceps 6½" (16.5 cm)

16-52 COLLIN Intestinal Forceps 8" (20.3 cm), with ½" (1.3 cm) wide jaws

16-54 COLLIN Intestinal Forceps 8" (20.3 cm), with 1" (2.5 cm) wide jaws

16-56 WILLIAMS Intestinal Forceps 6½" (16.5 cm), cross serrated jaws

16-62 PERCY Intestinal Forceps 6" (15.2 cm), serrated jaws









6-214

7-236









6-214 SINGLEY Tissue Forceps 9" (22.9 cm), with fenestrated serrated jaws 16-70 CHILD Intestinal Forceps 93/4" (24.8 cm), with interchangeable rubber laws Replacement Rubber Jaws for No. 16-70 forceps 2-62 COLLIN-DUVAL Tissue Forceps 8" (20.3 cm), with 1/2" (1.3 cm) wide jaws 16-72 16-80 WANGENSTEEN Tissue Forceps 9" (22.9 cm), rounded jaws with fine cross serrations BAINBRIDGE Intestinal Forceps 6" (15.2 cm), straight, longitudinal serrations with cross tip serrations 7-230 7-232 BAINBRIDGE Intestinal Forceps 6" (15.2 cm), curved, longitudinal serrations with cross tip serrations 7-234 BAINBRIDGE Intestinal Forceps 71/4" (18.4 cm), straight, longitudinal serrations with cross tip serrations

BAINBRIDGE Intestinal Forceps 71/4" (18.4 cm), curved, longitudinal serrations with cross tip serrations









16-99	DENNIS Anastomosis Clamp 9" (22.9 cm), longitudinal serrations
16-100	ALLEN Intestinal Clamp 8" (20.3 cm), longitudinal serrations with 1 x 2 teeth
16-102	Baby ALLEN Intestinal Clamp 61/2" (16.5 cm), delicate, longitudinal serrations with 1 x 2 teeth
16-108 16-110	Baby DOYEN Intestinal Forceps 6½" (16.5 cm), straight, flexible blades with diagonal serrations Baby DOYEN Intestinal Forceps 6½" (16.5 cm), curved, flexible blades with diagonal serrations













16-120 16-122

16-132




KOCHER Intestinal Forceps 101/2" (26.7 cm), straight, longitudinal serrations 16-120 16-122 KOCHER Intestinal Forceps 101/2" (26.7 cm), curved, longitudinal serrations

16-130 MAYO-ROBSON Intestinal Forceps 10" (25.4 cm), straight, Iongitudinal serrations


16-132 MAYO-ROBSON Intestinal Forceps 10" (25.4 cm), curved, longitudinal serrations

SCUDDER Intestinal Forceps 91/2" (24.1 cm), straight, smooth solid blades 16-140

30-1805 PEAN Intestinal Forceps 9" (22.9 cm), curved, with longitudinal serrations 30-1815 PEAN Intestinal Forceps 101/4" (26 cm), curved, with longitudinal serrations








16-170	FOSS Anterior Resection Clamp 111/2" (29.2 cm), S-shape with smooth jaws
16-172	BRUNNER Intestinal Forceps 91/2" (24.1 cm), longitudinal serrations
16-174	FEHLAND Intestinal Clamp 9% " (24.8 cm), right angle with 3% " (8.3 cm) long jaws
16-176	FOSS Intestinal Clamp 11¾" (29.8 cm), right angle with heavy jaws










STONE INTESTINAL CLAMPS

4" (10.2 cm)
ļ





16-229	Baby PAYR Pylorus Clamp 6" (15.2 cm)
16-230 16-232 16-234	PAYR Pylorus Clamp 8" (20.3 cm) PAYR Pylorus Clamp 11" (27.9 cm) PAYR Pylorus Clamp 13¾" (34.9 cm)
16-250	BEST Colon Clamp Set of 4, including 1 clamp 7.5 cm, 1 clamp 6 cm and 2 clamps 3.5 cm
16-252 16-254 16-256	BEST Colon Clamp 10½" (26.7 cm), with 3.5 cm jaws BEST Colon Clamp 10½" (26.7 cm), with 6 cm jaws BEST Colon Clamp 10½" (26.7 cm), with 7.5 cm jaws







17-100

17-115 17-215

17-500

PENCIL-TYPE MICRO INSTRUMENTS 71/4" (18.4 cm) overall length; extremely delicate working ends

Forceps, 1 mm wide serrated jaws, straight shaft 17-100

17-115 Scissors, 4 mm long blades, straight shaft 17-215 Scissors, 4 mm long blades, angled shaft

MILTEX Micro Dissecting Scissors 5%" (14.6 cm), double action, sharp points, spring handle 17-500











17-2150 17-2100

17-2160 17-2110

18-1558 18-1559

MILTEX MICRO SURGERY SCISSORS, sharp points

17-2060 17-2040	6" (15.2 cm) straight, 8 mm blades, flat handles
	6" (15.2 cm) curved, 8 mm blades, flat handles
17-2150	6" (15.2 cm) straight, 6 mm blades, round handles
17-2100	6" (15.2 cm) curved, 6 mm blades, round handles
17-2160	71/4" (18.1 cm) straight, 6 mm blades, round handles
17-2110	71/6" (18.1 cm) curved, 6 mm blades, round handles
18-1558	CASTROVIEJO Corneal Scissors 4" (10.2 cm), right, sub-miniature for micro surgery, 6 mm blades
18-1559	CASTROVIEJO Corneal Scissors 4" (10.2 cm), left, sub-miniature for micro surgery, 6 mm blades













18-1618 18-1619

18-1626

18-1631

18-1633 18-1634

19-2150-B

18-1618 Micro Iris Scissors 4" (10.2 cm), straight, sharp points 18-1619 Micro Iris Scissors 4" (10.2 cm), curved, sharp points 18-1625 McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), straight, sharp points McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), curved, sharp points 18-1626 Micro-VANNAS Scissors 31/4" (8.3 cm), straight, ultrafine blades 18-1630 18-1631 Micro-VANNAS Scissors 31/4" (8.3 cm), curved, ultrafine blades 18-1633 Micro Scissors 31/4" (8.3 cm), straight, superfine blades

18-1634 19-2150-B

Micro Scissors 31/4" (8.3 cm), curved, superfine blades BELLUCCI Micro Ear Scissors 2¾" (7 cm) shaft, 4 mm blades, ebony finish







17-302









SWISS JEWELER STYLE FORCEPS Non-Magnetic Stainless Steel

17-303

17-301	Style 1, fine 43/4" (12.1 cm)
17-302	Style 2, strong, fine, 4¾" (12.1 cm)
17-303	Style 3, narrow, fine, 43/4" (12.1 cm)
17-303C	Style 3C, narrow, fine, 4%" (11.2 cm
17-303X	Style 3, micro-fine, 4¾" (12.1 cm)
17-304	Style 4, fine, 4%" (11.2 cm)
17-305	Style 5, super fine, 4%" (11.2 cm)











17-306






SWISS JEWELER STYLE FORCEPS Non-Magnetic Stainless Steel


17-305X 17-306 17-307 17-307X	Style 5, micro-fine, 4%" (11.2 cm) Style 6, pickup, 4½" (11.4 cm) Style 7, fine curved, 4½" (11.4 cm) Style 7, micro-fine, curved. 4½" (11.4 cm)
17-1100 17-1150	RHOTON Micro Suture Forceps 7" (17.8 cm), straight, with 0.7 mm tying platform RHOTON Micro Forceps 7" (17.8 cm), straight,

RHOTON Micro Forceps 7" (17.8 cm), straight, 0.3 mm tips



219






17-2510	Micro ADSON Trissue Forceps 4% (12.1 cm), jaws 0.5 mm wide, 1 x 2 teem Micro ADSON Dressing Forceps 4% (12.1 cm), serrated jaws 0.5 mm wide
18-835 18-836	BONN Micro Suturing Forceps 294" (7 cm), 1 x 2 teeth 0.12 mm, with tying platform BONN Micro Suturing Forceps 334" (9.5 cm), 1 x 2 teeth 0.12 mm, with tying platform
18-837	McPHERSON Micro Suturing Forceps $3\frac{1}{2}$ " (8.9 cm), straight, 1 x 2 teeth 0.2 mm, with tying platform
18-838	McPHERSON Micro Suturing Forceps 3½" (8.9 cm), angled, 1 x 2 teeth 0.2 mm, with tying platform
18-948	McPHERSON Micro Iris Suturing Forceps 3½" (8.9 cm), straight, smooth platform 0.3 x 5 mm, no teeth
18-949	McPHERSON Micro Iris Suturing Forceps 3½" (8.9 cm), angled, smooth platform 0.3 x 5 mm, no teeth












18-867	BISHOP-HARMON Micro Tissue Forceps 3%" (8.6 cm), delicate 0.3 mm 1 x 2 teeth
18-953	CASTROVIEJO Micro Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.12 mm, with tying platform
19-2081-B 19-2100 19-2100-B 19-2090-B 19-2110 19-2110-B	Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, serrated jaws 0.6 mm wide, ebony finish Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, serrated jaws 0.8 mm wide Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, serrated jaws 0.8 mm wide, ebony finish Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, oval cup jaws 0.5 mm wide, ebony finish Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, oval cup jaws 0.8 mm wide Micro Alligator Ear Forceps 3¼" (8.3 cm) shaft, oval cup jaws 0.8 mm wide, ebony finish
19-2138	HOUSE Micro Strut Forceps 2¾" (7 cm) shaft, smooth jaws 0.8 mm wide
19-2150-B	BELLUCCI Micro Ear Scissors, 29/4" (7 cm) shaft, 4 mm blades, ebony finish





MILTEX MICRO SURGERY NEEDLE HOLDERS with round handle, 0.6 mm tips

17-1000	5¼" (13.3 cm) curved jaws
17-980	71/4" (18.1 cm) straight jaws
17-1010	7½" (18.1 cm) curved jaws
17-1020	51/4" (13.3 cm) curved jaws, with lock
17-1024	71/2" (18.1 cm) straight jaws, with lock
17-1025	71/e" (18.1 cm) curved jaws, with lock
17-1090TC	JACOBSON Microvascular Needle Holder 7" (17.8 cm), straight, with lock, with CARB-N-SERT jaws
17-1091TC	JACOBSON Microvascular Needle Holder 7" (17.8 cm), curved, with lock, with CARB-N-SERT jaws
17-1092TC	JACOBSON Microvascular Needle Holder 8½" (21.6 cm), straight, with lock, with CARB-N-SERT jaws




17-2600

17-2602

extremely delicate

Micro Surgery Instruments






	extremely delicate			// 1	
17-2603	JACOBSON Hemostatic Forceps 7" (17.8 cm), curved, extremely delicate				
	MICRO VES	SSEL CLIPS			
17-4000	straight 0.75 x 4 mm jaws	17-4050	slightly curved 1 x 6 mm jaws		
17-4005	straight 0.8 x 6 mm jaws	17-4060	fully curved 1 x 6 mm jaws	N H	
17-4010	straight 1 x 8 mm jaws	17-4070	45° angle 1 x 6 mm jaws		
17-4015	straight 1.5 x 10 mm jaws				I.
17-4030	Micro Clip Approximator Set,	2 x 6 mm jaw	s, 15 mm spread	V	
17-4100	Applying Forceps for Micro C	lips 5¾" (14.6	6 cm), with lock	17-4100	17-4102
17-4102	Applying Forceps for Micro C	lips 5½" (14 d	cm)		

JACOBSON Micro Mosquito Forceps 5" (12.7 cm), straight,


JACOBSON Micro Mosquito Forceps 5" (12.7 cm), curved,






7-310 7-315	DIETHRICH Micro Bulldog Clamp 11/4" (4.7 cm), straight jaws 8 mm x 1.2 mm DIETHRICH Micro Bulldog Clamp 13/4" (4.4 cm), angled jaws 8 mm x 1.2 mm
26-1040	SCHWARTZ Clip Applying Forceps 8" (20.3 cm)
26-1041 26-1042 26-1043	SCHWARTZ Vessel Clip 1" (2.5 cm), straight, jaws 1.7 x 8 mm SCHWARTZ Vessel Clip 1" (2.5 cm), slightly angled, jaws 1.7 x 8 mm SCHWARTZ Vessel Clip 1" (2.5 cm), strong angle, jaws 1.7 x 8 mm
17-4800 17-4801 17D-4802	Micro FRIEDMAN Rongeur 5½" (14 cm), straight 1.3 mm wide jaws, very delicate Micro FRIEDMAN Rongeur 5½" (14 cm), curved 1.3 mm wide jaws, very delicate Micro FRIEDMAN Rongeur 5½" (14 cm), 90° angled 1.3 mm wide jaws, very delicate






18-2 LANCASTER Eye Speculum 31/4" (8.3 cm), blades 15 x 5 mm 18-4 KNAPP Eye Speculum 31/4" (8.3 cm), blades 15 x 6 mm CASTROVIEJO Eye Speculum 4" (10.2 cm), large blades 18 x 5 mm O.D. 18-10 CASTROVIEJO Eye Speculum 3¾" (9.5 cm), medium blades 15 x 5 mm O.D. 18-12 CASTROVIEJO Eye Speculum 3¾" (9.5 cm), small blades 13 x 5 mm O.D. 18-14 18-20 LANGE Eye Speculum 35/8" (9.2 cm), blades 20 x 5 mm WILLIAMS Eye Speculum 31/4" (8.3 cm), large blades 15 x 5 mm 18-22 18-24 WILLIAMS Eye Speculum 23/4" (7 cm), small blades 11 x 5 mm










18-30	GUYTON-PARK Eye Speculum 3½" (8.9 cm), solid blades 14 mm wide, with suture posts
18-32	GUYTON-PARK Eye Speculum 3½" (8.9 cm), fenestrated blades 14 mm wide, with suture posts
18-33	MAUMENEE-PARK Eye Speculum 3½" (8.3 cm), solid blades 14 mm wide, with lock screws
18-34	Plain Wire Eye Speculum 21/2" (6.4 cm), blades 15 mm wide







18-35 18-36	BARRAQUER Eye Speculum 13/4" (4.4 cm), large blades 15 mm wide BARRAQUER Eye Speculum 13/4" (4.4 cm), small blades 10 mm wide
18-46	SMITH Eye Speculum 3" (7.6 cm), blades 15 mm wide, with concealed spring, chrome
18-47	SAUER Eye Speculum 11/4" (3.2 cm), infant size, 10 mm blades
18-50	MURDOCK Eye Speculum 21/4" (5.7 cm), blades 15 mm wide, chrome
18-54	WIENER Eye Speculum 2½" (5.7 cm), blades 15 mm, with extensions to cover eye lashes, chrome
18-62	McKINNEY Eye Speculum 21/6" (5.4 cm), blades 15 mm wide, chrome
18-63 18-64	COOK Eye Speculum 11%" (4.7 cm), infant size, 9 mm blades COOK Eye Speculum 11%" (4.7 cm), child size, 11 mm blades












18-66 18-68	MELLINGER Eye Speculum 2¾" (7 cm), solid blades 15 mm wide MELLINGER Eye Speculum 2¾" (7 cm), fenestrated blades 15 mm wide
18-80	AGRICOLA Lacrimal Sac Retractor 11/2" (3.8 cm), 3 x 3 pointed prongs 4 mm deep
18-84	MELLER Lacrimal Sac Retractor 27/8" (7.3 cm), 3 x 3 prongs 3 mm deep
18-88	GOLDSTEIN Lacrimal Sac Retractor 1½" (2.8 cm), 3 x 3 pointed prongs 5 mm deep
18-96	STEVENSON Lacrimal Sac Retractor 31/4" (8.3 cm), 3 x 3 prongs 3 mm deep










18-100 18-102	KNAPP (ROLLET) Retractor 51/4" (13.3 cm), with 4 sharp prongs KNAPP Retractor 51/4" (13.3 cm), with 4 blunt prongs
11D-27 11-30 11-32 11-34	KNAPP Lacrimal Sac Retractor 5" (12.7 cm), flexible shaft, 3 sharp prongs KNAPP Lacrimal Sac Retractor 6" (15.2 cm), flexible shaft, 1 sharp prong KNAPP Lacrimal Sac Retractor 6" (15.2 cm), flexible shaft, 2 sharp prongs KNAPP Lacrimal Sac Retractor 6" (15.2 cm), flexible shaft, 3 sharp prongs
18-110 18-112 18-114 18-116	DESMARRES Lid Retractor 5½" (14 cm), size 1, 11 mm DESMARRES Lid Retractor 5½" (14 cm), size 2, 13 mm DESMARRES Lid Retractor 5½" (14 cm), size 3, 15 mm DESMARRES Lid Retractor 5½" (14 cm), size 4, 17 mm
18-120	JAEGER Lid Plate 4" (10.2 cm), blades 20 mm and 23 mm





All units are constructed with jewelled hammers to minimize friction. Autoclavable. Scale markings are etched for accuracy. Scales are inclined, permitting reading from above.

Each unit supplied with certification to meet specification established by the Committee on Standardization of Tonometers of the A.A.O.O. Individually packaged in fitted case with weights and test block.


18-132 Original Model, with 3 weights (5.5 gr., 7.5 gr. and 10 gr.) Plunger, Footplate and Test Block are stainless steel; remainder made of rustless materials
 18-132A Same as No. 18-132, but with Allen Plunger Retractor*

18-134 Improved Model, with 4 weights (5.5 gr., 7.5 gr., 10 gr. and 15 gr.), mirror insert on scale reduces error of parallax; entirely made of stainless steel.

18-134A Same as No. 134, but with Allen Plunger Retractor*

18-139 Sterilizing Case for Schiotz Tonometers 6½" (16.5 cm) x 3½" (8.9 cm) x 1½" (3.2 cm) fitted to firmly hold tonometer, weights and test block


*With the Allen Plunger Retractor, the plunger is held in retracted position, flush with the surface of the footplate until tonometer is in place on cornea. This feature will prevent any possible damage to cornea tissue.






18-140

Ophthalmic Instruments




18-154



18-140	VON GRAEFE Cataract Knife No. 1, 1.5 x 25 mm
18-142	VON GRAEFE Cataract Knife No. 2, 1.7 x 27 mm
18-144	VON GRAEFE Cataract Knife No. 3, 2 x 30 mm
18-146	VON GRAEFE Cataract Knife No. 4, 2.3 x 33 mm
18-148	VON GRAEFE Cataract Knife No. 5, 2.5 x 35 mm
18-152	WHEELER Discission Knife 1 x 18 mm
18-154	KNAPP Cataract Knife 1 x 20 mm
18-170	BARRAQUER Blade Breaker and Holder 4½" (11.4 cm), with ratchet lock
18-171	CASTROVIEJO Blade Breaker and Holder 5¼" (13.3 cm), with push-pull lock
18-172	TROUTMAN-BARRAQUER Mini Blade Breaker and Holder 3%" (9.5 cm), angled laws 3 x 10 mm, with push-pull lock

18-148







18-372

18-399

Ophthalmic Instruments







18-374














DAVIS Foreign Body Spud 4¾" (12.1 cm), 0.8 mm blade 18-372 18-374 DIX Foreign Body Spud 4¾" (12.1 cm), 1.7 x 15 mm blade 18-376 WALTON Foreign Body Gouge 41/2" (11.4 cm), round end 18-380 ELLIS Foreign Body Spud 41/2" (11.4 cm), 1.2 mm blade, curved

18-384 ALVIS Foreign Body Spud 4%" (12.1 cm), oval curette 1 x 3 mm with fine point

HOSFORD Foreign Body Spud 4%" (12.1 cm), gouge blade 1 x 7 mm 18-388 rounded tip, angled to right 18-390

HOSFORD Foreign Body Spud 4¾" (12.1 cm), gouge blade 1 x 7 mm rounded tip, angled to left

18-398 Eye Magnet and Foreign Body Loop. The magnet probe and the Perlon loop are reversible in the pencil-shape handle. In plastic case.


18-398M Eye Magnet and Foreign Body Loop with adjustable Magnifier attached. Reversible plastic handle, metal parts are plated, in plastic case.

> MILTEX Eye Magnet; a powerful magnetic cylinder with one round and one pointed interchangeable tip. In velvet lined steel case.



SELTZER Hook 41/2" (11.4 cm), sharp point






18-400 18-402 18-402A 18-402B 18-404	DIX Foreign Body Needle and Spud, protected in sliding handle DIX Foreign Body Needle and Spud, protected in reversible screw handle DIX Foreign Body Needle with screw socket DIX Foreign Body Spud with screw socket DIX Foreign Body Needle and Gouge, protected in reversible screw handle		
18-412 18-414	TYRELL Hook 5" (12.7 cm), sharp TYRELL Hook 5" (12.7 cm), blunt		
18-416	GUTHRIE Fixation Hook 5" (12.7 cm), sharp double prongs 1.5 mm wide		2.6
18-418 18-420	WIENER Corneal Hook 5" (12.7 cm), sharp WIENER Corneal Hook 5" (12.7 cm), blunt	2009	
18-424	O'CONNOR Hook 4%" (12.1 cm), sharp		
18-428	NUGENT Hook 43/4" (12.1 cm), right angle, sharp		33
21-152	MILTEX Hook 4%" (12.1 cm), sharp prong 2 mm diameter	21-152	33-121

33-121





18-439	McGANNON Retractor 5" (12.7 cm), fenestrated blade 2.5 mm wide
18-440	STEVENS Tenotomy Hook 5" (12.7 cm), right angle 5 mm deep
18-442	GREEN Strabismus Hook 5" (12.7 cm), 9 mm deep
18-446	JAMESON Strabismus Hook 5" (12.7 cm), 10 mm deep
18-450 18-452 18-454	VON GRAEFE Strabismus Hook 5½" (14 cm), small size 8 mm long VON GRAEFE Strabismus Hook 5½" (14 cm), medium size 10 mm long VON GRAEFE Strabismus Hook 5½" (14 cm), large size 11 mm long
18-472	DAVIEL Lens Spoon 5" (12.7 cm), 3 mm wide
18-480 18-482	LEWIS Lens Loop 51/4" (13.3 cm), serrated loop 5.5 x 8 mm LEWIS Lens Loop 51/4" (13.3 cm), smooth loop 5.5 x 8 mm
18-484 18-486	WILDER Lens Loop 51/4" (13.3 cm), serrated loop 4 x 15 mm WILDER Lens Loop 51/4" (13.3 cm), smooth loop 4 x 15 mm











18-514




18-520

MEYHOEFER Chalazion Curette 5" (12.7 cm), size 00, 1 mm

HEBRA Chalazion Curette 5" (12.7 cm), oval cup 1 x 1.7 mm

HEBRA Chalazion Curette 5" (12.7 cm), oval cup 2 x 3.3 mm

HEBRA Chalazion Curette 5" (12.7 cm), oval cup 2.5 x 4 mm





18-498 18-499 18-500 18-502 18-504

MEYHOEFER Chalazion Curette 5" (12.7 cm), size 0, 1.5 mm MEYHOEFER Chalazion Curette 5" (12.7 cm), size 1, 1,7 mm MEYHOEFER Chalazion Curette 5" (12.7 cm), size 2, 2 mm MEYHOEFER Chalazion Curette 5" (12.7 cm), size 3, 2.5 mm MEYHOEFER Chalazion Curette 5" (12.7 cm), size 4, 3.5 mm

18-512 18-514

18-520

18-522

18-524

18-530

18-532

18-534

18-536

18-540

18-542

18-550

18-524

18-550 18-540 18-542

18-506

18-510

SKEELE Chalazion Curette 5" (12.7 cm), round cup 1 mm, serated edge

SKEELE Chalazion Curette 5" (12.7 cm), round cup 2 mm, serated edge

SKEELE Chalazion Curette 5" (12.7 cm), round cup 2.5 mm, serated edge

blunt edge, 2 mm

HEATH Chalazion Curette 4" (10.2 cm), size 0, one sharp and one blunt edge, 1/2 mm

HEATH Chalazion Curette 4" (10.2 cm), size 1, one sharp and one blunt edge, 1 mm)

HEATH Chalazion Curette 4" (10.2 cm), size 2, one sharp and one

HEATH Chalazion Curette 4" (10.2 cm), size 3, one sharp and one blunt edge, 3 mm


BUNGE Evisceration Spoon 51/2" (14 cm), small size 8 mm cup BUNGE Evisceration Spoon 51/2" (14 cm), large size 11 mm cup

WELLS Enucleation Spoon 6" (15.2 cm), cup 21 mm




18-584





18-562	LINDNER Cyclodialysis Spoon and Spatula 5% " (14 cm), malleable spatula end is graduated
18-564	CASTROVIEJO Spatula 5½" (14 cm), double end 3.5 mm wide angled and 2 mm wide slightly curved, sterling
18-566	SMITH-FISHER Iris Spatula 5" (12.7 cm), one flat end 2 mm wide, one spoon end 2 mm, sterling
18-568	FISHER Spoon and Needle 5% " (14 cm), spoon end 7 x 15 mm, needle 33 mm long
18-570	GREEN Spatula 6" (15.2 cm), one end 4 x 14 mm, other end 5 x 16 mm
18-573A	CASTROVIEJO Cyclodialysis Spatula $5\%^{\circ}$ (13.3 cm), blades 10 mm and 15 mm x 0.75 mm wide
18-574	CASTROVIEJO Cyclodialysis Spatula $5\%^{\circ}$ (13.3 cm), blades 10 mm and 15 mm x 0.5 mm wide
18-576	CASTROVIEJO Cyclodialysis Spatula 51/4" (13.3 cm), blade 10 x 1 mm
18-584	WECKER Iris Spatula 5" (12.7 cm), sterling blade 30 x 2 mm







18-613	KIRBY Muscle Hook and Expressor 6" (15.2 cm), double end
18-620 18-624	KIRBY Iris Forceps 4" (10.2 cm), straight, 1 x 2 teeth KIRBY Iris Forceps 4" (10.2 cm), curved, 1 x 2 teeth
18-628	KIRBY Fixation Forceps 4" (10.2 cm), straight 1 x 2 teeth, with slide lock
18-632	KIRBY Cornea-Scleral Forceps 4" (10.2 cm), angled tips with 1 x 2 teeth








CASTROVIEJO Caliper 31/4" (8.3 cm), a combination measuring and marking 18-650 caliper, graduated 0-20 mm THORPE Caliper 41/2" (11.4 cm), graduated in inches and mm, permits 18-657 measurements in deep areas JAMESON Caliper 39/4" (9.5 cm), graduated in inches and mm, with thin tips, chrome 18-658 Flexible Stainless Ruler 6" (15.2 cm) x ½" (1.3 cm), graduated in 1/32" (0.8 mm) fractions and mm 18-660 WADSWORTH-TODD Cautery 5¾" (14.6 cm), copperhead 9 mm dia. with platinum point 18-666










EYE DRESSING FORCEPS 4" (10.2 cm), serrated tips

Standard Pattern 0.8 mm wide tips		Extra Delicate Pattern 0.5 mm wide tips	
18-780	straight	18-779	straight
18-782	half curved	18-781	half curved
18-784	full curved	18-783	full curved






IRIS TISSUE FORCEPS 4" (10.2 cm), 1 x 2 teeth

	andard Pattern I mm wide tips		icate Pattern n wide tips
18-786 18-787 18-788 18-790	straight straight (Non Magnetic) half curved full curved	18-785 18-789 18-791	straight half curved full curved
18-793 18-794	O'CONNOR Iris Forceps 3% BRACKEN Iris Forceps 4" (1	, ,	ular jaws 0.5×9 mm with 1×2 teeth ine cross-serrated tips



















18-795

FOERSTER Eye Dressing Forceps 3¾" (9.5 cm), octagonal grip, straight, serrated tips 0.5 mm wide

18-798

18-795A

FOERSTER Eye Dressing Forceps 3¾" (9.5 cm), octagonal grip, full curved, serrated tips 0.5 mm wide

18-796

FOERSTER Iris Tissue Forceps 3¾" (9.5 cm), octagonal grip,

straight, 1 x 2 teeth, 0.5 mm wide

18-798 FOERSTER Iris Tissue Forceps 3¾" (9.5 cm), octagonal grip, full curved, 1 x 2 teeth, 0.5 mm wide

18-804

STEVENS Iris Forceps 4" (10.2 cm), straight, 1 x 2 teeth, 1 mm wide

18-814

GRAEFE Eve Dressing Forceps 2%" (7 cm), straight, serrated

18-816

GRAEFE Iris Forceps 2¾" (7 cm), straight, 1 x 2 teeth

18-818

GRAEFE Eye Dressing Forceps 23/4" (7 cm), curved, serrated GRAEFE Iris Forceps 23/4" (7 cm), curved, 1 x 2 teeth

18-820

GILL (CHANDLER) Iris Forceps 3" (7.6 cm), fine cross serrated tips. 0.6 mm wide

18-825 18-834

GILL Iris Forceps 27/6" (7.3 cm), tips 0.6 mm wide, with 1 x 2 teeth

18-827

TROUTMAN-BARRAQUER Colibri Type Corneal Utility Forceps 3" (7.6 cm), with 1 x 2 teeth, 0.4 mm, set at 45° angle, with tying platform

18-828


TROUTMAN-BARRAQUER Colibri Type Corneal Utility Forceps

3" (7.6 cm), with 1 x 2 teeth, 0.12 mm, set at 45° angle, with tying platform






18-828



18-829 18-830 18-835 18-836 18-837 18-838 18-866 18-867 18-868	BOTVIN Iris Forceps 2%" (7.3 cm), fine cross serrated tips, 0.5 mm widded and continuous services and			The state of the s
18-850	GRAEFE Fixation Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth			
18-852 18-854	GRAEFE Fixation Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth, Non-Magnetic GRAEFE Fixation Forceps 4%" (11.2 cm), standard jaws 4.5 mm wide, with fine teeth, with catch	18-840	18-850 18-852	18-854








18-860	WALDEAU Fixation Forceps 4% " (11.2 cm), concave jaws 5 mm wide, with fine teeth
18-870	GUYTON-NOYES Fixation Forceps 41/4" (10.8 cm), with slide catch, curved jaws, 4.5 mm wide, with fine teeth
18-876	CASTROVIEJO (ALVIN) Fixation Forceps 3¾" (9.5 cm), straight, 1 x 1 curved teeth extending beyond each other
18-878	CASTROVIEJO Fixation Forceps 33/4" (9.5 cm), straight, 2 x 2 curved teeth extending beyond each other
18-880 18-881	PIERSE Colibri Type Corneal Forceps 31/4" (8.3 cm), tips 0.12 mm PIERSE Colibri Type Corneal Forceps 31/4" (8.3 cm), tips 0.3 mm
18-883	CASTROVIEJO Colibri Type Corneal Forceps 4" (10.2 cm), tips 0.12 mm, with tying platform
18-896 18-897	GREEN Fixation Forceps 4" (10.2 cm), T-shape, 12 mm wide 19 x 20 teeth GREEN Fixation Forceps 4" (10.2 cm), T-shape, 5 mm wide 8 x 9 teeth






18-904	O'BRIEN Fixation Forceps 4" (10.2 cm), delicate 1 x 2 teeth, 0.9 mm wide, angled forward
18-906 18-908	ELSCHNIG Fixation Forceps 41/4" (10.8 cm), strong 1 x 2 teeth, 0.9 mm wide, angled upwards ELSCHNIG-O'CONNOR Fixation Forceps 41/4" (10.8 cm), with slide lock, strong 1 x 2 teeth, 1 mm wide, angled upwards
18-912 18-914	LESTER Fixation Forceps 3% " (9.5 cm), delicate 1 x 2 teeth, 0.6 mm wide LESTER Fixation Forceps 3% " (9.5 cm), delicate 2 x 3 teeth, 1 mm wide
18-916	HARMON Fixation Forceps 3" (7.6 cm), 4 x 5 teeth, delicate 2.5 mm wide jaws



18-946


Ophthalmic Instruments



18-947	HARMS Suturing Forceps 4%" (10.5 cm), curved, with 11 mm wide handles, tying platform, 0.6 mm wide at tip
18-948 18-949	McPHERSON Micro Iris Suturing Forceps $3\frac{1}{2}$ " (8.9 cm), straight, smooth platform 0.3 x 5 mm, no teeth McPHERSON Micro Iris Suturing Forceps $3\frac{1}{2}$ " (8.9 cm), angled, smooth platform 0.3 x 5 mm, no teeth
18-950 18-951 18-952	CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.9 mm, with tying platform CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.3 mm, with tying platform CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.6 mm, with tying platform
18-953	CASTROVIEJO Micro Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.12 mm, with tying platform
18-954	CASTROVIEJO Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.6 mm, with tying platform
18-955	CASTROVIEJO Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.9 mm, with tying platform
18-956 18-958	NUGENT Utility Forceps 4¼" (10.8 cm), smooth tips 1.2 mm wide NUGENT Utility Forceps 4¼" (10.8 cm), cross serrated tips 1.2 mm wide
18-962 18-964	McCullough Utility Forceps 4" (10.2 cm), smooth tips 1.5 mm wide McCullough Utility Forceps 4" (10.2 cm), cross serrated tips 1.5 mm wide

HARMS Suturing Forcess 41/4" (10.5 cm), straight, with 11 mm wide handles, twing platform, 0.6 mm wide at tip





18-966	McCULLOUGH Suturing Forceps 4" (10.2 cm), cross serrated jaws 1.5 mm wide with 1 x 2 teeth
18-970	WILLS HOSPITAL Utility Forceps 4" (10.2 cm), with cross serrated trying platform 1.2 mm wide
18-971	BONACCOLTO Utility Forceps 4" (10.2 cm), longitudinal serrations with cross serrations on tip, 1.2 mm wide with tapered tips
18-974	SAUER Suturing Forceps 3% " (8.9 cm), 1 x 2 curved teeth, overlapping each other, 0.6 mm wide
18-975	MANHATTAN EYE & EAR Suturing Forceps 3½" (8.9 cm), with smooth platform, 1 x 2 curved teeth, overlapping each other, 0.7 mm wide
18-978	SCHAAF Foreign Body Forceps 3¾" (9.5 cm), grooved tips















18-980

THORPE Corneal and Splinter Forceps 4" (10.2 cm), very sharp points, fine serrated tips

18-1000 18-1002

ARRUGA Capsule Forceps 4" (10.2 cm), standard curve ARRUGA Capsule Forceps 4" (10.2 cm), new curve

18-1034

MILTEX Extra-Capsular Forceps 4" (10.2 cm), improved pattern with 11 x 11 teeth, with long convex curve


18-1050

CASTROVIEJO Capsule Forceps 4" (10.2 cm), cross action, with stop to limit opening of blades

18-1058

SCHWEIGER Extra-Capsular Forceps 33/4" (9.5 cm), 4 x 5 teeth





18-1060 FUCHS Extra-Capsular Forceps 29/4" (7 cm), 5 x 5 teeth

18-1080 SHEPARD Intraocular Lens Forceps 41/2" (11.4 cm), 17 mm jaws, serrated tips

18-1081 CLAYMAN Intraocular Lens Forceps 41/2" (11.4 cm), 15 mm jaws, smooth tips

18-1100 LITTAUER Cilia Forceps 31/2" (8.9 cm), 4 mm wide jaws with fine horizontal serrations


18-1102 BEER Cilia Forceps 31/2" (8.9 cm), 4 mm wide smooth jaws

18-1104 ZIEGLER Cilia Forceps 31/2" (8.9 cm), 2 mm wide smooth jaws

18-1106 DOUGLAS Cilia Forceps 31/2" (8.9 cm)), thin tips 2.5 mm wide, with fine horizontal serrations







18-1107	SWISS Cilia and Suture Forceps 3%" (9.5 cm), 3 mm wide precision fitted diagonal jaws
18-1108 18-1110	BERGH Cilia Forceps 41/4" (10.8 cm), 10 mm wide jaws with horizontal serrations BERGH Cilia Forceps 31/2" (8.9 cm), 5 mm wide jaws with horizontal serrations
18-1111	BARRAQUER Cilia and Suture Forceps 41/2" (11.4 cm), with 5 mm smooth platform
18-1112	BENNETT Cilia Forceps 3" (7.6 cm), cup shaped jaws 3 mm diameter
18-1114	BEAUPRE Cilia Forceps 4½" (11.4 cm), smooth jaws 12 mm long, 45° angle
18-1116	GRADLE Cilia Forceps 4" (10.2 cm), precision ground jaws 1.4 mm wide















18-1136 BERKE Ptosis Forceps 4" (10.2 cm), blades 20 mm long with longitudinal serrations, with slide lock

18-1138 BERKE Ptosis Forceps 4" (10.2 cm), blades 27 mm long with longitudinal serrations, with slide lock

18-1150 ERHARDT Lid Forceps 31/2" (8.9 cm), serrated jaw 9.5 mm wide

18-1160 SNELLEN Entropium Forceps 3½" (8.9 cm), left, 35 mm wide plate 18-1162 SNELLEN Entropium Forceps 3½" (8.9 cm), right, 35 mm wide plate

18-1170 AYER Chalazion Forceps 3½" (8.9 cm), inside ring diameter 8.5 mm,

F1170 AYEH Chaiazion Forceps 3½" (8.9 cm), inside ring diameter 8.5 mm, solid blade 8 mm diameter

18-1174 LAMBERT Chalazion Forceps 31/2" (8.9 cm), inside ring diameter 8 mm

18-1178 HUNT Chalazion Forceps 3¾" (9.5 cm), inside ring diameter 12 mm

18-1182 HEATH Chalazion Forceps 4" (10.2 cm), inside ring measures 12 x 14 mm









18-1188



18-1208



18-1200 to 18-1204








18-1188 SPENCER Chalazion Forceps 4" (10.2 cm), inside ring diameter 10 mm, with slide lock 18-1190 WIES Chalazion Forceps 4" (10.2 cm), inside ring measures 9 x 12 mm 18-1200 DESMARRES Chalazion Forceps 31/2" (8.9 cm), small size, inside ring measures 11 x 17 mm 18-1202 DESMARRES Chalazion Forceps 31/2" (8.9 cm), medium size, inside ring measures 12 x 23 mm 18-1204 DESMARRES Chalazion Forceps 31/2" (8.9 cm), large size, inside ring measures 17 x 28 mm

18-1220	PRINCE Trachoma Forceps 31/2" (8.9 cm), rings measure 8 x 10 mm inside

FRANCIS Chalazion Forceps 3¾" (9.5 cm), inside ring measures 12 x 14 mm


18-1260 JAMESON Recession Forceps 3¾" (9.5 cm), left, with slide lock 18-1262 JAMESON Recession Forceps 3¾" (9.5 cm), right, with slide lock





18-1392 18-1394	BONN Miniature Iris Scissors $3\%"$ (8.9 cm), with 15 mm blades, straight BONN Miniature Iris Scissors $3\%"$ (8.9 cm), with 15 mm blades, curved		
18-1396 18-1398 5-300 5-302 5D-303 5-304 5-306 5-308	Iris Scissors 3½" (8.9 cm), with 20 mm blades, straight, delicate Iris Scissors 3½" (8.9 cm), with 20 mm blades, curved, delicate Iris Scissors 4" (10.2 cm), straight, standard pattern Iris Scissors 4" (10.2 cm), curved, standard pattern Iris Scissors 4½" (11.4 cm), curved sideways Iris Scissors 4½" (11.4 cm), straight, standard pattern Iris Scissors 4½" (11.4 cm), curved, standard pattern Iris Scissors 4½" (11.4 cm), angular	5-302 5-306	5-308





 18-1404
 MILTEX Iris Scissors 4" (10.2 cm), straight, heavy pattern

 18-1405
 MILTEX Iris Scissors 4" (10.2 cm), curved, heavy pattern

 18-1415*
 Iris Scissors 4" (10.2 cm), straight, standard pattern, ribbon type

 18-1416*
 Iris Scissors 4" (10.2 cm), curved, standard pattern, ribbon type

 18-1417*
 Iris Scissors 3%" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon type


 18-1418*
 Utility Scissors 3%" (9.5 cm), curved blunt blades, 15 mm long, ribbon type

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.





18-1420	KNAPP Iris Scissors 4" (10.2 cm), straight, sharp-blunt points
18-1422	KNAPP Iris Scissors 4" (10.2 cm), straight, blunt-blunt points
18-1426	KNAPP Iris Scissors 4" (10.2 cm), curved, sharp-blunt points
18-1428	KNAPP Iris Scissors 4" (10.2 cm), curved, blunt-blunt points
18-1430	Eye Scissors 4" (10.2 cm), with probe points, straight
18-1432	Eye Scissors 4" (10.2 cm), with probe points, curved
5-310	Eye Scissors 4½" (11.4 cm), with probe points, angled on side












18-1440





18-1460 to 18-1476



5D-324

LAGRANGE Scierectomy Scissors 41/4" (10.8 cm), curved, one serrated blade, angled shanks

18-1440

WILMER Conjunctival and Utility Scissors, 4" (10.2 cm), angled on flat

5-312 5-314 18-1453*

18-1454*

Strabismus Scissors 4" (10.2 cm), straight, standard pattern Strabismus Scissors 4" (10.2 cm), curved, standard pattern Strabismus Scissors 41/4" (10.8 cm), straight, ribbon type Strabismus Scissors 41/4" (10.8 cm), curved, ribbon type

STEVENS TENOTOMY SCISSORS

41/8" (10.5 cm)	4½" (11.4 cm)	
With Short Blades	With Long Blades	
18-1460	Straight Sharp Points	18-1470
18-1462	Straight Blunt Points	18-1472
18-1464	Curved Sharp Points	18-1474
18-1466	Curved Blunt Points	18-1476

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.



18-1460 to 18-1466



18-1470 to 18-1476













18-1477* STEVENS Tenotomy Scissors 3¾" (9.5 cm), straight, blunt points, ribbon type
18-1478* STEVENS Tenotomy Scissors 3¾" (9.5 cm), curved, blunt points, ribbon type
18-1480 WESTCOTT Utility Scissors 4½" (11.4 cm), curved, blunt tips
18-1482 SCHEIE-WESTCOTT Corneal Section Scissors 5" (12.7 cm), right, curved, blunt tips
18-1484 SCHEIE-WESTCOTT Corneal Section Scissors 5" (12.7 cm), left, curved, blunt tips

WESTCOTT Stitch Scissors 4½" (11.4 cm), curved, sharp points WESTCOTT Tenotomy Scissors 5½" (13.3 cm), right, wide handles

WESTCOTT Tenotomy Scissors 51/4" (13.3 cm), left, wide handles

*Rilbbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.




18-1487 18-1488

18-1486

18-1487 18-1488











Enucleation Scissors 5" (12.7 cm), light curve 18-1494 Enucleation Scissors 5" (12.7 cm), medium curve 18-1496 18-1498 Enucleation Scissors 5" (12.7 cm), full curve 18-1510 NOYES Iris Scissors 41/2" (11.4 cm), straight, sharp-sharp points NOYES Iris Scissors 41/2" (11.4 cm), straight, sharp-blunt points 18-1512

NOYES Iris Scissors 41/2" (11.4 cm), straight, blunt-blunt points 18-1514

18-1516 NOYES Iris Scissors 41/2" (11.4 cm), curved, sharp-sharp points

18-1520 BARRAQUER Iris Scissors 21/4" (5.7 cm), angled 7 mm blades, sharp points BARRAQUER Iris Scissors 21/4" (5.7 cm), angled 7 mm blades, blunt points 18-1521

18-1524 DeWECKER Iris Scissors 41/4" (10.8 cm), angled 7 mm blades, sharp-sharp points













18-1540

18-1550 18-1552

18-1554

18-1558 18-1559

18-1540	NOYES Iridectomy Scissors 5½" (14 cm), straight 7 mm blades, sharp-sharp points
18-1550	McCLURE Iris Scissors 3¾" (9.5 cm), angled over edge, 4 mm blades with very delicate sharp points
18-1552	McCLURE Iris Scissors 3¾" (9.5 cm), angled on flat, 4 mm blades with very delicate sharp points
18-1554	McCLURE Iris Scissors 3¾" (9.5 cm), angled on flat, regular 10 mm blades with sharp points
18-1558	CASTROVIEJO Corneal Scissors 4" (10.2 cm), right, sub-miniature for micro surgery, 6 mm blades
18-1559	CASTROVIEJO Corneal Scissors 4" (10.2 cm), left, sub-miniature for micro surgery, 6 mm blades
18-1560	CASTROVIEJO Corneal Section Scissors 4½" (11.4 cm), right, inner blade 1 mm longer, blunt tips
18-1562	CASTROVIEJO Corneal Section Scissors 4½" (11.4 cm), left, inner blade 1 mm longer, blunt tips



18-1562













18-1570 18-1572

18-1570

CASTROVIEJO Keratoplasty Scissors 3¾" (9.5 cm), angled blades


11 mm, sharp points CASTROVIEJO Corneal Scissors 3¾" (9.5 cm), angled blades 18-1572 11 mm, blunt points

CASTROVIEJO Corneal Scissors 3¾" (9.5 cm), curved, sharp points 18-1576 CASTROVIEJO Corneal Scissors 3¾" (9.5 cm), curved, blunt points 18-1578

18-1594 McGUIRE Corneal Scissors 41/8" (10.5 cm), right McGUIRE Corneal Scissors 41/8" (10.5 cm), left 18-1596

18-1600 AEBLI Corneal Scissors 4" (10.2 cm), straight 18-1602 AEBLI Corneal Scissors 4" (10.2 cm), angled on flat, right

AEBLI Corneal Scissors 4" (10.2 cm), angled on flat, left 18-1604















18-1619

18-1622

18-1624

18-1626

18-1630 18-1631

18-1618 18-1619	Micro Iris Scissors 4" (10.2 cm), straight, sharp points Micro Iris Scissors 4" (10.2 cm), curved, sharp points
18-1620	VANNAS Capsulotomy Scissors 3¼" (8.3 cm), straight, very delicate sharp points
18-1622	VANNAS Capsulotomy Scissors 3¼" (8.3 cm), curved, very delicate sharp points
18-1623	VANNAS Capsulotomy Scissors 3¼" (8.3 cm), angled to side, very delicate sharp points
18-1624	VANNAS Capsulotomy Scissors 31/4" (8.3 cm), angled on flat, very delicate sharp points
18-1625	McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), straight, sharp points
18-1626	McPHERSON-VANNAS Micro Iris Scissors 31/4" (8.3 cm), curved, sharp points
18-1630	Micro-VANNAS Scissors 31/4" (8.3 cm), straight ultra fine blades


Micro-VANNAS Scissors 31/4" (8.3 cm), curved ultra fine blades

Micro Scissors 31/4" (8.3 cm), straight super fine blades

Micro Scissors 31/4" (8.3 cm), curved super fine blades

HOLTH Corneoscleral Punch 31/2" (8.9 cm), 1 mm bite HOLTH Corneoscleral Punch 31/2" (8.9 cm), 1.5 mm bite





18-1718

18-1631

18-1633


18-1634

18-1716

18-1718







SPENCER Stitch Scissors 31/2" (8.9 cm), straight, delicate

9-101 SHORTBENT Stitch Scissors 31/2" (8.9 cm), curved to side, delicate

O'BRIEN Stitch Scissors 3¾" (9.5 cm), angled to side, sharp points 9-110

Eye Suture (GRADLE) Scissors 3¾" (9.5 cm), slightly curved, 18-1652

sharp points

9-100


18-1653* Eye Suture Scissors 4" (10.2 cm), ribbon type, slightly curved,

sharp points



*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.





18-1800	KALT Needle Holder 51/4" (13.3 cm), standard model, delicate jaws
18-1810	GREEN Needle Holder and Suturing Forceps 41/4" (10.8 cm), angled jaws, double end
18-1820 18-1824	CASTROVIEJO Needle Holder 5½" (14 cm), straight, without lock CASTROVIEJO Needle Holder 5½" (14 cm), curved, without lock
18-1828 18-1832	CASTROVIEJO Needle Holder 5½" (14 cm), straight, with lock CASTROVIEJO Needle Holder 5½" (14 cm), curved, with lock
18-1831 18-1833	CASTROVIEJO Needle Holder 5½" (14 cm), straight, extra delicate jaws, with lock CASTROVIEJO Needle Holder 5½" (14 cm), curved, extra delicate jaws, with lock
18-1834	PATON Needle Holder 41/6" (11.4 cm) straight without lock







8-26







18-1837	McPHERSON Needle Holder for microsurgery 4" (10.2 cm), curved smooth tapered jaws
18-1838	BARRAQUER Needle Holder 5" (12.7 cm), curved, without lock, 10 mm wide hollow round handle
18-1839	BARRAQUER Needle Holder 5" (12.7 cm), curved, with lock, 10 mm wide hollow round handle
18-1840	BARRAQUER Needle Holder 51/8" (13 cm), curved, without lock delicate tapered jaws, 8 mm·wide round handle


BOYNTON Needle Holder 5" (12.7 cm), delicate 18-1858 STEVENS Needle Holder 61/2" (16.5 cm), with ratchet lock





CARN-N-SERT NEEDLE HOLDERS

with tungsten carbide inserts insure firm holding of needles and sutures



18-1800TC KALT Needle Holder 5½" (14 cm), serrated Carb-N-Sert laws 16000 teeth p. sq. inch

18-1816TC MILTEX-CASTROVIEJO Needle Holder 41/8" (10.5 cm), round handles, with lock,

delicate smooth Carb-N-Sert jaws

18-1828TC CASTROVIEJO Needle Holder 5½" (14 cm), straight, with lock, smooth Carb-N-Sert jaws **18-1828A-TC** CASTROVIEJO Needle Holder 5½" (14 cm), straight, with lock, serrated Carb-N-Sert jaws

16000 teeth p. sq. inch

18-1832TC CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock, smooth Carb-N-Sert jaws

18-1832A-TC CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock, serrated Carb-N-Sert jaws

16000 teeth p. sq. inch



Ophthalmic Instruments

















7-280 Serrefine 11/2" (3.8 cm), straight, with suture hole Serrefine 21/2" (6.4 cm), straight, with suture hole 7-284 7-290 Serrefine 11/2" (3.8 cm), curved, with suture hole Serrefine 21/2" (6.4 cm), curved, with suture hole 7-294

7-310 DIETHRICH Micro Bulldog clamp 11/6" (4.7 cm), straight jaws 8 mm x 1.2 mm 7-315 DIETHRICH Micro Bulldog clamp 13/4" (4.4 cm), angled jaws 8 mm x 1.2 mm

7-8 HALSTED Mosquito Forceps 5" (12.7 cm), straight, extra delicate HALSTED Mosquito Forceps 5" (12.7 cm), curved, extra delicate 7-10 18-1934 HALSTED Mosquito Forceps 5" (12.7 cm), straight, Non-Magnetic HALSTED Mosquito Forceps 5" (12.7 cm), curved, Non-Magnetic 18-1936 7-24

HARTMAN Mosquito Forceps 31/2" (8.9 cm), straight, delicate 7-26 HARTMAN Mosquito Forceps 31/2" (8.9 cm), curved, delicate


JACOBSON Micro Mosquito Forceps 5" (12.7 cm), straight, extremely delicate 17-2600 JACOBSON Micro Mosquito Forceps 5" (12.7 cm), curved, extremely delicate 17-2602



17-2602

Ophthalmic Instruments





ARRUGA LACRIMAL TREPHINES


	AHHOUA	LAGISMAL THE T	IIILO	
18-19	enter Point 9 50-7 9 50-8 —	Diameter 7 mm 8 mm 9 mm 10 mm	Center Point 18-1952-7 18-1952-8 18-1952-9 18-1952-10	
18-1960	FREER Chis	sel 6½" (16.5 cm), c	urved blade 4 mm wide	
18-1962 18-1964			(15.9 cm), straight blade 5 mm wide .9 cm), straight blade 4.5 mm wide	18-2200
18-1968	FREER Elev	/ator 7" (17.8 cm), d	louble end sharp and blunt blades 5 mi	m wide
18-1990 18-1992 18-1994 18-1996	KERRISON KERRISON	RONGEUR 3" (7.6 RONGEUR 3" (7.6	cm) shaft, size 0, 3.5 mm bite cm) shaft, size 1, 4 mm bite cm) shaft, size 2, 4.5 mm bite cm) shaft, size 3, 5 mm bite	
18-2200	FOX Eye Sh	nields, aluminum, pk	rg. 12	














19-10A

19-21 to 19-25

19-30 to 19-36


BOU	CHERON EAR SPE	CULA, Round, Chrome		GRUBER EAR SI	PECULA, Oval,	
19-2	Set of 4 sizes 1,	2, 3 and 4	Standard Adult Sizes, Long Pattern, Chrome			
19-2-1	Single size 1	4.5 mm	19-21	Set of 4 sizes 1,	2, 3 and 4	
19-2-2 19-2-3 19-2-4	Single size 2 Single size 3 Single size 4	5.5 mm 6.5 mm 7.5 mm	19-22 19-23	Single size 1 Single size 2	4x5 mm 5x6 mm	
13-2-4	Single Size 4	7.5 11811	19-24	Single size 3	6x7 mm	
	GRUBER EAR SI Child and Int		19-25	Single size 4	7x8 mm	
	Light Short Pat	tern, Chrome	TC	YNBEE EAR SPECI	JLA, Round, Chro	ome
19-10	Set of 4 sizes 1, top case	2, 3 and 4 with screw	19-30	Set of 3 sizes 1,	2 and 3	
19-12	•	2, 3 and 4 without case	19-32	Single size 1	4 mm	
*** ****	00.07.70.200.7,	-,, -, -, -, -, -, -, -, -, -, -, -,	19-34	Single size 2	5 mm	
19-14	Single size 1	2x3 mm	19-36	Single size 3	6 mm	
19-16	Single size 2	3x4 mm		ŭ		
19-18	Single size 3	4x5 mm				
19-20	Single size 4	5x6 mm				
19-10A	Case only for set	t of 4 Gruber Specula, es				














19-71

19-70

HARTMAN EAR SPECULA.	. Round	١.	JLA	U	C	Е	P	S	R	Α	E	۷	ΑI	М	łΤ	٩F	Н
----------------------	---------	----	-----	---	---	---	---	---	---	---	---	---	----	---	----	----	---

19-40	Set of 3 sizes 3, 4 and 5 mm
19-48	Set of 4 sizes 4, 5, 6 and 7 mn
19-48-A	Single size 3 mm
19-48-B	Single size 4 mm
19-48-C	Single size 5 mm
19-48-D	Single size 6 mm
19-48-E	Single size 7 mm


FARRIOR EAR SPECULA, Oval,	angled	end
nonglare stainless	_	

19-50-1	Single size 1	3 x 4 mm
19-50-2	Single size 2	4 x 5 mm
19-50-3	Single size 3	4.5 x 5.5 mm
19-50-4	Single size 4	5 x 6 mm
19-50-5	Single size 5	5.5 x 6.5 mm
19-50-6	Single size 6	6 x 7 mm
19-50-7	Single size 7	6.5 x 7.5 mm
19-50-8	Single size 8	7 x 8 mm
19-50-9	Single size 9	7.5 x 8.5 mm
19-70	BRUENING Mag	gnifier Otoscope, Chrom

19-70 BRUENING Magnifier Otoscope, Chrome, 3.5x magnif., with 3 specula, small, medium and large sizes; complete with inflation bulb and tubing, in fitted case
 19-71 Set of 3 specula, small, medium and large sizes

19-71 Set of 3 specula, small, medium and large sizes for Bruening Otoscope, chrome
19-71-1 Single speculum, small size, chrome
19-71-2 Single speculum, medium size, chrome
19-71-3 Single speculum, large size, chrome






19-74




19-74 SIEGEL Magnifier Otoscope, chrome, 3.15x magnif., with 3 specula, small, medium and large sizes; complete with inflation bulb and tubing
 19-75 Rubber Bulb with tubing for Bruening and Siegel Otoscopes
 19-78 BARANY Noise Apparatus, with hard rubber ear tip






MILTEX TUNING FORKS, Aluminum Alloy


	•
19-100	C-64 vibrations
19-102	C-128 vibrations
19-104	C-256 vibrations
19-106	C-512 vibrations
19-108	C-1024 vibrations
19-110	C-2048 vibrations
19-112	C-4096 vibrations
19-120	Set of 5 Tuning Forks, C-128, C-256, C-512, C-1024 and C-2048 in fitted molded vinyl zippered case
19-120A	Fitted molded vinyl zippered case only for tuning fork set



















HARTMAN EUSTAC	CHIAN CATHETE	RS, bulbous tip, si	lver plated
19-240	Size 00	19-246	Size 2
19-242	Size 0	19-248	Size 3
19-244	Size 1		

	BUCK EA	AR CURETTES, 61/2"(16.5 cm)	
	SHARP	,,		JNT
Straight	Angled	Size	Straight	Angled
19-260	19-270	00	19-280	19-290
19-262	19-272	0	19-282	19-292
19-264	19-274	1	19-284	19-294
19-266	19-276	2	19-286	19-296
19-268	19-278	3	19-288	19-298
19-300 19-302		urette 6" (15.2 cm), si urette 6" (15.2 cm), lä	· ·	,
19-314 19-316 19-318	BILLEAU Flexible E	ar Loop 6½" (16.5 cn ar Loop 6½" (16.5 cn ar Loop 6½" (16.5 cn	n), medium size No.	2







19-350 WILDE Ear Forceps 5" (12.7 cm), angular with serrated tips

19-356 WILDE Ear Forceps 51/4" (13.3 cm), cross action, angular with serrated tips

19-360 BLAKE Ear Forceps 41/2" (11.4 cm), very delicate, angular with serrated tips

19-370 LUCAE Ear Forceps 51/2" (14 cm), bayonet shape with serrated tips

19-371 JANSEN Ear Forceps 61/4" (15.9 cm), bayonet, extra delicate with

2 mm wide serrated tips


JANSEN (GRUENWALD) Ear Forceps 61/4" (15.9 cm), bayonet 6-197 shape with serrated tips

6-198 ADSON Bayonet Dressing Forceps 81/4" (21 cm), serrated tips

> EAR SYRINGES, All Brass, Chrome Plated, with precision ground plunger, shaft graduated in ½ (.015 liter) ounces. complete with one tapered and one short bulbous tip, easily disassembled for sterilization


Without Shield	Capacity	With Shield
19-380	2 oz. (.06 liter)	19-381
19-382	3 oz. (.09 liter)	19-383
19-384	4 oz. (.12 liter)	19-385

19-390	Short Bulbous Tip only, for Ear Syringes, chrome
19-394	Taper Tip only, for Ear Syringes, chrome
19-396	Shield only, for Ear Syringes, chrome



19-380 to 19-396






19-400	QUIRE Mechanical Finger 4" (10.2 cm), for removal of foreign bodies, spring handle
19-404	NOYES Alligator Forceps 4" (10.2 cm), spring handle
19-410	HARTMAN Ear Forceps 5" (12.7 cm), delicate serrated jaws 2.5 x 6 mm
19-414	LITTAUER Ear Forceps 494" (12.1 cm), delicate serrated jaws 2.5 x 8 m








19-429 19-430	HARTMAN-HERZFELD Cup-Shape Forceps, alligator type, 3" (7.6 cm) shaft, 2 mm cup HARTMAN-HERZFELD Cup-Shape Forceps, alligator type, 3" (7.6 cm) shaft, 3 mm cup
19-434	HARTMAN-NOYES Ear Forceps, alligator type, 3¼" (8.3 cm) shaft, serrated jaws
19-437	WEINGARTNER Ear Forceps, alligator type, 3" (7.6 cm) shaft,concave serrated jaws
19-438	STRUEMPEL Ear Forceps, alligator type, 3" (7.6 cm) shaft, fenestrated oval spoon jaws
19-439	Duckbill Ear Forceps, alligator type, 3" (7.6 cm) shaft, very flat serrated jaws, rounded tip
19-442	HOFFMAN Ear Forceps, alligator type, 3" (7.6 cm) shaft, through cutting round jaws
19-464 19-465	BRUENING Ear Snare, bayonet shape, 2% " (6.4 cm), cannula, delicate pattern Stylet only for Bruening Ear Snare
19-468	Cut Snare Wire, 3" (7.6 cm) long, fine size 1 for ear snares

For Snare Wire sizes 4 and 5 see page 294 For Snare Wire sizes 6,7,8 and 9 see page 361






19-500 19-502 19-504	SEXTON Ear Knife 7" (17.8 cm), straight SEXTON Ear Knife 7" (17.8 cm), bayonet SEXTON Ear Knife 7" (17.8 cm), angled		
19-510 19-512	POLITZER Ear Knife 7" (17.8 cm), spear blade, bayonet POLITZER Ear Knife 7" (17.8 cm), spear blade, angled		
19-520	SEXTON Ear Knife 7½" (19.1 cm), in reversible swivel handle (use in straight or angled position)	1	
19-524	POLITZER Ear Knife 71/2" (19.1 cm), in reversible swivel handle (use in straight or angled position)	19-52	0 19-524







HOUSE SUCTION-IRRIGATION TUBES, 9 cm working length


19-562	1.8 mm Suction Tube and 1.8 mm Irrigating Tube		
19-563	2.5 mm Suction Tube and 2 mm Irrigating Tube		
19-564	3 mm Suction Tube and 2 r	nm Irrigating Tube.	
19-565	4.5 mm Suction Tube and 2	mm Irrigating Tube	
	FRAZIER FERGUSON SU	•	3 , 3
19-570	6 French (2 mm)	19-575	10 French (3.3 mm)
19-572	7 French (2.3 mm)	19-576	11 French (3.6 mm)
19-573	8 French (2.6 mm)	19-577	12 French (4 mm)
19-574	9 French (3 mm)		
19-580 19-582 19-584	BARON Suction Tube, 3 French (1 mm) with finger cut-off, 7.5 cm tube BARON Suction Tube, 5 French (1.7 mm) with finger cut-off, 7.5 cm tube BARON Suction Tube, 7 French (2.3 mm) with finger cut-off, 7.5 cm tube		
19-589	ANTHONY Suction Tube 61	⁄2" (16.5 cm), with 3 m	m lumen
19-590	DAY Attic Cannula, 3" (7.6)		cate, tapered, with












11-2	JANSEN Mastoid Retractor 4¼" (10.8 cm), spread 3" (7.6 cm), 3 x 3 blunt prongs 10 mm deep, standard model
11-4	JANSEN Mastoid Retractor $4\frac{1}{2}$ " (11.4 cm), spread 3" (7.6 cm), 4 x 4 blunt prongs 20 mm deep
11-6	ALLPORT Mastoid Retractor 4" (10.2 cm), spread 2½" (6.4 cm), 4 x 4 sharp prongs
11-8	JANSEN-WAGNER Masterid Retractor 5" (12.7 cm), 5 v 5 sharp propos, self-adjusting












11-10	DAVIS Retractor 4½" (11.4 cm), 3 x 3 blunt prongs
11-12 11-13	ALM Retractor 2%" (7 cm), spread 21/4" (5.7 cm), prongs 1/6" (.32 cm) deep ALM Retractor 3%" (9.8 cm), spread 3" (7.6 cm), prongs 1/2" (1.3 cm) deep
11-14 11-16	MILTEX Cross Action Retractor 4" (10.2 cm), sharp prongs MILTEX Cross Action Retractor 4" (10.2 cm), blunt prongs






WEITLAN	IER	RETR	RACT	ORS
---------	------------	------	------	-----

Sharp		Blunt
11-600	4" (10.2 cm), 2 x 3 teeth	11-602
11-608	51/2" (14 cm), 3 x 4 teeth	11-610
11-614	61/2" (16.5 cm), 3 x 4 teeth	11-616


BECKMAN-WEITLANER RETRACTORS, with hinged blades Sharp Blunt

11-630-SH 5½" (14 cm), 3 x 4 teeth 11-630-BL 11-634-SH 6½" (16.5 cm), 3 x 4 teeth 11-634-BL









19-668 19-669

19-655 SCHUHKNECHT Postauricular Retractor 5" (12.7 cm), 3 x 3 sharp teeth 3" (7.6 cm) opening

19-660 WULLSTEIN Endaural Retractor 5" (12.7 cm), 3 x 3 sharp teeth, ebony non-glare finish

19-668 PERKINS Endaural Retractor 5¼" (13.3 cm), left serrated solid blade, right side 3 teeth

19-669 PERKINS Endaural Retractor 5¼" (13.3 cm), right serrated solid blade, left side 3 teeth













19-700 to 19-720

SPRATT (BRUN) MASTOID CURETTES 61/2" (16.5 cm), oval cups


19-700	Size 5/0	1 9- 712	Size 2	
19-702	Size 4/0	19-714	Size 3	
19-704	Size 3/0	19-716	Size 4	
19-706	Size 2/0	19-718	Size 5	1
19-708	Size 0	19-720	Size 6	,
19-710	Size 1			
19-750	ALEXANDER Mastoid Chisel 7" (17.8 cm),			
19-752	ALEXANDER Mastoid Chisel 7" (17.8 cm)			
19-754	ALEXANDER Mastoid Chisel 7" (17.8 cm),			
19-756	ALEXANDER Mastoid Chisel 7" (17.8 cm)			
19-758	ALEXANDER Mastoid Chisel 7" (17.8 cm)			
19-760	ALEXANDER Mastoid Chisel 7" (17.8 cm)	, width 14 mm		
19-770	ALEXANDER Mastoid Gouge 7" (17.8 cm),	width 4 mm		
19-772	ALEXANDER Mastoid Gouge 7" (17.8 cm)	, width 6 mm		
19-774	ALEXANDER Mastoid Gouge 7" (17.8 cm).	width 8 mm		
19-776	ALEXANDER Mastoid Gouge 7" (17.8 cm).	width 10 mm		
19-780	ALEXANDER Mastoid Gouge 7" (17.8 cm)			
19-784	BISHOP Mastoid Gouge 6%" (17.1 cm), wi	dth 3 mm		
19-785	BISHOP Mastoid Gouge 6¾" (17.1 cm), wi			
19-790	GERZOG Mallet 71/2" (19.1 cm), head 8 oz.	(227 g) lead	filled with stainless jacket	:
19-792	LUCAE Mallet 7¾" (19.7 cm), head 9 oz. (2	255 g) solid sta	ainless	






Ear Instruments Mass












19-794

19-808


MILTEX Mallet 71/2" (19.1 cm), lightweight 7.5 oz. (213 g), 19-793 stainless head 1" (2.5 cm) dia. with 2 replaceable nylon caps MILTEX Mallet 71/2" (19.1 cm), lightweight 9 oz. (255 g), stainless head 19-794 1" (2.5 cm) dia. with one stainless face and one replaceable nylon cap 19-795 Nylon Cap Replacement for Nos. 19-793 and 19-794 Mallets 19-796 MEAD MALLET 71/2" (19.1 cm), head 8 oz. (227 g), with 2 replaceable nylon caps 19-797 Nylon Cap Replacement for No. 19-796 Mallet HARTMAN Mastoid Rongeur 53/4" (14.6 cm), strong curve 19-804

BANE Mastoid Rongeur 7" (17.8 cm), curved 18-1990 KERRISON Rongeur 3" (7.6 cm) shaft, size 0, 3.5 mm bite 18-1992 KERRISON Rongeur 3" (7.6 cm) shaft, size 1, 4 mm bite 18-1994 KERRISON Rongeur 3" (7.6 cm) shaft, size 2, 4.5 mm bite 18-1996 KERRISON Rongeur 3" (7.6 cm) shaft, size 3, 5 mm bite




18-1990 to 18-1996











19-850 to 19-856

ZAUFEL-JANSEN Rongeur 7" (17.8 cm), double action, slightly curved jaws 5 mm wide 19-840

BEYER Rongeur 7" (17.8 cm), double action, slightly curved jaws 3.5 mm wide $\,$ 19-844

H	IUSKIN HONGEURS,	double action	
Straight Jaws	S	Curved Jaws	
19-848	6" (15.2 cm) long,	3 mm bite 19-849	
19-850	7¼" (18.4 cm) long,	5 mm bite 19-852	
19-854	7¼" (18.4 cm) long,	7 mm bite 19-856	
	,	cm), straight, delicate, 2.5 mm bite cm), curved, delicate, 2.5 mm bite	



19-1230 19-1232

19-1230

19-1232










19-1234



17-4801







very delicate 17-4801 Micro Friedman Rongeur 51/2" (14 cm), curved 1.3 mm wide jaws,

very delicate 19-1330 LEMPERT Elevator 71/4" (18.4 cm), straight 2 mm blade.

light pattern LEMPERT Elevator 71/4" (18.4 cm), curved 2 mm blade, 19-1334

light pattern 19-1336 LEMPERT Elevator 71/4" (18.4 cm), curved 3 mm blade,

light pattern

LEMPERT Elevator 71/4" (18.4 cm), curved 4 mm blade, 19-1340 heavy pattern

19-1350 LEMPERT Endaural Curette 81/4" (21 cm), oval cup size 5/0 19-1352 LEMPERT Endaural Curette 81/4" (21 cm), oval cup size 4/0

19-1354 LEMPERT Endaural Curette 81/4" (21 cm), oval cup size 3/0

19-1356 LEMPERT Endaural Curette 81/4" (21 cm), oval cup size 2/0

19-1820 Ear Suction Tube, with 2 interchangeable tips 20 ga, and 22 ga,




19-1350 to 19-1356



19-1820







		19-2130
19-2060	WULLSTEIN Ear Scissors, 3" (7.6 cm) shaft, straight, very delicate	
19-2120	WULLSTEIN Ear Forceps 3" (7.6 cm) shaft, oval cup jaws 1 mm wide	19-2081-B 19-2100 19-2138
19-2130	WULLSTEIN Ear Forceps 3" (7.6 cm) shaft, round cup jaws 1 mm dia.	19-2100-B
19-2081-B	Micro Alligator Ear Forceps 3½" (8.3 cm) shaft, serrated jaws 0.6 mm wide, ebony finish	19-2090-B 19-2110
19-2100	Micro Alligator Ear Forceps, 3¼" (8.3 cm) shaft, serrated jaws 0.8 mm wide	19-2110-B
19-2100-B	Micro Alligator Ear Forceps, 31/4" (8.3 cm) shaft, serrated jaws 0.8 mm wide, ebony finish	19-2150-B
19-2090-B	Micro Alligator Ear Forceps, 31/4" (8.3 cm) shaft, oval cup jaws 0.5 mm wide, ebony finish	
19-2110	Micro Alligator Ear Forceps, 31/4" (8.3 cm) shaft, oval cup jaws 0.8 mm wide	
19-2110-B	Micro Alligator Ear Forceps, 31/4" (8.3 cm) shaft, oval cup jaws 0.8 mm wide, ebony finish	
19-2138	HOUSE Strut Forceps, 2¾" (7 cm) shaft, smooth jaws 0.8 mm wide	
19-2150-B	BELLUCCI Micro Ear Scissors, 2¾" (7 cm) shaft, 4 mm blades, ebony finish	19-2081-B to 19-2150-B







20-2	vienna nasai Speculum 5¾" (14.6 cm), standard pattern, large size
20-3	VIENNA Nasai Speculum 5¾" (14.6 cm), standard pattern, medium size
20-4	VIENNA Nasal Speculum 5%" (14.6 cm), standard pattern, small size
20-6	VIENNA Lightweight Nasal Speculum 5¾" (14.6 cm), adult size
20-7	VIENNA Lightweight Nasal Speculum 5¾" (14.6 cm), child size
20-12	TIECK (HALLE) Infant Nasal Speculum 5¾" (14 cm)
20-20	SONNENSCHEIN Nasal Speculum 6" (15.2 cm)
20-24	BECKMAN Nasal Speculum 6" (15.2 cm), angled on side











20-28 to 20-32


20-36-35 to 20-36-90

20-37-45 to 20-37-75

20-28	KILLIAN Septum Speculum 6" (15.2 cm), with set screw, 2" (5.1 cm) long blades
20-29	KILLIAN Septum Speculum 6" (15.2 cm), with set screw, 21/2" (6.4 cm) long blades
20-30	KILLIAN Septum Speculum 6" (15.2 cm), with set screw, 3" (7.6 cm) long blades
20-32	KILLIAN Septum Speculum 6" (15.2 cm), with set screw, 31/2" (8.9 cm) long blades
20-36-35	COTTLE Septum Speculum 6" (15.2 cm), with set screw, narrow 35 mm long blades, ebony finish
20-36-50	COTTLE Septum Speculum 6" (15.2 cm), with set screw, narrow 50 mm long blades, ebony finish
20-36-75	COTTLE Septum Speculum 6" (15.2 cm), with set screw, narrow 75 mm long blades, ebony finish
20-36-90	COTTLE Septum Speculum 6" (15.2 cm), with set screw, narrow 90 mm long blades, ebony finish
20-37-45	COTTLE Septum Speculum 5½" (14 cm), with set screw, extra delicate, very thin 45 mm long blades, ebony finish
20-37-75	COTTLE Septum Speculum 5½" (14 cm), with set screw, extra delicate, very thin 75 mm long blades, ebony finish








19-170 19-172	UEBE Applicator 7" (17.8 cm), triangular tip UEBE Applicator 7" (17.8 cm), cross serrated tip
19-180 19-182	BROWN Applicator 6" (15.2 cm), triangular tip, copper BROWN Applicator 6" (15.2 cm), cross serrated tip, copper
19-190 19-192	LATHBURY Applicator 5½" (14 cm), triangular tip LATHBURY Applicator 5½" (14 cm), cross serrated tip
19-200 19-202	FARRELL Applicator $6\%''$ (16.5 cm), triangular tip FARRELL Applicator $6\%''$ (16.5 cm), cross serrated tip
20-100	THEOBALD (FISH) Sinus Probe 6" (15.2 cm), sterling
20-102	SKILLERN Sinus Probe 6½" (16.5 cm), graduated, sterling

WILDE Nasal Dressing Forceps 61/2" (16.5 cm), angular

20-120










6-192








6-190 6-192	CUSHING Dressing Forceps $7\%"$ (18.4 cm), bayonet, serrated tips, scraper end CUSHING Tissue Forceps $7\%"$ (18.4 cm), bayonet, 1 x 2 teeth, scraper end
19-371	JANSEN Bayonet Dressing Forceps 6¼" (15.9 cm), extra delicate with 2 mm wide serrated tips
6-197	JANSEN (GRUENWALD) Bayonet Dressing Forceps 6¼" (15.9 cm), serrated tips
6-198	ADSON Bayonet Dressing Forceps 81/4" (21 cm), serrated tips
20-124	HARTMAN Nasal Dressing Forceps 7" (17.8 cm), improved delicate pattern, serrated tips
20-126	HARTMAN Nasal Dressing Forceps 71/4" (18.4 cm), regular model, serrated tips
20-130	LITTAUER Nasal Dressing Forceps 6¾" (17.1 cm), slender mortised shanks, serrated tips










20-140 20-142 20-144	NOYES Alligator Forceps 5½" (14 cm) shaft, serrated jaws NOYES Alligator Forceps 6½" (16.5 cm) shaft, serrated jaws NOYES Alligator Forceps 5½" (14 cm) shaft, 1 x 2 teeth
20-148	KNIGHT Nasal Forceps 6¾" (17.1 cm), standard pattern, serrated tips
20-150	MUSEHOLDT Nasal Dressing Forceps 7" (17.8 cm), slender long serrated blades






20-160	KNIGHT Nasal Scissors 6%" (17.1 cm), angled on side, standard pattern
20-164	HEYMAN Nasal Scissors 6¾" (17.1 cm), angled on side, rounded blades
20-168	SEILER Turbinate Scissors 61/4" (15.9 cm), angled on side
20-177 20-178	BRUENING Nasal Snare 9" (22.9 cm), bayonet shape for unobstructed view Stylet only for Bruening Snare
20-180 20-182 20-184	KRAUSE Nasal Snare 10" (25.4 cm), standard pattern with slightly angled cannula Stylet only for Krause Snare Cannula only for Krause Snare
20-186-4* 20-186-5*	Cut Snare Wire, 4" (10.2 cm) long with formed ends, size 4, for nasal snares Cut Snare Wire, 4" (10.2 cm) long with formed ends, size 5, for nasal snares

^{*}For Snare Wire size 1 see page 277
*For Snare Wire sizes 6, 7, 8 and 9 see page 361







FREER Septum Knife 61/4" (15.9 cm), "D" model blade 20-200 FREER Submucous Knife 71/2" (19.1 cm), angled blade 20-202

BALLENGER SWIVEL KNIVES 7¾" (19.7 cm)


Straight Width Bayonet 20-216 3 mm 20-222 20-218 4 mm 20-224 20-220 5 mm

PIERCE Double End Elevator 81/2" (21.6 cm), with right and left blades 20-230

18-1968 FREER Elevator 7" (17.8 cm), double end sharp and blunt blades 5 mm wide

20-238 ROGER Septum Elevator 7¾" (19.7 cm), double ended with curved semi-sharp and blunt blades





20-242	DUNNING Elevator and Curette 81/4" (21 cm), blunt elevator and sharp small oval curette
20-246	PENNINGTON Septum Elevator 81/4" (21 cm), double end for right and left side
20-258	HAJEK-BALLENGER Septum Elevator 7% (18.4 cm), double end straight semi-sharp and blunt blades
20-270 18-1960 20-274	FREER Submucous Chisel 6½" (16.5 cm), straight, 4 mm wide FREER Submucous Chisel 6½" (16.5 cm), curved, 4 mm wide FREER Septum Gouge 6½" (15.9 cm), straight, 6 mm wide
20-278 20-280 20-282 20-284	BALLENGER V-Shaped Chisel 6¼" (15.9 cm), 2 mm wide BALLENGER V-Shaped Chisel 6¼" (15.9 cm), 4 mm wide BALLENGER V-Shaped Chisel 6¼" (15.9 cm), 6 mm wide BALLENGER V-Shaped Chisel 6¼" (15.9 cm), 8 mm wide














20-312





HALLE Ethmoid Curette 81/2" (21.6 cm), stainless oval cup 20-300 on malleable copper shaft

MOSHER Ethmoid Curette 8¾" (22.2 cm), double end oval 20-308 cups, shaft graduated in ½ inches (1.3 cm)

20-312 PRATT Ethmoid Curette 81/4" (21.6 cm), double end, fenestrated cup on angled shaft and solid cup on straight shaft

COAKLEY ANTRUM CURETTES 6%" (17.1 cm)

Set of 6, Nos. 1 through 6 20-330 20-332 No. 1 Large Size 30° angle 20-334 No. 2 Large Size 60° angle

20-336 No. 3 Large Size 100° angle 20-338 No. 4 Small Size 30° angle


No. 5 Small Size 60° angle 20-340 No. 6 Triangle 30° angle 20-342

FAULKNER Ethmoid Curette 9" (22.9 cm), double end 7 mm 20-360 and 10 mm rings, sharp on both sides

20-378 WIENER Antrum Rasp 8" (20.3 cm), No. 1 small with trocar tip WIENER Antrum Rasp 8" (20.3 cm), No. 2 medium with trocar tip 20-380 20-382 WIENER Antrum Rasp 8" (20,3 cm), No. 3 large with trocar tip








20-400	COAKLEY Antrum Trocar 7" (17.8 cm), standard curve with 3 mm dia. cannula and irrigating extension	11 %	
20-404	PIERCE Antrum Trocar 7" (17.8 cm), with 5 mm dia. cannula and tubing connection		
20-420	LICHTWICZ Antrum Trocar Needle 4½" (11.4 cm), 16 gauge with beveled stylet and crossbar, Luer lock		
20-424	VAN ALYEA Antrum Cannula 4½" (11.4 cm), trocar tip 1.8 mm, Luer lock		
20-428 20-430	RUSKIN Antrum Trocar Needle 4" (10.2 cm), straight, 16 gauge with cross bar and tubing connection, Luer lock RUSKIN Antrum Trocar Needle, 4" (10.2 cm), curved, 16 gauge with cross bar and tubing connection, Luer lock	O	d
20-436	VAN ALYEA Frontal Sinus Cannula 6½" (16.5 cm), light curve, with tubing connection, Luer lock, sterling		<u>H</u>
20-438	VAN ALYEA Frontal Sinus Cannula 6½" (16.5 cm), strong curve, with tubing connection, Luer lock, sterling	20-436	20-438



10-220






20-442	blunt tip, 1.8 mm, Luer lock VAN ALYEA Frontal Cannula 4½" (11.4 cm), short curve, blunt tip, 1.8 mm, Luer lock
20-450	BOWER Sphenoid and Antrum Cannula 6½" (16.5 cm), dia. at tip 2 mm, lateral opening, Luer lock
20-456	KILLIAN (VON EICKEN) Antrum Cannula 5¾" (14.6 cm), dia. 2.5 mm, bulbous tip, center opening, sterling
20-460	FRAZIER FERGUSON Suction Tube 9 Fr. (3 mm) straight, with finger cut-off
19-574	FRAZIER FERGUSON Suction Tube 9 Fr. (3 mm) angled, with finger cut-off
20-466	COUPLAND Suction Tube 10" (25.4 cm), with 2 interchangeable tips 1.5 and 2 mm, chrome

ANTHONY Suction Tube 10" (25.4 cm), lumen 2 mm, chrome

VAN ALYEA Antrum Cannula 41/2" (11.4 cm), long curve,












20-500	KILLIAN Septum Forceps 694 " (17.1 cm), flat serrated jaws 25×5.5 mm, shaft 314 " (8.3 cm)
20-504	WATSON-WILLIAMS Polyp Forceps 7%" (19.7 cm), small jaws 6 mm wide, shaft 41/4" (10.8 cm)
20-506	WATSON-WILLIAMS Polyp Forceps 7¾" (19.7 cm), large jaws 9 mm wide, shaft 4¼" (10.8 cm)
20-510	LEWIS Septum Forceps 8" (20.3 cm), flat serrated jaws, 15 x 7 mm, shaft 43/4" (12.1 cm)
20-522	KNIGHT Septum Forceps $6\%"$ (17.1 cm), cup jaws 15 x 5 mm, shaft $3\%"$ (8.3 cm)












20-530	
20-532	

20-530 LUC Septum Forceps 7¾" (19.7 cm), cup jaws 9 x 5.5 mm
20-532 LUC Septum Forceps 7¾" (19.7 cm), cup jaws 11 x 6.5 mm
20-544 BRUENING Septum Forceps 7½" (19.1 cm), mortised shanks, small fenestrated cup jaws 7 mm wide
20-546 BRUENING Septum Forceps 7½" (19.1 cm), mortised shanks, large fenestrated cup jaws 8 mm wide
20-552 RUBIN Septal Morselizer 8" (20.3 cm), double action, deeply serrated jaws 17 x 5 mm, supplied with one slip-on guard
20-553 Slip-on guard only












20-554 20-556	JANSEN-MIDDLETON Septum Forceps 7½" (19.1 cm), spoon-shaped cups 15 x 5 mm JANSEN-MIDDLETON Septum Forceps 7½" (19.1 cm), through-cutting 19 x 4 mm
20-570 20-572	TAKAHASHI Ethmoid Forceps, 4½" (11.4 cm) shaft, 2.5 x 10 mm bite TAKAHASHI Ethmoid Forceps, 4½" (11.4 cm) shaft, 4 x 10 mm bite
20-580	WILDE (BLAKESLEY) Ethmoid Forceps 5" (12.7 cm) shaft, pointed fenestrated cups 3 x 11 mm
20-582	WILDE (BLAKESLEY) Ethmoid Forceps 5" (12.7 cm) shaft, pointed fenestrated cups 4 x 12 mm
20-584	WILDE (BLAKESLEY) Ethmoid Forceps 5" (12.7 cm) shaft, pointed fenestrated cups 5 x 13 mm
20-594	GRUENWALD Nasal Forceps, 41/2" (11.4 cm) shaft, through-cutting 2 x 10 mm bite
20-596	GRUENWALD Nasal Forceps, 4½" (11.4 cm) shaft, through-cutting 2.5 x 12 mm bite
20-598	GRUENWALD Nasal Forceps, 4½" (11.4 cm) shaft, through-cutting 3 x 12 mm bite
20-630	STRUYKEN Turbinate Forceps, $4\%^{\prime\prime}$ (11.4 cm) shaft, through-cutting 1 x 15 mm bite












18-1990	KERRISON Rongeur, 3" (7.6 cm) shaft, size 0, 3.5 mm bite
18-1992	KERRISON Rongeur, 3" (7.6 cm) shaft, size 1, 4 mm bite
18-1994	KERRISON Rongeur, 3" (7.6 cm) shaft, size 2, 4.5 mm bite
18-1996	KERRISON Rongeur, 3" (7.6 cm) shaft, size 3, 5 mm bite
20-660 20-662	HAJEK Antrum Rongeur, 5" (12.7 cm) shaft, up-bite HAJEK Antrum Rongeur, 5" (12.7 cm) shaft, down-bite
20-002	HAVEN ARRIGIN HUNGER, 5 (12.7 Cm) Shan, down-ble










20-800	BARNHILL Adenoid Curette 8½" (21.6 cm), size 0, 11 mm blade
20-802	BARNHILL Adenoid Curette 8½" (21.6 cm), size 1, 13 mm blade
20-804	BARNHILL Adenoid Curette 8½" (21.6 cm), size 2, 15 mm blade
20-806	BARNHILL Adenoid Curette 8½" (21.6 cm), size 3, 17 mm blade
20-808	BARNHILL Adenoid Curette 8½" (21.6 cm), size 4, 19 mm blade
20-809	VOGEL Adenoid Curette 8½" (21.6 cm), infant size, 7 mm blade
20-835	LOTHROP Uvula Retractor 8" (20.3 cm) shaft, grooved blade 40 x 18 mm with 5 mm lip
20-836	LOVE Uvula Retractor, 18 mm wide blade
20-838	LOVE Uvula Retractor, 20 mm wide blade









20-890







20-850 LAFORCE Adenotome 81/2" (21.6 cm), infant size 0, cuts 7 mm width 20-852 LAFORCE Adenotome 91/2" (24.1 cm), small size 1, cuts 11 mm width 20-854 LAFORCE Adenotome 91/2" (24.1 cm), medium size 2, cuts 13 mm width 20-856 LAFORCE Adenotome 91/2" (24.1 cm), large size 3, cuts 16 mm width 20-858-0 Blade 10 mm wide, solid, for size 0 LaForce Adenotome 20-858-1 Blade 13 mm wide, perforated, for size 1 LaForce Adenotome 20-858-2 Blade 16 mm wide, perforated, for size 2 LaForce Adenotome 20-858-3 Blade 19 mm wide, perforated, for size 3 LaForce Adenotome 20-872 SHAMBAUGH Adenotome 91/2" (24.1 cm), reverse curve, small size 1, cuts 11 mm width 20-874 SHAMBAUGH Adenotome 91/2" (24.1 cm), reverse curve, medium size 2, cuts 13 mm width 20-876 SHAMBAUGH Adenotome 91/2" (24.1 cm), reverse curve, large size 3, cuts 15 mm width 20-878-1 Blade 13 mm wide, perforated, for size 1 Shambaugh Adenotome 20-878-2 Blade 16 mm wide, perforated, for size 2 Shambaugh Adenotome 20-878-3 Blade 19 mm wide, perforated, for size 3 Shambaugh Adenotome

YANKAUER Pharyngeal Speculum







20-910 20-912	RONIS Adenoid Punch, 5" (12.7 cm) shaft curved upwards, small triangular jaws RONIS Adenoid Punch, 5" (12.7 cm) shaft curved upwards, large triangular jaws
20-922	MELTZER Adenoid Punch 5" (12.7 cm) shaft, small round jaws with basket, size 1
20-924	MELTZER Adenoid Punch 5" (12.7 cm) shaft, large round jaws with basket, size 2
20-926	MELTZER Adenoid Punch 5" (12.7 cm) shaft, small triangular jaws with basket, size 1
20-928	MELTZER Adenoid Punch 5" (12.7 cm) shaft, large triangular jaws with basket, size 2


















JOSEPH NASAL KNIVES 6" (15.2 cm), Double Edge

	Pointed Tips		Blunt Tips
21-2	Small blade, straight	21-16	Small blade, curved
21-4	Large blade, straight	21-18	Large blade, curved
21-8	Large blade, curved		
21-20	JOSEPH Button End Knife 6" (1	5.2 cm), strai	ght blade 4 x 30 mm
21-40	JOSEPH Button End Knife 61/4" (1	5.9 cm), angl	ed blade 4 x 12 mm
21-24	CONVERSE Bistouri 6" (15.2 cm)	, with button e	end, curved blade 4 x 35 mm
21-28	FOMON Knife 51/2" (14 cm), doub	le edge, blunt	tip, curved
21-32	COTTLE Nasal Knife 51/2" (14 cm), straight blac	de, 4 mm wide







18-1968	FREER Septum Elevator 7" (17.8 cm), double end, sharp and blunt blades 5 mm wide	
21-60	COTTLE Septum Elevator 9" (22.9 cm), double end, graduated in cm.	DELIS 23
21-62 21-63	McKENTY Septum Elevator 51/2" (14 cm), curved at tip, 4 mm wide McKENTY Septum Elevator 51/2" (14 cm), curved at tip, 5 mm wide	š
21-64	JOSEPH Periosteal Elevator 7" (17.8 cm), blade 6.5 mm wide, strong curve	
21-68	FOMON Periosteal Elevator 61/4" (15.9 cm), slightly curved blade 4.5 mm wide, standard pattern	
21-72	COTTLE Skin Elevator 7½" (19.1 cm), slightly curved blade 9 mm wide	11
73-62	MILTEX Periosteal Elevator 6" (15.2 cm), double end, slightly curved blade 6 mm wide and straight spear point rounded blade 3mm wide	
DELS23	SELDIN Periosteal Elevator 73/4" (19.7 cm), double end, slightly curved blades	DELS

10 mm and 13 mm wide










18-428 21-102 21-103

For additional Hooks see pages 188 - 194






21-107	MILLARD Thimble Retractor, %" (1.6 cm) working length, sharp single	e hook	
21-108	COTTLE Thumb Retractor, $1\%"\ (3.8\ cm)$ working length, sharp double 8 mm wide	le hook	
18-439	McGANNON Retractor 5" (12.7 cm), fenestrated blade 2.5 mm wide		
11D-27 11-30 11-32 11-34 11-36 11-38 11-40	Flexible Neck Rake Retractor 5" (12.7 cm), 3 sharp prongs Flexible Neck Rake Retractor 6" (15.2 cm), 1 sharp prong Flexible Neck Rake Retractor 6" (15.2 cm), 2 sharp prongs Flexible Neck Rake Retractor 6" (15.2 cm), 3 sharp prongs Flexible Neck Rake Retractor 6" (15.2 cm), 1 blunt prong Flexible Neck Rake Retractor 6" (15.2 cm), 2 blunt prongs Flexible Neck Rake Retractor 6" (15.2 cm), 3 blunt prongs		
11-50 11-52 11-54 11-56 11-58 11-60	Rigid Neck Rake Retractor 6" (15.2 cm), 1 sharp prong Rigid Neck Rake Retractor 6" (15.2 cm), 2 sharp prongs Rigid Neck Rake Retractor 6" (15.2 cm), 3 sharp prongs Rigid Neck Rake Retractor 6" (15.2 cm), 1 blunt prong Rigid Neck Rake Retractor 6" (15.2 cm), 2 blunt prongs Rigid Neck Rake Retractor 6" (15.2 cm), 3 blunt prongs		950
18-100 18-102	KNAPP (ROLLET) Retractor 5¼" (13.3 cm) with 4 sharp prongs KNAPP Retractor 5¼" (13.3 cm), with 4 blunt prongs	11-50 to 11-60	18-100 18-102









- 21-144 COTTLE Double Hook 51/2" (14 cm), two sharp prongs 12 mm wide
- 21-148 FOMON Retractor 61/2" (16.5 cm), two ball end prongs 11 mm wide
- 21-150 COTTLE Retractor 51/2" (14 cm), four blunt prongs, 2 center prongs recessed
- 21-152 MILTEX Skin Hook 43/4" (12.1 cm), sharp prong 2 mm dia.
- 21-153 JOSEPH Single Hook 61/4" (15.9 cm), one sharp prong
- 21-154 JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 2 mm
- 21-156 JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 5 mm
- 21-158 JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 7 mm
- 21-160 JOSEPH Double Hook 61/4" (15.9 cm), two sharp prongs 10 mm
- 21-430 MILTEX Tissue Twist Hook 51/2" (14 cm)
- 30-950 MILTEX Tenaculum Hook 8" (20.3 cm), extra delicate pointed hook, handle graduated in cm
- 33-121 SELTZER Dermal Biopsy Hook 41/2" (11.4 cm), sharp point
- 21-163 CONVERSE Alar Retractor 4" (10.2 cm), double end, with blunt double hooks 10 mm and 14 mm wide

For additional Hooks see pages 188 - 194








18-110	DESMARRES Alar Retractor 51/2" (14 cm), 11 mm wide blade	
21-165	COTTLE Knife Guide and Retractor 8" (20.3 cm), double end, ball tip prongs 10 mm wide, and 5 mm wide slotted hook	
11-73	RAGNELL Retractor 6" (15.2 cm), double end 3 x 8 mm and 5 x 15 mm blades, delicate	
11-74 11-76	SENN Retractor 6%" (16.2 cm), double end, 3 sharp prongs SENN Retractor 6%" (16.2 cm), double end, 3 blunt prongs	
11-80 11-81	MATHIEU Retractor 61/2" (16.5 cm), double end, 3 sharp prongs MATHIEU Retractor 61/2" (16.5 cm), double end, 3 blunt prongs	
11-82	MEYERDING Skin Hook and Retractor 6¼" (15.9 cm), double end	
11-88 11-90 11-92 11-94 11-96	MEYERDING Finger Retractor 7" (17.8 cm), style A. Blade $\frac{1}{4}$ " (.64 cm) x $\frac{1}{4}$ 6" (1.7 cm) MEYERDING Finger Retractor 7" (17.8 cm), style B. Blade $\frac{1}{4}$ 6" (.48 cm) x $\frac{1}{4}$ 8" (.64 cm) MEYERDING Finger Retractor 7" (17.8 cm), style C. Blade $\frac{1}{4}$ 8" (.32 cm) x $\frac{1}{4}$ " (.64 cm) MEYERDING Finger Retractor 7" (17.8 cm), style D. Blade $\frac{1}{4}$ 2" (.24 cm) x $\frac{1}{4}$ " (1 cm) MEYERDING Finger Retractor 7" (17.8 cm), style E. Blade $\frac{1}{4}$ 4" (.64 cm) x $\frac{1}{4}$ 8" (2.2 cm)	
23-1047	LUKENS Retractor 6%" (17.1 cm), 10 mm wide curved and 7 mm wide angular blades	

For a full line of Retractors see pages 151 - 187



23-1049 JACKSON Retractor 7" (17.8 cm), double end, 2 blunt prongs %*" (2.4 cm) wide, and solid blade %" (1 cm) wide

21-176 NEIVERT Knife Guide and Retractor 6½" (16.5 cm), double end,

21-176 NEIVERT Knife Guide and Retractor 6½" (16.5 cm), double end, open blade 14 mm wide

21-180 LOVE Nerve Retractor 8½" (21.6 cm), blade 7 mm wide

21-182 CONVERSE Blade Retractor 3¾" (9.5 cm), size 1
21-183 CONVERSE Blade Retractor 3¾" (9.5 cm), size 2
21-184 CONVERSE Blade Retractor 3¾" (9.5 cm), size 3
21-185 CONVERSE Blade Retractor 3¾" (9.5 cm), size 4


21-188 AUFRICHT Nasal Retractor-Speculum 6½" (16.5 cm), solid blade 45 mm long

21-190 AUFRICHT Nasal Retractor 7" (17.8 cm), fenestrated blade 45 mm long

20-270 FREER Septum Chisel 6½" (16.5 cm), straight 4 mm wide 18-1960 FREER Septum Chisel 6½" (16.5 cm), curved 4 mm wide 20-274 FREER Septum Gouge 6½" (15.9 cm), straight 6 mm wide 20-270 18-1960 20-274

For a full line of Retractors see pages 151 - 187






SHEEHAN OSTEOTOMES 61/4" (15.9 cm), hexagon handle

21-205-2 21-205-4 21-205-6 21-205-8	2 mm wide 4 mm wide 6 mm wide 8 mm wide	21-205-10 21-205-12 21-205-15	10 mm wide 12 mm wide 15 mm wide
21-206-4 21-206-7 21-206-9 21-206-12	COTTLE Osteotome 7" (17.8 cm COTTLE Osteotome 7" (17.8 cm COTTLE Osteotome 7" (17.8 cm COTTLE Osteotome 7" (17.8 cm	i), 7 mm wide i), 9 mm wide	, rounded corners , rounded corners
21-212	COTTLE Chisel 71/4" (18.4 cm), (curved 5 mm,	with depth graduations in cm.
21-215 21-216 21-217	CINELLI Osteotome 6½" (16.5 o CINELLI Osteotome 6½" (16.5 o CINELLI Osteotome 6½" (16.5 o	m), with doub m), with doub	le guards, 12 mm cut le guards, 14 mm cut
21-218	CINELLI Osteotome 61/2" (16.5 c	m), with doub	le quards, 16 mm cut







21-227	SILVER Osteotome 7" (17.8 cm), single guard, straight
21-228	SILVER Osteotome 7" (17.8 cm), single guard, right
21-229	SILVER Osteotome 7" (17.8 cm), single guard, left
21-232	NEIVERT Osteotome 8" (20.3 cm), straight, with single guard
21-233	NEIVERT-ANDERSON Osteotome 8" (20.3 cm), curved right, with single guard
21-234	NEIVERT-ANDERSON Osteotome 8" (20.3 cm), curved left, with single guard
20-278	BALLENGER V-Shaped Chisel 61/4" (15.9 cm), 2 mm wide
20-280	BALLENGER V-Shaped Chisel 61/4" (15.9 cm), 4 mm wide
20-282	BALLENGER V-Shaped Chisel 61/4" (15.9 cm), 6 mm wide
20-284	BALLENGER V-Shaped Chisel 61/4" (15.9 cm), 8 mm wide

BALLENGER SWIVEL KNIVES 7%" (19.7 cm)

Straight	Width	Bayonet
20-216	3 mm	20-222
20-218	4 mm	20-224
20-220	5 mm	***********
















21-282 21-284

19-790 GERZOG Mallet 71/2" (19.1 cm), head 8 oz. (226 g) lead filled with stainless jacket

19-792 LUCAE Mallet 7½" (19.1 cm), head 9 oz. (255 g), solid stainless

19-793 MILTEX Mallet 7½" (19.1 cm), lightweight 7.5 oz. (213 g), stainless head 1" (2.5 cm) dia. with 2 replaceable nylon caps

19-794 MILTEX Mallet 7½" (19.1 cm), lightweight 9 oz. (255 g), stainless head 1" (2.5 cm) dia. with 1 stainless face and 1 replaceable nylon cap

19-795 Nylon Cap Replacement for Nos. 19-793 and 19-794 Mallets

19-796 MEAD Mallet 71/2" (19.1 cm), head 8 oz. (226 g) with 2 replaceable nylon caps

19-797 Nylon Cap Replacement for No. 19-796 Mallet

21-280 JOSEPH Nasal Saw 71/2" (19.1 cm), straight

21-282 JOSEPH Nasal Saw 71/2" (19.1 cm), bayonet, right


21-284 JOSEPH Nasal Saw 71/2" (19.1 cm), bayonet, left

40-2000 MILTEX Nasal Saw 6" (15.2 cm), blade 2" (5.1 cm)

















CASTROVIEJO Caliper 31/4" (8.3 cm), a combination measuring and marking 18-650 caliper, graduated 0-20 mm


THORPE Caliper 41/2" (11.4 cm), graduated 18-657 in inches and mm

18-658 JAMESON Caliper 33/4" (9.5 cm), graduated in inches and mm, with thin tips


18-660 Flexible Stainless Ruler 6" (15.2 cm) x 1/2" (1.3 cm), graduated in 1/32" (.8 mm) fractions and mm

COTTLE Cartilage Crusher, 16 mm wide 21-300 anvil, with removing forceps









21-321	MARTINI Bone Curette 51/2" (14 cm), double end, 4 mm and 5 mm round cups
21-322	WILLIGER Bone Curette 5½" (14 cm), double end, 3 mm and 4 mm oval cups
27-920 27-924	VOLKMAN Bone Curette $51\!\!/\!\!z''$ (14 cm), double end, oval cups 5 x 10 mm and 6-x 20 mm VOLKMAN Bone Curette $51\!\!/\!z''$ (14 cm), double end, oval cups 6 x 10 mm and round cup 7 mm dia.
21-330	JOSEPH Rasp 61/4" (15.9 cm), fine cross serrations
21-338	COTTLE Rasp 8" (20.3 cm), backward cutting
21-342	MALTZ Rasp 7" (17.8 cm), backward cutting

















21-343 MALTZ Rasp 81/4" (21 cm), improved double end pattern, one end forward other end backward cutting 21-346 AUFRICHT Glabella Rasp 81/4" (21 cm), curved blade, forward cutting

21-354

21-348 AUFRICHT Glabella Rasp 81/4" (21 cm), curved blade, backward cutting 21-352

LEWIS Rasp 7½" (19.1 cm), fine teeth, backward cutting 21-354 LEWIS Rasp 71/2" (19.1 cm), coarse teeth, backward cutting

21-360 GALLAHER Rasp 71/2" (19.1 cm), coarse teeth, convex 21-368

FOMON Rasp 81/4" (21 cm), double end, 4-sided, fine and coarse teeth. both in flat and convex


21-369 FOMON Rasp 81/4" (21 cm), double end, coarse teeth flat and convex 21-370 FOMON Rasp 81/4" (21 cm), double end, fine teeth flat and convex 21-379 Stainless Wire Brush 61/4", for cleaning rasps and files, Autoclavable

3-1001 MILTEX Instrument Cleaning Brushes 71/4" (18.4 cm) Stainless steel bristles. Designed to clean instrument serrations, rasps, files and burs (set of 3), Autoclavable



3-1001





6-24	Dressing Forceps 41/2" (11.4 cm), delicate pattern
6-26	Dressing Forceps 5" (12.7 cm), delicate pattern
6-28	Dressing Forceps 51/2" (14 cm), delicate pattern
6-28XL	Dressing Forceps 51/2" (14 cm), delicate pattern, lightweight, fenestrated handles
6-30	Dressing Forceps 6" (15.2 cm), delicate pattern


EYE DRESSING FORCEPS 4" (10.2 cm), SERRATED TIPS

Standard Pattern 0.8 mm wide tips			licate Pattern n wide tips
18-780	Straight	18-779	Straight
18-782	Half curved	18-781	Half curved
18-784	Full curved	18-783	Full curved

For additional Dressing & Tissue Forceps see pages 64 - 82








6-90	Tissue Forceps, 1 x 2 teeth, 41/2" (11.4 cm), delicate pattern
6-92	Tissue Forceps, 1 x 2 teeth, 5" (12.7 cm), delicate pattern
6-94	Tissue Forceps, 1 x 2 teeth, 51/2" (14 cm), delicate pattern
6-96	Tissue Forceps, 1 x 2 teeth, 6" (15.2 cm), delicate pattern

IRIS TISSUE FORCEPS 4" (10.2 cm), 1 x 2 TEETH

Standard Pattern 0.8 mm wide tips			licate Pattern n wide tips
18-786 18-787 18-788 18-790	Straight Straight (Non Magnetic) Half curved Full curved	18-785 18-789 18-791	Straight Half curved Full curved
6-100	Dressing Forceps 5" (12.7 cm),	curved fine tips	

For additional Dressing & Tissue Forceps see pages 64 - 82





6-104 6-106 6-107 6-108 6-109 6-109A 18-866	SEMKEN Dressing Forceps 5" (12.7 cm), serrated SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth, curved tips SEMKEN Tissue Forceps 5" (12.7 cm), 2 x 3 teeth SEMKEN Dressing Forceps 6" (15.2 cm), serrated SEMKEN Tissue Forceps 6" (15.2 cm), 1 x 2 teeth BISHOP-HARMON Dressing Forceps 3%" (8.6 cm), fine cross-serrated tips 0.5 mm wide BISHOP-HARMON Micro Tissue Forceps 3%" (8.6 cm), delicate		
18-868	0.3 mm 1 x 2 teeth BISHOP-HARMON Tissue Forceps 3%" (8.6 cm), standard 0.7 mm 1 x 2 teeth		
6-112 6-114	HUDSON (EWALD) Dressing Forceps 4¾" (12.1 cm), serrated HUDSON (EWALD) Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth	71,675	
6-118 6-118A 6-118XL	ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate, angled ADSON Dressing Forceps 4¾" (12.1 cm), serrated, lightweight, fenestrated handles	6-118XL	17-2510
17-2510	Micro ADSON Dressing Forceps 4¾" (12.1 cm), serrated tips, 0.5 mm wide		

For additional Dressing & Tissue Forceps see pages 64 - 82









6-120A









6-120 ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, delicate

6-120A ADSON Tissue Forceps 43/4" (12.1 cm), 1 x 2 teeth, delicate, angled

6-120XL ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, lightweight, fenestrated handles

6-121 ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth and serrated, delicate 17-2500 Micro ADSON Tissue Forceps 49/4" (12.1 cm), tips 0.5 mm wide, 1 x 2 teeth

6-122 ADSON Tissue Forceps 43/4" (12.1 cm), 2 x 3 teeth, delicate

6-123 ADSON Tissue and Suture Forceps 4¾" (12.1 cm), 1 x 2 teeth with tying platform

6-123XL ADSON Tissue Forceps 4¾" (12.1 cm), with platform 1 x 2 teeth, lightweight,


fenestrated handles



For additional Dressing & Tissue Forceps see pages 64 - 82

Carb-N-Sert Forceps with tungsten carbide inserts are shown on pages 88, 89






6-124 6-124XL	BROWN-ADSON Tissue Forceps 49/4" (12.1 cm), 7×7 side grasping teeth BROWN-ADSON Tissue Forceps 49/4" (12.1 cm), 7×7 side grasping teeth, lightweight, fenestrated handles
6-125	DeBAKEY-ADSON Tissue Forceps 4¾" (12.1 cm), atraumatic tips 1.5 mm wide
6-126	BROWN-ADSON Tissue Forceps 4%" (12.1 cm), 7 x 7 side grasping teeth, angular
6-134	MARTIN Tissue Forceps 6" (15.2 cm), 7 x 8 teeth
6-138 6-140	BROWN Tissue Forceps 6" (15.2 cm), 8 x 8 side grasping teeth BROWN Tissue Forceps 8" (20.3 cm), 8 x 8 side grasping teeth
18-850 18-854	GRAEFE Fixation Forceps 4% " (11.2 cm), standard jaws 4.5 mm wide, with fine teeth GRAEFE Fixation Forceps 4% " (11.2 cm), standard jaws 4.5 mm wide, with fine teeth, with catch
24-570XL	Debakey Tissue Forceps 6" (15.2 cm), atraumatic teeth, 1.5 mm wide tips, lightweight, fenestrated handles











18-1107	Swiss Cilia and Suture Forceps 33/4" (9.5 cm), precision fitted diagonal jaws 3 mm wide
18-1136	BERKE Ptosis Forceps 4" (10.2 cm), blades 20 mm with longitudinal serrations, with slide lock
18-1138	BERKE Ptosis Forceps 4" (10.2 cm), blades 27 mm with longitudinal serrations, with slide lock
18-1150	ERHARDT Lid Forceps 31/2" (8.9 cm), serrated jaws 9.5 mm wide
18-1160 18-1162	SNELLEN Entropium Forceps $31/2^{\circ}$ (8.9 cm), left, 35 mm wide plate SNELLEN Entropium Forceps $31/2^{\circ}$ (8.9 cm), right, 35 mm wide plate
18-1202	DESMARRES Chalazion Forceps 3% " (8.9 cm), inside ring measures 12 x 23 mm
21-440	COTTLE Columella Forceps 41/2" (11.4 cm), jaws 10 mm wide
18-1252	PRINCE Muscle Forceps 3%" (9.5 cm)







HARTMAN Mosquito Forceps 31/2" (8.9 cm), curved

PROVIDENCE Hospital Forceps 51/2" (14 cm), straight

PROVIDENCE Hospital Forceps 51/2" (14 cm), curved

7-26

7-30

7-32














7-50 7-52	Baby CRILE Forceps 5½" (14 cm), straight, extra delicate Baby CRILE Forceps 5½" (14 cm), curved, extra delicate
7-70	Baby OCHSNER Forceps 5½" (14 cm), straight, 1 x 2 teeth, extra delicate
7-72	Baby OCHSNER Forceps 51/2" (14 cm), curved, 1 x 2 teeth,

Baby OCHSNER	Forceps 51/2" ((14 cm), curve	d, 1 x 2 teeth,
extra delicate			

7-110	Baby PEAN Forceps 51/2" (14 cm), straight, extra delicate
7-112	Baby PEAN Forceps 5½" (14 cm), curved, extra delicate

7-210	Baby MIXTER Forceps 51/4" (13.3 cm), fully curved, extra delicate
-------	---

7-310	DIETHRICH Micro 8 mm x 1.2 mm	Bulldog	Clamp	1%"	(4.7	cm), s	traight	jaws

DIETHRICH Micro Bulldog Clamp 1%" (4.4 cm), angled jaws 8 mm x 1.2 mm 7-315



7-304



7-302


7-305











MII TEV DI	VOLUE OF IDUCEDA	SCISSORS 43/4" (12.1 cm)


1711127	L. 10110 00110LITT 001000110 174 (14.1 01117		AVI
Straight		Curved		///
5-270	Sharp-Sharp Points	5-276		
5-271	Sharp-Sharp Points, serrated blade	5-277	d	
5-272	Sharp-Blunt Points	5-278	8.1	M
5-273	Sharp-Blunt Points, serrated blade	5-279		M
5-274	Blunt-Blunt Points	5-280	MII	
5-275	Blunt-Blunt Points, serrated blade	5-281		1111
	Angled on Side			
5-282	Sharp-Sharp Points			
5-282A	Sharp-Sharp Points, serrated blade			
			5-280	5-282
			5-281	5-282A

For Carb-N-Sert Scissors with tungsten carbide blades see pages 55 - 59

For SuperCut Scissors see pages 60 - 63







18-1396 18-1398	Iris Scissors 3½" (8.9 cm), with 20 mm blades, straight, very delicate Iris Scissors 3½" (8.9 cm), with 20 mm blades, curved, very delicate	
5-300 5-302	Iris Scissors 4" (10.2 cm), straight, delicate Iris Scissors 4" (10.2 cm), curved, delicate	
5-304 5-306 5-308 5D-303	Iris Scissors 4½" (11.4 cm), straight Iris Scissors 4½" (11.4 cm), curved Iris Scissors 4½" (11.4 cm), angled on side Iris Scissors 4½" (11.4 cm), curved sideways	
18-1415* 18-1416*	Iris Scissors 4" (10.2 cm), straight, standard pattern, ribbon type Iris Scissors 4" (10.2 cm), curved, standard pattern, ribbon type	1
18-1417*	Iris Scissors 31/2" (8.9 cm), straight, miniature 14 mm blades,	

18-1415* 18-1417* 18-1418* 18-1416*

For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59

For SuperCut Scissors see pages 60 - 63

ribbon type


ribbon type

18-1418*


*Ribbon type scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.

Iris Scissors 31/2" (8.9 cm), curved, miniature 14 mm blades,





18-1392 18-1394	BONN Miniature Iris Scissors 3% " (8.9 cm), with 15 mm blades, straight BONN Miniature Iris Scissors 3% " (8.9 cm), with 15 mm blades, curved
18-1404 18-1405	MILTEX Iris Scissors 4" (10.2 cm), straight, heavy pattern MILTEX Iris Scissors 4" (10.2 cm), curved, heavy pattern
18-1420 18-1422 18-1426 18-1428	KNAPP Iris Scissors 4" (10.2 cm), straight, sharp-blunt points KNAPP Iris Scissors 4" (10.2 cm), straight, blunt-blunt points KNAPP Iris Scissors 4" (10.2 cm), curved, sharp-blunt points KNAPP Iris Scissors 4" (10.2 cm), curved, blunt-blunt points
5-310	Eye Scissors with Probe Points, $4\frac{1}{2}$ " (11.4 cm), angled on side

















5-312 Strabismus Scissors 4" (10.2 cm), straight 5-314 Strabismus Scissors 4" (10.2 cm), curved


Strabismus Scissors 41/2"(11.4 cm), straight, ribbon type 18-1453* Strabismus Scissors 41/2"(11.4 cm), curved, ribbon type 18-1454*

BROPHY Gum and Suture Scissors 51/2" (14 cm), straight 5-316 5-318 BROPHY Gum and Suture Scissors 51/2" (14 cm), curved BROPHY Gum and Suture Scissors 51/2" (14 cm), curved, 5D-326

one serrated blade


STEVENS TENOTOMY SCISSORS

11/4" (10.5 cm) With Short Blades		4½" (11.4 cm) With Long Blades
18-1460	straight, sharp points	18-1470
18-1462	straight, blunt points	18-1472
18-1464	curved, sharp points	18-1474
18-1466	curved, blunt points	18-1476









18-1477*	STEVENS Tenotomy Scissors 3¾" (9.5 cm), straight, blunt points
	ribbon type

18-1478* STEVENS Tenotomy Scissors 3%" (9.5 cm), curved, blunt points, ribbon type

18-1419* Utility Scissors 3¾" (9.5 cm), curved blunt blades 15 mm long, ribbon type


18-1440 CONVERSE (WILMER) Nasal Scissors 4" (10.2 cm), angled on flat

33-240 BIRO Dermal Naevus Scissors 4" (10.2 cm), angled on flat, with sharp outside edges

9-100 SPENCER Stitch Scissors 31/2" (8.9 cm), straight

9-101 SHORTBENT Stitch Scissors 3½" (8.9 cm), curved to side

9-102 LITTAUER JR. Stitch Scissors 4½" (11.4 cm), light pattern




For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59 For SuperCut Scissors see pages 60 - 63

*Ribbon type scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.







9-103	NORTHBENT Stitch Scissors 43/4" (12.1 cm),
	curved, light pattern

9-104 LITTAUER Stitch Scissors 51/2" (14 cm), standard pattern

9-108 MILTEX Stitch Scissors 41/2" (11.4 cm), angled, delicate

9-110 O'BRIEN Stitch Scissors 33/4" (9.5 cm), angled to side, sharp points

WESTCOTT Stitch Scissors 41/2" (11.4 cm), curved, 18-1486 sharp points

Eye Suture (GRADLE) Scissors 3¾" (9.5 cm), 18-1652 slightly curved, sharp points


Eye Suture Scissors 4" (10.4 cm), slightly curved, 18-1653* sharp points, ribbon type

21-536 LITTLER Suture Carrying Scissors 4%" (11.8 cm), with suture hole in blades, curved

^{18-1653*} 21-536

^{*}Ribbon type scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished.





5-180	METZENBAUM Scissors 51/2" (14 cm), curved, standard	1
5-182	METZENBAUM Scissors 7" (17.8 cm), curved, standard	
5-284	METZENBAUM Scissors 5" (12.7 cm), curved, delicate	8
5-286	METZENBAUM Scissors 6" (15.2 cm), curved, delicate	Å
5-288	METZENBAUM Scissors 7" (17.8 cm), curved, delicate	
5-290 5-291	RAGNELL Undermining Scissors 5" (12.7 cm), curved with flattened dissecting tips RAGNELL Undermining Scissors 7" (17.8 cm), curved with flattened dissecting tips	
33-242	MILTEX Undermining Scissors 4% " (11.4 cm), curved blades, with flat pointed tip	
5-320	GOLDMAN-FOX Wound Debridement Scissors 5" (12.7 cm), curved fine tips, with one serrated blade	
5D-322	GOLDMAN-FOX Wound Debridement Scissors 5" (12.7 cm), double curved, with one serrated blade	5D-

For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59
For SuperCut Scissors see pages 60 - 63

LaGRANGE Scissors 41/4" (10.8 cm), curved, one serrated blade, angled shanks

5D-324













MILTEX Blepharoplasty Scissors 41/2" (11.4 cm), curved flat tips, one serrated 21-537TC and one CARB-N-SERT blade

KLEINERT-KUTZ Tenotomy Scissors 41/6" (12.4 cm), curved, rounded 21-539 dissecting blades

5-172 REYNOLDS Dissecting Scissors 6" (15.2 cm), curved tenotomy type dissecting tips

REYNOLDS Dissecting Scissors 7" (17.8 cm), curved tenotomy type 5-176 dissecting tips

21-558 QUIMBY Scissors 5" (12.7 cm), blades 15 mm, strongly curved and pointed QUIMBY Scissors 5" (12.7 cm), blades 15 mm, strongly curved and pointed, 21D-559 one serrated blade


TAYLOR Dural Scissors 51/2" (14 cm), with probe tip 21-560



For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59

For SuperCut Scissors see pages 60 - 63






5-152 5-154	SISTRUNK Scissors 5½" (14 cm), heavy pattern, straight SISTRUNK Scissors 5½" (14 cm), heavy pattern, slightly curved		5
21-600	FOMON Saber-Back Scissors 5" (12.7 cm), slightly curved blades with semi-sharp outer edges used for dissection		
21-602	FOMON Lower Lateral Scissors 5" (12.7 cm), strong curve, dissecting blades		
21-604	FOMON Dorsal Scissors 51/2" (14 cm), light pattern, angular shanks		
21-604TC	FOMON Dorsal Scissors 51/2" (14 cm), angular shanks, CARB-N-SERT blades, one serrated blade		
21-605	MILTEX Saber-Back Rhytidectomy Scissors 6½" (16.5 cm), curved with semi-sharp outer edges used for dissection		
21-606	CASTANARES Facelift Scissors 6¼" (15.9 cm), straight, short strong blades	21-605	21-606

For Carb-N-Sert Scissors with tungsten carbide blades see pages 55 - 59 For SuperCut Scissors see pages 60 - 63







21-607 GORNEY Facelift Scissors 73/4" (19.7 cm), straight, one serrated blade

21-608 GORNEY Facelift Scissors 73/4" (19.7 cm), slightly curved, one serrated blade

21-610 COTTLE Dorsal Scissors 61/2" (16.5 cm), heavy pattern with rounded blades, angular shanks

COTTLE Bulldog Nasal Scissors 41/2" (11.4 cm), slightly curved 21-612 blunt blades

21-630 JOSEPH Nasal Scissors 5¾" (14.6 cm), straight, sharp points 21-632 JOSEPH Nasal Scissors 53/4" (14.6 cm), curved, sharp points

21-621 GORNEY Septal Scissors 73/4" (19.7 cm), double action, angled handle, serrated blades with 134" (4.4 cm) cut


21-622 CAPLAN Nasal Bone Scissors 8" (20.3 cm), double action, serrated blades

21-622 21-621

For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59

For SuperCut Scissors see pages 60 - 63





21-624	KAZANJIAN Nasal Cutting Forceps 71/2" (19.1 cm), light model		1 1
21-626 21-627	KLEINERT-KUTZ (RUSKIN) Bone Cutting Forceps 6" (15.2 cm), straight KLEINERT-KUTZ (RUSKIN) Bone Cutting Forceps 6" (15.2 cm), curved		V
21-628	RUSKIN (ROWLAND) Nasal Hump Forceps 7" (17.8 cm), double action, narrow jaws		
21-657	WALSHAM Septum Straightener 91/4" (23.5 cm), Universal, with two convex blades		
21-658	WALSHAM Septum Straightener 91/4" (23.5 cm), right, with one convex and one concave blade		ĝ
21-659	WALSHAM Septum Straightener 9½" (23.5 cm), left, with one convex and one concave blade)
21-660	ASH Septum Straightener 9" (22.9 cm), angled blades flat inside		
21-900 21-901	OBWEGESER Mandibular Awl for wire suturing 5½" (14 cm), slightly curved OBWEGESER Mandibular Awl for wire suturing 7" (17.8 cm), strong curve	U	
21-902	OBWEGESER Zygomatic Arch Awl for wire suturing 9" (22.9 cm), slightly curved	21-901	21-902




9.96

9D-135

Curved

Plastic Surgery Instruments






0-00	TIENTITY THE COLLING COLOROTO OF THE COLLING	W	a A
	WIRE CUTTING SCISSORS 4¾" (12.1 cm)		
9-120 9D-121 9-122 9D-123 9-124 9-124TC	Straight, one serrated blade Straight Curved, one serrated blade Curved Angled to side, one serrated blade Angled to side, one serrated blade Angled to side, one serrated blade, Carb-N-Sert inserts	0 /=	
	WIRE CUTTING SCISSORS 41/4" (10.8 cm)		
9D-132	Straight, one serrated blade		
9D-133	Straight	9D-132	9D-134
9D-134	Curved, one serrated blade	9D-133	9D-135

HEATH Wire Cutting Scissors 6¼" (15.9 cm), one serrated blade





For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130

8-8

8D-36

WEBSTER Needle Holder 5" (12.7 cm), delicate,

HALSEY Needle Holder 5" (12.7 cm), smooth jaws

HALSEY Needle Holder 5" (12.7 cm), serrated jaws

smooth jaws, standard pattern WEBSTER Needle Holder 4¾" (12.1 cm),

extra delicate, smooth jaws

8-6


8-7

8-8

8D-36











For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130








For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130







CRILE-WOOD Needle Holder 6" (15.2 cm) 8-50 8-51 CRILE-WOOD Needle Holder 7" (17.8 cm) GILLIES-SHEEHAN Needle Holder, with Suture Scissors 61/2" (16.5 cm), straight 8-58 GILLIES-SHEEHAN Needle Holder, with Suture Scissors 61/2" (16.5 cm), curved 8-59 CRILE Needle Holder 6" (15.2 cm) 8-94 8-96 SAROT Needle Holder 71/8" (18.1 cm), light model












MILTEX Micro Surgery Needle Holders with round handles, 0.6 mm tips

17-1000	5¼" (13.3 cm), curved jaws
17-980	71/6" (18.1 cm), straight jaws
17-1010	71/a" (18.1 cm), curved jaws
17-1020	51/4" (13.3 cm), curved jaws, with lock
17-1024	71/8" (18.1 cm), straight jaws, with loci
17-1025	71/a" (18.1 cm), curved jaws, with lock








18-1800 KALT Needle Holder 51/4" (13.3 cm), standard model, delicate jaws GREEN Needle Holder and Suture Forceps 41/4" (10.8 cm), 18-1810 angled jaws, double end 18-1820 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, without lock 18-1824 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, without lock 18-1828 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock 18-1832 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock

> For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130






18-1831 18-1833	CASTROVIEJO Needle Holder 5½" (14 cm), extra delicate jaws, straight, with lock CASTROVIEJO Needle Holder 5½" (14 cm), extra delicate jaws, curved, with lock
18-1834	PATON Needle Holder 41/2" (11.4 cm), straight, without lock
18-1837	McPHERSON Needle Holder for microsurgery 4" (10.2 cm), curved smooth tapered jaws
18-1838 18-1839 18-1840	BARRAQUER Needle Holder 5" (12.7 cm), curved, without lock, 10 mm wide hollow round handle BARRAQUER Needle Holder 5" (12.7 cm), curved, with lock, 10 mm wide hollow round handle BARRAQUER Needle Holder 434" (12.1 cm), curved, without lock, delicate tapered jaws, 8 mm wide round handle

For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130








22-2 22-4	WIEDER Tongue Depressor, small blade 28 mm wide, serrated WIEDER Tongue Depressor, large blade 36 mm wide, serrated
22-6	BOSWORTH Tongue Depressor, fenestrated blade 25 mm wide
22-10	ANDREWS Tongue Depressor, serrated blade 21 mm wide
22-16	Tongue Depressor 6¾" (17.1 cm), double end, 20 mm and 25 mm wide, cross serrated
22-200	MILTEX Lip Retractor 6" (15.2 cm), double end 25 mm and 36 mm wide
22D-202	COLUMBIA Lip Retractor 5% (14 cm), double end 35 mm and 50 mm wide
22D-209	UNIVERSITY OF MINNESOTA Cheek Retractor 5% (14 cm), double end 14 mm and 20 mm wide









22-210	Cheek Retrac	ctor 7¾" (19.7 cm), with swivel l	blade 35 mm wide
2-120 2-122	DAVIS Mouth Gag 61/4" (15.9 cm), double bite, left, without blades DAVIS Mouth Gag 61/4" (15.9 cm), double bite, right, without blades		
	BLADES FOR	R DAVIS OR McIVOR MOUTH	GAGS
W/o Ether Tube	Size No.	Size	W/Ether Tube

W/o Ether	Size	•	W/Ether
Tube	No.	Size	Tube
2-130	1	3/4" x 1" (1.9 x 2.5 cm) Tooth Plate	
2-132	2	1" x 2¾" (2.5 x 7 cm) Tongue Blade	2-132A
2-134	3	1" x 31/4" (2.5 x 8.3 cm) Tongue Blade	2-134A
2-136	4	1" x 3¾" (2.5 x 9.5 cm) Tongue Blade	2-136A
2-138	5	1" x 4" (2.5 x 10.2 cm) Tongue Blade	2-138A

RING TONGUE RETRACTOR BLADES

For use with Davis or McIvor Mouth Gags. Grooved blade holds endotracheal tube firmly in place.


	CHOOVER	nade iloido endotracinear ta	de minny in piace.	
Left	Size No.	Size		Right
2-140	2	7/8" x 21/4" (2.2 x 5.7 cm)	Tooth Plate	2-144
2-141	3	11/8" x 3" (2.9 x 7.6 cm)	Tongue Blade	2-145
2-142	4	11/4" x 31/2" (3.2 x 8.9 cm)	Tongue Blade	2-146
2-143	5	11/2" x 33/4" (3.8 x 9.5 cm)	Tongue Blade	2-147

2-150 DOYEN-JANSEN Mouth Gag 51/2" (14 cm), with grip lock









2-180	McIVOR Mouth Gag 5¾" (14.6 cm), complete with 3 tongue blades,
	one each Nos. 2-132, 2-134 and 2-136.

2-182 McIVOR Mouth Gag 53/4" (14.6 cm), frame only

2-160 JENNINGS Mouth Gag, adult size, 6" (15.2 cm) wide 2-162 JENNINGS Mouth Gag, child size, 43/4" (12.1 cm) wide 2-164 JENNINGS Mouth Gag, infant size, 3¾" (9.5 cm) wide

2-154 DENHARDT Mouth Gag 5¾" (14.6 cm), improved model, with grip lock


> MOLT Mouth Gags with silicone guards

63-10 Pediatric 4" (10.2 cm) 63-15 Child 41/2" (11.4 cm) 63-20

Adult 51/2" (14 cm)









2-20	MAUNDER Oral Screw, white plastic		
2-50	COLLIN Tongue Seizing Forceps 6¾" (17.1 cm), jaws 25 mm wide		//
2-60 2-62	YOUNG Tongue Seizing Forceps 61/2" (16.5 cm), with soft rubber jaws Rubber Jaws for YOUNG Tongue Forceps		
2-100	PYNCHON Suction Tube 11" (27.9 cm), tip with large openings, with removable tubing connection, chrome		To the last of the
2-104 2-104SS	YANKAUER Suction Tube 11%" (29.8 cm), with removable tip and tubing connection, chrome YANKAUER Suction Tube 11%" (29.8 cm), with removable tip and tubing connection		
2-106 2-106SS	Tip for YANKAUER Suction Tube, chrome Tip for YANKAUER Suction Tube	2-106 2-106SS	
2-109	YANKAUER Pediatric Suction Tube 8½" (21.6 cm), with removable tip, delicate pattern		2-109











18-1968







2-110 ANDREWS-PYNCHON Suction Tube 10" (25.4 cm), delicate pattern, chrome 2-110SS ANDREWS-PYNCHON Suction Tube 10" (25.4 cm), delicate pattern 22-300 BLAIR Cleft Palate Elevator 7¾" (19.7 cm), L-Shape, blunt BLAIR Cleft Palate Elevator 7¾" (19.7 cm), right angle, semi-sharp tip 22-302 22-304 BLAIR Cleft Palate Elevator 8" (20.3 cm), curved blade, blunt FREER Elevator 7" (17.8 cm), double end, sharp and blunt blades 5 mm wide 18-1968 BROPHY Scissors 51/2" (14 cm), straight, sharp points 5-316 BROPHY Scissors 51/2" (14 cm), curved, sharp points 5-318 5D-326 BROPHY Scissors 51/2" (14 cm), curved, one serrated blade 9-106 LONG Oral Surgery Stitch Scissors 6" (15.2 cm)














22D-420

22-484

22D-422

22D-424

22-480 to 22-482

22D-420 BLUMENTHAL Oral Surgery Rongeur 6" (15.2 cm), beaks at 45° angle
22D-422 Mini-BLUMENTHAL Oral Surgery Rongeur 4½" (11.4 cm), beaks at 30° angle


22D-424 Mini-BLUMENTHAL Oral Surgery Rongeur 41/2" (11.4 cm), beaks at 45° angle

 22-480
 BLUMENTHAL Oral Surgery Rongeur 6" (15.2 cm), beaks at 30° angle

 22-481
 BLUMENTHAL Oral Surgery Rongeur 7" (17.8 cm), beaks at 30° angle

 22-482
 BLUMENTHAL Oral Surgery Rongeur 6" (15.2 cm), beaks at 90° angle


MEAD Oral Surgery Rongeur 61/2" (16.5 cm), No. 1A pattern, slightly angled beaks



22-484














17-4800 17-4801

17-4800

Micro-FRIEDMAN Rongeur 5½" (14 cm), straight 1.3 mm wide jaws, very delicate

Micro-FRIEDMAN Rongeur 51/2" (14 cm), 17-4801 curved 1.3 mm wide jaws, very delicate

Micro-FRIEDMAN Rongeur 5½" (14 cm), 90° angled 1.3 mm wide jaws, very delicate 17D-4802

HARTMAN Rongeur 5¾" (14.6 cm), 19-804

strong curve

BANE Rongeur 7" (17.8 cm), slight curve 19-808


BEYER Rongeur, double action, 7" (17.8 cm), slight curve 19-844

FRIEDMAN Oral Surgery Rongeur 5½" (14 cm), very delicate, 22-494 slightly curved beaks









22-486 22-488 22-490	Oral Surgery Rongeur 5½" (14 cm), No. 4 pattern, slightly curved beaks Oral Surgery Rongeur 6½" (16.5 cm), No. 4A pattern, slightly curved beaks Oral Surgery Rongeur 6¾" (17.1 cm), No. 5 pattern, side cutting, slightly curved		
22D-426	Oral Surgery Rongeur 6½" (16.5 cm), No. 5A pattern, side cutting, slightly curved		'
22D-428	Oral Surgery Rongeur 6" (15.2 cm), No. 5S pattern, side cutting, slightly curved		
22D-432	Oral Surgery Rongeur 6½" (16.5 cm), No. 1 pattern, strong curved beaks		
22-492	Oral Surgery Rongeur $6\%^{\prime\prime}$ (16.5 cm), No. 6A pattern, end cutting, curved beaks		
22-500TC 22-502TC	Wire Twister 6" (15.2 cm), with Carb-N-Sert jaws Wire Twister 8" (20.3 cm), with Carb-N-Sert jaws	22-492	22-500TC 22-502TC



Tonsil Instruments












WHITE Tonsil Seizing Forceps 71/2" (19.1 cm), medium curve, open ring 22-550 22-552 WHITE Tonsil Seizing Forceps 91/2" (24.1 cm), medium curve, open ring COLVER Tonsil Seizing Forceps 71/2" (19.1 cm), straight, open ring 22-556 COLVER Tonsil Seizing Forceps 71/2" (19.1 cm), medium curve, open ring 22-558 22-560 BALLENGER Tonsil Seizing Forceps 81/2" (21.6 cm), slight curve, 3 x 3 prongs, open ring 22-586 TIVNEN Tonsil Seizing Forceps 8" (20.3 cm), angled 3 x 3 prongs, open ring

Tonsil Instruments











7-604 7-606

7-614 7-616

7-604	BALLENGER Sponge Forceps 7" (17.8 cm), straight, serrated loop jaws
7-606	BALLENGER Sponge Forceps 7" (17.8 cm), straight, smooth loop jaws
7-614	BALLENGER Sponge Forceps 7" (17.8 cm), curved, serrated loop jaws
7-616	BALLENGER Sponge Forceps 7" (17.8 cm), curved, smooth loop jaws
7-190	MICULICZ Tonsil Hemostatic Forceps $7 \ensuremath{\mathcal{V}}\xspace$ (19.1 cm), slight curve, 1 x 2 teeth
22-660	SCHNIDT Tonsil Hemostatic Forceps 7½" (19.1 cm), slight curve, open ring
22-662	SCHNIDT Tonsil Hemostatic Forceps 7½" (19.1 cm), slight curve, closed rings







22-666

Tonsil Instruments








2-110

2-110SS

2-109

Tonsil Instruments






9-154 9-156	DUPUYS-WEISS Ligature Needle 8½" (21.6 cm), right DUPUYS-WEISS Ligature Needle 8½" (21.6 cm), left
22-802	CANFIELD Tonsil Knife 81/2" (21.6 cm), curved
22-808	FISHER Tonsil Knife and Dissector 8½" (21.6 cm), sharp tip
22-810	DEAN Tonsil Knife 71/2" (19.1 cm), straight, double edge
22-854	HURD Tonsil Dissector and Pillar Retractor 9" (22.9 cm), standard pattern, 9 x 25 mm blade
22-856	MORRISON-HURD Tonsil Dissector and Pillar Retractor 9" (22.9 cm), flat handle, 8 x 25 mm serrated blade
22-890	COLVER Tonsil Examining Retractor 8" (20.3 cm), blunt tip
22-900	LEWIS Tonsil Screw 7¾" (19.7 cm), two sharp points, standard small size



Tonsil Instruments





5-182 METZENBAUM Tonsil Scissors 7" (17.8 cm), curved, standard pattern 5-288 METZENBAUM Tonsil Scissors 7" (17.8 cm), curved, delicate pattern 22-962 LILLY Tonsil Scissors 71/2" (19.1 cm), strong curve BOETTCHER Tonsil Scissors 71/4" (18.4 cm), curved 5-256 DEAN Tonsil Scissors 6¾" (17.1 cm), blades angled on flat, one serrated blade 5-264

For Carb-N-Sert and SuperCut Metzenbaum Scissors see pages 57 and 61.

Tonsil Instruments I









22-1000 EVE Tonsil Snare 11" (27.9 cm), standard model,

without ratchet

22-1016

22-1002 EVE Tonsil Snare 11" (27.9 cm), with ratchet

22-1010 BECK-SCHENCK Tonsil Snare 10" (25.4 cm), complete with

one straight and one vedder tip Medium Vedder Tip, opening 23 x 19 mm


22-1020 Straight Tip 21/2" (6.4 cm) for BECK-SCHENCK Tonsil Snare

22-1030 SAGE Automatic Ratchet Tonsil Snare, improved pattern, with removable straight tip

22-1030



Tonsil Instruments



TYDING Tonsil Snare, complete with straight tip 22-1041

22-1042 Straight Tip for TYDING Snare

HARTMAN Tonsil Punch, 41/2" (11.4 cm) shaft, 9 mm 22-1052 round cut, with basket

*TONSIL SNARE WIRE, pre-cut 4" (10.2 cm), bent ends


22-1060-8 22-1060-6 Size 6 size 8 22-1060-9 22-1060-7 Size 7 size 9

BROWN Tonsil Snare, complete with one straight and 22-1070

one vedder tip

Straight tip for BROWN Snare 22-1072


Vedder Tip 23 x 19 mm opening, for BROWN Snare 22-1074



22-1070 to 22-1074

Laryngeal Instruments





LARYNGEAL MIRRORS, Boilable, Thin Model


Complete with octagon			Mirror only with
threaded handle	Size No.	Diameter	threaded stem
23-1-000	000	10 mm	23-29-000
23-2-00	00	12 mm	23-30-00
23-4-0	0	14 mm	23-320
23-6-1	1	16 mm	23-34-1
23-8-2	2	18 mm	23-36-2
23-10-3	3	20 mm	23-38-3
23-12-4	4	22 mm	23-40-4
23-14-5	5	24 mm	23-42-5
23-16-6	6	26 mm	23-44-6
23-18-7	7	28 mm	23-46-7

23-60 Octagon Threaded Handle for Laryngeal Mirrors
23-62 UNIVERSAL Chuck Handle for Laryngeal Mirrors, round, chrome




Laryngeal Instruments







23-232










Laryngeal Instruments






	UNIVERSAL LARYNGEAL TIPS for use with Nos. 23-399 to 23-410 Cannula and No. 23-420 Handle
23-282	BRUENING Cup Jaws, 3 mm round, straight
23-286	Cup Jaws, 3 mm round, angled-up-
23-320 23-322 23-324 23-326	KRAUSE Square Basket Jaws, through-cutting, small, straight KRAUSE Square Basket Jaws, through-cutting, large, straight KRAUSE Oval Basket Jaws, through-cutting, small, straight KRAUSE Oval Basket Jaws, through-cutting, large, straight
23-399 23-400 23-402 23-410	UNIVERSAL CANNULA for use with Laryngeal Tips and Universal Handles 15 cm straight 20 cm straight 30 cm straight 20 cm curved
23-420	UNIVERSAL Handle for use with all Cannula and Tips



Bronchial and Esophageal Instruments






BRONCHOSCOPIC FORCEPS Double Articulated Jaws

23-600	Serrated Grasping Jaws, 30 cm shaft, 4 mm pass diameter
23-622	KAHLER Alligator Grasping Jaws, 40 cm shaft, 4 mm pass diameter
23-640 23-642 23-644	Serrated Jaws 3 x 12 mm, 30 cm shaft, 4 mm pass diameter Serrated Jaws 3 x 12 mm, 40 cm shaft, 4 mm pass diameter Serrated Jaws 3 x 12 mm, 50 cm shaft, 4 mm pass diameter
23-650	Papilloma Grasping Jaws 7 mm wide, with fine serrations, 30 cm shaft
23-680 23-690 23-720	JACKSON Cup Jaws 2 mm dia. straight, 30 cm shaft JACKSON Cup Jaws 4 mm dia. straight, 30 cm shaft JACKSON Cup Jaws 4 mm dia. angular, 30 cm shaft
23-804	Serrated Jaws 2 x 8 mm, Miniature size, 20 cm shaft, 3 mm pass diameter

Bronchial and Esophageal Instruments






JACKSON ENDOSCOPIC FORCEPS, Standard Pattern

23-755 Serrated Grasping Jaws 4 x 20 mm, 30 cm shaft, passes 5 mm tube 23-756 Serrated Grasping Jaws 4 x 20 mm, 30 cm shaft, passes 5 mm tube








23-1047	LUKENS Trachea Retractor 6% " (17.1 cm), 10 mm wide curved and 7 mm wide angular blades
23-1048	HUPP Trachea Hook 61/2" (16.5 cm), sharp point
23-1049	JACKSON Trachea Retractor 7" (17.8 cm), double end, 2 blunt prongs $\%$ (2.4 cm) wide, and solid blade $\%$ (1 cm) wide
23-1050	JACKSON Trachea Tenaculum 5¾" (14.6 cm), sharp point
23-1051	NEWS Tracheotomy Hook 6" (15.2 cm), sharp point






11D-27	Trachea Re	tractor 5" (12.7 cm), w	ith Flexible Nec	k, 3 sharp prongs		
TRA	TRACHEA RETRACTORS 6" (15.2 cm), with Flexible Neck					0
11-30	Sharp	1 Prong	11-36	Blunt		
11-32	Sharp	2 Prongs	11-38	Blunt	ľ	
11-34	Sharp	3 Prongs	11-40	Blunt		
ŢI	RACHEA RETF	RACTORS 6" (15.2 cm), with Rigid Ne	ck		
11-50	Sharp	1 Prong	11-56	Blunt		
11-52	Sharp	2 Prongs	11-58	Blunt		
11-54	Sharp	3 Prongs	11-60	Blunt		
11-73	RAGNELL blades, deli	Retractor 6" (15.2 cm) cate	1			
11-74 11-76	SENN Retractor 6%" (16.2 cm), sharp, double end SENN Retractor 6%" (16.2 cm), blunt, double end					e de la companya de l
11-80 11-81	, , , , , , , , , , , , , , , , , , , ,				11-80	11-82
11-82	MEYERDIN	IG Skin Hook and Reti	11-81			













23-1060 JACKSON Trachea Bistouri 61/2" (16.5 cm), curved blade 2" (5.1 cm), long, with probe tip 23-1402 MOSHER Life Saving Tube 8" (20.3 cm) size 20 Fr. (6.6 mm) 23-1420 TROUSSEAU Trachea Dilator 51/4" (13.3 cm), standard size, jaws 3 cm long x 5 mm wide at tip TROUSSEAU Trachea Dilator 41/2" (11.4 cm), neo natal size, jaws 2 cm long x 3 mm wide at tip 23-1421 LaBORDE Trachea Dilator 51/2" (14 cm), three bladed 23-1430



			Cat. No. 23-1139SS	Cat. No. 23-1243SS	Cat. No. 23-1350SS	Cat. No. 23-1360ST
	Actual	O.D.	to	to	to	and
Size	Size	mm	23-1219SS	23-1251SS	23-1352SS	23-1361ST
00	0	4	1-5/16" (33mm)			
0	0	5	1-1/2" (38mm)			
1		5.5	1-3/4" (44mm)			
2		6	2" (51mm)	1-3/4" (44mm)		
3		7	2-9/32" (58mm)	1-7/8" (48mm)		
4		8	2-7/16" (62mm)	2-1/16" (52mm)		
_5		9	2-5/8" (68mm)	2-3/16" (56mm)		
6		10	2-11/16" (69mm)	2-3/16" (56mm)		
.7		11	2-11/16" (69mm)	2-3/16" (56mm)		2-5/32" (55mm)
8		12	2-11/16" (69mm)	2-3/16" (56mm)	1-5/8" (42mm)	2-5/32" (55mm)
0		40	0.4414011 (00)	0.0402 (50)		
9		13	2-11/16" (69mm)	2-3/16" (56mm)		
10		14	2-11/16" (69mm)	2-3/16" (56mm)	1-5/8" (42mm)	
12		16			1-5/8" (42mm)	
				anath in inches as	A 1919 1	ured on inside curve

370

Length in inches and millimeters measured on inside curve







23-1139SS




JACKSON TRACHEA TUB	S, Improved Model with Rotating Lock
---------------------	--------------------------------------


23-114055	Size U	23-115255	Size 6				
23-1142SS	Size 1	23-1154SS	Size 7				
23-1144SS	Size 2	23-1156SS	Size 8				
23-1146SS	Size 3	23-1158SS	Size 9				
23-1148SS	Size 4	23-1159SS	Size 10				
JACKSON TRACHEA TUBES, Regular Length, Original Model with Wing Lock							
23-1199SS	Size 00	23-1210SS	Size 5				
23-1200SS	Size 0	23-1212SS	Size 6				
23-1202SS	Size 1	23-1214SS	Size 7				
23-1204SS	Size 2	23-1216SS	Size 8				
23-1206SS	Size 3	23-1218SS	Size 9				
23-1208SS	Size 4	23-1219SS	Size 10				
IACKCON TE	MOUTE A TUREO	Chart Langth Original Madal with	Mina Look				

JACKSON	INAUNEA IU	DEO	, anon Lengin, Ongina Model with	wing Lo	CK
23-124355	Size	2	23-1248SS	Size	7
23-1244SS	Size	3	23-1249SS	Size	8
23-1245SS	Size	4	23-1250SS	Size	9
23-1246SS	Size	5	23-1251SS	Size	10
23-1247SS	Size	6			









23-1350SS MARTIN Laryngectomy Tube, Size 8

23-1351SS MARTIN Laryngectomy Tube, Size 10

23-1352SS MARTIN Laryngectomy Tube, Size 12

23-1360ST LEWIS Laryngectomy Tube, Size 7, tapered, sterling

23-1361ST LEWIS Laryngectomy Tube, Size 8, tapered, sterling

CLEANING BRUSHES for Trachea and Laryngectomy Tubes Nylon Bristles with Stainless Wire Handles


23-1600 Small Brush for tube sizes 2 to 5, overall length 6" (15.2 cm)
23-1602 Large Brush for tube sizes 6 to 12, overall length 8" (20.3 cm)

23-1600 23-1602



Atraumatic Forceps


Enlarged cross sections of our basic patterns:



Manufactured in the Miltex tradition, these high quality instruments feature:

- Individual hand finishing assures soft spring action of the handles and perfect balance
- A scientific milling process produces an exclusive design on the atraumatic clamping surfaces, providing fine serrations that are of non-traumatic quality, but that will grip the tissue firmly
- Basic jaw types include the DeBakey type and the Cooley type





DeBAKEY THORACIC TISSUE FORCEPS,

Lio min moo apo				
24-550	6" (15.2 cm)	24-554	91/2" (24.1 cm)	
24-552	7½" (19.7 cm)	24-556	12" (30.5 cm)	
	DeBAKEY VASCULA	B TISSUE FORCE	FPS	

DeBAKEY VASCULAR TISSUE FORCEPS, 1.5 mm wide tips

 24-570
 6" (15.2 cm)

 24-570XL
 6" (15.2 cm), lightweight with fenestrated handles

 24-572
 8" (20.3 cm)

 24-574
 9½" (24.1 cm)


 6-125
 DeBAKEY-ADSON Tissue Forceps 4¾" (12.1 cm), 1.5 mm wide tips

24-580 COOLEY Vascular Tissue Forceps 6" (15.2 cm), tips 2 mm wide


24-582 COOLEY Vascular Tissue Forceps 8" (20.3 cm), tips 2 mm wide

DeBAKEY BULLDOG CLAMPS, Cross Action


Size	Tips	Curved
7.5 cm	2 cm	24-650
8.5 cm	3 cm	24-652
10.5 cm	4.5 cm	24-654
12 cm	6 cm	24-656
	7.5 cm 8.5 cm 10.5 cm	7.5 cm 2 cm 8.5 cm 3 cm 10.5 cm 4.5 cm








24-700 to 24-704




GLOVER BULLDOG CLAMPS, with Spring Tension Adjusting Screw

Straight	Size	Jaws	Curved
24-700	7 cm	2 cm	24-750
24-702	8 cm	2.5 cm	24-752
24-704	9 cm	4 cm	24-754
	11 cm	6 cm	24-756










DeBAKEY RING HANDLED BULLDOG CLAMPS

24-800	Straight 5" (12.7 cm)
24-802	Curved 4¾" (12.1 cm)
24-804	90° Angle 4" (10.2 cm)
24-806	45° Angle 4%" (12.4 cm)
24-808	S-Curved 43/4" (12.1 cm)







LEE Right Angle Bronchus and General Purpose Vascular Clamp 9" (22.9 cm) 24-900

GLOVER Patent Ductus Clamp 8" (20.3 cm), straight,

3 cm jaws

24-1000

24-1002


GLOVER Patent Ductus Clamp 8" (20.3 cm), angular,

3 cm jaws


GLOVER Coarctation Clamp 8¾" (22.2 cm), straight, 24-1100

GLOVER Coarctation Clamp 8¾" (22.2 cm), angular, 24-1102


6 cm jaws








24-1150	DeBAKEY Patent Ductus and Peripheral Vascular Clamp 8" (20.3 cm) straight, 3 cm jaws		
24-1152	DeBAKEY Patent Ductus and Peripheral Vascular Clamp 8" (20.3 cm) angular, 3 cm jaws		
24-1180	DeBAKEY Peripheral Vascular Clamp 7" (17.8 cm), angled jaws 6.5 cm		
24-1190	DeBAKEY Peripheral Vascular Clamp 8" (20.3 cm), S-Shaped		
24-1200	DeBAKEY Peripheral Vascular Clamp 6" (15.2 cm), curved jaws		
		24-1190	24-1200








24-1400DeBAKEY Vascular Clamp for Aortic Aneurism 10"(25.4 cm), slightly curved jaws 7.5 cm24-1404DeBAKEY Vascular Clamp for Aortic Aneurism 10½" (26.7 cm), slightly curved jaws 8.5 cm24-1408DeBAKEY Vascular Clamp for Aortic Aneurism 12¾" (31 cm), slightly curved jaws 10 cm24-1412DeBAKEY Vascular Clamp for Aortic Aneurism 12¾" (32.4 cm), slightly curved jaws 12 cm





24-1500

DeBAKEY Vascular Clamp for Aortic Aneurism 9½" (24.1 cm), strongly curved jaws


DeBAKEY MULTIPURPOSE CLAMPS, OBTUSE ANGLE, 60°

24-1600 9½" (24.1 cm) long, 4 cm jaws

24-1602 8¼" (21 cm) long, 5 cm jaws


24-1604 91/2" (24.1 cm) long, 6.5 cm jaws

24-1606 12" (30.5 cm) long, 10 cm jaws



Atraumatic Forceps Atraumatic






DeBAKEY MULTIPURPOSE CLAMPS, RIGHT ANGLE, 90°


 24-1650
 9" (22.9 cm) long,
 4 cm jaws

 24-1654
 8" (20.3 cm) long,
 6.5 cm jaws


 24-1656
 10½" (26.7 cm) long,
 10 cm jaws











 24-2100
 SATINSKY-DeBAKEY Vascular Clamp for Vena Cava 9%" (24.8 cm), small


 24-2110
 SATINSKY-DeBAKEY Vascular Clamp for Vena Cava 9%" (24.8 cm), medium

 24-2120
 SATINSKY-DeBAKEY Vascular Clamp for Vena Cava 9%" (24.8 cm), large


 24-2500
 COOLEY Pediatric Vascular Clamp for Anastomosis 6½" (16.5 cm), with 5mm graduations on jaws









24-2700







DeBAKEY Pediatric Multi-Purpose Clamp 6½" (16.5 cm), angled 30° DeBAKEY Pediatric Multi-Purpose Clamp 6" (15.2 cm), angled 60° 24-2710 DeBAKEY Pediatric Multi-Purpose Clamp 24-2720 51/2" (14 cm), angled 90° 24-2960 DeBAKEY-DERRA Vascular Clamp for Anastomosis 61/4" (15.9 cm), small jaws DeBAKEY-DERRA Vascular Clamp for Anastomosis 6%" (17.1 cm), medium jaws 24-2961 DeBAKEY-DERRA Vascular Clamp for Anastomosis 7" (17.8 cm), large jaws 24-2962












COOLEY-DERRA Pediatric Vascular Clamps for Anastomosis 61/2" (16.5 cm)

24-3000 small 24-3010 medium 24-3020 large


COOLEY Pediatric Vascular Clamps 5%" (14 cm) 5 mm calibrations on outer sides of jaws

24-3200 angled jaws 25 mm, straight handles 24-3210 straight jaws 25 mm, angled handles 24-3220 straight jaws 25 mm, straight handles












COOLEY CLAMPS Pediatric Vascular Clamps, 5 mm calibrations on outer sides of jaws 4¾" (12.1 cm)

24-3230 24-3240 24-3250	right angle jaws 17 mm, straight handles right angle jaws 25 mm, straight handles right angle jaws 35 mm, straight handles
24-3260	COOLEY-SATINSKY Clamp 51/2" (14 cm)
24-3270	COOLEY Spoon Shape Clamp 51/2" (14 cm







25-14







DOYEN RIB ELEVATORS AND RASPATORIES 7" (17.8 cm)

25-2 25-4 25-6	Adult Size, right Adult Size, left Child Size, right
25-14	MATSON-ALEXANDER Rib Elevator and Stripper 8¾" (22.2 cm
25-18	MATSON Rib Elevator and Stripper 8¾" (22.2 cm)
25-22	ALEXANDER-FARABEUF Periosteotome 81/4" (21 cm), standard size, 12 mm and 24 mm wide ends
25-24	ALEXANDER-FARABEUF Periosteotome 6" (15.2 cm), pediatric size, 10 mm and 14 mm wide ends








25-106

25-130

25-132






25-100 to 25-106

25-28	SEDILLOT Periosteal Elevator 7½" (19.1 cm), blade 17 mm wide
25-98	FINOCHIETTO-DeBAKEY Infant Rib Spreader, blade 19 mm deep x 21 mm wide, 85 mm spread, light construction
	FINOCHIETTO RIB SPREADERS, Original Model
25-100	Child Size, blades 38 mm deep x 45 mm wide, 6" (15.2 cm) spread, straight arms
25-102	Standard Size, blades 48 mm deep x 65 mm wide, 8" (20.3 cm) spread, straight arms
25-104	Standard Size, blades 48 mm deep x 65 mm wide, 10" (25.4 cm) spread, straight arms

Standard Size, blades 48 mm deep x 65 mm wide,

Standard Size, blades 48 mm deep x 65 mm wide,

Standard Size, blades 48 mm deep x 65 mm wide,




10" (25.4 cm) spread, curved arms Standard Size, blades 48 mm deep x 65 mm wide, 25-134


12" (30.5 cm) spread, straight arms


8" (20.3 cm) spread, curved arms

12" (30.5 cm) spread, curved arms









BURFORD-FINOCHIETTO RIB SPREADERS

with 2 pairs interchangeable blades 48 mm deep x 65 mm wide, and 72 mm deep x 65 mm wide

25-140 8" (20.3 cm) Spread

25-142 10" (25.4 cm) Spread

25-144 12" (30.5 cm) Spread

BURFORD-FINOCHIETTO RIB SPREADERS, LIGHT WEIGHT ALUMINUM with 2 pairs interchangeable blades 43 mm deep x 63 mm wide. and 65 mm deep x 63 mm wide

25-150 8" (20.3 cm) Spread

25-152 10" (25.4 cm) Spread

25-154 12" (30.5 cm) Spread

25-158


HAIGHT Baby Rib Spreader, 6.5 cm arms, 32 x 32 mm blades, 11 cm spread


25-160

HAIGHT Baby Rib Spreader, 9 cm arms, 32 x 32 mm blades, 11 cm spread









11-676

HARKEN RIB SPREADERS and SCAPULA RETRACTORS, blades 45 mm deep x 65 mm wide, with detachable scapula blades 43 mm deep x 50 mm wide

25-170

8" (20.3 cm) Spread

25-172

10" (25.4 cm) Spread

11-672


BADGLEY Rib Spreader, 4 x 4 blunt teeth, 2" (5.1 cm) x 13/4" (4.4 cm), 61/2" (16.5 cm) spread

11-676

TUFFIER Rib Spreader, 2" (5.1 cm) x 13/4" (4.4 cm), grooved blades, 61/2" (16.5 cm) spread

391






25-200	BAILEY Rib Contractor, Standard model
25-202	BAILEY-GIBBON Rib Contractor, with 9 mm longer arms than standard model
25-204	BAILEY Pediatric Rib Contractor, light-weight construction
11-314	HARRINGTON Splanchnik Retractor 12" (30.5 cm), large blade 21/2" (6.4 cm) wide
11-316	HARRINGTON Splanchnik Retractor 12" (30.5 cm), small blade 11/2" (3.8 cm) wide
11-318	HARRINGTON Splanchnik Retractor 91/2" (24.1 cm), extra small blade 1" (2.5 cm) wide









11-246 to 11-249







25-272

25-274

25-274

SAUERBRUCH Retractor 9" (22.9 cm), Size 2, blade 7/8" (2.2 cm) x 21/4" (5.7 cm) SAUERBRUCH Retractor 9" (22.9 cm), Size 3, blade 3/4" (1.9 cm) x 27/4" (7.3 cm)

RIC	HARDSON RETHACTORS 91/2" (24	i.1 cm)
With Loop Handle	Blade Size	With Hollow Grip Handle
11-246	1" (2.5 cm) x ¾" (1.9 cm)	11-250
11-247	11/4" (3.2 cm) x 1" (2.5 cm)	11-252
11-248	1½" (3.8 cm) x 1½" (3.8 cm)	11-254
11-249	2" (5.1 cm) x ¾" (1.9 cm)	11-256
11-349	OLLIER Retractor 81/2" (21.6 cm), b	lunt, 4 prongs 21/4" (5.7 cm) x 13/4" (4.4 cm)
11-350	ISRAEL Retractor 9" (22.9 cm), blu	nt, 4 prongs 11/2" (3.8 cm) x 15/8" (4.1 cm)
11-352	ISRAEL Retractor 9" (22.9 cm), blu	nt, 6 prongs 21/4" (5.7 cm) x 11/4" (4.4 cm)
11-590	HIBBS Retractor 91/4" (23.5 cm), bla	ade 1" (2.5 cm) x 3" (9.6 cm) deep

For additional Retractors see pages 151 - 187














25-302 ALLISON Lung Retractor 1294" (32.4 cm), standard size 2" (5.1 cm) wide blade
25-306 DAVIDSON Scapula Retractor, blade 31½" (8.9 cm) wide x 314" (8.3 cm) deep
25-350 BETHUNE Rib Shears 131¾" (34.9 cm), double curved handles
25-360 GLUCK Rib Shears 81½" (21.6 cm)
25-362 STILLE Pattern Rib Shears 81½" (21.6 cm)










25-378	STILLE-GIERTZ Rib Shears 10" (25.4 cm)
25-384 25-385 21-628	RUSKIN Bone Cutting Forceps 7% " (19.1 cm), straight standard blades RUSKIN Bone Cutting Forceps 7% " (19.1 cm), slightly angled standard blades RUSKIN Bone Cutting Forceps 7% " (19.1 cm), straight delicate blades
25-390 25-392	STILLE-LISTON Bone Cutting Forceps 109/4" (27.3 cm), straight STILLE-LISTON Bone Cutting Forceps 109/4" (27.3 cm), curved
25-396	STILLE-HORSLEY Bone Cutting Forceps 101/2" (26.7 cm), angled blades







25-440LEBSCHE Sternum Knife 10" (25.4 cm)25-454STILLE-LUER Rongeur, 8¾" (22.2 cm), straight jaws 9 x 15 mm25-456STILLE-LUER Rongeur, 8¾" (22.2 cm), curved jaws 9 x 15 mm25-458STILLE-LUER Duckbill Rongeur, 9" (22.9 cm), angled on side jaws 6 x 12 mm












19-850 to 19-856

25-511

25-480	SAUERBRUCH Rongeur 12" (30.5 cm), straight, square jaws 17 mm wide
19-850 19-852 19-854 19-856	RUSKIN Rongeur 7¼" (18.4 cm), straight, 5 mm bite RUSKIN Rongeur 7¼" (18.4 cm), curved, 5 mm bite RUSKIN Rongeur 7¼" (18.4 cm), straight, 7 mm bite RUSKIN Rongeur 7¼" (18.4 cm), curved, 7 mm bite
25-510 25-511	LEKSELL Rongeur 91/2" (24.1 cm), standard curve, 3 mm bite LEKSELL Rongeur 91/2" (24.1 cm), standard curve, 8 mm bite





5-128	MAYO-HARRINGTON Dissecting Scissors 9" (22.9 cm), straight
5-130	MAYO-HARRINGTON Dissecting Scissors 9" (22.9 cm), curved

NELSON LOBECTOMY SCISSORS


		~~	1 8
Straight	Length	Curved	
5-186	9" (22.9 cm)	5-188	
5-190	11" (27.9 cm)	5-192	WA I
5-194	14½" (36.8 cm)	5-196	///
5-200	HARRINGTON (WILLAUER) Lober straight	ctomy Scissors 111/2" (29.2 cm),	
5-202	HARRINGTON (WILLAUER) Lober curved	ctomy Scissors 11½" (29.2 cm),	
5-243	JORGENSON Scissors 9" (22.9 cm	n), strongly curved blades	
5-244 5-245	THOREK Scissors 71/4" (18.4 cm) THOREK Scissors 10" (25.4 cm)		5-244
5-260	SWEET Esophageal Scissors 81/2"	(21.6 cm), angled on flat	5-245

For Carb-N-Sert Scissors with tungsten carbide inserts see pages 55 - 59

For SuperCut Scissors see pages 60 - 63








6-140 6-141	BROWN Side Grasping Tissue Forceps 8" (20.3 cm), 8 \times 8 teeth BROWN Side Grasping Tissue Forceps 10" (25.4 cm), 8 \times 8 teeth
6-142 6-144 6-145	RUSSIAN Tissue Forceps 6" (15.2 cm) RUSSIAN Tissue Forceps 8" (20.3 cm) RUSSIAN Tissue Forceps 10" (25.4 cm)
6-200 6-202 6-204	KELLY Tissue Forceps 9" (22.9 cm), 1 x 2 teeth, slotted Gutch handle KELLY Tissue Forceps 9" (22.9 cm), 2 x 3 teeth, slotted Gutch handle KELLY Tissue Forceps 9" (22.9 cm), 3 x 4 teeth, slotted Gutch handle
6-210	MAYO-RUSSIAN Tissue Forceps 9" (22.9 cm), fenestrated jaws, slotted grooved handles
6-212	DAVIS Tissue Forceps 10" (25.4 cm), Russian type jaws
6-214	SINGLEY (TUTTLE) Tissue Forceps 9" (22.9 cm), fenestrated serrated jaws
16-80	WANGENSTEEN Tissue Forceps 9" (22.9 cm), rounded jaws with fine cross serrations






16-12	ALLIS Tissue Forceps 10" (25.4 cm), 5 x 6 teeth
16-40	LOCKWOOD Tissue Forceps 7¾" (19.7 cm)
16-46 16-48	BABCOCK Tissue Forceps 81/4" (21 cm) BABCOCK Tissue Forceps 91/2" (24.1 cm)
16-56	WILLIAMS Tissue Forceps 61/2" (16.5 cm), fenestrated jaws with cross serrated fins













COLLIN-DUVAL Lung Grasping Forceps 8" (20.3 cm), with 1/2" (1.3 cm) wide jaws

COLLIN-DUVAL Lung Grasping Forceps 8" (20.3 cm), with 16-52 1/2" (1.3 cm) wide jaws

COLLIN-DUVAL Lung Grasping Forceps 8" (20.3 cm), with 16-54 1" (2.5 cm) wide jaws


COLLIN-DUVAL Lung Grasping Forceps 9" (22.9 cm), with 1" (2.5 cm) wide jaws 25-754

LOVELACE Lung Grasping Forceps 8" (20.3-cm), straight, serrated jaws 1" (2.5 cm) wide 25-760

LOVELACE Lung Grasping Forceps 7%" (18.4 cm), angled shanks, serrated jaws 1" (2.5 cm) wide 14-66







25-782 25-784 25-786	RUMEL Thoracic Forceps 9" (22.9 cm), medium curve, 2%" (6.1 cm) long RUMEL Thoracic Forceps 9" (22.9 cm), full curve, 2%" (6.1 cm) long jaws RUMEL Thoracic Forceps 9" (22.9 cm), angular, 2%" (6.1 cm) long jaws	•
7-250 7-252 7-254	MOSQUITO Clamp 8½" (21 cm), straight MOSQUITO Clamp 8½" (21 cm), curved MOSQUITO Clamp 8½" (21 cm), angular	
17-2603	JACOBSON Forceps 7" (17.8 cm), curved, extremely delicate	
25-800	SAROT Intra-Thoracic Artery Forceps 9½" (24.1 cm), with 2½" (6.4 cm) long jaws	
7-255	GEMINI-MIXTER Thoracic Forceps 5½" (14 cm), delicate full curved jaws	
7-256	GEMINI-MIXTER Thoracic Forceps 7" (17.8 cm), delicate full curved jaws	
7-257	GEMINI-MIXTER Thoracic Forceps 8" (20.3 cm), delicate full curved jaws	()
7-258	GEMINI-MIXTER Thoracic Forceps 9" (22.9 cm), delicate full curved jaws	7-255 to 7-259


GEMINI-MIXTER Thoracic Forceps 11" (27.9 cm), delicate


7-259

full curved jaws













7-206		9" (22.9 cm), full curved jaws, horizontal
	serrations	

25-812

25-814

25-812 MIXTER Thoracic Forceps 9" (22.9 cm), right angle jaws, longitudinal serrations

25-814 MIXTER Thoracic Forceps 11" (27.9 cm), right angle jaws, longitudinal serrations

25-821 MILTEX Artery and Undermining Forceps 11" (27.9 cm), a delicate lightweight instrument with fine, serrated 45 mm long jaws


25-834 BRIDGE Deep Surgery Forceps 11" (27.9 cm), straight, very delicate BRIDGE Deep Surgery Forceps 11" (27.9 cm), curved, very delicate 25-835

25-836 KANTROWITZ Thoracic Forceps 71/2" (19.1 cm), delicate right angle jaws


25-837 KANTROWITZ Thoracic Forceps 91/2" (24.1 cm), delicate right angle jaws

KANTROWITZ Thoracic Forceps 11" (27.9 cm), delicate 25-838

right angle jaws







25-840	FINOCHIETTO Thoracic Forceps 9" (22.9 cm), slightly curved jaws
25-848	FINOCHIETTO Thoracic Forceps and Ligature Carrier 9½" (24.1 cm), angular jaws with longitudinal serrations and suture holes at tip

25-880 SAROT Bronchus Clamp 9" (22.9 cm), Left, adult size, 35 mm long jaws


25-882 SAROT Bronchus Clamp 9" (22.9 cm), Right, adult size, 35 mm long jaws

25-1100 SATINSKY Vascular Clamp for Vena Cava 10¾" (27.3 cm), inside jaw length 32 mm x 9 mm depth

25-1104 SATINSKY Vascular Clamp for Vena Cava 10%" (27.3 cm), inside jaw length 50 mm x 13 mm depth


25-1108 SATINSKY Vascular Clamp for Vena Cava 6½" (27.3 cm), pediatric size, inside jaw length 32 mm x 7 mm depth

25-1160 CRAFOORD Coarctation Clamp 9½" (24.1 cm), slightly curved jaws 65 mm with longitudinal serrations









5-246	POTTS-SMITH Scissors 71/2" (19.1 cm), angled on side 25°,
-------	---

20 mm blades

5-247 POTTS-SMITH Scissors 71/2" (19.1 cm), angled on side 45°,

20 mm blades

POTTS-SMITH Scissors 71/2" (19.1 cm), angled on side 60°, 5-248

20 mm blades

POTTS-SMITH Scissors 7" (17.8 cm), delicate, 5-249A

angled on side 25°, 13 mm blades

POTTS-SMITH Scissors 7" (17.8 cm), delicate, angled on side 45°, 13 mm blades 5-249B

5-249C POTTS-SMITH Scissors 7" (17.8 cm), delicate,

angled on side 60°, 13 mm blades

25-1240B DeBAKEY-POTTS Scissors 6%" (17.1 cm), blunt tips,

angled on side 45°

25-1240C DeBAKEY-POTTS Scissors 6%" (17.1 cm), blunt tips,

angled on side 60°

5-176 REYNOLDS Dissecting Scissors 7" (17.8 cm), curved tenotomy

type dissecting tips

REYNOLDS Dissecting Scissors 7" (17.8 cm), curved, tenotomy 5-178

type dissecting tips, one serrated blade
















7-310 7-315

25-1280	SATINSKY	Scissors	10"	(25.4	cm),	S-curved	

6-154	POTTS-SMITH Dressing Forceps 7" (17.8 cm), serrated
6-156	POTTS-SMITH Dressing Forceps 81/4" (21 cm), serrated
6-158	POTTS-SMITH Dressing Forceps 91/2" (24.1 cm), serrated
6-159	POTTS-SMITH Dressing Forceps 12" (30.5 cm), serrated
6-160	POTTS-SMITH Tissue Forceps 7" (17.8 cm), 1 x 2 teeth and serrated tips
6-162	POTTS-SMITH Tissue Forceps 81/4" (21 cm), 1 x 2 teeth and serrated tips
6-164	POTTS-SMITH Tissue Forceps 91/2" (24.1 cm), 1 x 2 teeth and serrated tips
6-164A	POTTS-SMITH Tissue Forceps 12" (30.5 cm), 1 x 2 teeth and serrated tips

JOHNS HOPKINS BULLDOG CLAMPS

	BOTTLE TIOT KIND DOLLDOG OLD INI	
Straight	Size	Curved
7-280	1½" (3.8 cm)	7-290
7-282	2" (5.1 cm)	7-292
7-284	21/2" (6.4 cm)	7-294
7-286	3" (7.6 cm)	7-296
7-288	3½" (8.9 cm)	7-298

7-310	DIETHRICH Micro Bulldog Clamp 1%" (4.7 cm), straight jaws 8 mm x 1.2 mm
7-315	DIETHRICH Micro Bulldog Clamp 1%" (4.7 cm), angled jaws 8 mm x 1.2 mm

25-1420 Bulldog Clamp Applying Forceps 91/2" (24.1 cm)














8-51TC CRILE-WOOD (DeBAKEY) Needle Holder 7" 8-52TC CRILE-WOOD (DeBAKEY) Needle Holder 8" CRILE-WOOD (DeBAKEY) Needle Holder 9" 8-53TC

CRILE-WOOD (DeBAKEY) Needle Holder 12"

(17.8 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch (20.3 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch (22.9 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch CRILE-WOOD (DeBAKEY) Needle Holder 101/2" (26.7 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch

(30.5 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch

8-57TC 8-62TC

8-55TC

COOLEY Needle Holder 6" (15.2 cm), extra delicate, serrated Carb-N-Sert jaws 4000 teeth p. sq. inch COOLEY Needle Holder 61/4" (17.1 cm), extra delicate, serrated Carb-N-Sert jaws 4000 teeth p. sq. inch

8-63TC

RYDER-Type Intracardiac Needle Holder 5" (12.7 cm), extra delicate, serrated Carb-N-Sert jaws

8-65TC 4000 teeth p. sq. inch

8-66TC RYDER-Type Intracardiac Needle Holder 6" (15.2 cm), extra delicate, serrated Carb-N-Sert jaws

4000 teeth p. sq. inch

8-67TC RYDER-Type Intracardiac Needle Holder 7" (17.8 cm), extra delicate, serrated Carb-N-Sert jaws


4000 teeth p. sq. inch

RYDER-Type Intracardiac Needle Holder 9" (22.9 cm), extra delicate, serrated Carb-N-Sert jaws 8-69TC

4000 teeth p. sq. inch



17-1090TC to 17-1092TC



8-961C	SAROT Needle Holder 7/4" (18.4 cm), serrated Carb-N-Seri Jaws 2500 teeth p. sq. inch
8-114TC	SAROT Needle Holder 101/2" (26.7 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch
8-101TC	RYDER Needle Holder 5" (12.7 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-103TC	RYDER Needle Holder 6" (15.2 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-104TC	RYDER Needle Holder 7" (17.8 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-105TC	RYDER Needle Holder 8" (20.3 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-106TC	RYDER Needle Holder 9" (22.9 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-107TC	RYDER Needle Holder 101/4" (26 cm), serrated Carb-N-Sert jaws 4000 teeth p. sq. inch
8-134TC	WANGENSTEEN Needle Holder 10%" (27.3 cm), serrated Carb-N-Sert jaws 2600 teeth p. sq. inch
17-1090TC	JACOBSON Microvascular Needle Holder 7" (17.8 cm), straight, with lock, serrated Carb-N-Sert jaws 16000 teeth p. sq. inch
17-1091TC	JACOBSON Microvascular Needle Holder 7" (17.8 cm), curved, with lock, serrated Carb-N-Sert jaws 16000 teeth p. sq. inch
17-1092TC	JACOBSON Microvascular Needle Holder 8½" (21.6 cm), straight, with lock, serrated Carb-N-Sert laws 16000 teeth p. sq. inch



Neurosurgical Instruments















CUSHING Perforator Drill 1/2" (1.3 cm) diameter 26-2 26-16 McKENZIE Perforator Drill, 13 mm diameter 26-18 McKENZIE Enlarging Bur, 16 mm diameter D'ERRICO Perforator and Enlarging Bur, 16 mm diameter 26-20 Above Burs and Drills will fit Hudson Braces HUDSON Cerebellar Extension 4" (10.2 cm), for use with Hudson Braces 26-24



26-33-1 26-33-2 26-33-3 26-33-4	HUDSON Bur, 9 mm diameter HUDSON Bur, 14 mm diameter HUDSON Bur, 16 mm diameter HUDSON Bur, 22 mm diameter
26-30	HUDSON Brace Set, standard thumbscrew chuck, complete with 4 Burs (1 ea. Nos. 26-33-1 through 4)
26-32 26-35	HUDSON Brace 10½" (26.7 cm), with standard thumbscrew chuck HUDSON Brace 11" (27.9 cm), with improved snap-lock chuck
26-40	STILLE Pattern Cranial Drill Set, complete with 4 twist drills (2.5 to 4 mm), 4 burs (4 to 16 mm) and 1 trephine (10 mm). Chrome plated drill, stainless points.
26-42	STILLE Pattern Cranial Drill 10" (25.4 cm), chrome plated
	Penlacement Points for Stille Pattern Drille

Replacement Points for Stille Pattern Drills


	•		
Twist	Drills	Roun	d Burs
26-44-1	2.5 mm	26-46-1	4 mm
26-44-2	3 mm	26-46-2	8 mm
26-44-3	3.5 mm	26-46-3	11 mm
26-44-4	4 mm	26-46-4	16 mm
26-48	Trephine 10 mm dia		













26-28A to 26-28E

26-60 26-62

GALT SKULL TREPHINES, with Removable Handle For Use With Hudson Brace (26-30), Chrome Plated. Cutters and Adjustable Centering Points are Stainless

26-28A 1/2" (1.3 cm) dia. 26-28D 7/8" (2.2 cm) dia. 26-28B %" (1.6 cm) dia, 26-28E 1" (2.5 cm) dia. 26-28C 3/4" (1.9 cm) dia,

RANEY-CRUTCHFIELD Skull Traction Tong 120 mm 26-60 spread, incl. 2 drill points


26-62 RANEY-CRUTCHFIELD Drill Point

JANSEN Retractor 41/4" (10.8 cm), spread 3" (7.6 cm), 11-2 3 x 3 blunt prongs, 10 mm deep


JANSEN Scalp Retractor 41/4" (10.8 cm), spread 3" (7.6 cm), 11-4 4 x 4 blunt prongs, 20 mm deep

11-8 JANSEN-WAGNER Retractor 5" (12.7 cm), 5 x 5 sharp

prongs, self adjusting











11-10	DAVIS Scalp Retractor 41/2" (11.4 cm), 3 x 3 blunt prongs 10 mm deep

26-91 TAYLOR Spinal Retractor 7¼" (18.4 cm), blade 1¼" (3.2 cm) x 3" (7.6 cm) TAYLOR Spinal Retractor 71/4" (18.4 cm), blade 11/4" (3.2 cm) x 4" (10.2 cm) 26-92


ADSON Cerebellar Retractor 71/2" (19.1 cm), 4 x 4 teeth 3/4" (1.9 cm) 11-700

x 34" (1.9 cm), straight arms

11-702 ADSON Cerebellar Retractor 71/2" (19.1 cm), 4 x 4 teeth 3/4" (1.9 cm)


x 3/4" (1.9 cm), angled arms

BECKMAN-WEITLAN	ER RETRACTORS	6, with hinged blades, 3 x 4 teeth	
Sharp	Size	Blunt	a
11-630-SH	5½" (14 cm)	11-630-BL (Cone Scalp Retr.)	((
11-634-SH (Cone Lamin. Retr.)	6½" (16.5 cm)	11-634-BL	//
11-642-SH (Cone Lamin. Retr.)	9" (22.9 cm)	11-642-BL	



11-630 to 11-642







26-119 26-120 26-122 26-124	GIGLI Saw Handle, Standard Pattern, T-Shaped, turnable GIGLI-STRULLY Saw Handle, Standard Pattern, T-Shape, rigid, with snap-lock GIGLI-STRULLY Saw Handle, Standard Pattern, T-Shape, turnable, with snap-lock GIGLI-STRULLY Saw Handle, Standard Pattern, Loop Grip, rigid, with snap-lock
26-130 26-132 26-134	GIGLI Saw 12" (30.5 cm), Standard twisted wire type GIGLI Saw 20" (50.8 cm), Standard twisted wire type GIGLI Saw 30" (76.2 cm), Standard twisted wire type
26-138	BAILEY GIGLI Saw Guide 12" (30.5 cm), flexible
26-140 26-141 26-142 26-143	MICHELE Trephine 6¼" (15.9 cm) shaft, ¾6" (.47 cm) O.D., graduated ¼" (.64 cm) to 1¼" (3.2 cm) MICHELE Trephine 6¼" (15.9 cm) shaft, ½2" (.56 cm) O.D., graduated ¼" (.64 cm) to 1¼" (3.2 cm) MICHELE Trephine 6¼" (15.9 cm) shaft, ¾6" (.79 cm) O.D., graduated ¼" (.64 cm) to 1¼" (3.2 cm) Marker for No. 26-142 Michele Trephine











19-808

26-190 26-192

26-198

26-200

19-808	BANE Rongeur 7" (17.8 cm), curved
26-190 26-192	LUER Rongeur 7" (17.8 cm), straight jaws 8 mm bite LUER Rongeur 7" (17.8 cm), curved jaws 8 mm bite
26-198	BACON Ronguer 81/4" (21 cm), angled on side, 5 mm bite
26-200	STOOKEY Rongeur 7%" (19.7 cm), angled on side, 3 mm bite


26-220 ADSON Cranial Rongeur 8" (20.3 cm), straight jaws 7 x 13 mm 26-222 ADSON Cranial Rongeur 8" (20.3 cm), curved jaws 7 x 13 mm



26-220 26-222







ZAUFEL-JANSEN Rongeur 7" (17.8 cm), slightly curved jaws 5 mm wide 19-840

19-844 BEYER Rongeur 7" (17.8 cm), slightly curved jaws 3.5 mm wide

25-510 LEKSELL Laminectomy Rongeur 91/2" (24.1 cm), standard curve, 3 mm bite

LEKSELL Laminectomy Rongeur 91/2" (24.1 cm), standard curve, 25-511 8 mm bite


BUSKIN BONGEURS

	110011111111111111111111111111111111111	
Straight Jaws	Size	Curved Jaws
19-848	6" (15.2 cm) long, 3 mm bite	19-849
19-850	71/4" (18.4 cm) long, 5 mm bite	19-852
19-854	71/4" (18.4 cm) long, 7 mm bite	19-856




19-850 to 19-856









26-255 to 26-258


26-261 26-262

26-264 to 26-266

26-255 26-256 26-258	SMITH-PETERSEN Laminectomy Rongeur 7½" (19.1 cm), straight jaws, 3 mm bite SMITH-PETERSEN Laminectomy Rongeur 7½" (19.1 cm), slightly curved jaws, 3 mm bite SMITH-PETERSEN Laminectomy Rongeur 9½" (24.1 cm), straight jaws, 3 mm bite
26-261 26-262	FULTON Laminectomy Rongeur 9½" (24.1 cm), straight jaws 7 x 18 mm FULTON Laminectomy Rongeur 9½" (24.1 cm), slightly curved jaws 7 x 18 mm
26-264 26-265 26-266	ECHLIN Duckbill Rongeur 9" (22.9 cm), 2 mm bite ECHLIN Duckbill Rongeur 9" (22.9 cm), 3 mm bite ECHLIN Duckbill Rongeur 9" (22.9 cm), 4 mm bite








STILLE-LUER Rongeur 8¾" (22.2 cm), straight jaws 9 x 15 mm 25-454 STILLE-LUER Rongeur 8¾" (22.2 cm), curved jaws 9 x 15 mm 25-456 STILLE-LUER Duckbill Rongeur, 9" (22.9 cm), angled on side jaws 6 x 12 mm 25-458 DeVILBISS Cranial Rongeur 81/4" (21 cm), with 2 blades 26-270








20-660 HAJEK-KOFLER Laminectomy Rongeur 5" (12.7 cm) shaft, up-bite, 4 x 5 mm
4 HAJEK-KOFLER Laminectomy Rongeur 5" (12.7 cm) shaft, down-bite, 4 x 5 mm

SCHLESINGER CERVICAL RONGEURS, THIN TIPS, 6" (15.2 cm) SHAFT

Up-Bite		Down	-Bite
26-326	3 mm	26-328	3 mm
26-327	5 mm	26-329	5 mm








SPURLING-KERRISON CERVICAL RONGEURS

Up-Bite		Down-Bite 6" (15.2 cm) Shaft		40° Forward Angle	
26-330	3 mm	26-332	3 mm	26-334	3 mm
26-331	5 mm	****	*****	26-335	5 mm
		8" (20.3 cr	n) Shaft		
26-336	3 mm	26-338	3 mm	26-339A	3 mm
26-337	5 mm			26-339B	5 mm

LOVE-KERRISON CERVICAL RONGEURS

Up-I	3ite		Dow	n-Bite
		6" (15.2 cm) Shaft		
26-340	3 mm		26-370	3 mm






MILTEX-KERRISON CERVICAL RONGEURS 7" (17.8 cm) SHAFT, IMPROVED HANDLE

	2 mm	3 mm	5 mm
Up-Bite	26-392A	26-392B	26-392
Down-Bite	26-393A	26-393B	26-393
40° Forward Angle	26-388	26-389	26-390







CUSHING PITUITARY RONGEURS 2 x 10 mm CUP JAWS

	Straight	Curved-Up	Curved-Dowr
5" (12.7 cm) Shank	26-394	26-395	26-396
7" (17.8 cm) Shank	26-397	26-398	26-399


LOVE-GRUENWALD PITUITARY RONGEURS 3 x 10 mm CUP JAWS


	Straight	Curved-Up	Curved-Down
5" (12.7 cm) Shank	26-400	26-410	26-420
7" (17.8 cm) Shank	26-402	26-412	26-422

SPURLING PITUITARY RONGEURS 4 x 10 mm CUP JAWS

	Straignt	Curvea-up	Carvea-Dow
5" (12.7 cm) Shank	26-430	****	
7" (17.8 cm) Shank	26-435	26-436	26-437







26-440 OLDBERG Pituitary Rongeur 7" (17.8 cm) Shank, round cup jaws 6 mm dia.

CLOWARD PITUITARY RONGEURS 6 x 10 mm CUP JAWS, 6" (15.2 cm) SHANKS

Straight 26-480 Curved-Up 26-481 Curved-Down 26-482

SCHLESINGER PITUITARY RONGEURS, SERRATED CUP JAWS, STRAIGHT 6" (15.2 cm) SHANKS

26-483


2 x 10 mm bite

26-484

3 x 10 mm bite

26-485

4 x 10 mm bite






21-560 5-250 5D-251 5-252 5D-253

5-262

probe tips

Neurosurgical Instruments






FERRIS-SMITH PITUITARY CUP JAW RONGEURS, 7" (17.8 cm) SHANK

	Straight	Curved-Up	Curved-Down
2 x 10 mm Jaws			26-470
3 x 10 mm Jaws	26-461	26-466	26-471
4 x 10 mm Jaws	26-462		26-472
6 x 12 mm Jaws	26-474 (Clow	ard Patt.)	
KELLY (ADSON) Ga	nglion Scissors 61/4 nglion Scissors 61/4	ith probe tip ' (15.9 cm), straight, sh ' (15.9 cm), straight, sh	, ,
with one serrated bla			1
. 1	•	' (15.9 cm), curved, sha	1 4
KELLY (ADSON) Ga with one serrated bla		' (15.9 cm), curved, sha	arp points,
STRULLY Neurologic	eal Scissors 8" (20 :	3 cm) slightly curved bl	ades with

5-262






9-32	McKENZIE Brain Clip Holding Rack, for 20 clips	
9-74	SMITHWICK Clip Applying Forceps 9" (22.9 cm), curved jaws, straight handles	A.
9-80	McKENZIE Clip Applying Forceps 6" (15.2 cm), curved jaws, straight handles	Of the second
9-90	SCOVILLE Clip Applying Forceps, 6" (15.2 cm) shaft	
26-590 26-591	RANEY Scalp Clip Applying Forceps 61/4" (15.9 cm) RANEY Reusable Scalp Clips, black finish	








26-660	ADSON Aneurysm Needle 9" (22.9 cm)
26-700	SCOVILLE Spatula 8" (20.3 cm), malleable ends $\%$ " (1.6 cm) and $\%$ " (2.2 cm) wide
26-710	CHILDREN'S HOSPITAL Spatula 8% (21.6 cm), malleable, ends are $\%$ (1.3 wide), silver plated
26-724 26-726	DAVIS Spatula 7" (17.8 cm) x ½" (1.3 cm), nickel silver, malleable DAVIS Spatula 7" (17.8 cm) x ¾" (1.9 cm), nickel silver, malleable
26-730 26-732	CUSHING Spatula Spoon 6%" (17.1 cm), spoon ends ½" (1.3 cm) and %" (1.9 cm) wide CUSHING Spatula Spoon 7%" (19.7 cm), spoon ends %" (1.9 cm) and 1" (2.5 cm) wide
26-760 26-762	ADSON Cannula 9 Fr. (3 mm), 10 cm graduated in cm, open end, sterling ADSON Cannula 9 Fr. (3 mm), 14 cm graduated in cm, open end, sterling
26-770 26-772 26-774	ELSBERG Cannula 11 cm, 3-hole closed end, 5 Fr., (1.7 mm), Luer lock, sterling ELSBERG Cannula 11 cm, 3-hole closed end, 7 Fr., (2.3 mm), Luer lock, sterling ELSBERG Cannula 11 cm, 3-hole closed end, 9 Fr., (3 mm), Luer lock, sterling
26-778	FRAZIER Needle 9 cm, 3-hole closed blunt trocar point, 9 Fr. (3 mm)






26-780 26-782	ADSON Suction Tube 6" (15.2 cm), with finger valve and obturator, 11 Fr. (3.6 mm)			
20-702	ADSON Suction Tube 8" (20.3 cm), with finger valve and obturator, 15 Fr. (5 mm)			
	HOUSE SU	CTION - IRRIGATING TUBES	, 9 cm working leng	gth
	19-562	1.8 mm Suction Tube and 1	1.8 mm Irrigating Tu	ıbe
	19-563	2.5 mm Suction Tube and	2 mm Irrigating Tu	ube
	19-564	3 mm Suction Tube and	2 mm Irrigating Tu	ube
	19-565	4.5 mm Suction Tube and	2 mm Irrigating To	ube
	FRAZIER SUCTION TUBES, with finger cut-off, angled			
	19-570	6 French (2 mm)	19-575	10 French (3.3 mm)
	19-572	7 French (2.3 mm)	19-576	11 French (3.6 mm)
	19-573	8 French (2.6 mm)	19-577	12 French (4 mm)
	19-574	9 French (3 mm)		
19-580	BARON Suction Tube, 3 French (1 mm), with finger cut-off			
19-582	BARON Su	ction Tube, 5 French (1.7 mm)	, with finger cut-off	
19-584	BARON Su	ction Tube, 7 French (2.3 mm)	, with finger cut-off	



6-118

Neurosurgical Instruments





6-118A	ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate, angled
6-118XL	ADSON Dressing Forceps 4¾" (12.1 cm), serrated, lightweight, fenestrated handles
17-2510	Micro ADSON Dressing Forceps 4%" (12.1 cm), serrated tips, 0.5 mm wide
6-120	ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, delicate
6-120A	ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, delicate, angled
6-120XL	ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, lightweight,

ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate

fenestrated handles 6-121 ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth and cross serrated tips, delicate


Micro ADSON Tissue Forceps 4¾" (12.1 cm), tips 0.5 mm wide, 17-2500 1 x 2 teeth

6-122 ADSON Tissue Forceps 43/4" (12.1 cm), 2 x 3 teeth, delicate ADSON Tissue and Suture Forceps 4¾" (12.1 cm), 1 x 2 teeth with 6-123 tying platform ADSON Tissue Forceps 43/4" (12.1 cm), w/platform, 1 x 2 teeth, 6-123XL

lightweight, fenestrated handles


17-2500 6-123 6-123XL






6-124 6-124XL	BROWN-ADSON Tissue Forceps 494" (12.1 cm), 7 x 7 side grasping teeth BROWN-ADSON Tissue Forceps 494" (12.1 cm), 7 x 7 side grasping teeth, lightweight, fenestrated handles		
6-126	BROWN-ADSON Tissue Forceps 4%" (12.1 cm), 7 x 7 side grasping teeth, angular	77.	
6-104 6-106 6-107 6-108 6-109 6-109A	SEMKEN Dressing Forceps 5" (12.7 cm), serrated SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth, curved tips SEMKEN Tissue Forceps 5" (12.7 cm), 2 x 3 teeth SEMKEN Dressing Forceps 6" (15.2 cm), serrated SEMKEN Tissue Forceps 6" (15.2 cm), 1 x 2 teeth		
6-166	CUSHING Dressing Forceps 7" (17.8 cm), serrated, Gutch handle, dissecting end		
6-168	CUSHING Tissue Forceps 7" (17.8 cm), 1 x 2 teeth, Gutch handle, dissecting end		
6-188	CUSHING Dressing Forceps 7" (17.8 cm), serrated handle, serrated jaws, dissecting end	6-188	€ 6-189
6-189	CUSHING Tissue Forceps 7" (17.8 cm), serrated handle, 1 x 2 teeth, dissecting end		








6-138 6-140	BROWN Tissue Forceps 6" (15.2 cm), double rows of 8 x 8 teeth BROWN Tissue Forceps 8" (20.3 cm), double rows of 8 x 8 teeth
6-180 6-182 6-184 6-186	GERALD Dressing Forceps 7" (17.8 cm), straight, serrated GERALD Tissue Forceps 7" (17.8 cm), straight, 1 x 2 teeth GERALD Dressing Forceps 7" (17.8 cm), curved, serrated GERALD Tissue Forceps 7" (17.8 cm), curved, 1 x 2 teeth
6-190 6-192	CUSHING Bayonet Dressing Forceps 71/4" (18.4 cm), serrated, scraper end CUSHING Bayonet Tissue Forceps 71/4" (18.4 cm), 1 x 2 teeth, scraper end
6-197	GRUENWALD (JANSEN) Bayonet Dressing Forceps 61/2" (16.5 cm), serrated
6-198	ADSON Bayonet Dressing Forceps 81/4" (21 cm), serrated
26-956 26-958	FERRIS SMITH Tissue Forceps 7" (17.8 cm), serrated platform tips with 1 x 2 teeth FERRIS SMITH Tissue Forceps 7" (17.8 cm), serrated platform tips with 2 x 3 teeth









7-240 7-242



26-1040

7-10	extra delicate HALSTED Mosquito Forceps 5" (12.7 cm), curved, extra delicate
17-2600	JACOBSON Micro Mosquito Forceps 5" (12.7 cm), straight, extremely delicate
17-2602	JACOBSON Micro Mosquito Forceps 5" (12.7 cm), curved, extremely delicate
17-2603	JACOBSON Hemostatic Forceps 7" (17.8 cm), curved, extremely delicate
7-102	DANDY Scalp Hemostatic Forceps 5½" (14 cm), curved sideways, serrated
7-240 7-242	ADSON Hemostatic Forceps $74\mbox{\ensuremath{''}}$ (18.4 cm), straight ADSON Hemostatic Forceps $74\mbox{\ensuremath{''}}$ (18.4 cm), curved
26-1040 26-1041	SCHWARTZ Clip Applying Forceps 8" (20.3 cm) SCHWARTZ Vessel Clip 1" (2.5 cm), straight jaws 1.7 x 8 mm
26-1042	SCHWARTZ Vessel Clip 1" (2.5 cm), slightly angled jaws 1.7 x 8 mm
26-1043	SCHWARTZ Vessel Clip 1" (2.5 cm), strong angled jaws 1.7 x 8 mm


HALSTED Mosquito Forceps 5" (12.7 cm), straight,

For Micro Vessel Clips see page 223

7-8








FRAZIER Dura Separator 61/2" (16.5 cm), angled blade 1/2" (1.3 cm) long, 26-1054 from 1/8" (.32 cm) to 1/16" (.16 cm) wide CUSHING Dura Hook 61/2" (16.5 cm), blunt 26-1080 26-1084 FRAZIER Dura Hook 5" (12.7 cm), sharp point FRAZIER Dura Hook 5" (12.7 cm), blunt point 26-1086 ADSON Dura and Dissecting Hook 8" (20.3 cm), sharp blade 5 mm at right angle 26-1090 26-1092 ADSON Dura Hook 8" (20.3 cm), blunt 5 mm hook at right angle 21-430 MILTEX Tissue Twist Hook 51/2" (14 cm) DANDY Nerve Hook 9" (22.9 cm), straight 26-1150






26-1156	GRAHAM Nerve Hook 7" (17.8 cm), blunt, 6 mm deep
26-1170 26-1174	SMITHWICK Nerve Hook 7" (17.8 cm), blunt, 10 mm deep SMITHWICK Sympathectomy Hook and Dissector 12" (30.5 cm)
11-100 11-102	CUSHING Decompression Retractor 81/2" (21.6 cm), fenestrated blade CUSHING Vein and Nerve Retractor 9" (22.9 cm), blade 11 x 13 mm
26-1326	SACHS Nerve Separator-Spatula 71/2" (19.1 cm), blunt blade 3 mm wide, slightly curved








LOVE Nerve Root Retractor 81/2" (21.6 cm), straight, 7 mm blade 21-180 LOVE Nerve Root Retractor 81/2" (21.6 cm), angled 45°, 7 mm blade 26-1352 LOVE Nerve Root Retractor 81/2" (21.6 cm), angled 90°, 7 mm blade 26-1354 SCOVILLE Nerve Root Retractor 9" (22.9 cm), straight, 8 mm wide blade 26-1360 SCOVILLE Nerve Root Retractor 9" (22.9 cm), angular, 8 mm wide blade 26-1362






26-1403	CUSHING "LITTLE JOKER" Periosteal Elevator 6" (15.2 cm)
26-1410 26-1412	LANGENBECK Periosteal Elevator 7½" (19.1 cm), narrow model, 7 mm wide blunt blade LANGENBECK Periosteal Elevator 7½" (19.1 cm), wide model, 17 mm wide sharp blade
18-1968	FREER Elevator 7" (17.8 cm), double end sharp and blunt blades 5 mm wide
26-1430 26-1432	ADSON Periosteal Elevator 61/2" (16.5 cm), chisel edge, straight, 8 mm wide ADSON Periosteal Elevator 61/2" (16.5 cm), blunt, curved, 6 mm wide
26-1434	ADSON Periosteal Elevator 6½" (16.5 cm), semi-sharp, curved, 7 mm wide
26-1436	ADSON Periosteal Elevator 61/8" (16.5 cm), semi-sharp, straight, 6 mm wide








PENFIELD DISSECTORS

26-1450	Style No. 1	7¼" (18.4 cm), double end, 6 mm sharp cup and 6 mm blunt dissector
26-1451	Style No. 2	7%" (19.7 cm), double end, wax packer and 6 mm blunt dissector, slightly curved
26-1452	Style No. 3	7% (19.1 cm), double end, wax packer and 6 mm blunt dissector, fully curved
26-1453	Style No. 4	8¾" (22.2 cm), single end, light pattern, 3 mm blunt dissector, slightly curved
26-1454	Style No. 5	11% (29.2 cm), double end, slightly curved, blunt dissector blades 7 mm and 8 mm wide






26-1620 26-1624



26-1627 to 26-1640



19-700 to 19-720



26-1650 to 26-1660

26-1620 26-1624 SCOVILLE Ruptured Disc Curette 10" (25.4 cm), angled up, cup 4 x 10 mm SCOVILLE Ruptured Disc Curette 10" (25.4 cm), angled down, cup 4 x 10 mm

HIBBS SPRATT SPINAL FUSION CURETTES 9" (22.9 cm), Oval Cup

26-1627	Size 3/0	26-1634	Size 3
26-1628	Size 2/0	26-1636	Size 4
26-1629	Size 0	26-1638	Size 5
26-1630	Size 1	26-1640	Size 6
26-1632	Size 2		


SPRATT (BRUN) BONE CURETTES 61/2" (16.5 cm), Oval Cups


19-700	Size 5/0	19-712	Size 2
19-702	Size 4/0	19-714	Size 3
19-704	Size 3/0	19-716	Size 4
19-706	Size 2/0	19-718	Size 5
19-708	Size 0	19-720	Size 6
10.710	Size 1		

BRUN BONE CURETTES 6¾" (17.1 cm), Round Cups

26-1650	Size 1	7 mm	26-1656	Size 4	11 mm
26-1652	Size 2	8 mm	26-1658	Size 5	12 mm
26-1654	Size 3	9 mm	26-1660	Size 6	14 mm









11-672 BADGLEY Laminectomy Retractor, 6½" (16.5 cm) spread, 4 x 4 blunt teeth 2" (5.1 cm) x 1¾" (4.4 cm)

11-660 FRAZIER Laminectomy Retractor, 3¼" (8.3 cm) spread, 3 x 4 blunt teeth

11-676 TUFFIER Retractor, 6½" (16.5 cm) spread, 2" (5.1 cm) x 1¾" (4.4 cm) grooved blades







11-240 MEYERDING Laminectomy Retractor 9½" (24.1 cm), large blade 3½" (8.9 cm) x 2" (5.1 cm)

11-242 MEYERDING Laminectomy Retractor 9½" (24.1 cm), medium blade 3" (7.6 cm) x 1" (2.5 cm)

11-244 MEYERDING Laminectomy Retractor 9½" (24.1 cm), small blade 2" (5.1 cm) x ¾" (1.6 cm)

11-642-SH BECKMAN-WEITLANER Retractor, hinged, 9" (22.9 cm), sharp 3 x 4 teeth

11-642-SH BECKMAN-WEITLANER Retractor, hinged, 9" (22.9 cm), sharp 3 x 4 teeth 11-642-BL BECKMAN-WEITLANER Retractor, hinged, 9" (22.9 cm), blunt 3 x 4 teeth

MEYERDING SELF-RETAINING LAMINECTOMY RETRACTORS, 3%" (9.5 cm) spread

 11-650
 1" (2.5 cm) x 2%" (6.1 cm) blades

 11-652
 1" (2.5 cm) x 2%" (5.4 cm) blades

 11-654
 1" (2.5 cm) x 1%" (4.4 cm) blades

 11-656
 %" (1.9 cm) x 2%" (6.1 cm) blades


11-684 BECKMAN-ADSON Laminectomy Retractor 12½" (31.8 cm), hinged, sharp 4 x 4 teeth 1" (2.5 cm) wide x ¾" (1.9 cm, deep


11-686 BECKMAN-ADSON Laminectomy Retractor 12½" (31.8 cm), hinged, sharp 4 x 4 teeth 1½" (3.2 cm) wide x 1¾" (4.4 cm) deep













26-1973 BECKMAN-EATON Laminectomy Retractor 123/4" (32.4 cm), hinged, sharp 7 x 7 teeth 13/4" (4.4 cm) wide x 2" (5.1 cm) deep 11-365 Pediatric GELPI Retractor 51/2" (14 cm), with grip lock GELPI Retractor 63/4" (17.1 cm), with grip lock 11-366 INGE Lamina Spreader 6" (15.2 cm), jaws open to 11/6" (2.8 cm) 26-2000 26-2001 INGE Lamina Spreader 6" (15.2 cm), jaws open to 11/8" (2.8 cm), without teeth 26-2005 INGE Lamina Spreader 91/2" (24.1 cm), jaws open to 11/4" (4.4 cm)



27-2	KERN Bone Holding Forceps 91/2" (24.1 cm), standard size, plain
27-4	KERN Bone Holding Forceps 91/2" (24.1 cm), standard size, with ratchet
27-6	KERN Bone Holding Forceps 6" (15.2 cm), pediatric size, plain
27-7	KERN Bone Holding Forceps 6" (15.2 cm), pediatric size, with ratchet
27-20 27-21 27-22 27-23	LANE Bone Holding Forceps 13" (33 cm), plain LANE Bone Holding Forceps 13" (33 cm), with ratchet LANE Bone Holding Forceps 17" (43.2 cm), plain LANE Bone Holding Forceps 17" (43.2 cm), with ratchet





27-14 FARABEUF-LAMBOTTE Bone Holding Forceps, 10¼" (26 cm)
 27-16 LEWIN Bone Holding Forceps 7" (17.8 cm), serrated jaws

27-24 DINGMAN Small Bone and Cartilage Clamp 7½" (19.1 cm), serrated jaws, 2 x 2 teeth. 4 mm wide

441









27-47





27-30 to 27-34

27-30

27-32

27-34


LOWMAN Bone Clamp 5" (12.7 cm), calibrated in 1/4" (.32 cm), 1 x 2 prong jaws 3/4" (1.9 cm) wide
LOWMAN Bone Clamp 7" (17.8 cm), calibrated in 1/4" (.32 cm), 1 x 2 prong jaws 1" (2.5 cm) wide
LOWMAN Bone Clamp 8" (20.3 cm), calibrated in 1/4" (.32 cm), 1 x 2 prong jaws 11/4" (3.2 cm) wide

27-42 LOWMAN-HOGLUND Bone Clamp 7" (17.8 cm), calibrated in 1/4" (.32 cm), 2 x 2 prong swivel jaws 1" (2.5 cm) wide


27-44 LOWMAN-HOGLUND Bone Clamp 8" (20.3 cm), calibrated in 1/4" (.32 cm), 2 x 2 prong swivel jaws 13/4" (3.5 cm) wide


27-46 LAMBERT-LOWMAN Bone Clamp 7" (17.8 cm), 2 x 2 prong jaws 1" (2.5 cm) wide
27-47 LAMBERT-LOWMAN Bone Clamp 8" (20.3 cm), 2 x 2 prong jaws 1%" (3.5 cm) wide
27-60 VAN BUREN Sequestrum Forceps 9" (22.9 cm), serrated jaws angled on side
27-110 Bone Wire Guide 8%" (22.2 cm), U.S. Pattern

27-144 POLK Finger Goniometer 6" (15.2 cm), graduated in inches and centimeters











27-159TC

9-128 Double Action Parallel Pliers and Wire Cutter 41/2" (11.4 cm), nickel plated

22-500TC Wire Twister 6" (15.2 cm), with Carb-N-Sert jaws 22-502TC Wire Twister 8" (20.3 cm), with Carb-N-Sert jaws

MILTEX Cannulated Pin and Wire Cutter 71/2" (19.1 cm), end cutting up to 27-158

3.2 mm (1/6") soft wire or 1.6 mm (1/16") hard wire

C-N-S Pin and Wire Cutter 7" (17.8 cm), double action, CARB-N-SERT cutting blades, cuts up to 2 mm (% ") soft wire or 1.6 mm (1/16") hard wire 27-159TC

DIAMOND Pin and Wire Cutter 61/2" (16.5 cm), cuts up to 2 mm (5/4") diameter 27-160

hard wire, stainless



27-160

Orthopedic Instruments 27-16/ITC 27-16/STC 27-16/STC 27-176/ITC 27-176/ITC

27-161TC MILTEX Double Action Pin Cutter 15" (38.1 cm), heavy duty, with CARB-N-SERT cutting blades, will cut up to 4.8 mm (%") diameter, stainless

27-163TC C-N-S Front and Side Pin and Wire Cutter 7" (17.8 cm), Double Action Flush Cutting CARB-N-SERT Blades cuts up to 2 mm (%4") soft wire or 1.6 mm (%4") hard wire


27-164TC C-N-S Pin & Wire Puller 7" (17.8 cm), double action, with CARB-N-SERT 4 mm wide jaws
 27-165TC C-N-S Pin & Wire Puller 7" (17.8 cm), double action, with CARB-N-SERT 6 mm wide jaws


27-165TC

27-176TC MILTEX Pin Puller 5½" (14 cm) CARB-N-SERT serrated top jaw measures 1 cm in length. Grooved cross-serrated opposing jaw. (Groove measures 1 cm in length.)

For Orthopedic Suture Wire & Wire Cutting Scissors - see pages 135 & 141.







26-40 STILLE Pattern Drill Set, complete with 4 twist drills (2.5 to 4 mm), 4 burs (4 to 16 mm) and 1 trephine (10 mm). Chrome plated drill, stainless points

26-42 STILLE Pattern Drill, chrome plated

26-42

Replacement Points for Stille Pattern Drill

Twist Drills		Round E	Burs
26-44-1	2.5 mm	26-46-1	4 mm
26-44-2	3 mm	26-46-2	8 mm
26-44-3	3.5 mm	26-46-3	11 mm
26-44-4	4 mm	26-46-4	16 mm
	26-48	Trephine 10 mm dia.	

27-180 BUNNELL Hand Drill, Lightweight, chrome plated, with '%' (4 mm) capacity stainless Jacobs chuck. Cannulated for use with Kirschner wires and thin drills. Chuck key stores in handle.

KIRSCHNER WIRES 12" (30.5 cm), Trocar Point with Flat Chuck End


27-202	.035" (0.9 mm) diameter
27-204	.045" (1.1 mm) diameter
27-206	.0625" (1.6 mm) diameter

KIRSCHNER WIRES 71/2" (19.1 cm), Trocar Points at both ends

27-208	.035" (0.9 mm) diameter
27-210	.045" (1.1 mm) diameter
27-212	.0625" (1.6 mm) diameter

27-202 to 27-206 27-208 to 27-212





27-213 27-214	MINI-HOHMANN Retractor 6¼" (15.9 cm), blade 6 mm wide MINI-HOHMANN Retractor 6¼" (15.9 cm), blade 8 mm wide
27-215 27-216 27-217 27-218 27-219	HOHMANN Retractor 9¼" (23.5 cm), blade 10 mm wide HOHMANN Retractor 9¼" (23.5 cm), blade 18 mm wide HOHMANN Retractor 9¼" (23.5 cm), blade 25 mm wide HOHMANN Retractor 9¾" (23.5 cm), blade 43 mm wide HOHMANN Retractor 9¾" (24.8 cm), blade 70 mm wide
11-580 11-582	BENNETT Bone Elevator and Retractor 9½" (24.1 cm), small blade 1¾" (4.4 cm) wide BENNETT Bone Elevator and Retractor 9¾" (24.8 cm), large blade 2½" (6.4 cm) wide
27-223	BLOUNT Knee Retractor 7" (17.8 cm), blade 1¾" (4.4 cm) deep, %" (1.6 cm) wide tapering to ¾6" (.48 cm)

For additional Retractors see pages 151 - 187

3/4" (1.9 cm) x 1%" (3.5 cm)

SMILLIE Meniscus Hook Retractor 5%4" (14.6 cm), blade




27-224













HAYS Hand and Finger Retractor 6" (15.2 cm), angled blade 1% " (3.8 cm) long, %" (1 cm) wide tapering to 1% (.32 cm) 27-225

11-500	VOLKMAN Hook 81/2" (21.6 cm), sharp, hollow handle
11-502	VOLKMAN Hook 81/2" (21.6 cm), sharp, ring handle


VOLKMAN Hook Retractor 9" (22.9 cm), sharp point 15 mm deep 27-230 VOLKMAN Hook Retractor 9" (22.9 cm), blunt tip 15 mm deep 27-232

VOLKMAN RETRACTORS 81/2" (21.6 cm), hollow handle

Sharp		Blunt
11-504	2 prongs	11-512
11-506	3 prongs	11-514
11-508	4 prongs	11-516

For additional Retractors see pages 151 - 187





	VOLKMAN RETRACTORS 8" (20.3 cm), Ring Handle		
	Sharp		Blunt
	11-530	2 prongs	11-538
	11-532	3 prongs	11-540
	11-534	4 prongs	11-542
	11-536	6 prongs	11-544
11-349	OLLIER Rake Retract 21/4" (5.7 cm) x 13/4" (ctor 8½" (21.6 cm), blun 4.4 cm)	t, 4 prongs
11-350	ISRAEL Rake Retract 11/2" (3.8 cm) x 15/8" (ctor 9" (22.9 cm), blunt, 4.1 cm)	4 prongs
11-352	ISRAEL Rake Retract 21/4" (5.7 cm) x 13/4" (ctor 9" (22.9 cm), blunt, (4.4 cm)	6 prongs

For additional Retractors see pages 151 - 187















27-260 to 27-264

27-270 to 27-274


27-280 27-282

19-750 to 19-760

19-770 to 19-780

U.S. ARMY Pattern Chisel 61/2" (16.5 cm), 6 mm wide 27-260 27-262 U.S. ARMY Pattern Chisel 61/2" (16.5 cm), 12 mm wide 27-264 U.S. ARMY Pattern Chisel 61/2" (16.5 cm), 18 mm wide U.S. ARMY Pattern Osteotome 61/2" (16.5 cm), 6 mm wide 27-270 U.S. ARMY Pattern Osteotome 61/2" (16.5 cm), 12 mm wide 27-272 27-274 U.S. ARMY Pattern Osteotome 61/2" (16.5 cm), 18 mm wide 27-280 U.S. ARMY Pattern Gouge 61/2" (16.5 cm), 6 mm wide U.S. ARMY Pattern Gouge 61/2" (16.5 cm), 12 mm wide 27-282 19-750 ALEXANDER Chisel 7" (17.8 cm), 4 mm wide ALEXANDER Chisel 7" (17.8 cm), 6 mm wide 19-752 19-754 ALEXANDER Chisel 7" (17.8 cm), 8 mm wide ALEXANDER Chisel 7" (17.8 cm), 10 mm wide 19-756 ALEXANDER Chisel 7" (17.8 cm), 12 mm wide 19-758 19-760 ALEXANDER Chisel 7" (17.8 cm), 14 mm wide 19-770 ALEXANDER Gouge 7" (17.8 cm), 4 mm wide ALEXANDER Gouge 7" (17.8 cm), 6 mm wide 19-772 19-774 ALEXANDER Gouge 7" (17.8 cm), 8 mm wide ALEXANDER Gouge 7" (17.8 cm), 10 mm wide 19-776 19-778 ALEXANDER Gouge 7" (17.8 cm), 12 mm wide 19-780 ALEXANDER Gouge 7" (17.8 cm), 14 mm wide







27-331	HOKE Osteotome 51/2" (14 cm) straight 1/8" (.32 cm)
27-332	HOKE Osteotome 51/2" (14 cm) straight 3/16" (.48 cm)
27-333	HOKE Osteotome 51/2" (14 cm) straight 1/4" (.64 cm)
27-334	HOKE Osteotome 51/2" (14 cm) straight 5/16" (.79 cm)
27-335	HOKE Osteotome 51/2" (14 cm) straight 3/4" (.95 cm)
27-336	HOKE Osteotome 51/2" (14 cm) straight 1/2" (1.3 cm)

HIBBS CHISELS 9" (22.9 cm)		HIBBS OSTEOTOMES 9" (22.9 cm)			
Straight	Size	Straight	Size	Curved	
27-418	1/4" (.64 cm)	27-438	1/4" (.64 cm)	27-438C	
27-420	%" (.95 cm)	27-440	%" (.95 cm)	27-440C	
27-422	½" (1.3 cm)	27-442	½" (1.3 cm)	27-442C	
27-424	%" (1.6 cm)	27-444	%" (1.6 cm)	27-444C	
27-426	3/4" (1.9 cm)	27-446	¾" (1.9 cm)	27-446C	
27-427	7/8" (2.2 cm)	27-447	%" (2.2 cm)	27-447C	
27-428	1" (2.5 cm)	27-448	1" (2.5 cm)	27-448C	
27-430	11/4" (3.2 cm)	27-449	11/8" (2.8 cm)	27-449C	
27-432	11/2" (3.8 cm)	27-450	1¼" (3.2 cm)	27-450C	
		27-452	11/6" (3.8 cm)	27-452C	








27-458C to 27-466C

HIRRS COLIGES 0" (22.0 cm)







27-458 to 27-468

LAMBOTTE OSTEOTOMES 9" (22 9 cm)

HIDDS GOUGES 9 (22.9 GH)			LAMBOT	IL COLLOLOMICO 2 (22.0 0111
Straight	Size	Curved	Straight	Width	Curved
27-458	1/4" (.64 cm)	27-458C	27-488	6 mm (¼")	27-499
27-460	3/8" (.95 cm)	27-460C	27-490	13 mm (½")	27-500
27-462	½" (1.3 cm)	27-462C	27-492	19 mm (¾")	27-502
_	%" (1.6 cm)	27-464C	27-494	25 mm (1")	27-504
27-466	3/4" (1.9 cm)	27-466C	27-496	32 mm (11/4")	27-506
27-468	1" (2.5 cm)		27-498	38 mm (1½")	27-508

Straight	Width	Curved
27-510	1/4" (.64 cm)	27-530
27-512	%" (.95 cm)	27-532
27-514	1/2" (1.3 cm)	27-534
27-516	%" (1.6 cm)	27-536


SMITH-PETERSEN OSTEOTOMES 8" (20.3 cm)

3/4" (1.9 cm) 27-518 27-538 1" (2.5 cm) 27-520 27-540 27-522 11/4" (3.2 cm) 27-542









SMITH.	PFT	FRSEN	GOUGES	Ŗ۳	(20.3 cm)

SWISS PATTERN OSTEOTOMES 5" (12.7 cm)

Straight	Width		Width		Width
27-549	1/4" (.64 cm)	27-599	2 mm	26-602A	12 mm
27-550	%" (.95 cm)	27-600	4 mm	27-603	14 mm
27-551	½" (1.3 cm)	27-601	6 mm	27-604	18 mm
27-553	%" (1.6 cm)	27-601A	8 mm	27-605	20 mm
27-554	¾" (1.9 cm)	27-602	10 mm		
27-556	1" (2.5 cm)				
27-558	11/4" (3.2 cm)				


27-609

Sterilizer – Storage Case 5½" (14 cm) x 6½" (16.5 cm) x ½" (1.3 cm), fitted to hold 9 Swiss Pattern Osteotomes 2 mm to 20 mm

27-610 Set of 9 Swiss Pattern Osteotomes 2 mm to 20 mm, complete in fitted Sterilizer Case







19-790	GERZOG Mallet 71/2" (19.1 cm), head 8 oz. (227 g) lead filled with stainless jacket
19-792	LUCAE Mallet 7½" (19.1 cm), head 9 oz. (255 g) solid stainless
19-793	MiLTEX Mallet 7'/2" (19.1 cm), lightweight 7.5 oz. (213 g), stainless head 1" (2.5 cm) dia. with 2 replaceable nylon caps
19-794	MILTEX Mallet 7½" (19.1 cm), lightweight 9 oz. (255 g), stainless head 1" (2.5 cm) dia. with 1 stainless face and 1 replaceable nylon cap
19-795	Nylon Cap Replacement for Nos. 19-793 and 19-794 Mallets
19-796 19-797	MEAD Mallet 7½" (19.1 cm), head 8 oz. (227 g), with 2 replaceable nylon caps Nylon Cap Replacement for No. 19-796 Mallet
27-660 27-662	MILTEX Mallet 11" (27.9 cm), head 16 oz. (454 g) solid stainless MILTEX Mallet 11" (27.9 cm), head 32 oz. (908 g) solid stainless







27-700



27-702







27-680

27-680

MURPHY-LANE Bone Skid 12" (30.5 cm)

27-700 PUTTI Bone Rasp 10½" (26.7 cm), flat ends, one side curved on flat, other side curved sideways

27-702 PUTTI Bone Rasp 12" (30.5 cm), round and half-round tapered curved ends

27-703 Bone File 9%" (24.8 cm), bayonet shank, half-round blade 3%" (8.9 cm) x %" (1.3 cm)

26-1410 LANGENBECK Periosteal Elevator 7½" (19.1 cm), narrow model, 7 mm wide blunt blade

26-1412 LANGENBECK Periosteal Elevator 7½" (19.1 cm), wide model, 17 mm wide sharo blade

25-28 SEDILLOT Periosteal Elevator 71/2" (19.1 cm), blade 17 mm wide



















27-746

27-752


27-75

27-75

27-820 to 27-826


27-744	SAYRE Periosteal Elevator 6¾" (17.1 cm), one straight blunt and one curved sharp end
27-746	SAYRE Periosteal Elevator 6%" (17.1 cm), one straight blunt and one curved blunt end
27-747	CHANDLER Elevator 71/2" (19.1 cm), blade 1/2" (1.3 cm) tapering to 5/16" (.79 cm)
27-748	CHANDLER Elevator 8" (20.3 cm), blade ¾" (1.9 cm) tapering to ½" (1.3 cm)
27-749	CHANDLER Elevator 91/4" (23.5 cm), blade 1" (2.5 cm) tapering to 11/16" (1.7 cm)
27-750	CHANDLER Elevator 10" (25.4 cm), blade 114" (3.2 cm) tapering to 34" (1.9 cm)
27-752	CREGO Periosteal Elevator 8" (20.3 cm), slightly curved blade 7 mm
27-753	CREGO Periosteal Elevator 8" (20.3 cm), fully curved blade 7 mm
27-754	CREGO Periosteal Elevator 8" (20.3 cm), fully curved blade 14 mm
27-820	KEY Periosteal Elevator 7" (17.8 cm), ¼" (.64 cm) wide
27-821	KEY Periosteal Elevator 7½" (19.1 cm), ¾" (.95 cm) wide
27-822	KEY Periosteal Elevator 7½" (19.1 cm), ½" (1.3 cm) wide
27-824	KEY Periosteal Elevator 8" (20.3 cm), 3/4" (1.9 cm) wide
27-826	KEY Periosteal Elevator 81/4" (21 cm), 1" (2.5 cm) wide











27-924 27-926 27-928

SPRATT (BRUN) BONE CURETTES, Oval Cups

61/2" (16.5 cm) Long		9" Long (22.9 cm)	
19-700	Size 5/0	26-1627	Size 3/0
19-702	Size 4/0	26-1628	Size 2/0
19-704	Size 3/0	26-1629	Size 0
19-706	Size 2/0	26-1630	Size 1
19-708	Size 0	26-1632	Size 2
19-710	Size 1.	26-1634	Size 3
19-712	Size 2	26-1636	Size 4
19-714	Size 3	26-1638	Size 5
19-716	Size 4	26-1640	Size 6
19-718	Size 5		
19-720	Size 6		

VOLKMAN DOUBLE END CURETTES


27-920	51/2" (14 cm) long, oval cups 5 x 10 mm and 6 x 20 mm
27-924	5" (12.7 cm) long, oval cup 6 x 10 mm and round cup 7 mm dia.
27-926	61/2" (16.5 cm) long, oval cup 8 x 14 mm and round cup 10 mm dia.
27-928	8" (20.3 cm) long, oval cup 9 x 18 mm and round cup 11 mm dia.



27-972

27-974








MARTINI Bone Curette 51/2" (14 cm), double end, 4 mm and 21-321 5 mm round cups WILLIGER Bone Curette 51/2" (14 cm), double end, 3 mm and 21-322 4 mm oval cups 27-970 SMILLIE Meniscus Knife 71/4" (18.4 cm), straight chisel blade SMILLIE Meniscus Knife 7" (17.8 cm), curved left 27-972 SMILLIE Meniscus Knife 7" (17.8 cm), curved right 27-974 27-980 LOWE-BRECK Cartilage Knife 101/2" (26.7 cm), straight edge with guards DOWNING Cartilage Knife 10" (25.4 cm), concave edge with guards 27-984









27-990 27-992	MARTIN Cartilage Clamp 7" (17.8 cm) MARTIN Cartilage Clamp 8" (20.3 cm)
27-1000	MARTIN Cartilage Scissors 8" (20.3 cm), serrated non-slip cutting edges
27-1016	MILTEX Tendon Pulling Forceps, 5" (12.7 cm) curved shaft, alligator type
27-1020 27-1022 27-1024 27-1026 27-1028	BUNNELL Tendon Stripper 6" (15.2 cm), inside dia. 3 mm BUNNELL Tendon Stripper 6" (15.2 cm), inside dia. 4 mm BUNNELL Tendon Stripper 6" (15.2 cm), inside dia. 5 mm BUNNELL Tendon Stripper 6" (15.2 cm), inside dia. 6 mm BUNNELL Tendon Stripper 6" (15.2 cm), inside dia. 7 mm














22-480 to 22-482

19-804

19-808

26-198

DEAN (BLUMENTHAL) RONGEURS 6" (15.2 cm), 3 x 10 mm bite tapered to 4 mm

22-480

slightly curved jaws

22-482

jaws 90° angled

22-481

CICHERELLI (REINER) Rongeur 7" (17.8 cm), slightly curved jaws 3 x 8 mm bite tapered to 4 mm

19-804

HARTMAN Rongeur 53/4" (14.6 cm), strong curve, 4 x 12 mm bite

19-808

BANE Rongeur 7" (17.8 cm), curved jaws, 5 x 15 mm bite tapered to 7 mm

26-198

BACON Rongeur 81/4" (21 cm), angled on side, 5 mm bite


26-190

LUER Rongeur 7" (17.8 cm), straight jaws 8 mm bite LUER Rongeur 7" (17.8 cm), curved jaws 8 mm bite


26-192 22D-420

BLUMENTHAL Rongeur 6" (15.2 cm), 45° angled jaws

3 x 12 mm bite








22D-420











22-484

27-1295 27-1297

19-848 19-849

19-850 to 19-856

22-484	MEAD Rongeur 6½" (16.5 cm), slightly angled jaws 4 x 12 mm bite
--------	---

KLEINERT-KUTZ Bone Rongeur 6" (15.2 cm), slightly curved jaws 2 x 13 mm 27-1295 KLEINERT-KUTZ Synovectomy Rongeur 6" (15.2 cm), full curved jaws 3 x 8 mm 27-1297

RUSKIN RONGEURS

Straight	Sizes	Curved
19-848	6" (15.2 cm) long, 3 mm bite	19-849
19-850	71/4" (18.4 cm) long, 5 mm bite	19-852
19-854	71/4" (18.4 cm) long, 7 mm bite	19-856

ZAUFEL-JANSEN Rongeur 7" (17.8 cm), slightly curved jaws 5 mm wide 19-840




19-840













25-511

25-456

BEYER Rongeur 7" (17.8 cm), slightly curved jaws 3.5 mm wide 19-844


LEKSELL Laminectomy Rongeur 91/2" (24.1 cm), standard curve, 3 mm bite 25-510 LEKSELL Laminectomy Rongeur 91/2" (24.1 cm), standard curve, 8 mm bite 25-511


STILLE-LUER Rongeur 8¾" (22.2 cm), straight jaws 9 x 15 mm 25-454 STILLE-LUER Rongeur 8¾" (22.2 cm), curved jaws 9 x 15 mm 25-456 STILLE-LUER Duckbill Rongeur 9" (22.9 cm), angled on side 25-458

jaws 6 x 12 mm














26-264 to 26-266

27-1390 to 27-1394

27-1400 to 27-1404


26-264	ECHLIN Duckbill Rongeur 9" (22.9 cm), 2 mm bite
26-265	ECHLIN Duckbill Rongeur 9" (22.9 cm), 3 mm bite
26-266	ECHLIN Duckbill Rongeur 9" (22.9 cm), 4 mm bite


40-227-SS LITTAUER-LISTON Bone Cutting Forceps 6" (15.2 cm), straight 20 mm blades

LISTON BONE CUTTING FORCEPS

Straight	Size	Angled on Flat
27-1390	5½" (14 cm)	27-1400
27-1392	7½" (19.1 cm)	27-1402
27-1394	81/2" (21.6 cm)	27-1404














21-626 21-627


25-385

RUSKIN BONE CUTTING FORCEPS

Straight Size Angled 21-626 standard blades 51/2" (14 cm) 21-627 standard blades 71/2" (19.1 cm) 25-384 25-385 21-628 delicate blades 71/2" (19.1 cm)

25-390 STILLE-LISTON Bone Cutting Forceps 103/4" (27.3 cm), straight STILLE-LISTON Bone Cutting Forceps 10¾" (27.3 cm), curved 25-392 25-396

STILLE-HORSLEY Bone Cutting Forceps 101/2" (26.7 cm), angled blades




463



26-138


65-28



screws, burs etc., during cleaning, rinsing or sterilization








27-1550	LANGENBECK Metacarpal Saw 91/4" (23.5 cm), blade 41/2" (11.4 cm) long
40-2000	MILTEX Bone Saw 6" (15.2 cm), blade 2" (5.1 cm)
27-1560 27-1561	SATTERLEE Bone Saw 12" (30.5 cm), light metal handle with stainless 9" (22.9 cm) blade Blade for No. 27-1560 Satterlee Saw
27-1564 27-1565-12 27-1565-18	CHARRIERE Bone Saw 14" (35.6 cm), with one each 12 mm and 18 mm stainless blades 12 mm wide blade for No. 27-1564 saw 18 mm wide blade for No. 27-1564 saw
27-1600 27-1602	LISTON Amputating Knife 11%" (29.8 cm), with 6½" (16.5 cm) blade LISTON Amputating Knife 13%" (34.9 cm), with 8" (20.3 cm) blade

Cast Instruments














27-3002 U.S. Pattern Cast Knife 5½" (14 cm), wood handle, 1½" (3.8 cm) blade
27-3004 REINER Cast Knife 7" (17.8 cm), metal handle, 1½" (3.8 cm) blade
27-3020 ENGEL Cast Saw 6" (15.2 cm)
27-3022 WIGMORE Cast Saw 7¾" (19.7 cm)

5-800 Bandage and Utility Scissors 8½" (21.6 cm), with fluoride coated blades to help prevent adhesive or fiberglass casting material from sticking to the blade surfaces. Neon green plastic finger

rings, serrated blade. Autoclavable.


5-998 Utility and Light Cast Shears 6½" (16.5 cm), black plastic finger rings, serrated blade. Autoclavable.

5-1000 Utility and Light Cast Shears 7½" (19.1 cm), black plastic finger rings, serrated blade. Autoclavable.





Cast Instruments



5-568 ESMARCH Cast Shears 8" (20.3 cm)

5-547 HERCULES Cast Shears 7½" (19.1 cm), serrated blade

27-3070 BRUNS Cast Shears 91/4" (23.5 cm), serrated blade

5-600 MILTEX Cast and Utility Shears 81/2" (21.6 cm), with wire

cutter and serrated blade

27-3059 STILLE Cast Shears 8" (20.3 cm), Pediatric Size,

serrated blades


27-3060

STILLE Cast Shears 11" (27.9 cm), Medium Size,

serrated blades

27-3062 STILLE Cast Shears 15" (38.1 cm), Standard Size,

serrated blades



Cast Instruments






27-3080	WOLFE-BOEHLER Cast Breaker 10" (25.4 cm), jaws 15 mm wide with deep serrations
27-3090	BEESON Cast Spreader 12" (30.5 cm), blades 2" (5.1 cm) wide
27-3094	HENNIG Cast Spreader 11" (27.9 cm), blades 2" (5.1 cm) wide
27-3100	MILTEX 3-Prong Cast Spreader 9" (22.9 cm), with spring action for one-hand operation








28-1	PRATT Rectal Speculum 8½" (21.6 cm) long, with set screw, blades 3½" (8.9 cm) long x 1" (2.5 cm) wide at distal end
28-2	SIMS Rectal Speculum 6" (15.2 cm) long, with set screw, fenestrated blades $3\frac{1}{2}$ " (8.9 cm) long x $\%$ " (1.6 cm) wide at distal end
28-4	BODENHAMMER Rectal Speculum 6" (15.2 cm) long, with set screw, blades 3½" (8.9 cm) long x ¾" (1.9 cm) wide at distal end
28-6	BARR Anal Retractor 8" (20.3 cm), self-retaining, blades 234" (7 cm) x 3/8" (2.2 cm)










28-65

28-70 to 28-

28-64	SMITH (BUIE) Anal Retractor, self-retaining, swivel blades 2½" (6.4 cm) x ½" (2.2 cm) open to 4" (10.2 cm)
28-65	SMITH (BUIE) Anal Retractor, self-retaining, extra long swivel blades 4" (10.2 cm) x 1/8" (2.2 cm) open to 4" (10.2 cm)
28-70	SAWYER Rectal Retractor 11" (27.9 cm), blade %" (2.2 cm) x 21/2" (6.4 cm)
28-72	SAWYER Rectal Retractor 11" (27.9 cm), blade 11/8" (2.8 cm) x 3" (7.6 cm)
28-74	SAWYER Rectal Retractor 11" (27.9 cm), blade 11/2" (3.8 cm) x 31/2" (8.9 cm)

28-8 COOK Rectal Retractor 8" (20.3 cm), triple bladed, self-retaining, 4" (10.2 cm) deep blades, lower blade ½" (1.3 cm) wide





28-10	PENNINGTON Rectal Speculum 3%" (8.6 cm) long x 1" (2.5 cm) O.D.,
	with boilable obturator, 43/4" (12.1 cm) long handle, chrome

28-11 Obturator for No. 28-10 Pennington Speculum

28-12 FANSLER-IVES Rectal Speculum 31/4" (8.3 cm) long x 11/6" (2.8 cm) O.D.,

with boilable obturator, 23/4" (7 cm) long handle, chrome

28-14 Obturator for No. 28-12 Fansler-Ives Speculum

28-20 HIRSCHMAN Anoscope, Small size 21/4" (5.7 cm) long x 1/16" (1.4 cm) I.D.


28-22 HIRSCHMAN Anoscope, Medium size 25/6" (6.7 cm) long x 11/46" (1.7 cm) l.D.


28-24 HIRSCHMAN Anoscope, Large size $3\frac{1}{2}$ " (8.9 cm) long x $\frac{7}{6}$ " (2.2 cm) l.D.




28-20 to 28-24











KELLY Anoscope 2" (5.1 cm) long x 1/8" (2.2 cm) O.D.


KELLY Proctoscope 51/2" (14 cm) long x 1/4" (2.2 cm) O.D.

28-40 28-44






28-48	HIRSCHMAN Proctoscope 5" (12.7 cm) long x $^{15}\!\!/\!\! s^{e}$ " (2.4 cm) O.D.
28-50	FANSLER Operating Speculum 2%" (6.1 cm) long x 1% " (3.5 cm) O.D., slotted tube
28-52	BARR-SHUFORD Rectal Speculum 4% " (11.4 cm) tapered from %" (1.6 cm) to 2" (5.1 cm), chrome
28-68	HILL-FERGUSON Retractor 8½" (21.6 cm), small blade, 2.2 x 6.4 cm
28-69	HILL-FERGUSON Retractor 8¾" (22.2 cm), large blade, 3 x 7.6 cm
28-79	CHELSEA EATON Anal Speculum 2%" (6.7 cm) long x %" (2.2 cm) O.D.
28-80	CHELSEA EATON Anal Speculum 2%" (6.7 cm) long x 11/4" (3.2 cm) O.D.




Rectal Instruments 28-100 28-104 28-106 28-108 28-109 28-120 28-128 28-122 28-130 28-124 28-100 LARRY Probe and Grooved Director 7" (17.8 cm) 28-104 BARR Fistula Probe, with 8" (20.3 cm) sterling shaft, eye at tip, 111/2" (29.2 cm) overall BUIE Fistula Probe, with 51/2" (14 cm) sterling shaft, 63/4" (17.1 cm) overall 28-106 PRATT Rectal Probe, with 7¾" (19.7 cm) sterling shaft, 11" (27.9 cm) overall 28-108 28-109 BARR Double Ended Probe 7" (17.8 cm), rigid 28-120 PRATT Crypt Hook 81/2" (21.6 cm), narrow style 28-122 PRATT Crypt Hook 81/2" (21.6 cm), short style 28-124 PRATT Crypt Hook 81/2" (21.6 cm), long style 28-128 STEWART Crypt Hook 81/2" (21.6 cm), short style 28-130 STEWART Crypt Hook 81/2" (21.6 cm), long style 28-131 ROSSER Crypt Hook 6¾" (17.1 cm) 28-140 BARR Crypt Hook 10" (25.4 cm), short and long hook ends 28-131 28-140



28-150 28-151	McGIVNEY Hemorrhoidal Ligator, improved model with offset handle for better vision and comfort. 7" (17.8 cm) working length. Can be disassembled entirely for cleaning. Supplied with drum loading cone and 100 Latex O-Rings. Latex O-Rings for No. 28-150 Ligator. Bag with 100 rings.
28-152	McGIVNEY Hemorrhoid Grasping Forceps 7% " (19.1 cm), angled shanks. Specially designed for use with McGivney Ligator.
28-153	Loading Cone for No. 28-150 Ligator







7-106	PRATT T-Shape Tissue Grasping Forceps 6" (15.2 cm)
7-220	PENNINGTON Tissue Grasping Forceps 6" (15.2 cm), triangular jaws
16-30	ADAIR-ALLIS Tissue Forceps 6% " (16.5 cm), 9 x 10 teeth
16-3R	THOMS, ALLIS Tissue Engage 8" (20.3 cm), 6 x 7 teeth












16-10 16-24	ALLIS Tissue Forceps 7½" (18.4 cm), 5 x 6 teeth, straight ALLIS Tissue Forceps 7½" (18.4 cm), 5 x 6 teeth, 90° angled jaws	
28-196	BUIE Pile Clamp 8½" (21.6 cm), curved angiotribe jaws	
28-200	FERGUSON Tenaculum Forceps 6¾" (17.1 cm), straight jaws	
7-260 7-262 7-264 7-266	FERGUSON Angiotribe Forceps 6½" (16.5 cm), straight FERGUSON Angiotribe Forceps 6½" (16.5 cm), curved FERGUSON Angiotribe Forceps 7½" (19.1 cm), straight FERGUSON Angiotribe Forceps 7½" (19.1 cm), curved	
















28-240	MILLER Rectal Scissors 634" (17.1 cm), bayonet shape, delicate serrated blades with blunt points
28-244	PRATT Rectal Scissors 81/4" (21 cm), S-shaped, blunt points
5-250 5D-251 5-252	KELLY Rectal Scissors 6½" (15.9 cm), sharp points, straight KELLY Rectal Scissors 6½" (15.9 cm), sharp points, straight, with one serrated blade KELLY Rectal Scissors 6½" (15.9 cm), sharp points, curved
5D-253	KELLY Rectal Scissors 614" (15.9 cm), sharp points, curved, with one serrated blade










28-294 to 28-298

28-280	BUIE Biopsy Forceps 13" (33 cm) shaft, 3.5 mm bite			
28-294 28-296	YEOMAN Biopsy Forceps 9" (22.9 cm) shaft, 3.5 x 10 mm bite, with basket YEOMAN Biopsy Forceps 13" (33 cm) shaft, 3.5 x 10 mm bite, with basket			
28-298	YEOMAN Biopsy Forceps 15" (38.1 cm) shaft, 3.5 x 10 mm bite, with basket			





YEOMAN BIOPSY FORCEPS with ROTATING SHAFTS, Straight 4 x 8 mm Bite, With Basket

Complete with Ring Handle	Shaft Length	Shaft Only	
28-300	10" (25.4 cm)	28-306	
28-302	14" (35.6 cm)	28-308	
28-303	14" (35.6 cm) with teeth	28-309	
28-304	16" (40.6 cm)	28-310	




TURRELL BIOPSY FORCEPS with ROTATING SHAFTS, Angled Jaws 4 x 8 mm Bite, With Basket

Complete with Ring Handle	Shaft Length	Shaft Only
28-320	10" (25.4 cm)	28-326
28-322	14" (35.6 cm)	28-328
28-323	14" (35.6 cm) with teeth	28-329
28-324	16" (40.6 cm)	28-330

28-340 Ring Handle fitting Yeoman and Turrell Biopsy Forceps Nos. 28-300 to 28-330

28-341 Pistol Grip Handle fitting Yeoman and Turrell Biopsy Forceps Nos. 28-300 to 28-330

NOTE: All parts of Nos. 28-300 through 28-341 are interchangeable



MILTEX IMPROVED BIOPSY FORCEPS with ROTATING SHAFTS in Pistol Grip Spring Handle. Shafts easily disassembled from handle. Constructed without small assembly parts. All shafts and handles are interchangeable.


Jaws

Shaft

Complete with

	Pistol Grip Handle	Length	3.5 x 10 mm Bite	Shaft only	
	28-380	25.4 cm (10")	straight	28-410	
	28-382	35.6 cm (14")	straight	28-412	
	28-384	25.4 cm (10")	angled	28-414	
	28-386	35.6 cm (14")	angled	28-416	
28-400	Pistol Grip Spring 28-416. This han	Pistol Grip Spring Handle for Miltex Improved Biopsy Forceps Nos. 28-380 through 28-416. This handle also fits Uterine Biopsy Forceps Nos. 30-1490 to 30-1495.			
28-360	Rectal Cotton Ca	Rectal Cotton Carrier, 15" (38.1 cm) shaft, with spring action			
28-370	BUIE Suction Tub	BUIE Suction Tube 16" (40.6 cm), with finger valve, 15 Fr. (5 mm)			






١	/ A Ki	DUDEN	HOSTLIDAL	COLINDO	11" (27.9 cm)	
١	IΛΩN	RIIHEN		SCH INI IS	11" (2) U.cmi	

29-2 29-4	Set of 6 Sounds in Khaki Roll, even sizes 14 Fr. (4.6 mm) to 24 Fr. (8 mm) Set of 12 Sounds in Khaki Roll, even sizes 8 Fr. (2.6 mm) to 30 Fr. (10 mm)					
29-6-8	8 Fr. (2.6 mm)	29-6-20	20 Fr. (6.6 mm)	29-6-32	32 Fr. (10.6 mm)	
29-6-10	10 Fr. (3.3 mm)	29-6-22	22 Fr. (7.3 mm)	29-6-34	34 Fr. (11.2 mm)	
29-6-12	12 Fr. (4 mm)	29-6-24	24 Fr. (8 mm)	29-6-36	36 Fr. (11.9 mm)	
29-6-14	14 Fr. (4.6 mm)	29-6-26	26 Fr. (8.6 mm)	29-6-38	38 Fr. (12.5 mm)	
29-6-16	16 Fr. (5.3 mm)	29-6-28	28 Fr. (9.2 mm)	29-6-40	40 Fr. (13.2 mm)	
29-6-18	18 Fr. (6 mm)	29-6-30	30 Fr. (10 mm)			
	DITTEL URETHRAL SOUNDS 111/2" (29.2 cm), STRAIGHT					

29-8-8	8 Fr. (2.6 mm)	29-8-18	18 Fr. (6 mm)	29-8-28	28 Fr. (9.2 mm)
29-8-10	10 Fr. (3.3 mm)	29-8-20	20 Fr. (6.6 cm)	29-8-30	30 Fr. (9.9 mm)
29-8-12	12 Fr. (4 mm)	29-8-22	22 Fr. (7.3 mm)	29-8-32	32 Fr. (10.6 mm)
29-8-14	14 Fr. (4.6 mm)	29-8-24	24 Fr. (8 mm)	29-8-34	34 Fr. (11.2 mm)
29-8-16	16 Fr. (5.3 mm)	29-8-26	26 Fr. (8.6 mm)		, ,











LEFORT URETHRAL SOUNDS 11" (27.9 cm), With Standard 1-56 Thread for Filiform (not included)

29-10-8	8 Fr. (2.6 mm)	29-10-20	20 Fr. (6.6 mm
29-10-10	10 Fr. (3.3 mm)	29-10-22	22 Fr. (7.3 mm
29-10-12	12 Fr. (4 mm)	29-10-24	24 Fr. (8 mm)
29-10-14	14 Fr. (4.6 mm)	29-10-26	26 Fr. (8.6 mm
29-10-16	16 Fr. (5.3 mm)	29-10-28	28 Fr. (9.2 mm
29-10-18	18 Fr. (6 mm)	29-10-30	30 Fr. (10 mm)

McCREA INFANT URETHRAL SOUNDS 7" (17.8 cm) CURVED, PLAIN TIP Also used as Female Sounds

29-14-8	8 Fr. (2.6 mm)	29-14-24	24 Fr. (8 mm)
29-14-10	10 Fr. (3.3 mm)	29-14-26	26 Fr. (8.6 mm)
29-14-12	12 Fr. (4 mm)	29-14-28	28 Fr. (9.2 mm)
29-14-14	14 Fr. (4.6 mm)	29-14-30	30 Fr. (10 mm)
29-14-16	16 Fr. (5.3 mm)	29-14-32	32 Fr. (10.6 mm)
29-14-18	18 Fr. (6 mm)	29-14-34	34 Fr. (11.2 mm)
29-14-20	20 Fr. (6.6 mm)	29-14-36	36 Fr. (11.9 mm)
29-14-22	22 Fr. (7.3 mm)		

WALTHER URETHRAL SOUNDS 11" (27.9 cm) Shaft Tapered Downward toward Handle

29-18-16	16 Fr. (5.3 mm)	29-18-28	28 Fr. (9.2 mm)
29-18-18	18 Fr. (6 mm)	29-18-30	30 Fr. (9.9 mm)
29-18-20	20 Fr. (6.6 mm)	29-18-32	32 Fr. (10.6 mm)
29-18-22	22 Fr. (7.3 mm)	29-18-34	34 Fr. (11.2 mm)
29-18-24	24 Fr. (8 mm)	29-18-36	36 Fr. (11.9 mm)
20-18-26	26 Fr /8.6 mm\		








29-32-8 29-32-10	Female Catheter 5¾" (14.6 cm), Size 8 Fr. (2.6 mm) Female Catheter 5¾" (14.6 cm), Size 10 Fr. (3.3 mm)
29-32-12	Female Catheter 5¾" (14.6 cm), Size 12 Fr. (4 mm)
29-32-14	Female Catheter 5¾" (14.6 cm), Size 14 Fr. (4.6 mm)
29-32-16	Female Catheter 5%" (14.6 cm), Size 16 Fr. (5.3 mm)
29-32-18	Female Catheter 5¾" (14.6 cm), Size 18 Fr. (6 mm)

WALTHER FEMALE DILATOR-CATHETERS 51/4" (13.3 cm)

29-33-12	12 Fr. (4 mm)	29-33-22	22 Fr. (7.3 mm)	29-33-32	32 Fr. (10.6 mm)
29-33-14	14 Fr. (4.6 mm)	29-33-24	24 Fr. (8 mm)	29-33-34	34 Fr. (11.2 mm)
29-33-16	16 Fr. (5.3 mm)	29-33-26	26 Fr. (8.6 mm)	29-33-36	36 Fr. (11.9 mm)
29-33-18	18 Fr. (6 mm)	29-33-28	28 Fr. (9.2 mm)	29-33-38	38 Fr. (12.5 mm)
29-33-20	20 Fr (6.6 mm)	29-33-30	30 Fr. (10 mm)		







VAN BUREN Catheter Guide 17" (43.2 cm), size 8 Fr. (2.6 mm) 29-48


29-72 KOLLMANN Urethral Dilator 10" (25.4 cm), curved, dilates from 20 Fr. (6.6 mm) to 45 Fr. (14.9 mm), uses Handle No. 29-78

KOLLMANN Dilator Handle, adjustable from 20 Fr. (6.6 mm) to 45 Fr. (14.9 mm), 29-78 for use with Dilator No. 29-72

OTIS Urethrotome, straight, with 2 knife blades, adjusts from 20 Fr. (6.6 mm) to 45 Fr. (14.9 mm) 29-80 11" (27.9 cm) working length. In fitted wooden case

OTIS Urethrotome Replacement Blade

29-81



		()		I	
29-102	MATHIEU Urethral Forceps 11½" (29.2 cm) shaft, serrated jaws 3.5 mm x 17 mm			ŭ	
28-280	BRAASCH (BUIE) Bladder Biopsy Forceps 13" (33 cm) shaft, 3.5 mm cup jaws				
11-390	MASSON-JUDD Bladder Retractor, with 2 swivel blades 1" (2.5 cm) x 2¾" (7 cm)		11-392		

MASSON-JUDD Bladder Retractor, with 2 swivel blades and center blade 1" (2.5 cm) x 2¾" (7 cm)

486

11-392






8-102	STRATTE Needle Holder 9" (22.9 cm), double curved jaws and bent shanks
29-232	MILLIN Capsule Holding Forceps $9\frac{1}{2}$ " (24.1 cm), jaws 11 mm wide, 9 x 10 teeth
29-280	WERTHEIM Pedicle Clamp 9% " (24.8 cm), curved jaws 2" (5.1 cm) long
29-282	WERTHEIM-CULLEN Pedicle Clamp 8½" (21.6 cm), right angle jaws 2" (5.1 cm) long


29-284 RANDALL Kidney Stone Forceps 9¼" (23.5 cm), quarter curved 29-286 RANDALL Kidney Stone Forceps 8¼" (21 cm), half curved

29-288 RANDALL Kidney Stone Forceps 7%" (19.7 cm), three quarter curved

29-290 RANDALL Kidney Stone Forceps 7½" (19.1 cm), full curved









29-300	YOUNG Renal Pedicle Clamp 9" (22.9 cm), curved, longitudinal serrations with cross serrations at tip
29-302	STILLE Pattern Vessel Clamp 9" (22.9 cm)
29-304	HERRICK Kidney Pedicle Clamp 9½" (24.1 cm), double angle jaws, serrated
29-308	MAYO (GUYON) Vessel Clamp 9¼" (23.5 cm), double curved jaws, serrated
29-310TC	BIHRLE Dorsal Clamp — Tungsten Carbide Needle Holder 11" (27.9 cm). Designed to aid in the fixation and suturing of the dorsal vein complex during retropubic prostatectomy for prostate cancer during cystectomy.















GRAVES VAGINAL SPECULA, ORIGINAL MODEL

30-10 30-15 30-20	Small Size, 3/4" (1.9 cm) x 3" (7.6 cm) Medium Size 13/4" (3.5 cm) x 4" (10.2 cm) Large Size 11/2" (3.8 cm) x 4/2" (11.4 cm)
30-10B 30-15B 30-20B	GRAVES Vaginal Speculum, Small Size ¾" (1.9 cm) x 3" (7.6 cm), ebony finish GRAVES Vaginal Speculum, Medium Size 1¾" (3.5 cm) x 4" (10.2 cm), ebony finish GRAVES Vaginal Speculum, Large Size 1½" (3.8 cm) x 4½" (11.4 cm), ebony finish
30-25	GRAVES Improved Vaginal Speculum, Medium Size 1%" (3.5 cm) x 4" (10.2 cm), wide angle blades
30-30	GRAVES Open Side Vaginal Speculum, Medium Size 1%" (3.5 cm) x 4" (10.2 cm), wide angle blades





PEDERSON VAGINAL SPECULA, NARROW BLADES


30-50	Small Size, ½" (1.3 cm) x 3" (7.6 cm)
30-55	Medium Size %" (2.2 cm) x 4" (10.2 cm)
30-60	Large Size 1" (2.5 cm) x 4%" (12.1 cm)
30-62	PEDERSON Extra Narrow Vaginal Speculum %" (1.6 cm) x 4" (10.2 cm) blades

CUSCO VAGINAL SPECULA, with Folding Handles


With Center Screw		Blade Sizes	With Side Screw
_	Extra Narrow	3/4" (1.9 cm) x 31/4" (8.3 cm)	30-70
30-75	Small	11/4" (3.2 cm) x 31/2" (8.9 cm)	30-90
30-80	Medium	1%" (3.5 cm) x 4" (10.2 cm)	30-95
30-85	Large	1%" (3.5 cm) x 41/2" (11.4 cm)	30-100







30-120 to 30-130



30-140 to 30-142


30-120	COLLIN Vaginal Speculum, Small Size	11/4" (3.2 cm) x 33/4" (9.5 cm)
30-125	COLLIN Vaginal Speculum, Medium Size	e 1%" (3.5 cm) x 4" (10.2 cm)

30-130 COLLIN Vaginal Speculum, Large Size 1%" (4.1 cm) x 41/2" (11.4 cm)


DeVILBISS Vaginal Specula Two sets of screws permit parallel,or bivalve opening of the blades, for examination or operative use


30-140	Small Size, 3½" (8.9 cm) x ¾" (2.2 cm) blades
30-141	Medium Size, 4" (10.2 cm) x 11/4" (3.2 cm) blades
30-142	Large Size, 5" (12.7 cm) x 1%" (3.5 cm) blades
30-180	GUTMANN Vaginal Speculum, Examining Type,

1%" (3.5 cm) x 4" (10.2 cm), self-retaining









		W/I 7
30-174	SIMS Vaginal Speculum, Single End, Small Size 11/16" (2.7 cm) x 21/2" (6.4 cm)	W
30-176	SIMS Vaginal Speculum, Single End, Medium Size 11/4" (3.2 cm) x 3" (7.6 cm)	
30-178	SIMS Vaginal Speculum, Single End, Large Size 1½" (3.8 cm) x 3½" (8.9 cm)	S. COLUMN
30-210	SIMS Vaginal Speculum, Double End, Small Size 1" (2.5 cm) x 2½" (6.4 cm) and 1½" (3.2 cm) x 3" (7.6 cm)	
30-215	SIMS Vaginal Speculum, Double End, Medium Size 11/6" (2.8 cm) x 3" (7.6 cm) and 11/4" (3.2 cm) x 31/6" (8.9 cm)	
30-220	SIMS Vaginal Speculum, Double End, Large Size 1%" (3.5 cm) x 3½" (8.9 cm) and 1½" (3.8 cm) x 4" (10.2 cm)	
	AUVARD WEIGHTED VAGINAL SPECULA 9" (22.9 cm)	
30-186	2 lbs. (909 g), straight blade 11/2" (3.8 cm) x 31/4" (8.3 cm), slightly angled (75°)	
30-187	2.5 lbs. (1136 g), straight blade 11/2" (3.8 cm) x 31/4" (8.3 cm), slightly angled (75°)	
30-188	3 lbs. (1363 g), straight blade 1½" (3.8 cm) x 3¼" (8.3 cm), slightly angled (75°)	
20.102	2.5 lbs (1136 a) blade ("(10.2 cm) long v 13/4" (2.9 cm) wide taparing to	30-195

2.5 lbs. (1136 g), blade $4^{\rm m}$ (10.2 cm) long x 1%" (3.8 cm) wide, tapering to 1%" (4.4 cm) width at tip, blade at acute angle (45°)


STEINER-AUVARD weighted Vaginal Speculum, 2.5 lbs. (1136 g), slightly curved blade 5% " (14 cm) x 11% (3.2 cm), at 90° angle


30-193

30-195













30-300 JACKSON Vaginal Retractor 7" (17.8 cm), Blade 11/2" (3.8 cm) x 3" (7.6 cm) 30-305 JACKSON Vaginal Retractor 7" (17.8 cm), Blade 11/2" (3.8 cm) x 31/2" (8.9 cm) 30-310 JACKSON Vaginal Retractor 7" (17.8 cm), Blade 11/2" (3.8 cm) x 4" (10.2 cm) 30-320 EASTMAN Vaginal Retractor 7" (17.8 cm), Blade 1½" (3.8 cm) x 3" (7.6 cm) 30-325 EASTMAN Vaginal Retractor 7" (17.8 cm), Blade 11/2" (3.8 cm) x 31/2" (8.9 cm) 30-330 EASTMAN Vaginal Retractor 7" (17.8 cm), Blade 11/2" (3.8 cm) x 4" (10.2 cm) 30-340 DE LEE Universal Retractor 91/2" (24.1 cm), Blade 21/4" (7 cm) x 2" (5.1 cm) 30-395 DOYEN Vaginal Retractor 9" (22.9 cm), Blade 1¾" (4.4 cm) x 2½" (6.4 cm) 30-406 DOYEN Vaginal Retractor 9" (22.9 cm), Blade 21/6" (5.4 cm) x 3" (7.6 cm) 30-415 DOYEN Vaginal Retractor 9" (22.9 cm), Blade 21/4" (5.7 cm) x 31/2" (8.9 cm)



■ Gynecological and Obstetrical Instruments ■ Market Mark







30-370	O'SULLIVAN-O'CONNOR Self-Retaining Vaginal Retractor, complete with 3 interchangeable blades ½" (1.3 cm) x 1½" (4.7 cm), 1½" (3.8 cm) x 3½" (8.9 cm) and 1¾" (4.4 cm) x 3½" (8.9 cm)			
11-180	HEANEY Hysterectomy Retractor 101/2" (26.7 cm), Blade 1" (2.5 cm) x 13/4" (4.4 cm)			
11-182	HEANEY Hysterectomy Retractor 101/2" (26.7 cm), Blade 11/6" (2.8 cm) x 21/6" (7.3 cm)			
11-184	HEANEY Hysterectomy Retractor 10½" (26.7 cm), Blade 11%" (2.8 cm) x 31½" (8.3 cm)			
11-185	HEANEY-SIMON Vaginal Retractor 11%" (29.2 cm). Blade 1" (2.5 cm) x 4%" (11.4 cm)			








11-374






11-365 Pediatric GELPI Perineal Retractor 51/2" (14 cm), with grip-lock 11-366 GELPI Perineal Retractor 6¾" (17.1 cm), with grip-lock 11-368 GELPI Perineal Retractor 6¾" (17.1 cm), with ball stops and grip-lock

RIGBY Vaginal Retractor 63/4" (17.1 cm), with grip-lock







HANK UTERINE DILATORS 11" (27.9 CM), DOUBLE END

30-500 Set of 6 Dilators in Khaki Roll, Sizes 11-12, (3.6-4 mm), 13-14 (4.3-4.6 mm), 15-16 (5-5.3 mm), 17-18 (5.6-6 mm), 19-20 (6.3-6.6 mm), and 21-22 (7-7.3 mm) French


30-505-78	7-8 Fr. (2.3-2.6 mm)	30-505-1516	15-16 Fr. (5-5.3 mm)	30-505-2324	23-24 Fr. (7.6-8 mm)
30-505-910	9-10 Fr. (3-3.3 mm)	30-505-1718	17-18 Fr. (5.6-6 mm)	30-505-2526	25-26 Fr. (8.3-8.6 mm)
30-505-1112	11-12 Fr. (3.6-4 mm)	30-505-1920	19-20 Fr. (6.3-6.6 mm)	30-505-2728	27-28 Fr. (8.9-9.2 mm)
30-505-1314	13-14 Fr. (4.3-4.6 mm)	30-505-2122	21-22 Fr. (7-7.3 mm)	30-505-2930	29-30 Fr. (9.6-10 mm)

HEGAR UTERINE DILATORS 7" (17.8 cm), SINGLE END


Set of 14 Dilators in Khaki Holl, Sizes 4 to 17 mm				
Set of 26 Dilators in Khaki Roll, Sizes 1 to 26 mm				
1 mm	30-530-85	8.5 mm	30-530-16	16 mm
1.5 mm	30-530-9	9 mm	30-530-165	16.5 mm
2 mm	30-530-95	9.5 mm	30-530-17	17 mm
2.5 mm	30-530-10	10 mm	30-530-175	17.5 mm
3 mm	30-530-105	10.5 mm	30-530-18	18 mm
3.5 mm	30-530-11	11 mm	30-530-185	18.5 mm
4 mm	30-530-115	11.5 mm	30-530-19	19 mm
4.5 mm	30-530-12	12 mm	30-530-20	20 mm
5 mm	30-530-125	12.5 mm	30-530-21	21 mm
5.5 mm	30-530-13	13 mm	30-530-22	22 mm
6 mm	30-530-135	13.5 mm	30-530-23	23 mm
6.5 mm	30-530-14	14 mm	30-530-24	24 mm
7 mm	30-530-145	14.5 mm	30-530-250	25 mm
7.5 mm	30-530-150	15 mm	30-530-26	26 mm
8 mm	30-530-155	15.5 mm		
	Set of 26 Dilators in Kha 1 mm 1.5 mm 2 mm 2.5 mm 3 mm 3.5 mm 4 mm 4.5 mm 5 mm 5 mm 6 mm 6.5 mm 7 mm 7.5 mm	Set of 26 Dilators in Khaki Roll, Sizes 1 mm 30-530-85 1.5 mm 30-530-95 2 mm 30-530-95 2.5 mm 30-530-10 3.5 mm 30-530-115 4.5 mm 30-530-115 4.5 mm 30-530-12 5 mm 30-530-13 6 mm 30-530-13 6.5 mm 30-530-14 7 mm 30-530-145 7.5 mm 30-530-150	Set of 26 Dilators in Khaki Roll, Sizes 1 to 26 mm 1 mm 30-530-85 8.5 mm 1.5 mm 30-530-9 9 mm 2 mm 30-530-95 9.5 mm 2.5 mm 30-530-10 10 mm 3 mm 30-530-105 10.5 mm 3.5 mm 30-530-11 11 mm 4 mm 30-530-115 11.5 mm 4.5 mm 30-530-12 12 mm 5 mm 30-530-125 12.5 mm 5.5 mm 30-530-13 13 mm 6 mm 30-530-14 14 mm 7 mm 30-530-145 14.5 mm 7.5 mm 30-530-150 15 mm	Set of 26 Dilators in Khaki Roll, Sizes 1 to 26 mm 1 mm 30-530-85 8.5 mm 30-530-16 1.5 mm 30-530-9 9 mm 30-530-165 2 mm 30-530-95 9.5 mm 30-530-17 2.5 mm 30-530-10 10 mm 30-530-175 3 mm 30-530-105 10.5 mm 30-530-18 3.5 mm 30-530-11 11 mm 30-530-185 4 mm 30-530-115 11.5 mm 30-530-19 4.5 mm 30-530-12 12 mm 30-530-20 5 mm 30-530-125 12.5 mm 30-530-21 5.5 mm 30-530-13 13 mm 30-530-22 6 mm 30-530-13 13.5 mm 30-530-23 6.5 mm 30-530-14 14 mm 30-530-24 7 mm 30-530-150 15 mm 30-530-26







30-555



HEGAR UTERINE DILATORS 71/2" (19.1 cm), DOUBLE END

Set of 8 Dilators in Khaki Roll, sizes 3-5 mm to 17-18 mm 30-550 30-555-1112 11-12 mm 3-4 mm 30-555-34 30-555-1314 13-14 mm 30-555-56 5-6 mm

15-16 mm 30-555-1516 7-8 mm 30-555-78 17-18 mm 9-10 mm 30-555-1718 30-555-910

PRATT UTERINE DILATORS 111/2" (29.2 cm), DOUBLE END

Set of 8 Dilators in Khaki Roll, sizes 13-15 Fr. (4.3-5 mm) 30-561 to 41-43 Fr. (13.5-14.2 mm)

13-15 Fr. (4.3-4.6 mm) 30-560-2931 29-31 Fr. (9.6-10.2 mm) 30-560-1315 30-560-3335 33-35 Fr. (10.9-11.6 mm) 17-19 Fr. (5.6-6.3 mm) 30-560-1719 30-560-3739 37-39 Fr. (12,2-12.9 mm) 30-560-2123 21-23 Fr. (7-7.6 mm) 30-560-4143 41-43 Fr. (13.5-14.2 mm) 30-560-2527 25-27 Fr. (8.3-9 mm)














30-620 GOODELL Uterine Dilator 13" (33 cm), Large Model, with gauge and corrugated blades GOODELL Uterine Dilator 11" (27.9 cm), Small Model, with gauge and corrugated blades 30-625 30-650 SIMS Uterine Sound 13" (33 cm), graduated in cm, malleable, silver plated 30-655 SIMPSON Uterine Sound 13" (33 cm), graduated in inches, malleable, silver plated 7-620 BOZEMAN Uterine Dressing Forceps 101/2" (26.7 cm), straight 7-622 BOZEMAN Uterine Dressing Forceps 101/2" (26.7 cm), double curved, standard model 7-624 BOZEMAN Uterine Dressing Forceps 101/2" (26.7 cm), double curved, with one large finger ring


CHERON Uterine Dressing Forceps 10½" (26.7 cm), rings angled sideways, serrated jaws 8 x 25 mm



30-835







MAIER UTERINE DRESSING FORCEPS 10" (25.4 cm)

Serrated Jaws 5 x 30 mm	

Without Ratchet With Ratchet 30-850 30-840 Straight 30-845 Curved 30-855 FOERSTER Sponge Forceps 91/2" (24.1 cm), straight, serrated jaws

7-600 7-602 FOERSTER Sponge Forceps 91/2" (24.1 cm), straight, smooth jaws 7-610 FOERSTER Sponge Forceps 91/2" (24.1 cm), curved, serrated jaws 7-612 FOERSTER Sponge Forceps 91/2" (24.1 cm), curved, smooth jaws

7-630 FLETCHER Cervix Holding Forceps 91/2" (24.1 cm), straight, serrated jaws 7-632 FLETCHER Cervix Holding Forceps 91/2" (24.1 cm), curved, serrated jaws (This instrument can also be used as a polyp forceps.)

V930-996 VANTAGE Polyp Forceps 101/2" (26.7 cm), angled shanks, delicate narrow jaws




■ Gynecological and Obstetrical Instruments ■ 30-940 30-950 30-951 to 30-955 30-965 30-970 30-940 MILTEX I.U.D. Extractor Hook 101/4" (26 cm) 30-950 MILTEX Uterine Tenaculum Hook 8" (20.3 cm), extra delicate pointed hook, handle graduated in cm 29-102 EMMETT UTERINE TENACULUM HOOKS 9" (22.9 cm), Sharp Points 30-951 Style 1 Half-Curved, small 30-954 Style 4 Half-Curved, large 30-952 Style 2 Acute Angle 30-955 Style 5 Double Angle 30-953 Style 3 Right Angle 30-965 SCHROEDER-BRAUN Uterine Tenaculum Forceps 994" (24.8 cm), straight 30-970 KAHN Uterine Tenaculum Forceps 91/2" (24.1 cm), angled jaws and bent shanks 29-102 MATHIEU I.U.D. Removal Forceps 111/2" (29.2 cm) shaft, serrated jaws 3.5 mm x 17 mm 30-975


DUPLAY Uterine Tenaculum Forceps 11" (27.9 cm), double curved

30-975





















16-38 THOMS-ALLIS Vulsellum Forceps 8" (20.3 cm), 6 x 7 teeth TEALE Vulsellum Forceps 10" (25.4 cm), curved sideways, 3 x 4 teeth, 5 mm wide 30-997 KELLY Vulselium Forceps 9%" (23.5 cm), curved sideways, 4 x 5 teeth, 7.7 mm wide 30-998 30-1000 SCHROEDER Uterine Vulsellum Forceps 10" (25.4 cm), straight SCHROEDER Uterine Vulsellum Forceps 30-1005 10" (25.4 cm), curved on flat 30-1010 SCHROEDER Uterine Vulsellum Forceps 10" (25.4 cm), curved sideways









SKENE Uterine Vulsellum Forceps 91/2" (24.1 cm), curved sideways 30-1020

30-1025 JACOBS Uterine Vulsellum Forceps 81/2" (21.6 cm), straight, serrated jaws 2 x 2 teeth

JACOBS Uterine Vulsellum Forceps 81/2" (21.6 cm), curved sideways, 30-1030 serrated jaws 2 x 2 teeth

> MUSEUX UTERINE VULSELLUM FORCEPS 91/2" (24.1 cm), 2 x 2 teeth


Straight Jaw Width 30-1050-6 6 mm 30-1050-8 8 mm 30-1050-10 10 mm

12-80 M.G.H. (Lahey) Forceps 61/4" (15.9 cm), 3 x 3 teeth, straight

30-1077 GORDON Uterine Vulsellum Forceps 9" (22.9 cm), straight, heavy pattern, 3 x 3 teeth








6-128	HEANEY Hissue Forceps 51/2" (14 cm), broad handles, 2 x 3 teeth and cross serrated tips
6-134 6-136	MARTIN Tissue Forceps $6''$ (15.2 cm), 7×8 teeth MARTIN Tissue Forceps $7\frac{1}{2}''$ (19.1 cm), 7×8 teeth
6-204	KELLY Tissue Forceps 9" (22.9 cm), 3 x 4 teeth, slotted serrated handles
6-142 6-144 6-145	RUSSIAN Tissue Forceps 6" (15.2 cm) RUSSIAN Tissue Forceps 8" (20.3 cm) RUSSIAN Tissue Forceps 10" (25.4 cm)
6-210	MAYO-RUSSIAN Tissue Forceps 9" (22.9 cm), fenestrated jaws, slotted grooved handles








30-1150 BILLROTH Tumor Forceps 11" (27.9 cm), 4 x 4 prongs 30-1160 DOYEN Tumor Screw 6¾" (17.1 cm), screw diameter 1" (2.5 cm)

30-1185 COLLIN (BUXTON) Uterine Forceps 101/2" (26.7 cm), swivel jaws

30-1190 SOMER Uterine Elevating Forceps 9" (22.9 cm), curved jaws





SIMS UTERINE CURETTES 11" (27.9 cm), Sharp Blades on Malleable Shanks

30-1205-00	Size 00	30-1205-3	Size 3
30-1205-0	Size 0	30-1205-4	Size 4
30-1205-1	Size 1	30-1205-5	Size 5
30-1205-2	Size 2	30-1205-6	Size 6

THOMAS UTERINE CURETTES 11" (27.9 cm), Blunt Blades on Malleable Shanks

30-1225-1	Size 1	30-1225-4	Size 4
30-1225-2	Size 2	30-1225-5	Size 5
30-1225-3	Size 3	30-1225-6	Size 6


30-1255-2 GREENE Uterine Curette 11" (27.9 cm), Serrated Blunt Blade on Malleable Shank, Size 2

30-1255-4 GREENE Uterine Curette 11" (27.9 cm), Serrated Blunt Blade on Malleable Shank, Size 4

30-1255-6 GREENE Uterine Curette 11" (27.9 cm), Serrated Blunt Blade on Malleable Shank, Size 6







BUMM UTERINE CURETTES 12" (30.5 cm), Extra Large Sizes on Malleable Shanks

	Sharp	3126	DIUIT	
	30-1270-1	1	30-1275-1	
	30-1270-2	2		
	30-1270-3	3		
	30-1270-4	4	30-1275-4	
30-1350		rical Speculum 9½" (24 strated jaws 7 x 25 mm	1 cm), with gauge and set s	screw, offset
30-1351	MILTEX-KOGAN Endocervical Speculum 11" (27.9 cm), improved pattern with ratchet and gauge. Extra delicate fenestrated jaws 3 mm wide at tip x 25 mm long			
30-1355		etractor 8¼" (21 cm) wit m) for an extra wide vie	h extra long thumb ratchet r wing area.	etracts up





30-1	380	SIMS Uterine Biopsy Curette 10" (25.4 cm), extra small sharp loop, outside meas. 2.5 x 5 mm
30-1 30-1		KEVORKIAN-YOUNGE Endocervical Biopsy Curette 12" (30.5 cm), loop size 3 x 12 mm, with basket KEVORKIAN-YOUNGE Endocervical Biopsy Curette 12" (30.5 cm), loop size 3 x 12 mm, without basket
30-1	385	HEANEY Uterine Biopsy Curette 9" (22.9 cm), serrated loop
30-1	390	NOVAK Endometrial Suction Biopsy Curette 9¾" (24.8 cm), with Luer hub, 4 mm dia., standard 12 Fr. size
30-1	391	NOVAK Endometrial Suction Biopsy Curette 9¾" (24.8 cm), with Luer hub, 2 mm dia., extra delicate
30-1	392	NOVAK Endometrial Suction Biopsy Curette 9¾" (24.8 cm), with Luer hub, 3 mm dia., delicate '
30-1	395	RANDALL Endometrial Suction Biopsy Curette 9¾" (24.8 cm), with Luer hub, 4 mm dia.



30-1430


30-1435

30-1440

angled shanks

S-shape, triangular bite

Gynecological and Obstetrical Instruments 📖



30-1400	GAYLOR Uterine Biopsy Forceps 8½" (21.6 cm), straight shanks, curved jaws, 5 mm bite, interlitted cups
30-1405	THOMS-GAYLOR Uterine Biopsy Forceps 8½" (21.6 cm), angled shanks, curved jaws, 5 mm bite, interfitted cups
30-1410	THOMS-GAYLOR Uterine Biopsy Forceps 91/2" (24.1 cm), angled shanks, curved jaws, 5 mm bite, interfitted cups
30-1420	SCHUBERT Uterine Biopsy Forceps 8½" (21.6 cm), straight
30-1425	SCHUBERT Uterine Biopsy Forceps 8½" (21.6 cm),


SCHUBERT Uterine Biopsy Forceps 11" (27.9 cm),

SCHUBERT Uterine Biopsy Forceps 11" (27.9 cm),

GELLHORN Uterine Biopsy Forceps 91/2" (24.1 cm), curved, oval bite, 5 x 8 mm

30-1435 30-1440







30-1442	TISCHLEH Cervical Biopsy Punch Forceps 9%" (24.8 cm) shart, 6 x 3 x 1.5 mm bite
30-1443	BURKE (Modified Baby TISCHLER) Cervical Biopsy Punch Forceps $7\%4^\circ$ (19.7 cm) shaft, $4.5\times3\times1.5$ mm bite
30-1445	TOWNSEND 'Mini-bite' Cervical Biopsy Punch Forceps 7% " (19.1 cm) shaft, $4.2\times2.3\times1$ mm bite with locking mechanism
30-1480	MILTEX Cup Jaw Uterine Biopsy Forceps 9" (22.9 cm) shaft, plain not interfitted cups. 6 mm bite









30-1482

KEVORKIAN-YOUNGE Uterine Biopsy Forceps 9½" (24.1 cm) shaft, 3.5 x 8 mm bite, teeth on lower jaw


EPPENDORFER Uterine Biopsy Forceps, 81/2" (21.6 cm) shaft 30-1485

WITTNER Uterine Biopsy Forceps 8% (21.6 cm) shaft, with teeth on lower jaw, standard model with angled tip 30-1486

WITTNER Uterine Biopsy Forceps 81/2" (21.6 cm) shaft, with teeth 30-1487

on lower jaw, straight tip






MILTEX ROTATING SHAFT UTERINE BIOPSY FORCEPS 9" (22.9 cm) shafts are interchangeable on Pistol Grip Spring Handles

	- () 3
30-1490 30-1491	KEVORKIAN-YOUNGE Biopsy Forceps, 3.5 x 8 mm bite, teeth on lower jaw, complete with Pistol Grip Handle KEVORKIAN-YOUNGE Biopsy Forceps, Shaft only
30-1492 30-1493	EPPENDORFER-KRAUSE Biopsy Forceps, oval 4.5 x 6.5 mm bite, complete with Pistol Grip Handle EPPENDORFER-KRAUSE Biopsy Forceps, Shaft only
30-1494 30-1495	TISCHLER Biopsy Forceps, 3 x 7 mm bite, complete with Pistol Grip Handle TISCHLER Biopsy Forceps, Shaft only
28-400	Pistol Grip Spring Handle fitting all Miltex Rotating Shaft Uterine Biopsy Forceps (this handle also fits Rectal Biopsy Punches Nos. 28-380 to 28-416)








30-1710	HEANEY Hysterectomy Forceps 81/4" (21 cm), heavy pattern, single tooth, curved HEANEY Hysterectomy Forceps 81/4" (21 cm), heavy pattern, double tooth, curved HEANEY Hysterectomy Forceps 93/4" (24.8 cm), extra heavy pattern, double tooth, curved
30-1716	HEANEY-BALLENTINE Hysterectomy Forceps 8½" (21 cm), longitudinal serrations, single tooth, curved
30-1718	BALLENTINE Hysterectomy Forceps 8½" (21.6 cm), longitudinal serrations, single tooth, straight
30-1720	GARLAND Hysterectomy Forceps $8\frac{1}{2}$ " (21 cm), longitudinal serrations, single tooth, angled on flat







LONG Hysterectomy Forceps 71/2" (19.1 cm), longitudinal serrations, 30-1730 1 x 2 teeth, curved


PHANEUF Uterine Artery Forceps 8" (20.3 cm), horizontal serrated jaws, 30-1735 1 x 2 teeth, straight

PHANEUF Uterine Artery Forceps 8" (20.3 cm), horizontal serrated jaws, 30-1740 1 x 2 teeth, angled on flat

> PEAN HYSTERECTOMY FORCEPS, with Horizontal Serrations Straight Size Curved 7-124 7-142 (20.3 cm) 7-144 7-126 (22.9 cm) 7-128 101/4" (26 cm) 7-146 7-130 12" (30.5 cm) 7-148

PEAN Hysterectomy Forceps 9" (22.9 cm), curved, with 30-1805 longitudinal serrations


30-1815 PEAN Hysterectomy Forceps 101/4" (26 cm), curved, with longitudinal serrations



30-1815








7-190	MICULICZ Peritoneal Clamp 8" (20.3 cm), curved, 1 x 2 teeth
29-280	WERTHEIM Vaginal Forceps 9%" (24.8 cm), curved, 2" (5.1 cm) jaws
29-282	WERTHEIM-CULLEN Pedicle Forceps 8½" (21.6 cm), right angle 2" (5.1 cm) jaws
5-160 5-162	DOYEN Scissors 7" (17.8 cm), straight, heavy pattern DOYEN Scissors 7" (17.8 cm), curved, heavy pattern
5-166	FERGUSON Scissors 7" (17.8 cm), angled on flat








SIMS UTERINE SCISSORS 8" (20.3 cm)

	Straight 5-220	Points Sharp-Sharp	Curved 5-226	1	
	5-222	Sharp-Blunt	5-228		Ø
	5-224	Blunt-Blunt	5-230	SI.	
5-250		e Scissors 6¼" (15.9 cm), s			
5D-251	KELLY Uterine Scissors 61/4" (15.9 cm), sharp points, straight, with one serrated blade				
5-252	KELLY Uterin	e Scissors 61/4" (15.9 cm), s	harp points, curved	// \}	6
5D-253	KELLY Utering with one serra	e Scissors 61/4" (15.9 cm), s ited blade	harp points, curved,		
5-128	MAYO Uterine	e Scissors 9" (22.9 cm), stra	ight		
5-130	MAYO Uterine	e Scissors 9" (22.9 cm), cun	/ed		
5-186	NELSON Uter	ine Scissors 9" (22.9 cm), s	traight		
5-188	NELSON Uter	rine Scissors 9" (22.9 cm), c	urved		
30-2190	BRAUN Epision guarded lower	otomy Scissors 5½" (14 cm) r blade	, angled to side,	00	
30-2195	BRAUN Epision guarded lower	otomy Scissors 81⁄24" (21.6 o r blade	m), angled to side,	5-186 5-188	30-2190 30-2195

For Carb-N-Sert Scissors with Tungsten Carbide inserts and Super Cut Scissors – see pages 55-63.







30-2200	PINARD Fetal Heart Stethoscope, aluminum
30-2205	DE LEE-HILLIS Obstetrical Stethoscope with padded head support, rigid metal arm. Complete with binaurals, chrome

DE LEE Obstetrical Stethoscope, same as No. 30-2205 30-2210 but with flexible arm, chrome

MARTIN Pelvimeter, graduated in centimeters and inches 30-2230 30-2235 COLLYER Pelvimeter, graduated in centimeters and inches



Gynecological and Obstetrical Instruments











30-2250 IOWA Membrane Puncturing Forceps 10¼" (26 cm), double curved, 6 x 6 teeth
 30-2260 BEACHAM Anniotome 10¼" (26 cm)
 30-2300 SIMPSON Obstetrical Forceps, Long Model, 14" (35.6 cm)
 30-2305 SIMPSON Obstetrical Forceps, Short Model, 12" (30.5 cm)



Gynecological and Obstetrical Instruments















DE LEE Obstetrical Forceps, Short Model, 30-2320 12" (30.5 cm)

ELLIOTT Obstetrical Forceps, Short Model, 30-2330 121/2" (31.8 cm)

PIPER Obstetrical Forceps 171/2" (44.5 cm) 30-2335

30-2345 BILL Traction Handle, for use with any O.B. Forceps

having standard type handle



■ Gynecological and Obstetrical Instruments





30-2400







30-2350 KIELLAND Obstetrical Forceps 15½" (39.4 cm)
 30-2355 KIELLAND-LUIKART Obstetrical Forceps 15½" (39.4 cm), solid blades
 30-2390 LUIKART Obstetrical Forceps 15½" (39.4 cm), solid blades

SIMPSON-LUIKART Obstetrical Forceps 14" (35.6 cm), solid blades




■ Gynecological and Obstetrical Instruments ■









McLEAN-LUIKART Obstetrical Forceps 15¾" (40 cm), solid blades, long shanks 30-2410 30-2415 McLEAN-TUCKER-LUIKART Obstetrical Forceps 15¾" (40 cm), solid blades,

long shanks

DEWEY Obstetrical Forceps 14" (35.6 cm), with Axis Traction 30-2430

■ Gynecological and Obstetrical Instruments







30-2445





30-2435 HAWKS-DENNEN Obstetrical Forceps 15" (38.1 cm)
 30-2445 BAILEY-WILLIAMSON Obstetrical Forceps 15½" (39.4 cm)
 30-2450 MURLESS Head Extractor 10" (25.4 cm), with hinged shaft and sliding collar lock


30-2505 KELLY Placenta Forceps 121/2" (31.8 cm), serrated ring jaws, slightly curved







Gynecological and Obstetrical Instruments













30-2530 WINTER Placenta Forceps 11" (27.9 cm), cup jaws, straight 30-2535 WINTER Placenta Forceps 11" (27.9 cm), cup jaws, curved 30-2560 KANE Umbilical Clamp 31/2" (8.9 cm) 30-2575 Umbilical Scissors 41/8" (10.5 cm), American Pattern 30-2580 BUSCH Umbilical Scissors 61/2" (16.5 cm), curved on side

30-2605 DUBOIS Decapitation Scissors 101/2" (26.7 cm), curved, extra heavy pattern, blunt points

30-2632 IOWA Trumpet Pudendal Needle Guide 51/2" (14 cm)







	FOX DERMAL CURETTES 51/2" (14 cm), Round		
33-2	6 mm dia.	33-5	3 mm dia.
33-3	5 mm dia.	33-5A	2 mm dia.
33-4	4 mm dia.	33-5B	1 mm dia.
	PIFFARD DERMAL CURETTES 61	/4" (15.9 cm),	Oval
33-6	Size 1 (3 mm dia.)	33-11	Size 4 (8 mm dia)
33-8	Size 2 (5 mm dia.)	33-11A	Size 5 (10 mm dia.)
33-10	Size 3 (6 mm dia.)		, ,
	PIFFARD DERMAL CURETTES 51/2" (14 cm), oval, narrow handles		
33-16-00	Size 00 (1 mm dia.)	33-16-2	Size 2 (5 mm dia)
33-16-0	Size 0 (1.5 mm dia.)	33-16-3	Size 3 (6 mm dia.)
33-16-1	Size 1 (3 mm dia.)		
33-12	MILTEX Dermal Curette 61/4" (15.9	cm), oval spo	on size 0 (3 mm dia.)
33-13	MILTEX Dermal Curette 61/4" (15.9 cm), oval spoon size 1 (4 mm dia.)		
33-14	MILTEX Dermal Curette 61/4" (15.9 cm), oval spoon size 2 (5 mm dia.)		
33-15	MILTEX Dermal Curette 61/4" (15.9	cm), oval spo	on size 3 (6 mm dia.)
40-54	MILTEX Verucca Curette 5" (12.7 o	cm), cup 4 mm	ı dia.




to 40-57/3

to 40-58/3

18-520

Dermal Instruments









MILTEX EXCAVATOR CURETTES, Single End 5" (12.7 cm)

W/O Hole	Diameter	W/Hole
40-57/0	1 mm	40-58/0
40-57/1	1.5 mm	40-58/1
40-57/2	2 mm	40-58/2
40-57/3	2.5 mm	40-58/3

MILTEX EXCAVATOR CURETTES, Double End 51/2" (14 cm)

40-58/1-2 Small 1.5 mm and Medium 2 mm dia. Small 1.5 mm and Large 2.5 mm dia. 40-58/1-3

40-58/1-3


BUCK CURETTES 61/2" (16.5 cm), Sharp

Straight	Size	Angled
19-260	00	19-270
19-262	0	19-272
19-264	1	19-274
19-266	2	19-276
19-268	3	19-278

18-522 18-524	SKEELE Curette 5" (12.7 cm), round 2 mm cup, serrated edge SKEELE Curette 5" (12.7 cm), round 2.5 mm cup, serrated edge	
18-530	HEATH Curette 4" (10.2 cm), Size 0, one sharp and one blunt edge, 0	0.5 mm
18-532	HEATH Curette 4" (10.2 cm), Size 1, one sharp and one blunt edge,	1 mm
18-534	HEATH Curette 4" (10.2 cm), Size 2, one sharp and one blunt edge,	2 mm
18-536	HEATH Curette 4" (10.2 cm), Size 3, one sharp and one blunt edge,	3 mm

SKEELE Curette 5" (12.7 cm), round 1 mm cup, serrated edge





	KEYES DERMAL F	PUNCHES 4" (10.2	cm)
33-20	2 mm	33-25	5 mm
33-22	3 mm	33-26	6 mm
33-24	4 mm	33-28	8 mm
33-30	KEYES Dermal Punch Set, 3½" (8.9 cm) Handle, with screw-in Keyes Knives 2, 3, 4, 5, 6 and 7 mm mounted on storage and sterilization rack		
	Replacement Parts	for No. 33-30 Pun	ich Set
33-30-1	Handle		
33-30-2	2 mm Head	33-30-5	5 mm Head
33-30-3	3 mm Head	33-30-6	6 mm Head
33-30-4	4 mm Head	33-30-7	7 mm Head

MILTEX Sterile Disposable Biopsy Punches – made by Kai Each Sterile Biopsy Punch is sealed in a peel-apart see-through package. Packed 50 Biopsy Punches to a Dispenser Box. Sterilized by EO.

33-31A	1.5 mm	33-34	4 mm
33-31	2 mm	33-35	5 mm
33-32	3 mm	33-36	6 mm
33-33	3.5 mm	33-37	8 mm



FEATURES:

- EALUMES.
 Seamless stainless steet razor sharp cutting edge
 Ribbed handle provides maximum control and comfort
 Size marked on every punch for easy identification
 Sterile ready to use every time

- Convenient use once and discard Eliminates need for re-sharpening Clean, precise specimens for accurate diagnosis Wide range of sizes









5-300








VAUGHAN Abscess Knife Set 5" (12.7 cm), one straight 20 mm lancet and one curved 12 mm knife, with protective 33-40 screw-on caps, and pocket clip

Iris Scissors 4" (10.2 cm), straight 5-300 5-302 Iris Scissors 4" (10.2 cm), curved Iris Scissors 41/2" (11.4 cm), curved sideways 5D-303 5-304 Iris Scissors 41/2" (11.4 cm), straight Iris Scissors 41/2" (11.4 cm), curved 5-306 5-308 Iris Scissors 41/2" (11.4 cm), curved sideways

5-308 5-306

For Carb-N-Sert Scissors with tungsten carbide inserts and SuperCut Scissors see pages 55 - 63






18-1392 18-1394	BONN Miniature Iris Scissors 3½" (8.9 cm), with 15 mm blades, straight BONN Miniature Iris Scissors 3½" (8.9 cm), with 15 mm blades, curved		
18-1396 18-1398	Iris Scissors 3½" (8.9 cm), with 20 mm blades, straight, delicate Iris Scissors 3½" (8.9 cm), with 20 mm blades, curved, delicate		
18-1404 18-1405	MILTEX Iris Scissors 4" (10.2 cm), straight, heavy pattern MILTEX Iris Scissors 4" (10.2 cm), curved, heavy pattern		
18-1415* 18-1416*	Iris Scissors 4" (10.2 cm), straight, standard pattern, ribbon type Iris Scissors 4" (10.2 cm), curved, standard pattern, ribbon type	18-1415*	18-1417*
18-1417* 18-1418*	Iris Scissors 3%" (9.5 cm), straight, miniature pattern with 14 mm blades, ribbon type Iris Scissors 3%" (9.5 cm), curved, miniature pattern with 14 mm blades, ribbon type	18-1416*	18-1418*
	biddes, sepon type		


*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished







18-1419*	Utility Scissors 3¾" (9.5 cm), curved blunt blades, 15 mm long, ribbon type
18-1420 18-1422 18-1426 18-1428	KNAPP Iris Scissors 4" (10.2 cm), straight, sharp-blunt points KNAPP Iris Scissors 4" (10.2 cm), straight, blunt-blunt points KNAPP Iris Scissors 4" (10.2 cm), curved, sharp-blunt points KNAPP Iris Scissors 4" (10.2 cm), curved, blunt-blunt points
18-1440	WILMER Utility Scissors 4" (10.2 cm), angled on flat
5-312 5-314 18-1453* 18-1454*	Strabismus Scissors 4" (10.2 cm), straight Strabismus Scissors 4" (10.2 cm), curved Strabismus Scissors 4¼"(10.8 cm), straight, ribbon type Strabismus Scissors 4¼"(10.8 cm), curved, ribbon type




*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished











18-1478*

18-1510 to 18-1514

STEVENS TENOTOMY SCISSORS

18-1460 18-1462 18-1464 18-1466	Straight sharp points 41/6" (10.5 cm), short blades Straight blunt points 41/6" (10.5 cm), short blades Curved sharp points 41/6" (10.5 cm), short blades Curved blunt points 41/6" (10.5 cm), short blades	M	M	M
18-1470 18-1472 18-1474 18-1476	Straight sharp points 4½" (11.4 cm), long blades Straight blunt points 4½" (11.4 cm), long blades Curved sharp points 4½" (11.4 cm), long blades Curved blunt points 4½" (11.4 cm), long blades			
18-1477* 18-1478*	STEVENS Tenotomy Scissors 3¾" (9.5 cm), straight, blunt points, ribbon type STEVENS Tenotomy Scissors 3¾" (9.5 cm), curved, blunt points, ribbon type			
21-539	KLEINERT-KUTZ Tenotomy Scissors 4% " (12.4 cm), curved, rounded dissecting blades			
18-1510 18-1512 18-1514 18-1516	NOYES Iris Scissors 4½" (11.4 cm), straight, sharp-sharp points NOYES Iris Scissors 4½" (11.4 cm), straight, sharp-blunt points NOYES Iris Scissors 4½" (11.4 cm), straight, blunt-blunt points NOYES Iris Scissors 4½" (11.4 cm), curved, sharp-sharp points	18-1516	18-1480	18-1486

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished


WESTCOTT Eye Scissors 41/2" (11.4 cm), curved, blunt points

WESTCOTT Eye Scissors 41/2" (11.4 cm), curved, sharp points

18-1480

18-1486






33-240	BIRO Dermal Naevus Scissors 4" (10.2 cm), angled on flat, w/sharp outside edge
33-242	MILTEX Undermining Scissors 41/2" (11.4 cm), curved blades, with flat pointed tip
18-1652 18-1653*	GRADLE Scissors 3%" (9.5 cm), slightly curved, sharp points GRADLE Scissors 4" (10.2 cm), ribbon type, slightly curved, sharp points
21-537TC	MILTEX Blepharoplasty Scissors 4% (11.4 cm), curved flat tips, one serrated and one Carb-N-Sert blade
21-630 21-632	JOSEPH Scissors 5%" (14.6 cm), straight, sharp points JOSEPH Scissors 5%" (14.6 cm), curved, sharp points

21-630 21-632

*Ribbon Type Scissors have extra large finger rings; shanks and rings have satin finish; blades are highly polished






5-180 5-284	METZENBAUM Scissors 5½" (14 cm), curved, standard METZENBAUM Scissors 5" (12.7 cm), curved, delicate	The second secon
5-290	RAGNELL Undermining Scissors 5" (12.7 cm), curved with flattened dissecting tips	The same of the sa
5-320	GOLDMAN-FOX Wound Debridement Scissors 5" (12.7 cm), curved fine tips with one serrated blade	
5D-322	GOLDMAN-FOX Wound Debridement Scissors 5" (12.7 cm), curved, one serrated blade, angled shanks	$\cup \cup$
5D-324	LA GRANGE Scissors 41/4" (10.8 cm), curved, one serrated blade, angled shanks	5D-324

For Carb-N-Sert Scissors with tungsten carbide inserts and SuperCut Scissors see pages 55 - 63








9-100 9-101 9-102 9-103	SPENCER Stitch Scissors 3½" (8.9 cm), delicate SHORTBENT Stitch Scissors 3½" (8.9 cm), curved, delicate LITTAUER JR. Stitch Scissors 4½" (11.4 cm), light pattern NORTHBENT Stitch Scissors 4¾" (12.1 cm), curved, light pattern		
9-104	LITTAUER Stitch Scissors 5½" (14 cm), standard pattern		
9-108	MILTEX Stitch Scissors 4½" (11.4 cm), angled, delicate		
9-110	O'BRIEN Stitch Scissors 3¾" (9.5 cm), angled, sharp points	9-108	9-110






8-4	DERF Needle Holder	43/4"	(19.1	cm\
Ų~ ~	DELLI LAGGRIG LIGITIGE	474	12.1	UIII)

8-6 WEBSTER Needle Holder 5" (12.7 cm), delicate smooth jaws,

standard pattern

8-7 WEBSTER Needle Holder 4¾" (12.1 cm), extra delicate smooth jaws

8-8 HALSEY Needle Holder 5" (12.7 cm), smooth jaws 8D-36 HALSEY Needle Holder 5" (12.7 cm), serrated jaws



For additional Needle Holders see pages 112 - 121

For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130













8-14 OLSEN-HEGAR Needle Holder-Scissor Combination 43/4" (12.1 cm), serrated jaws, extra delicate 8-14A

OLSEN-HEGAR Needle Holder-Scissor Combination 4%" (12.1 cm), smooth jaws, extra delicate

8-26 BOYNTON Needle Holder 5" (12.7 cm), delicate

8-28 MATHIEU Needle Holder 51/2" (14 cm), with spring handles

8-40 BAUMGARTNER Needle Holder 51/2" (14 cm)


18-1800 KALT Needle Holder 51/4" (13.3 cm), standard model, delicate jaws

GREEN Needle Holder and Suturing Forceps 41/4" (10.8 cm), 18-1810 angled jaws, double end

For additional Needle Holders see pages 112 - 121

For Carb-N-Sert Needle Holders with tungsten

carbide inserts see pages 123 - 130
















18-1828 18-1831 18-1834 18-1832 18-1833

18-1820 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, without lock 18-1824 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, without lock 18-1828 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, with lock 18-1832 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, with lock 18-1831 CASTROVIEJO Needle Holder 51/2" (14 cm), straight, extra delicate jaws, with lock 18-1833 CASTROVIEJO Needle Holder 51/2" (14 cm), curved, extra delicate jaws, with lock 18-1834 PATON Needle Holder 41/2" (11.4 cm), straight, without lock BARRAQUER Needle Holder 5" (12.7 cm), curved, without lock, 18-1838 10 mm wide hollow round handles 18-1839 BARRAQUER Needle Holder 5" (12.7 cm), curved, with lock, 10 mm wide hollow round handles 18-1840 BARRAQUER Needle Holder 51/8" (13 cm), curved, delicate tapered jaws, 8 mm wide hollow round handles

For additional Needle Holders see pages 112 - 121

For Carb-N-Sert Needle Holders with tungsten carbide inserts see pages 123 - 130



18-1840













6-118







6-104 SEMKEN Dressing Forceps 5" (12.7 cm), serrated

6-106 SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth

6-107 SEMKEN Tissue Forceps 5" (12.7 cm), 1 x 2 teeth, curved tips

SEMKEN Tissue Forceps 5" (12.7 cm), 2 x 3 teeth 6-108

6-118 ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate 6-118A ADSON Dressing Forceps 4¾" (12.1 cm), serrated, delicate, angled ADSON Dressing Forceps 43/4" (12.1 cm), serrated, lightweight, 6-118XL

fenestrated handles

17-2510 Micro ADSON Dressing Forceps 43/4" (12.1 cm), serrated tips, 0.5 mm wide

6-120 ADSON Tissue Forceps 43/4" (12.1 cm), 1 x 2 teeth, delicate ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, delicate, angled 6-120A

6-120XL ADSON Tissue Forceps 43/4" (12.1 cm), 1 x 2 teeth, lightweight. fenestrated handles


For additional Dressing and Tissue Forceps see pages 64 - 82











6-121	ADSON Tissue Forceps 4¾" (12.1 cm), 1 x 2 teeth, and cross
	serrated tips, delicate

17-2500	Micro ADSON Tissue Forceps 4¾" (12.1 cm), tips
	0.5 mm wide 1 x 2 teeth

6-122 ADSON Ti	ssue Forceps 4¾" (12.1 cm), 2 x 3 teeth, delicate	
6-123 ADSON Tis	ssue and Suture Forceps 43/4" (12.1 cm), 1 x 2 tee	t

ADSON Tissue Forceps 4¾" (12.1 cm), w/platform, 1 x 2 teeth, lightweight, fenestrated handles

6-124 BROWN-ADSON Tissue Forceps 4%4" (12.1 cm), 7 x 7 side grasping teeth


6-124XL BROWN-ADSON Tissue Forceps 4¾" (12.1 cm), 7 x 7 side grasping teeth, lightweight, fenestrated handles

6-125 DeBAKEY-ADSON Tissue Forceps 4¾" (12.1 cm), atraumatic tips 1.5 mm wide

24-570XL DeBAKEY Tissue Forceps 6' (15.2 cm), atraumatic tips, 1.5 mm wide, lightweight, fenestrated handles

6-126 BROWN-ADSON Tissue Forceps 4¾" (12.1 cm), 7 x 7 side grasping teeth, angular

For additional Dressing and Tissue Forceps see pages 64 - 82



24-570XL

6-125



6-126

6-123XL


















EYE DRESSING FORCEPS 4" (10.2 cm), SERRATED TIPS

Standard Pattern 0.8 mm wide tips			licate Pattern n wide tips
18-780	Straight	18-779	Straight
18-782	Half curved	18-781	Half curved
18-784	Full curved	18-783	Full curved

For additional Dressing and Tissue Forceps see pages 64 - 82






IRIS TISSUE FORCEPS 4" (10.2 cm), 1 x 2 TEETH

	1110 110001 1011011 0 4 (1	U.Z. OHIJ, 1 X Z. 1 L	.L. [11	B / I	
	ard Pattern		elicate Pattern	M	
0.8 111	m wide tips	0.5 111	m wide tips		700
18-786	Straight	18-785	Straight		112
18-787	Straight (Non Magnetic)	18-789	Half curved		
18-788	Half curved	18-791	Full curved		
18-790	Full curved				
18-850	GRAEFE Tissue Forceps 4%" with fine teeth	(11.2 cm), standa	ard jaws 4.5 mm wide,		on an
18-854	GRAEFE Tissue Forceps 4%" with fine teeth, with catch	(11.2 cm), standa	ard jaws 4.5 mm wide,		
				18-850	18-854

For additional Dressing & Tissue Forceps see pages 64 - 82








18-866	BISHOP-HARMON Dressing Forceps 3%" (8.6 cm), fine cross serrated tips 0.5 mm wide
18-868	BISHOP-HARMON Tissue Forceps 3%" (8.6 cm), standard 0.7 mm 1 x 2 teeth
18-1107	Swiss Cilia and Suture Forceps 3%" (9.5 cm), precision fitted diagonal jaws 3 mm wide
18-876	CASTROVIEJO (ALVIN) Fixation Forceps 3¾" (9.5 cm), straight, 1 x 1 curved teeth extending beyond each other
18-878	CASTROVIEJO Fixation Forceps $3\%''$ (9.5 cm), straight, 2 x 2 curved teeth extending beyond each other
18-950 18-951 18-952 18-953	CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.9 mm, with tying platform CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.3 mm, with tying platform CASTROVIEJO Suturing Forceps 4" (10.2 cm), 1 x 2 teeth 0.6 mm, with tying platform CASTROVIEJO Micro Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.12 mm, with tying platform
18-954	CASTROVIEJO Suturing Forceps 4* (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.6 mm, with tying platform
18-955	CASTROVIEJO Suturing Forceps 4" (10.2 cm), with 11 mm wide handles, 1 x 2 teeth 0.9 mm, with tying platform

For additional Fixation and Suturing Forceps see pages 238, 244 - 248 & 251






21-430	MILTEX Tissue Twist Hook 5½" (14 cm)	Î	7
26-1084 26-1086	FRAZIER Skin Hook 5" (12.7 cm), sharp point FRAZIER Skin Hook 5" (12.7 cm), blunt point		
33-121	SELTZER Dermal Biopsy Hook 41/2" (11.4 cm), sharp point	l	
18-412 18-414	TYRELL Skin Hook 5" (12.7 cm), sharp point TYRELL Skin Hook 5" (12.7 cm), blunt point		
18-416	GUTHRIE Hook 5" (12.7 cm), sharp double prongs 1.5 mm wide	milita	
21-85 21-86 21-87	KLEINERT-KUTZ Hook 5" (12.7 cm), small size 3 mm dia. KLEINERT-KUTZ Hook 5" (12.7 cm), medium size 5 mm dia. KLEINERT-KUTZ Hook 5" (12.7 cm), large size 7 mm dia.	The state of the s	
21-88 21-89	GILLIES (CONVERSE) Skin Hook 7" (17.8 cm), large size GILLIES (CONVERSE) Skin Hook 7" (17.8 cm), small size	21-88 21-89	

For additional Hooks see pages 188 - 194








21-90 21-92 21-94	COTTLE Skin Hook 5½" (14 cm), small, shallow curve COTTLE Skin Hook 5½" (14 cm), small, deep curve COTTLE Skin Hook 5½" (14 cm), large, deep curve
21-102 21-103	FREER Skin Hook 6" (15.2 cm), sharp double prongs 2.5 mm wide FREER Skin Hook 6" (15.2 cm), blunt double prongs 2.5 mm wide
21-144	COTTLE Double Hook 5% " (14 cm), two sharp prongs 12 mm wide
21-152	MILTEX Skin Hook 4¾" (12.1 cm), sharp prong 2 mm dia.
21-153 21-154 21-156 21-158 21-160	JOSEPH Single Hook 6¼" (15.9 cm), one sharp prong JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs JOSEPH Double Hook 6¼" (15.9 cm), two sharp prongs






18-100	ROLLET Retractor 51/4" (13.3 cm), with 4 sharp prongs
18-102	KNAPP Retractor 51/4" (13.3 cm), with 4 blunt prongs
18-439	McGANNON Retractor 5" (12.7 cm), fenestrated blade 2.5 mm wide
11D-27	Flexible Neck Rake Retractor 5" (12.7 cm), 3 sharp prongs
	RAKE RETRACTORS 6" (15.2 cm)

	HAKE HETHAUTURS 6" (15.2 cm)	
Sharp	With Flexible Neck	Blunt
11-30	1 prong	11-36
11-32	2 prongs	11-38
11-34	3 prongs	11-40
	RAKE RETRACTORS 6" (15.2 cm)	
Sharp	With Rigid Neck	Blunt
11-50	1 prong	11-56
11-52	2 prongs	11-58
11-54	3 prongs	11-60













11-12	ALM Retractor 2¾" (7 cm), spread 2¼" (5.7 cm), prongs %" (.32 cm) deep
11-13	ALM Retractor 3%" (9.8 cm), spread 3" (7.6 cm), prongs 1/2" (1.3 cm) deep
11-14 11-16	MILTEX Cross Action Retractor 4" (10.2 cm), sharp prongs MILTEX Cross Action Retractor 4" (10.2 cm), blunt prongs
11-74 11-76	SENN Retractor 6%" (16.2 cm), sharp, double end SENN Retractor 6%" (16.2 cm), blunt, double end
11-80 11-81	MATHIEU Retractor 6½" (16.5 cm), sharp, double end MATHIEU Retractor 6½" (16.5 cm), blunt, double end







18-570





40-94 40-96


DESMARRES Eye Lid Retractor 5½" (14 cm), size 1, 11 mm
DESMARRES Eye Lid Retractor 51/2" (14 cm), size 2, 13 mm
DESMARRES Eye Lid Retractor 51/2" (14 cm), size 3, 15 mm
DESMARRES Eye Lid Retractor 51/2" (14 cm), size 4, 17 mm
GREEN Spatula 6" (15.2 cm), slightly curved ends 4 x 14 mm and 5 x 16 mm
Finger Ring Cutter $61\!\!/2"$ (16.5 cm), with circular saw blade, chrome Replacement Blade for Finger Ring Cutter, chrome
LOCKE Elevator 4½" (11.4 cm), small blade 3 x 15 mm LOCKE Elevator 5" (12.1 cm), large blade 5 x 20 mm
Flexible Stainless Ruler 6" (15.2 cm) x ½" (1.3 cm), graduated in ½2" (.8 mm) fractions & mm

18-660

For Measuring Calipers see page 239







33-200 33-201 33-202	SCHAMBERG Comedone Extractor 3¾" (9.5 cm), with square loop ends SCHAMBERG Comedone Extractor 3¾" (9.5 cm), standard pattern with crimped small loop SCHAMBERG Comedone Extractor 3¾" (9.5 cm), standard pattern with extra fine loops, small loop crimped
33-203	Comedone Extractor 4" (10.2 cm), double end, flat and wire end expressors
33-204 33-206	SAALFIELD Comedone Extractor 5%" (14.6 cm), with acne lancet, professional pattern SAALFIELD Comedone Extractor 4%" (10.8 cm), with acne lancet, professional pattern
33-210	UNNA Comedone Extractor 5%" (14.6 cm), spoon ends












33-214	MILTEX Comedone Extractor 6" (15.2 cm), one fenestrated cup and one delicate curved lancet with protective cap
33-215	GREG Comedone Extractor 5" (12.7 cm), one fenestrated cup and chuck with interchangeable No. 65 lancet blade, protected by cap
33-216	Replacement Blade No. 65 for No. 33-215 Comedone Extractor
33-230 33-232 33-234	Finger Nail Drill 4" (10.2 cm), standard pattern, complete with 3 drill points stored in handle Finger Nail Drill 4½" (11.4 cm), with Chuck Head, complete with 3 drill points stored in handle Set of 3 Drill Points, small, medium and large size, for Finger Nail Drills, chrome
33-500	MILTEX Depilatory Forceps 51/4" (13.3 cm), flat square tips with deep serrations



40-200

Dermal Instruments


40-215



40-210 40-212


40-200	Nail Nipper 41/2" (11.4 cm), concave jaws, single spring, chrome
40-205 40-207	Nail Nipper $4\frac{1}{2}$ " (11.4 cm), concave jaws, barrel spring, chrome Nail Nipper $5\frac{1}{2}$ " (14 cm), concave jaws, barrel spring, chrome
40-210-CH 40-210-SS	Nail Nipper $5\frac{1}{2}$ " (14 cm), concave jaws, double spring, chrome Nail Nipper $5\frac{1}{2}$ " (14 cm), concave jaws, double spring, stainless
40-212-CH 40-212-SS	Nail Nipper $5\frac{1}{2}$ " (14 cm), straight jaws, double spring, chrome Nail Nipper $5\frac{1}{2}$ " (14 cm), straight jaws, double spring, stainless
40-215-CH 40-215-SS	Nail Nipper $51\!\!/\!\!z''$ (14 cm), angled concave jaws, double spring, chrome Nail Nipper $51\!\!/\!\!z'''$ (14 cm), angled concave jaws, double spring, stainless
40-217	Nail Nipper 5" (12.7 cm), concave jaws, double spring, stainless
40-219	Double Action Nail Nipper 6" (15.2 cm), sharp concave jaws, barrel spring, stainless

40-205 40**-**207



40-217






40-220	Tissue and Cuticle Nipper 4" (10.2 cm), straight jaws, ring handles, chrome	
40-224	Nail Splitter 4" (10.2 cm), straight jaws, double spring, stainless	
40-225A	Nail Splitter 41/2" (11.4 cm), delicate straight jaws, double spring, stainless	
40-225-CH 40-225-SS	Nail Splitter 4½" (11.4 cm), straight jaws, single spring, chrome Nail Splitter 4½" (11.4 cm), straight jaws, double spring, stainless	0
40-226A	Nail Splitter 5" (12.7 cm), straight jaws, extra narrow, double spring, stainless	
40-226-CH 40-226-SS	Nail Splitter 5" (12.7 cm), straight jaws, single spring, chrome Nail Splitter 5" (12.7 cm), straight jaws, double spring, stainless	
40-227-CH 40-227-SS	Nail Splitter 6" (15.2 cm), straight jaws, single spring, chrome Nail Splitter 6" (15.2 cm), straight jaws, double spring, stainless	\
40-228	Nail Splitter 6" (15.2 cm), straight jaws, extra narrow, double spring, stainless	
40-230	Nail Splitter 5" (12.7 cm), English Anvil Pattern, double spring, stainless	₩ 40-234
40-234 40-235	PLATYPUS Nail Pulling Forceps 5½" (14 cm), standard wide jaws, stainless PLATYPUS Nail Pulling Forceps 5½" (14 cm), narrow jaws, stainless	40-235 40-235



40-260

Dermal Instruments



40-245-CH 40-245-SS	Tissue and Cuticle Nipper 4" (10.2 cm), convex jaws, chrome Tissue and Cuticle Nipper 4" (10.2 cm), convex jaws, stainless
40-250-CH 40-250-SS	Tissue and Cuticle Nipper 41/2" (11.4 cm), convex jaws, chrome Tissue and Cuticle Nipper 41/2" (11.4 cm), convex jaws, stainless
40-251-CH 40-251-SS	Tissue and Cuticle Nipper 5" (12.7 cm), convex jaws, chrome Tissue and Cuticle Nipper 5" (12.7 cm), convex jaws, stainless
40-252	Tissue and Cuticle Nipper 5" (12.7 cm), convex jaws, heavy pattern, stainless
40-425	Nail Scissors 3½" (8.9 cm), curved blades, chrome
40-435	Cuticle Scissors 3½" (8.9 cm), curved blades, standard pattern, chrome
40-445	Cuticle Scissors 31/2" (8.9 cm), curved blades, extra delicate, chrome


Ingrown Nail Splitting Scissors 6" (15.2 cm), stainless



Post Mortem Instruments













34-52


34-50 34-52	VIRCHOW Brain Sectioning Knife, blade 180 mm x 20 mm, double edge VIRCHOW Brain Sectioning Knife, blade 250 mm x 40 mm, double edge
34-80	Dissecting Scissors 4% (11.4 cm), straight, two sharp points, economy grade
34-92	RICHTER Dissecting Scissors 5" (12.7 cm), angular blades, one sharp and one blunt point $$
34-98	Enterotomy Scissors 8" (20.3 cm), with hook blade
30-2605	DUBOIS Scissors 101/2" (26.7 cm), curved, blunt points, extra heavy pattern
25-360	GLUCK Rib Shears 81/2" (21.6 cm)



For a complete selection of Scissors see pages 18 - 52



Post Mortem Instruments



			ana a
34-256	Aneurysm Needle 6½" (16.5 cm), blunt tip with eye		
34-250	Dissecting Tenaculum 6½" (16.5 cm), sharp point		
34-220	Post Mortem Hammer 91/2" (24.1 cm), with hook handle		
34-210	VIRCHOW Skull Breaker 3" (7.6 cm), blade ½" (1.3 cm) wide		Λ
40-227-SS	LITTAUER-LISTON Bone Cutting Forceps 6" (15.2 cm), straight 20 mm blades	The Process of the Pr	
25-362	STILLE Pattern Rib Shears 81/2" (21.6 cm)	4	





Suggested Surgical Instrument Sets

The following lists have been compiled as a guide in determining basic instrumentation requirements.

	Page Nos.
General Surgery	. 555-561
Genito-Urinary Surgery	. 562-564
Rectal Surgery	. 565-566
Obstetrical and Gynecological Surgery	. 567-570
Neurosurgery	. 571-575
Orthopedic Surgery	. 576-582
Thoracic and Cardiac Surgery	. 583-588
Ophthalmic Surgery	. 589-591
Ear, Nose and Throat Surgery	. 592-596
Suture	597



Surgical Instrument Sets GENERAL SURGERY

BASIC LAPAROTOMY SET

Quantity	Catalog No.	Description
1 ea.	2-104-SS	Yankauer Suction Tube
1 ea.	3-70	Forceps Holder 8" (20.3 cm)
1 ea.	3-71	Forceps Holder 10" (25.4 cm)
1 ea.	3-72	Forceps Holder 12" (30.5 cm)
2 ea.	4-7	Knife Handle No. 3
2 ea.	4-8	Knife Handle No. 4
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-16	Operating Scissors 51/2" (14 cm), str. S/B
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str., Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	5-186TC	Metzenbaum Scissors 9" (22.9 cm), str., Carb-N-Sert blades
1 ea.	5-188TC	Metzenbaum Scissors 9" (22.9 cm), cvd., Carb-N-Sert blades
1 ea.	6-8	Thumb Dressing Forceps 5½" (14 cm)
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
1 ea.	6-118	Adson Dressing Forceps 4%" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-144	Russian Tissue Forceps 8" (20.3 cm)
1 ea.	6-166	Cushing Dressing Forceps 7" (17.8 cm), serr.
1 ea.	6-168	Cushing Tissue Forceps 7" (17.8 cm), 1 x 2 teeth
6 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
6 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
6 ea.	7-36	Kelly Forceps 5½" (14 cm), str.
12 ea.	7-38	Kelly Forceps 51/2" (14 cm), cvd.
6 ea.	7-138	Rochester-Pean Forceps 6¼" (15.9 cm), cvd.
2 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm), cvd.
4 ea.	7-150	Rochester-Ochsner Forceps 61/4" (15.9 cm), str., 1 x 2 teeth
1 ea.	7-202	Mixter Forceps 71/4" (18.4 cm), fully cvd.
1 ea.	7-206	Mixter Forceps 9" (22.9 cm)
1 ea.	7-210	Baby Mixter Forceps 51/4" (13.3 cm), fully cvd.
8 ea.	7-506	Backhaus Towel Clamp 5¼" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str. serr.
1 ea.	8-44TC	Mayo-Hegar Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws

Surgical Instrument Sets **GENERAL SURGERY**



BASIC I APAROTOMY SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	8-48TC	Mayo-Hegar Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	8-130TC	Mayo-Hegar Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
1 ea.	10-312	Poole Suction Tube, str.
2 ea.	11-120	Goelet Retractor
1 ea.	11-122	U.S. Army Retractor, Set of 2
1 ea.	11-130	Ribbon Retractor ¾" (1.9 cm) x 13" (33 cm)
1 ea.	11-134	Ribbon Retractor 11/4" (3.2 cm) x 13" (33 cm)
1 ea.	11-140	Ribbon Retractor 2" (5.1 cm) x 13' (33 cm)
1 ea.	11-202	Deaver Retractor 1" (2.5 cm) x 12" (30.5 cm)
1 ea.	11-208	Deaver Retractor 2" (5.1 cm) x 12" (30.5 cm)
1 ea.	11-246	Richardson Retractor, small
1 ea.	11-248	Richardson Retractor, large
1 ea.	11-267	Kelly Retractor 21/2" (6.4 cm) x 2" (5.1 cm)
1 ea.	11-432	Balfour Retractor w/interchangeable blades
2 ea.	14-48	Lahey Gall Duct Forceps 7½" (19.1 cm)
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
2 ea.	16-12	Allis Tissue Forceps 10" (25.4 cm)
2 ea.	16 -4 4	Babcock Tissue Forceps 61/4" (15.9 cm)
2 ea.	16-48	Babcock Tissue Forceps 91/4" (23.5 cm)
1 ea.	18-660	Stainless Steel Ruler
4 ea.	22-662	Schnidt Tonsil Forceps 71/2" (19.1 cm)
1 ea.	24-572	DeBakey Tissue Forceps 8" (20.3 cm)

MINOR BASIC INSTRUMENT SET

Quantity	Catalog No.	Description
1 ea.	2-104-SS	Yankauer Suction Tube 11¾" (29.8 cm)
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-16	Operating Scissors 5½" (14 cm), str. S/B
1 ea.	5-120TC	Mayo Dissecting Scissors 5½" (14 cm), str. Carb-N-Sert blades
1 ea.	5-122TC	Mayo Dissecting Scissors 5½" (14 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	6-8	Dressing Forceps 5½" (14 cm)
4 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.



GENERAL SURGERY

MINOR BASIC INSTRU	JMENT SET	(Cont'd)	
--------------------	-----------	----------	--

Quantity	Catalog No.	Description
4 ea.	7-4	Halsted Mosquito Forceps, 5" (12.7 cm), cvd.
4 ea.	7-38	Kelly Forceps, 51/2" (14 cm), cvd.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str. serr.
1 ea.	8-44TC	Mayo-Hegar Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	8-50TC	Crile Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	10-26-SS	Probe with Eye 51/2" (14 cm)
1 ea.	10-82	Grooved Director with Probe Tip and Tongue Tie 51/2" (14 cm)
2 ea.	11-74	Senn Retractor, double end, sharp
2 ea.	11-552	Volkman Retractor, 3 prong, blunt
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr. (3 mm)
2 ea.	22-662	Schnidt Hemostatic Forceps 7½" (19.1 cm), slight curve

PLASTIC SET

PLASTIC SET		
Quantity	Catalog No.	Description
1 ea.	5-304	Iris Scissors 41/2" (11.4 cm), str.
1 ea.	5-306	Iris Scissors 41/2" (11.4 cm), cvd.
1 ea.	5-180	Baby Metzenbaum Scissors 51/2" (14 cm), cvd.
1 ea.	5-270	Plastic Surgery Scissors 4¾" (12.1 cm), str. S/S
1 ea.	5-290	Ragnell Scissors 5" (12.7 cm)
1 ea.	6-118	Adson Dressing Forceps, 4¾" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps, 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-124	Brown-Adson Tissue Forceps
2 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
2 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
1 ea.	7-38	Kelly Forceps 5½" (14 cm), cvd.
1 ea.	7-48	Crile Forceps 61/4" (15.9 cm), cvd.
2 ea.	7-255	Gemini-Mixter Forceps 51/2" (14 cm), full curved
8 ea.	7-504	Backhaus Towel Clamp 31/2" (8.9 cm)



GENERAL SURGERY

PLASTIC	SET ((Cont'd)
---------	-------	----------

Quantity	Catalog No.	Description
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm)
1 ea.	8-6TC	Webster Needle Holder 5' (12.7 cm), Carb-N-Sert jaws
1 ea.	8-50TC	Crile Wood N.H. 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	11-2	Jansen Retractor
2 ea.	11-74	Senn Retractor, sharp
1 ea.	11-602	Weitlaner Retractor 4" (10.2 cm), blunt
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr.
2 ea.	21-153	Joseph Skin Hook, single
2 ea.	21-156	Joseph Skin Hook, double 5 mm

HAND SET

IIMID OLI		
Quantity	Catalog No.	Description
2 ea.	4-7	Knife Handle No. 3
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-182A-TC	Metz Scissors 7" (17.8 cm), cvd. del. Carb-N ₇ Sert blades
1 ea.	5-182TC	Metz Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	5-304TC	Iris Scs 4¾" (11.4 cm), str. Carb-N-Sert blades
2 ea.	6-8	Dressing Forceps 5½" (14 cm)
2 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
2 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
2 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
6 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
6 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
6 ea.	7-38	Kelly Forceps 51/2" (14 cm), cvd.
4 ea.	7-504	Backhaus Towel Clamp 31/2" (8.9 cm)
2 ea.	7-604	Foerster Sponge Forceps 7" (17.8 cm), str. serr.
2 ea.	8-42TC	Mayo-Hegar NH 5" (12.7 cm), Carb-N-Sert jaws
2 ea.	11-50	Rigid Shaft Retractor, 1 sharp prong
2 ea.	11-102	Cushing Vein Retractor
2 ea.	11-538	Volkman Retractor, ring handle, 2 blunt prongs
2 ea.	16-8	Allis Tissue Forceps 6" (15.2 cm), 5 x 6 T



GENERAL SURGERY

PEDIATRIC	LAPAROTOMY	SET (add Plastic Set)

Quantity	Catalog No.	Description
1 ea.	11-112	Baby Farabeuf Retractor 4¾" (12.1 cm)
1 ea.	11-144	Pediatric Ribbon Retractor %" (1.6 cm)
2 ea.	11-219	Pediatric Deaver Retractor ¾" (1.9 cm)
2 ea.	11-246	Richardson Retractor 1" (2.5 cm) x ¾" (1.9 cm)
1 ea.	11-430	Pediatric Balfour Retractor
6 ea.	16-4	Baby Allis Tissue Forceps 51/2" (14 cm), 4 x 5 T
1 ea.	18-660	Stainless Steel Ruler

APPENDECTOMY AND HERNIA SET

Quantity	Catalog No.	Description
2 ea.	4-7	Knife Handle No. 3
2 ea.	4-8	Knife Handle No. 4
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Dissecting Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	6-8	Dressing Forceps 51/2" (14 cm)
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-144	Russian Tissue Forceps 8" (20.3 cm)
2 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
6 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
4 ea.	7-38	Kelly Forceps 51/2" (14 cm), cvd.
6 ea.	7-48	Crile Hemostatic Forceps 61/4" (15.9 cm), cvd.
6 ea.	7-154	Rochester Ochsner Forceps 8" (20.3 cm), 1 x 2 T, str.
8 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 9" (22.9 cm), str. serr.
1 ea.	8-44TC	Mayo-Hegar Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	10-310	Poole Suction Tube, cvd.
1 ea.	11-122	U.S. Army Retractor, Set of 2
1 ea.	11-130	Ribbon Retractor 3/4" (1.9 cm), x 13" (33 cm)
1 ea.	11-136	Ribbon Retractor 11/2" (3.8 cm) x 13" (33 cm)



GENERAL SURGERY

APPENDECTOMY AND HERNIA SET (Cont'd)

Quantity	Catalog No.	Description
2 ea.	11-247	Richardson Retractor 11/4" (3.2 cm) x 1" (2.5 cm)
2 ea.	11-267	Kelly Retractor 21/2" (6.4 cm) x 2" (5.1 cm)
4 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
2 ea.	16-44	Babcock Tissue Forceps 61/4" (15.9 cm)

COMMON DUCT INSTRUMENTS (add Basic Laparotomy Set)

COMMON DOCT INSTROMENTS (add basic caparotomy set)		
Catalog No.	Description	
7-258	Gemini-Mixter Forceps 9" (22.9 cm)	
14-2	Ochsner (Fenger) Gall Duct Probe	
14-8	Bakes Common Duct Dilator, Set of 9	
14-26-L	Ferguson Gall Scoop, large	
14-26-M	Ferguson Gall Scoop, medium	
14-26-S	Ferguson Gall Scoop, small	
14-49	Lahey Gall Duct Forceps 9" (22.9 cm)	
29-284	Randall Stone Forceps, 1/4 curve	
29-286	Randall Stone Forceps, 1/2 curve	
29-288	Randall Stone Forceps, ¾ curve	
29-290	Randall Stone Forceps, full curve	
	Catalog No. 7-258 14-2 14-8 14-26-L 14-26-M 14-26-S 14-49 29-284 29-286	

BOWEL RESECTION SET (and Basic Langrotomy Set)

DOWNER THE DECITION OF I (and basic capacitomy only		
Quantity	Catalog No.	Description
8 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
2 ea.	7-180	Rochester Carmalt Forceps 8" (20.3 cm), str.
1 ea.	16-56	Williams Tissue Forceps 61/2" (16.5 cm)
2 ea.	16-100	Allen Intestinal Clamp 8" (20.3 cm)
2 ea.	16-160	Doyen Intestinal Forceps 9" (22.9 cm), str.
2 ea.	16-230	Payr Pylorus Clamp 8" (20.3 cm)
1 ea.	24-552	DeBakey Tissue Forcens 73/4" (19.7 cm)

ABDOMINOPERINEAL RESECTION SET (add Basic Laparotomy Set and Bowel Resection Set)

Quantity	Catalog No.	Description
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str., Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd., Carb-N-Sert blades



Surgical Instrument Sets GENERAL SURGERY

ABDOMINOPERINEAL RESECTION SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	5-130TC	Mayo Dissecting Scissors 9" (22.9 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	6-210	Mayo-Russian Tissue Forceps 9" (22.9 cm)
2 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm), cvd.
2 ea.	7-154	Rochester-Ochsner Forceps 8" (20.3 cm), str.
1 ea.	7-259	Gemini-Mixter Forceps 11" (27.9 cm)
1 ea.	8-52TC	Crile-Wood Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	8-134TC	Wangensteen Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
1 ea.	11-350	Israel Retractor, 4 prongs
1 ea.	11-352	Israel Retractor, 6 prongs
2 ea.	11-542	Volkman Retractor, 4 prongs, blunt
4 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T

THYROIDECTOMY SET (add Plastic Set)

Quantity	Catalog No.	Description
1 ea.	5-182B-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., sharp points, Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
2 ea.	7-98	Lovelace Forceps 61/4" (15.9 cm), str., 1 x 2 T
1 ea.	11-20	Lahey Retractor
1 ea.	11-22	Green Retractor 81/2" (21.6 cm)
1 ea.	11-23	Little Retractor 7¾" (19.7 cm)
1 ea.	11-102	Cushing Vein Retractor 9" (22.9 cm)
6 ea.	12-80	Lahey Traction Forceps 61/4" (15.9 cm)



GENITO-URINARY SURGERY

Catalog No.	Description
5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), delicate, Carb-N-Sert blades
7-206	Mixter Gallbladder Forceps 9" (22.9 cm)
7-252	Mosquito Forceps 81/4" (21 cm), cvd.
7-258	Gemini Mixter Forceps 9" (22.9 cm)
7-259	Gemini Mixter Forceps 11" (27.9 cm)
7-622	Bozeman Sponge Forceps 101/2" (26.7 cm), double curve
8-100TC	Heaney Needle Holder 81/4" (21.6 cm), Carb-N-Sert jaws
8-114TC	Sarot Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
10-26-ST	Probe w/eye 5½" (14 cm), sterling
10-82	Grooved Director 51/2" (14 cm), w/probe tip
11-102	Cushing Vein Retractor
13-68	Ochsner Trocar, impr. patt., 20 Fr. (6.6 mm)
14-56	Grey Cystic Duct Forceps 9" (22.9 cm), Set of 2
16-48	Babcock Tissue Forceps 91/4" (23.5 cm)
21-180	Love Nerve Retractor, str.
24-2110	Satinsky-DeBakey Forceps 91/2" (24.1 cm)
26-1352	Love Nerve Retractor, 45° angle
29-284	Randall Stone Forceps, ¼ cvd.
29-286	Randall Stone Forceps, ½ cvd.
29-288	Randall Stone Forceps, ¾ cvd.
29-290	Randall Stone Forceps, full cvd.
29-304	Herrick Pedicle Clamp 91/2" (24.1 cm)
29-308	Mayo Guyon Vessel Clamp 91/4" (23.5 cm)
	5-182A-TC 7-206 7-252 7-258 7-259 7-622 8-100TC 8-114TC 10-26-ST 10-82 11-102 13-68 14-56 16-48 21-180 24-2110 26-1352 29-284 29-286 29-288 29-290 29-304

PYELOPLASTY AND URETEROPLASTY SET (add Basic Laparotomy Set and Nephrectomy Set)

Quantity	Catalog No.	Description
1 ea.	5-300	Iris Scissors 4" (10.2 cm)
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
2 ea.	6-124	Brown-Adson Tissue Forceps 4¾" (12.1 cm), 7 x 7 T
12 ea.	7-10	Halsted Mosquito Forceps 5" (12.7 cm), cvd. delicate
1 ea.	8-7A-TC	Webster Needle Holder 4¾" (12.1 cm), Carb-N-Sert jaws
1 ea.	14-27	Mayo Common Duct Scoop, small
1 ea.	14-28	Mayo Common Duct Scoop, medium
1 ea.	18-660	Stainless Steel Ruler
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr. (3 mm)



GENITO-URINARY SURGERY



SUPRAPUBIC AND RETROPUBIC PROSTATECTOMY SET (add Basic Laparotomy Set)

Quantity	Catalog No.	Description
1 ea.	5-128TC	Mayo Dissecting Scissors 9" (22.9 cm), str., Carb-N-Sert blades
1 ea.	5-130TC	Mayo Dissecting Scissors 9" (22.9 cm), cvd., Carb-N-Sert blades
6 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm), cvd.
2 ea.	11-247	Richardson Retractor 11/4" (3.2 cm) x 1" (2.5 cm)
2 ea.	11-266	Kelly Retractor 2" (5.1 cm) x 11/2" (3.8 cm)
1 ea.	11-392	Masson-Judd Bladder Retractor, with center blade
1 ea.	29-2	Van Buren Sounds, Set of 6
1 ea.	29-310TC	Bihrle-Dorsal Clamp/Needle Holder, 11" (27.9 cm), Carb-N-Sert jaws

VASECTOMY - MEATOTOMY SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-304	Iris Scissors 41/2" (11.4 cm), str.
1 ea.	5-306	Iris Scissors 41/2" (11.4 cm), cvd.
1 ea.	5-180TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
2 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
2 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
4 ea.	7-504	Backhaus Towel Clamp 31/2" (8.9 cm)
1 ea.	8-42TC	Mayo-Hegar Needle Holder 5" (12.7 cm), Carb-N-Sert jaws
1 ea.	16-3	Baby Allis Tissue Forceps 5" (12.7 cm)

TUROPI ASTY SET (add Basic Lanarotomy Set)

Quantity	Catalog No.	Description
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-300	Iris Scissors 4" (10.2 cm), str.
1 ea.	5-302	Iris Scissors 4" (10.2 cm), cvd.
1 ea.	6-124	Brown-Adson Tissue Forceps 4¾" (12.1 cm), 7 x 7 T
1 ea.	6-154	Potts Smith Dressing Forceps 7" (17.8 cm), serr.
4 ea.	7-10	Halsted Mosquito Forceps 5" (12.7 cm), delicate, cvd.
1 ea.	8-7A-TC	Webster Needle Holder 4¾" (12.1 cm), Carb-N-Sert jaws
1 ea.	10-26-ST	Probe w/eye 51/2" (14 cm), sterling
1 ea.	17-303C	Swiss Jewelers Style Forceps No. 3C
		Cont'd on next page





GENITO-URINARY SURGERY

TUBOPLASTY SET (Cont'd)

Quantity	Catalog No.	Description
2 ea.	18-416	Guthrie Hook, double prongs
1 ea.	18-740 to 18-750	Bowman Sterling Probes, Set of 6
1 ea.	18-950	Castroviejo Suturing Forceps
1 ea.	18-1625	McPherson-Vannas (Iris) Scissors 31/4" (8.3 cm), str.
1 ea.	18-1831	Castroviejo Needle Holder, w/lock, delicate
1 ea.	20-460	Frazier Ferguson Suction Tube, str. 9 Fr. (3 mm)
2 ea.	21-152	Miltex Skin Hook, 2 mm



RECTAL SURGERY

RECTAL INSTRUMENT SE	Т
----------------------	---

Quantity	Catalog No.	Description
1 ea.	5-16	Operating Scissors 5½" (14 cm), str. S/B
1 ea.	5-122TC	Mayo Operating Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str., Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd., Carb-N-Sert blade
1 ea.	6-10	Dressing Forceps 6" (15.2 cm)
4 ea.	7-46	Crile Forceps 61/4" (15.9 cm), str.
2 ea.	7-138	Rochester-Pean Forceps 61/4" (15.9 cm), cvd.
4 ea.	7-220	Pennington Forceps 6" (15.2 cm)
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
2 ea.	11-553	Volkman Retractor, 4 prong, sharp
1 ea.	16-10	Allis Tissue Forceps, 71/4" (18.4 cm), 5 x 6 T
1 ea.	28-2	Sims Rectal Speculum
1 ea.	28-20	Hirschman Anoscope, small
1 ea.	28-22	Hirschman Anoscope, medium
1 ea.	28-24	Hirschman Anoscope, large
1 ea.	28-70	Sawyer Rectal Retractor %" (2.2 cm) x 21/2" (6.4 cm)
1 ea.	28-74	Sawyer Rectal Retractor 1½" (3.8 cm) x 3½" (8.9 cm)
1 ea.	28-100	Larry Probe and Grooved Director 7" (17.8 cm)
1 ea.	28-108	Pratt Rectal Probe 7¾" (19.7 cm) shaft
1 ea.	28-131	Rosser Crypt Hook 6¾" (17.1 cm)
1 ea.	28-370	Buie Suction Tube 15 Fr. (5 mm)
1 ea.	28-386	Miltex Rotating Biopsy Forceps 14" (35.6 cm)

PILONIDAL CYST EXCISION SET (add Rectal Instrument Set)

Quantity	Catalog No.	Description
1 ea.	10-28-ST	Probe w/eye 6" (15.2 cm), sterling
1 ea.	11-553	Volkman Retractor 4 prong, blunt
1 ea.	19-714	Spratt Brun Curette 6½" (16.5 cm), oval, size 3

28-12



RECTAL SURGERY

HEMORRHOIDE	CTOMY	SET
-------------	-------	-----

Quantity	Catalog No.	Description	
2 ea.	4-7	Knife Handle No. 3	
1 ea.	5-180TC	Metzenbaum Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades	
1 ea.	6-8	Dressing Forceps 51/2" (14 cm)	
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T	
4 ea.	7 - 2	Halsted Mosquito Forceps 5" (12.7 cm), str.	
6 ea.	7-506	Backhaus Towel Clamp 51/2" (14 cm)	
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str. serr.	
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T	
1 ea.	28-22	Hirschman Anoscope, medium	
1 ea.	28-32	Brinkerhoff Rectal Speculum, medium	
1 ea.	28-70	Sawyer Rectal Retractor 7/8" (2.2 cm) x 21/2" (6.4 cm)	
1 ea	28-74	Sawyer Rectal Retractor 11/2" (3.8 cm) x 31/2" (8.9 cm)	
FOR OFFICE PRO	CEDURE:		
1 ea.	28-150	McGivney Hemorrhoidal Ligator. Supplied with loading cone and 100 latex O-rings	
1 ea.	28-152	McGivney Hemorrhoid Grasping Forceps 7" (19.1 cm)	

(2.8 cm), chrome

Fansier-ives Rectal Speculum 31/4" (8.3 cm) long x 11/6"

1 ea.



OBSTETRICAL AND GYNECOLOGICAL SURGERY

DILATATION AND CURETTAGE SET			
Quantity	Catalog No.	Description	
1 ea.	4-7	Knife Handle No. 3	
1 ea.	4-11	Knife Handle No. 4L	
1 ea.	6-14	Dressing Forceps 8" (20.3 cm)	
1 ea.	6-128	Heaney Tissue Forceps 51/2" (14 cm), 2 x 3 T	
4 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)	
1 ea.	7-600	Foerster Sponge Forceps 9½" (24.1 cm), str. serr.	
1 ea.	7-622	Bozeman Dressing Forceps 101/2" (26.7 cm), double cvd.	
1 ea.	8-100TC	Heaney Needle Holder 8" (20.3 cm), Carb-N-Sert jaws	
1 ea.	30-15	Graves Vaginal Speculum, medium size	
1 ea.	30-193	Auvard Weighted Speculum 21/2 lbs. (1136 g)	
1 ea.	30-300	Jackson Vaginal Retractor, 11/2" (3.8 cm) x 3" (7.6 cm)	
1 ea.	30-305	Jackson Vaginal Retractor, 11/2" (3.8 cm) x 31/2" (8.9 cm)	
1 ea.	30-310	Jackson Vaginal Retractor, 11/2" (3.8 cm) x 4" (10.2 cm)	
1 ea.	30-550	Hegar Uterine Dilators, Set of 8	
1 ea.	30-650	Sims Uterine Sound	
1 ea.	30-965	Schroeder Braun Tenaculum Forceps 9¾" (24.8 cm), str.	
1 ea.	30-1025	Jacobs Vulsellum Forceps 81/2" (21.6 cm), str.	
1 ea.	30-1205-1	Sims Uterine Curette, sharp, size 1	
1 ea.	30-1205-3	Sims Uterine Curette, sharp, size 3	
1 ea.	30-1205-5	Sims Uterine Curette, sharp, size 5	
1 ea.	30-1380	Sims Uterine Biopsy Curette	

BASIC VAGINAL INSTRUMENT SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-11	Knife Handle No. 4L
1 ea.	5-16	Operating Scissors 5½" (14 cm), str. S/B
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6%4" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	6-8	Dressing Forceps 5½" (14 cm)
1 ea.	6-14	Dressing Forceps 8" (20.3 cm)
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
1 ea.	6-48	Tissue Forceps 8" (20.3 cm), 1 x 2 T





OBSTETRICAL AND GYNECOLOGICAL SURGERY

BASIC	VAGINAL	INSTRUMENT	SET	(Cont'd)	
-------	---------	------------	-----	----------	--

Quantity	Catalog No.	Description	
4 ea.	7-36	Kelly Forceps 51/2" (14 cm), str.	
4 ea.	7-38	Kelly Forceps 5½" (14 cm), cvd.	
4 ea.	7-138	Rochester-Pean Forceps 61/4" (15.9 cm), cvd.	
2 ea.	7-150	Rochester-Ochsner Forceps 61/4" (15.9 cm), str.	
8 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)	
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str. serr.	
1 ea.	7-622	Bozeman Dressing Forceps 101/2" (26.7 cm), double curved	
1 ea.	11-116	Parker Retractor, Set of 2	
1 ea.	11-200	Deaver Retractor 1" (2.5 cm) x 9" (22.9 cm)	
2 ea.	16-10	Allis Tissue Forceps 71/4" (18.4 cm), 5 x 6 T	
1 ea.	30-193	Auvard Weighted Speculum, 21/2 lbs. (1136 g)	
1 ea.	30-215	Sims Retractor, medium	
2 ea.	30-965	Schroeder-Braun Tenaculum Forceps 91/2" (24.1 cm)	
1 ea.	30-1020	Skene Vulsellum Forcens 91/2" (24.1 cm), str.	

ABDOMINAL GYNECOLOGICAL SET (add Basic Laparotomy Set)

Quantity	Catalog No.	Description	
1 ea.	4-11	Knife Handle No. 4L	
1 ea.	5-130TC	Mayo Dissecting Scissors 9" (22.9 cm), cvd., Carb-N-Sert blades	
1 ea.	6-52	Tissue Forceps 12" (30.5 cm), 1 x 2 T	
1 ea.	6-128	Heaney Tissue Forceps, 2 x 3 T	
2 ea.	7-124	Rochester-Ochsner Forceps 8" (20.3 cm), str.	
1 ea.	7-600	Foerster Sponge Forceps 9" (22.9 cm), str. serr.	
2 ea.	8-100TC	Heaney Needle Holder 8" (20.3 cm), curved, Carb-N-Sert jaws	
1 ea.	11-410	O'Sullivan - O'Connor Retractor	
2 ea.	16-10	Allis Tissue Forceps 71/4" (18.4 cm), 5 x 6 T	
1 ea.	30-15	Graves Vaginal Speculum, medium	
2 ea.	30-965	Schroeder-Braun Tenaculum Forceps	
1 ea.	30-1010	Schroeder Vulsellum Forceps, cvd. sideways	
1 ea.	30-1025	Jacobs Vulsellum Forceps, str.	
1 ea.	30-1710	Heaney Hysterectomy Forceps 81/4" (21 cm), double tooth	
1 ea.	30-1716	Heaney-Ballantine Hysterectomy Forceps, cvd.	



OBSTETRICAL AND GYNECOLOGICAL SURGERY

LAPAROSCOPY SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
2 ea.	7-91	Coller Forceps 61/4" (15.9 cm), cvd.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 9" (22.9 cm), str. serr.
1 ea.	8-51TC	Crile-Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	11-185	Heaney-Simon Retractor
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
1 ea.	30-15	Graves Vaginal Speculum, medium
1 ea.	30-193	Auvard Weighted Speculum 2½ lbs. (1136 g)
1 ea.	30-650	Sims Uterine Sound, grad. in cm
1 ea.	30-965	Schroeder-Braun Tenaculum Forceps 9¾" (24.8 cm)

BARTHOLIN CYST EXCISION SET

DANTHOLINGIS	LYCIDION OF I	
Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str., Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	6-8	Dressing Forceps 51/2" (14 cm)
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
2 ea.	7-4	Halstead Mosquito Forceps 5" (12.7 cm), cvd.
2 ea.	7-138	Rochester Pean Forceps 61/4" (15.9 cm), cvd.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str. serr.
1 ea.	8-46TC	Mayo Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
2 ea.	16-10	Allis Tissue Forceps 71/2" (19.1 cm), 5 x 6 T
1 ea.	30-193	Auvard Weighted Speculum, 21/2 lbs. (1136 g)
1 ea.	30-215	Sims Vaginal Speculum, double end, medium



OBSTETRICAL AND GYNECOLOGICAL SURGERY

VAGINAL TUBAL LIGATION S

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-11	Knife Handle No. 4L
1 ea.	5-124TC	Mayo Dissecting Scissors 6%" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
1 ea.	6-145	Russian Tissue Forceps 10" (25.4 cm)
6 ea.	7-38	Kelly Forceps 51/2" (14 cm), cvd.
2 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm), cvd.
2 ea.	7-259	Gemini-Mixter Forceps 11" (27.9 cm)-
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 9½" (24.1 cm), str. serr.
1 ea.	8-100TC	Heaney Needle Holder 8" (20.3 cm), cvd., Carb-N-Sert jaws
1 ea.	11-185	Heaney-Simon Retractor 111/2" (29.2 cm)
1 ea.	11-202	Deaver Retractor 1" (2.5 cm) x 12" (30.5 cm)
4 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
2 ea.	16-10	Allis Tissue Forceps 71/4" (18.4 cm), 5 x 6 T
2 ea.	16-46	Babcock Tissue Forceps 81/4" (21 cm)
1 ea.	30-15	Graves Vaginal Speculum, medium
1 ea.	30-193	Auvard Weighted Speculum, 21/2 lbs. (1136 g)
1 ea.	30-650	Sims Uterine Sound, grad. in cm
2 ea.	30-965	Schroeder-Braun Tenaculum Forceps 9½" (24.1 cm), str.
2 ea.	30-1710	Heaney Hysterectomy Forceps 81/4" (21 cm), double tooth



NEUROSURGERY

BASIC CRANIOTOMY SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-120TC	Mayo Dissecting Scissors 51/2" (14 cm), str., Carb-N-Sert blades
1 ea.	5-122TC	Mayo Dissecting Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	5-246	Potts-Smith Scissors 71/2" (19.1 cm), 25° angle
1 ea.	5-250	Adson Scissors 61/4" (15.9 cm), str., sharp points
1 ea.	6-109A	Semkin Tissue Forceps 6" (15.2 cm), 1 x 2 T
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm)
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-168	Cushing Tissue Forceps 7" (17.8 cm), 1 x 2 T
1 ea.	6-198	Adson Bayonet Forceps 81/4" (21 cm), serr.
6 ea.	7-8	Halsted Mosquito Forceps 5" (12.7 cm), str., delicate
6 ea.	7-10	Halsted Mosquito Forceps 5" (12.7 cm), cvd., delicate
12 ea.	7-102	Dandy Scalp Hemostat 5½" (14 cm)
2 ea.	7-120	Rochester-Ochsner Forceps 61/4" (15.9 cm), str.
2 ea.	7-122	Rochester-Ochsner Forceps 71/4" (18.4 cm), str.
2 ea.	7-240	Adson Hemostatic Forceps 7¼" (18.4 cm), str.
2 ea.	7-242	Adson Hemostatic Forceps 71/4" (18.4 cm), cvd.
8 ea.	7-504	Backhaus Towel Clamp 3½" (8.9 cm)
8 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str., serr.
1 ea.	8-50TC	Crile-Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	8-51TC	Crile-Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
	9-2	Michel Wound Clips - 11 mm, one rack w/25 clips
1 ea.	9-124	Wire Cutting Scissors 49/4" (12.1 cm), angular
1 ea.	11-2	Jansen Retractor 41/4" (10.8 cm), blunt
1 ea.	11-8	Jansen-Wagner Retractor 5" (12.7 cm), sharp
2 ea.	11-74	Senn Retractor 6%" (16.2 cm), double end, sharp
2 ea.	11-564	Murphy Retractor 71/2" (19.1 cm), 4 prong, sharp
1 ea.	11-684	Beckman Retractor 121/2" (31.8 cm), sharp
1 ea.	11-700	Adson Retractor 7½" (19.1 cm), sharp
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm)
1 ea.	18-1968	Freer Elevator 7" (17.8 cm), double end
1 ea.	19-572	Frazier Ferguson Suction Tube, 7 Fr. (2.3 mm)

Surgical Instrument Sets NEUROSURGERY



BASIC CRANIOTOMY SET (Cont'd)

DAGIO CHANGIO	Juli OE: (Oom a)	
Quantity	Catalog No.	Description
1 ea.	19-574	Frazier Ferguson Suction Tube, 9 Fr. (3 mm)
1 ea.	19-840	Zaufel-Jansen Rongeur 7" (17.8 cm)
1 ea.	21-180	Love Nerve Retractor 8½" (21.6 cm), str.
1 ea.	25-458	Stille-Luer Rongeur 9" (22.9 cm), angled
1 ea.	26-2	Cushing Perforator Drill
1 ea.	26-24	Cerebellar Extension
1 ea.	26-32	Hudson Brace
1 pr.	26-124	Gigli-Strully Saw Handle, Set of 2
2 ea.	26-132	Gigli Saw 20" (50.8 cm)
1 ea.	26-138	Bailey Gigli Saw Guide 12" (30.5 cm)
1 ea.	26-190	Luer Rongeur 7" (17.8 cm), str.
1 ea.	26-326	Schlesinger Cervical Rongeur 6" (15.2 cm), 3 mm upwards bite
1 ea.	26-389	Miltex-Kerrison Cervical Rongeur 7" (17.8 cm), 3 mm bite, 40° forward
1 ea.	26-397	Cushing Pituitary Rongeur
1 ea.	26-700	Scoville Brain Spatula 8" (20.3 cm)
1 ea.	26-730	Cushing Spatula Spoon 6¾" (17.1 cm)
1 ea.	26-1054	Frazier Dura Separator 6½" (16.5 cm)
2 ea.	26-1084	Frazier Dura Hook 5" (12.7 cm), sharp
1 ea.	26-1090	Adson Dura Hook 8" (20.3 cm)
1 ea.	26-1150	Dandy Nerve Hook 9" (22.9 cm), str.
1 ea.	26-1352	Love Nerve Retractor 81/2" (21.6 cm), angled 45°
1 ea.	26-1354	Love Nerve Retractor 8½" (21.6 cm), angled 90°
1 ea.	26-1403	Cushing Little Joker Periosteal Elevator 6" (15.2 cm)
1 ea.	26-1412	Langenbeck Periosteal Elevator 71/2" (19.1 cm)
1 ea.	26-1432	Adson Periosteal Elevator 61/2" (16.5 cm), blunt

VENTRICULOGRAPHY SET

Quantity	Catalog No.	Description
2 ea.	4-7	Knife Handle No. 3
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-250	Kelly Scissors 61/4" (15.9 cm), str., sharp points



NEUROSURGERY

VENTRICULOGRAPHY SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	6-109A	Semkin Tissue Forceps 6" (15.2 cm), 1 x 2 T
1 ea.	6-186	Cushing Tissue Forceps 7" (17.8 cm), 1 x 2 T
1 ea.	6-190	Cushing Bayonet Forceps 71/4" (18.4 cm), serrated
6 ea.	7-102	Dandy Scalp Hemostat 51/2" (14 cm)
4 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str., serr.
1 ea.	8-51TC	Crile-Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	11-2	Jansen Retractor 41/4" (10.8 cm), blunt
1 ea.	11-8	Jansen-Wagner Retractor 5" (12.7 cm), sharp
1 ea.	11-74	Senn Retractor 6%" (16.2 cm), double end, sharp
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr. (3 mm)
1 ea.	19-704	Spratt Curette 61/2" (16.5 cm), size 3/0
1 ea.	19-708	Spratt Curette 61/2" (16.5 cm), size 0
1 ea.	26-2	Cushing Perforator Drill
1 ea.	26-24	Cerebellar Extension
1 ea.	26-32	Hudson Brace
1 ea.	26-772	Elsberg Cannula 7 Fr. (2.3 mm)
1 ea.	26-778	Frazier Needle 9 Fr. (3 mm)
1 ea.	26-1084	Frazier Dura Hook 5" (12.7 cm), sharp
1 ea.	26-1403	Cushing Little Joker Elevator 6" (15.2 cm)
1 ea.	26-1434	Adson Periosteal Elevator 6½" (16.5 cm), cvd.

NEUROSURGICAL-LAMINECTOMY SET (add Basic Onhopedic Soft Tissue Set)

talog No.	Description
247	Potts-Smith Scissors 7" (17.8 cm), 45° angle
262	Strully Scissors 8" (20.3 cm), with probe tips
118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
190	Cushing Dressing Forceps 71/4" (18.4 cm), bayonet, serr.
192	Cushing Tissue Forceps 71/4" (18.4 cm), bayonet, 1 x 2 T
240	Adson Hemostatic Forceps 71/4" (18.4 cm), str.
242	Adson Hemostatic Forceps 71/4" (18.4 cm), cvd.
600	Foerster Sponge Forceps 9" (22.9 cm), str. serr.
B2TC	Julian Needle Holder 81/4" (21 cm), Carb-N-Sert jaws
	talog No. 247 262 118 120 190 192 240 242 600 82TC



NEUROSURGERY

NEUROSURGICAL-LAMINECTOMY	SET (Cont'd)
---------------------------	--------------

Quantity	Catalog No.	Description
1 ea.	11-366	Gelpi Retractor 6%" (17.1 cm)
1 ea.	11-684	Beckman-Adson Retractor 121/2" (31.8 cm), blade 3/4" (1.9 cm) deep
1 ea.	19-702	Spratt Bone Curette 61/2" (16.5 cm), size 4/0
1 ea.	19-706	Spratt Bone Curette 61/2" (16.5 cm), size 2/0
1 ea.	19-708	Spratt Bone Curette 61/2" (16.5 cm), size 0
1 ea.	20-180	Love Nerve Retractor 8½" (21.6 cm), str.
2 ea.	24-550	DeBakey Tissue Forceps 6" (15.2 cm)
1 ea.	26-24	Hudson Cerebellar Attachment
1 ea.	26-30	Hudson Brace Set, with 4 burs
1 ea.	26-91	Taylor Spinal Retractor 71/4" (18.4 cm)
1 ea.	26-388	Miltex-Kerrison Rongeur 7" (17.8 cm), 2 mm 40° forward angle
1 ea.	26-389	Miltex-Kerrison Rongeur 7" (17.8 cm), 3 mm 40° forward angle
1 ea.	26-390	Miltex-Kerrison Rongeur 7" (17.8 cm), 5 mm 40° forward angle
1 ea.	26-392A	Miltex-Kerrison Rongeur 7" (17.8 cm), 2 mm up-bite
1 ea.	26-393A	Miltex-Kerrison Rongeur 7" (17.8 cm), 2 mm down-bite
1 ea.	26-400	Love-Gruenwald Rongeur 5" (12.7 cm), str., 3 x 10 mm
1 ea.	26-410	Love-Gruenwald Rongeur 5" (12.7 cm), curved up, 3 x 10 mm
1 ea.	26-422	Love-Gruenwald Rongeur 7" (17.8 cm), curved down, 3 x 10 mm
1 ea.	26-660	Adson Aneurysm Needle 9" (22.9 cm)
1 ea.	26-730	Cushing Spatula Spoon 6¾" (17.1 cm)
1 ea.	26-1150	Dandy Nerve Hook 9" (22.9 cm), str.
1 ea.	26-1174	Smithwick Hook and Dissector 12" (30.5 cm), double end
1 ea.	26-1360	Scoville Nerve Root Retractor 9" (22.4 cm), str.
1 ea.	26-1362	Scoville Nerve Root Retractor 9" (22.9 cm), angular
1 ea.	26-1434	Adson Periosteal Elevator, cvd., semi-sharp
1 ea.	26-1451	Penfield Dissector, style 2, double end
1 ea.	26-1973	Beckman-Eaton Retractor 12¾" (32.4 cm), blade 2" (5.1 cm) deep
1 ea.	27-331	Hoke Osteotome 51/2" (14 cm), 1/a" (.32 cm)
1 ea.	27-333	Hoke Osteotome 5½" (14 cm), ¼" (.64 cm)
1 ea.	27-335	Hoke Osteotome 51/2" (14 cm), %" (.95 cm)
1 ea.	27-336	Hoke Osteotome 51/2" (14 cm), 1/2" (1.3 cm)



NEUROSURGERY

SKULL TRACTION SET

Quantity	Catalog No.	Description
Quantity	Catalog 140.	
1 ea.	5-140	Mayo Dissecting Scissors 6¾" (17.1 cm), str.
1 ea.	6-96	Tissue Forceps 6" (15.2 cm), delicate, 1 x 2 T
2 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), straight
2 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), curved
8 ea.	7-504	Backhaus Towel Clamp 31/2" (8.9 cm)
1 ea.	7-600	Foerster Sponge Forceps 9" (22.9 cm), str. serr.
1 ea.	26-40	Stille Pattern Cranial Drill Set, w/set of drill points
1 ea.	26-60	Raney-Crutchfield Skull Tracing Tong w/2 points

Surgical Instrument Sets ORTHOPEDIC SURGERY



BASIC MAJOR ORTHOPEDIC SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
2 ea.	4-8	Knife Handle No. 4
1 ea.	5-120TC	Mayo Dissecting Scissors 51/2" (14 cm), str., Carb-N-Sert blades
1 ea.	5-122TC	Mayo Dissecting Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str., Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd., Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	5-568	Esmarch Bandage and Cast Shears 8" (20.3 cm)
1 ea.	6-46	Tissue Forceps 6" (15.2 cm), 1 x 2 T
1 ea.	6-86	Tissue Forceps 6" (15.2 cm), 4 x 5 T
2 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
6 ea.	7-48	Crile Forceps 61/4" (15.9 cm), cvd.
4 ea.	7-154	Rochester-Ochsner Forceps 8" (20.3 cm), str.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	8-48TC	Mayo-Hegar Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	8-50TC	Crile Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	11-122	U.S. Army Retractor 81/2" (21.6 cm), set of 2
1 ea.	11-240	Meyerding Retractor 91/2" (24.1 cm), 31/2" (8.9 cm) x 2" (5.1 cm)
1 ea.	11-244	Meyerding Retractor 9½" (24.1 cm), 2" (5.1 cm) x %" (1.6 cm)
1 ea.	11-350	Israel Retractor 9" (22.9 cm), 4 prong, blunt
1 ea.	11-548	Volkman Retractor 8¾" (22.2 cm), 4 prong, sharp
2 ea.	11-553	Volkman Retractor 8" (20.3 cm) x ¾" (1.9 cm), 4 prong, blunt
1 ea.	11-580	Bennett Bone Elevator 9½" (24.1 cm), small
1 ea.	11-582	Bennett Bone Elevator 9½" (24.1 cm), large
1 ea.	11-590	Hibbs Retractor 91/4" (23.5 cm)
1 ea.	11-684	Beckman Retractor 12" (30.5 cm), hinged blades
6 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
1 ea.	18-660	Stainless Steel Ruler 6" (15.2 cm)
1 ea.	19-573	Frazier Ferguson Suction Tube 8 Fr. (2.6 mm)
1 ea.	19-577	Frazier Ferguson Suction Tube 12 Fr. (4 mm)
1 ea.	19-844	Beyer Rongeur 7" (17.8 cm), cvd.
2 ea.	22-662	Schnidt Tonsil Forceps 7½" (19.1 cm)
1 ea.	25-385	Ruskin Bone Cutting Forceps 7½" (19.1 cm), angled
1 ea.	25-456	Stille-Luer Rongeur 8¾" (22.2 cm), cvd.



ORTHOPEDIC SURGERY

BASIC MAJOR ORTHOPEDIC SET (Cont'd)

Quantity	Catalog No.	Description
2 ea.	27-2	Kern Bone Holding Forceps 9" (22.9 cm)
1 ea.	27-12	Lambotte Bone Holding Forceps 10¾" (27.3 cm)
1 ea.	27-21	Lane Bone Holding Forceps 13" (33 cm), w/ratchet
1 ea.	27-30	Lowman Bone Clamp 5" (12.7 cm)
1 ea.	27-34	Lowman Bone Clamp 8" (20.3 cm)
·2 ea.	27-230	Bone Hook 9" (22.9 cm), sharp
1 ea.	27-438	Hibbs Osteotome 9" (22.9 cm), str. 1/4" (.64 cm)
1 ea.	27-442	Hibbs Osteotome 9" (22.9 cm), str. 1/2" (1.3 cm)
1 ea.	27-446	Hibbs Osteotome 9" (22.9 cm), str. ¾" (1.9 cm)
1 ea.	27-449	Hibbs Osteotome 9" (22.9 cm), str. 11/4" (2.8 cm)
1 ea.	27-458	Hibbs Gouge 9" (22.9 cm), str. 1/4" (.64 cm)
1 ea.	27-462	Hibbs Gouge 9" (22.9 cm), str. 1/2" (1.3 cm)
1 ea.	27-466	Hibbs Gouge 9" (22.9 cm), str. ¾" (1.9 cm)
1 ea.	27-440C	Hibbs Osteotome 9" (22.9 cm), cvd. %" (.96 cm)
1 ea.	27-444C	Hibbs Osteotome 9" (22.9 cm), cvd. %" (1.6 cm)
1 ea.	27-448C	Hibbs Osteotome 9" (22.9 cm), cvd. 1" (2.5 cm)
1 ea.	27-460C	Hibbs Gouge 9" (22.9 cm), cvd. %" (.96 cm)
1 ea.	27-464C	Hibbs Gouge 9" (22.9 cm), cvd. %" (1.6 cm)
1 ea.	27-660	Miltex Mallet 11" (27.9 cm), 16 oz. (454 g)
1 ea.	27-744	Sayre Elevator 6¾" (17.1 cm), sharp
1 ea.	27-820	Key Elevator 7" (17.8 cm), 1/4" (.64 cm) wide
1 ea.	27-822	Key Elevator 7½" (19.1 cm), ½" (1.3 cm) wide
1 ea.	27-826	Key Elevator 81/4" (21 cm), 1" (2.5 cm) wide

BASIC ORTHOPEDIC SOFT TISSUE SET

BASIC ORTHUPE	DIC SUFT HISSUE SET	
Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-8	Knife Handle No. 4
1 ea.	5-120TC	Mayo Dissecting Scissors 51/2" (14 cm), str. Carb-N-Sert blades
1 ea.	5-122TC	Mayo Dissecting Scissors 51/2" (14 cm), cvd. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	5-514	Lister Bandage Scissors 5½" (14 cm)
1 ea.	6-46	Tissue Forceps 6" (15.2 cm), 1 x 2 T
1 ea.	6-86	Tissue Forceps 6" (15.2 cm), 4 x 5 T

Surgical Instrument Sets ORTHOPEDIC SURGERY



BASIC ORTHOPEDIC SOFT TISSUE SET (Cont'd)

BASIC URITIO	OPEDIC SOFT TISSUE SET (CONT	a)
Quantity	Catalog No.	Description
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
4 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
4 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
6 ea.	7-38	Kelly Forceps, cvd.
4 ea.	7-66	Kocher Forceps 51/2" (14 cm), str.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	8-6TC	Webster Needle Holder 5" (12.7 cm), Carb-N-Sert jaws
1 ea.	8-50TC	Crile Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	<i>∂</i> 11-74	Senn Retractor 61/4" (15.9 cm), sharp
1 ea.	11-76	Senn Retractor 61/4" (15.9 cm), blunt
1 ea.	11-122	U.S. Army Retractor, Set of 2
1 ea.	11-160	Langenbeck Retractor 8" (20.3 cm), blade %" (.96 cm) x 11/4" (3.2 cm)
1 ea.	11-164	Langenbeck Retractor 8" (20.3 cm), blade ½" (1.3 cm) x 1½" (3.8 cm)
1 ea.	11-546	Volkman Retractor, 2 prongs, sharp
1 ea.	11-547	Volkman Retractor, 3 prongs, sharp
1 ea.	11-590	Hibbs Retractor 91/4" (23.5 cm), double end
1 ea.	11-608	Weitlaner Retractor 51/2" (14 cm), sharp
1 ea.	11-610	Weitlaner Retractor 5½" (14 cm), blunt
4 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm), 4 x 5 T
1 ea.	18-660	Stainless Steel Ruler
1 ea.	18-1968	Freer Elevator, double end
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr. (3 mm), angled
1 ea.	19-794	Miltex Mallet 71/2" (19.1 cm), nylon and stainless faces
1 ea.	19-848	Ruskin Rongeur 6" (15.2 cm), str.
1 ea.	20-460	Frazier Ferguson Suction Tube 9 Fr. (3 mm), str.
2 ea.	21-153	Joseph Hook 61/4" (15.9 cm), single prong
2 ea.	21-158	Joseph Hook 61/4" (15.9 cm), double prong 7 mm
1 ea.	25-384	Ruskin Bone Cutting Forceps 71/2" (19.1 cm), str.
1 ea.	26-1434	Adson Periosteal Elevator, curved, semi sharp
1 ea.	27-744	Sayre Elevator, double end, sharp
1 ea.	27-820	Key Elevator 7" (17.8 cm), 1/4" (.64 cm) wide
1 ea.	27-822	Key Elevator 7½" (19.1 cm), ½" (1.3 cm) wide



ORTHOPEDIC SURGERY

KNEE SET		
Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6%" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-180TC	Metzenbaum Scissors 51/2" (14 cm), cvd. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	5-514	Lister Bandage Scissors 51/2" (14 cm)
1 ea.	6-44	Tissue Forceps 5½" (14 cm), 1 x 2 T
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
2 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-134	Martin Tissue Forceps 6" (15.2 cm), 7 x 8 T
6 ea.	7-38	Kelly Forceps 51/2" (14 cm), cvd.
4 ea.	7-150	Rochester-Ochsner Forceps 61/4" (15.9 cm), str.
6 ea.	7-504	Backhaus Towel Clamp 3½" (8.9 cm)
1 ea.	8-6TC	Webster Needle Holder 5" (12.7 cm), Carb-N-Sert jaws
1 ea.	8-51TC	Crile Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	11-74	Senn Retractor 61/2" (16.5 cm), sharp
1 ea.	11-76	Senn Retractor 61/2" (16.5 cm), blunt
1 ea.	11-122	U.S. Army Retractor, Set of 2
1 ea.	11-546	Volkman Retractor, 2 prong, sharp
1 ea.	11-548	Volkman Retractor, 4 prong, sharp
1 ea.	11-610	Weitlaner Retractor 51/2" (14 cm), blunt
4 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm)
1 ea.	18-660	Stainless Steel Ruler
1 ea.	18-1968	Freer Elevator, double end
1 ea.	19-794	Miltex Mallet 7" (17.8 cm), nylon and stainless faces
1 ea.	26-1092	Adson Hook 8" (20.3 cm), blunt
1 ea.	26-1156	Graham Hook 7" (17.8 cm), blunt
1 ea.	26-1434	Adson Periosteal Elevator, curved, semi-sharp
2 ea.	27-223	Blount Knee Retractor 7" (17.8 cm)
1 ea.	27-224	Smillie Hook Retractor 5¾" (14.6 cm)
1 ea.	27-744	Sayre Elevator, double end, sharp
1 ea.	27-821	Key Elevator 7½" (17.8 cm), %" (.96 cm) wide
1 ea.	27-826	Key Elevator 81/4" (21 cm), 1" (2.5 cm) wide



ORTHOPEDIC SURGERY

KNEE SET (Cont'd)		
Quantity	Catalog No.	Description
1 ea.	27-970	Smillie Meniscus Knife 71/4" (18.4 cm), str.
1 ea.	27-972	Smillie Meniscus Knife 7" (17.8 cm), cvd. left
1 ea.	27-974	Smillie Meniscus Knife 7" (17.8 cm), cvd. right
1 ea.	27-980	Lowe-Breck Cartilage Knife 10½" (26.7 cm)
1 ea.	27-984	Downing Cartilage Knife 10" (25.4 cm)
1 ea.	27-990	Martin Cartilage Clamp 7" (17.8 cm)
1 ea.	27-1000	Martin Cartilage Scissors 8" (20.3 cm)

BASIC TOTAL HIP SET		
Quantity	Catalog No.	Description
3 ea.	2-104SS	Yankauer Suction Tube
2 ea.	4-8	Knife Handle No. 4
1 ea.	4-9	Knife Handle No. 7
1 ea.	4-11	Knife Handle No. 4L
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-130TC	Mayo Dissecting Scissors 9" (22.9 cm), cvd. Carb-N-Sert blades
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd. Carb-N-Sert blades
1 ea.	5-184TC	Metzenbaum Scissors 8" (20.3 cm), cvd. Carb-N-Sert blades
1 ea.	6-44	Tissue Forceps 51/2" (14 cm), 1 x 2 T
1 ea.	6- 84	Tissue Forceps 5½" (14 cm), 4 x 5 T
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-142	Russian Tissue Forceps 6" (15.2 cm)
1 ea.	6-144	Russian Tissue Forceps 8" (20.3 cm)
6 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
6 ea.	7-44	Crile Forceps 5½" (14 cm), cvd.
6 ea.	7-48	Crile Forceps 61/4" (15.9 cm), cvd.
6 ea.	7-91	Coller Forceps 61/4" (15.9 cm), cvd.
2 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm), cvd.
4 ea.	7-154	Rochester-Ochsner Forceps 8" (20.3 cm), str.
2 ea.	7-258	Gemini-Mixter Forceps 9" (22.9 cm)
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 9½" (24.1 cm), str. serr.
		Cont'd on next page



ORTHOPEDIC SURGERY

BASIC	TOTAL	HIP SET	(Cont'd)
-------	-------	---------	----------

DAGIC TOTAL III		
Quantity	Catalog No.	Description
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	8-48TC	Mayo-Hegar Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	8-51TC	Crile Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	11-260	Richardson-Eastman Retractor, Set of 2
1 ea.	11-366	Gelpi Retractor 6¾" (17.1 cm)
2 ea.	11-548	Volkman Retractor 4 prong, sharp
1 ea.	11-582	Bennett Bone Elevator 9¾" (24.8 cm), large blade
1 ea.	11-590	Hibbs Retractor 91/4" (23.5 cm)
1 ea.	11-614	Weitlaner Retractor 6¾" (16.5 cm), sharp
1 ea.	11-700	Adson Retractor 71/4" (19.1 cm)
6 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm)
4 ea.	16-10	Allis Tissue Forceps 71/4" (18.4 cm)
1 ea.	18-660	Stainless Steel Ruler
1 ea.	25-396	Stille-Horsley Bone Cutting Forceps 101/2" (26.7 cm)
1 ea.	25-456	Stille-Luer Rongeur 8¾" (22.2 cm), cvd.
1 ea.	25-458	Stille-Luer Rongeur 9" (22.9 cm), duckbill
2 ea.	26-122	Gigli Strully Saw Handle
1 ea.	26-132	Gigli Saw 20" (50.8 cm)
1 ea.	26-138	Bailey Saw Guide
1 ea.	26-1412	Langenbeck Periosteal Elevator, wide model
1 ea.	26-1628	Spratt Curette 9" (22.9 cm), size 00
1 ea.	26-1632	Spratt Curette 9" (22.9 cm), size 2
1 ea.	26-1636	Spratt Curette 9" (22.9 cm), size 4
1 ea.	26-1640	Spratt Curette 9" (22.9 cm), size 6
1 ea.	26-1973	Beckman-Eaton Retractor 12¾" (32.4 cm)
1 ea.	27-215	Hohman Retractor 91/4" (23.5 cm), 10 mm wide
1 ea.	27-217	Hohman Retractor 91/4" (23.5 cm), 25 mm wide
1 ea.	27-219	Hohman Retractor 91/4" (23.5 cm), 70 mm wide
1 ea.	27-230	Volkman Bone Hook 9" (22.9 cm), sharp
1 ea.	27-438	Hibbs Osteotome 9" (22.9 cm), 1/4" (.64 cm)
1 ea.	27-440	Hibbs Osteotome 9" (22.9 cm), %" (.96 cm)
1 ea.	27-442	Hibbs Osteotome 9" (22.9 cm), 1/2" (1.3 cm)
1 ea.	27-444	Hibbs Osteotome 9" (22.9 cm), %" (1.6 cm)
1 ea.	27-446	Hibbs Osteotome 9" (22.9 cm), 3/4" (1.9 cm)



ORTHOPEDIC SURGERY

BASIC TOTAL HIP SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	27-660	Miltex Mallet, 16 oz. (454 g)
1 ea.	27-680	Murphy Bone Skid 12" (30.5 cm)
1 ea.	27-702	Putti Bone Rasp 12" (30.5 cm), round and half-round
1 ea.	27-826	Key Periosteal Elevator 1" (2.5 cm)

AMPUTATION SET (add Orthopedic Soft Tissue Set)

Quantity	Catalog No.	Description
2 ea.	26-122	Gigli Strully Saw Handle
1 ea.	26-132	Gigli Saw 20" (50.8 cm)
1 ea.	26-138	Bailey Saw Guide
1 ea.	26-190	Luer Bone Ronguer 7" (17.8 cm), straight
1 ea.	26-1412	Langenbeck Periosteal Elevator 7½" (19.1 cm), wide, sharp
1 ea.	27-438	Hibbs Osteotome 9" (22.9 cm), 1/4" (.64 cm)
1 ea.	27-442	Hibbs Osteotome 9" (22.9 cm), 1/2" (1.3 cm)
1 ea.	27-446	Hibbs Osteotome 9" (22.9 cm), ¾" (1.9 cm)
1 ea.	27-448	Hibbs Osteotome 9" (22.9 cm), 1" (2.5 cm)
1 ea.	27-450	Hibbs Osteotome 9" (22.9 cm), 11/4" (3.2 cm)
1 ea.	27-700	Putti Bone Rasp 10½" (26.7 cm), flat ends
1 ea.	27-1392	Liston Bone Cutting Forceps 71/2" (19.1 cm), straight
1 ea.	27-1560	Satterlee Saw, 9" (22.9 cm) blade
1 ea.	27-1602	Liston Amputation Knife, 8" (20.3 cm) blade



Surgical Instrument Sets **THORACIC AND CARDIAC SURGERY**

BASIC RIB INSTRUMENT SET

Quantity	Catalog No.	Description
1 ea.	8-55TC	DeBakey Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
6 ea.	16-12	Allis Tissue Forceps 10" (25.4 cm)
2 ea.	24-1602	DeBakey Multipurpose Clamp 81/4" (21 cm), 5 cm jaws
1 ea.	25-2	Doyen Rib Elevator and Raspatory, right
1 ea.	25-4	Doyen Rib Elevator and Raspatory, left
1 ea.	25-14	Matson-Alexander Periosteotome 8¾" (22.2 cm)
1 ea.	25-18	Matson Rib Elevator and Stripper 8¾" (22.2 cm)
1 ea.	25-142	Burford-Finochietto Rib Spreader, w/2 pr. blades
1 ea.	25-200	Bailey Rib Contractor
1 ea.	25-302	Allison Lung Retractor 123/4" (32.4 cm)
1 ea.	25-350	Bethune Rib Shears 131/2" (34.9 cm)
1 ea.	25-396	Stille-Horsley Bone Cutting Forceps 101/2" (26.7 cm)
1 ea.	25-480	Sauerbruch Rongeur 12" (30.5 cm), square jaws 17 mm
1 ea.	25-510	Leksell Rongeur 91/2" (24.1 cm), 3 mm bite
2 ea.	25-754	Duval Tissue Grasping Forceps 9" (22.9 cm), 1" (2.5 cm) wide

THORACOTOMY SET (add Basic Laparotomy Set)

Quantity	Catalog No.	Description
1 ea.	2-10455	Yankauer Suction Tube
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-182B-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., del., Carb-N-Sert blades
1 ea.	5-189TC	Metzenbaum Scissors 9" (22.9 cm), cvd., del., Carb-N-Sert blades
1 ea.	6-136	Martin Tissue Forceps 71/2" (19.1 cm), 7 x 8 T
1 ea.	6-158TC	Potts-Smith Tissue Forceps 9¾" (24.8 cm), Carb-N-Sert jaws
12 ea.	7-38	Kelly Forceps 5½" (14 cm), cvd.
12 ea.	7-91	Coller Forceps 61/4" (15.9 cm), cvd.
12 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	8-52TC	Crile Wood Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	8-55TC	Crile Wood Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
1 ea.	8-81TC	New Orleans Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	8-130TC	Mayo-Hegar Needle Holder 101/2" (26.7 cm), Carb-N-Sert jaws
1 ea.	11-130	Ribbon Retractor ¾" (1.9 cm), x 13" (33 cm)
1 ea.	11-134	Ribbon Retractor 11/4" (3.2 cm) x 13" (33 cm)
1 ea.	11-140	Ribbon Retractor 2" (5.1 cm) x 13" (33 cm)





THORACIC AND CARDIAC SURGERY

THOR	ACO:	TOMY	SET	(Cont'd)

111011100101111	JET (COINT)	
Quantity	Catalog No.	Description
1 ea.	11-246	Richardson Retractor 1" (2.5 cm) x ¾" (1.9 cm)
1 ea.	11-248	Richardson Retractor 11/2" (3.8 cm) x 11/2" (3.8 cm)
1 ea.	11-266	Kelly Retractor 2" (5.1 cm) x 11/2" (3.8 cm)
1 ea.	11-268	Kelly Retractor 3" (7.6 cm) x 21/2" (6.4 cm)
1 ea.	11-314	Harrington Retractor 12" (30.5 cm), 21/2" (6.4 cm) wide
1 ea.	11-316	Harrington Retractor 12" (30.5 cm), 11/2" (3.8 cm) wide
2 ea.	11-549	Volkman Retractor 6 prongs, sharp
2 ea.	11-553	Volkman Retractor 4 prongs, blunt
1 ea.	18-1968	Freer Elevator, double end
2 ea.	24-572	DeBakey Tissue Forceps 8" (20.3 cm)
2 ea.	24-1602	DeBakey Multipurpose Clamp
1 ea.	25-2	Doyen Rib Raspatory, right
1 ea.	25-4	Doyen Rib Raspatory, left
1 ea.	25-18	Matson Rib Elevator and Stripper 8%" (22.2 cm)
1 ea.	25-22	Alexander-Farabeauf Periosteotome 81/4" (21 cm)
1 ea.	25-142	Burford-Finochietto Rib Spreader, 2 pr. blades
1 ea.	25-200	Bailey Rib Contractor
1 ea.	25-302	Allison Lung Retractor 1234" (32.4 cm)
1 ea.	25-456	Stille-Luer Bone Rongeur 8¾" (22.2 cm), cvd.
1 ea.	25-480	Sauerbruch Rib Rongeur 12" (30.5 cm), str. 17 mm sq. jaws
2 ea.	25-754	Duval Tissue Grasping Forcep 9" (22.9 cm), 1" (2.5 cm) wide
2 ea.	25-814	Mixter Thoracic Forceps 11" (27.9 cm)
= 00,		······································

CLOSED THORACOSTOMY SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-120TC	Mayo Dissecting Scissors 51/2" (14 cm), str. Carb-N-Sert blades
2 ea.	7-46	Crile Forceps 61/4" (15.9 cm), str.
6 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str., serrated
1 ea.	13-9	Patterson Trocar, Size 2, 7.5 mm
1 ea.	13-11	Patterson Trocar, Size 4, 12 mm



THORACIC AND CARDIAC SURGERY

OPEN THORACOSTOMY SET (add Thoracotomy Instruments)

Quantity	Catalog No.	Description
1 ea.	25-350	Bethune Rib Shears 13¾" (34.9 cm)
1 ea.	25-360	Gluck Rib Shears 81/2" (21.6 cm)

LOBECTOMY AND SEGMENTAL LUNG RESECTION SET (add Thoracotomy Instruments)

Quantity	Catalog No.	Description
1 ea.	5-192	Nelson Lobectomy Scissors 11" (27.9 cm), cvd.
6 ea.	7-142	Rochester-Pean Forceps 8" (20.3 cm)
2 ea.	7-202	Mixter Forceps 71/4" (18.4 cm)
4 ea.	16-10	Allis Tissue Forceps 71/4" (18.4 cm)
2 ea.	16-12	Allis Tissue Forceps 10" (25.4 cm)
2 ea.	16-54	Duval Lung Grasping Forceps 8" (20.3 cm), 1" (2.5 cm) wide
1 ea.	25-142	Burford-Finochietto Rib Spreader, 2 pr. blades
2 ea.	25-800	Sarot Intra-Thoracic Artery Forceps 9½" (24.1 cm)
2 ea.	25-814	Mixter Thoracic Forceps 11" (27.9 cm)
1 ea.	25-880	Sarot Bronchus Clamp 9" (22.9 cm), left
i ea.	25-882	Sarot Bronchus Clamp 9" (22.9 cm), right

BASIC OPEN HEART SET

Quantity	Catalog No.	Description
2 ea.	2-104SS	Yankauer Suction Tube
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-8	Knife Handle No. 4
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-124TC	Mayo Dissecting Scissors 6¾" (17.1 cm), str. Carb-N-Sert blades
1 ea.	5-126TC	Mayo Dissecting Scissors 6¾" (17.1 cm), cvd. Carb-N-Sert blades
1 ea.	5-180TC	Metzenbaum Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades
2 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-184TC	Metzenbaum Scissors 8" (20.3 cm), cvd., Carb-N-Sert blades
1 ea.	5-192TC	Metzenbaum Scissors 11" (27.9 cm), cvd., Carb-N-Sert blades
1 ea.	5-193A-TC	Metzenbaum Scissors 12" (30.5 cm), cvd., Carb-N-Sert blades
2 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
1 ea.	6-144	Russian Tissue Forceps 8" (20.3 cm)
1 ea.	6-168	Cushing Tissue Forceps 7" (17.8 cm), 1 x 2 T



THORACIC AND CARDIAC SURGERY

BASIC OPEN HEA	RT SET (Cont'd)	
Quantity	Catalog No.	Description
6 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
18 ea.	7-4	Halsted Mosquito Forceps 5" (12.7 cm), cvd.
6 ea.	7-44	Crile Forceps 51/2" (14 cm), cvd.
4 ea.	7-91	Coller. Forceps 61/4" (15.9 cm), cvd.
4 ea.	7-144	Rochester-Pean Forceps 9" (22.9 cm), cvd.
2 ea.	7-154	Rochester-Ochsner Forceps 8" (20.3 cm), str.
2 ea.	7-202	Mixter Forceps 71/4" (18.4 cm)
2 ea.	7-252	Mosquito Forceps 81/4" (21 cm), cvd.
1 ea.	7-258	Gemini-Mixter Forceps 9" (22.9 cm)
8 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
2 ea.	7-600	Foerster Sponge Forceps 91/2" (24.1 cm), str., serr.
1 ea.	8-46TC	Mayo-Hegar Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	8-51TC	DeBakey Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	8-53TC	DeBakey Needle Holder 9" (22.9 cm), Carb-N-Sert jaws
1 ea.	11-260	Richardson-Eastman Retractor, double end, Set of 2
4 ea.	16-12	Allis Tissue Forceps 10" (25.4 cm)
1 ea.	24-552	DeBakey Thoracic Tissue Forceps 7¾" (19.7 cm)
1 ea.	24-570	DeBakey Tissue Forceps 6" (15.2 cm)
1 ea.	24-574	DeBakey Tissue Forceps 91/2" (24.1 cm)
2 ea.	24-1000	Glover Patent Ductus Clamp, str.
1 ea.	24-1600	DeBakey Multipurpose Clamp 91/2" (24.1 cm), 4 cm jaws
1 ea.	24-1604	DeBakey Multipurpose Clamp 91/2" (24.1 cm), 6.5 cm jaws
1 ea.	24-2110	Satinsky-DeBakey Vascular Clamp 91/2" (24.1 cm), medium
4 ea.	V97-591	U.S. Pattern Tube Occluding Forceps 71/2" (19.1 cm)

MAJOR CARDIOVASCULAR SET (add Basic Laparotomy Set)

Quantity	Catalog No.	Description
1 ea.	4-10	Knife Handle No. 3L
1 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd. delicate, Carb-N-Sert blades
1 ea.	5-182B-TC	Metzenbaum Scissors 7" (17.8 cm), cvd. sharp points, Carb-N-Sert blades
1 ea.	5-246	Potts-Smith Scissors 71/2" (19.1 cm), 25° angle



2 ea.

Surgical Instrument Sets

THORACIC AND CARDIAC SURGERY

MAJOR CARDIOVASCU	ILAR SET (Cont'd)	
Quantity	Catalog No.	Description
1 ea.	5-247	Potts-Smith Scissors 71/2" (19.1 cm), 45° angle
1 ea.	5-248	Potts-Smith Scissors 71/2" (19.1 cm), 60° angle
1 ea.	5-249A	Potts-Smith Scissors 7" (17.8 cm), 25° angle, delicate
1 ea.	6-158TC	Potts-Smith Tissue Forceps 9¾" (24.8 cm), Carb-N-Sert jaws
12 ea.	7-32	Providence Hospital Forceps 5½" (14 cm), cvd.
1 ea.	8-53TC	DeBakey Needle Holder 9" (22.9 cm), Carb-N-Sert jaws
2 ea.	11-23	Little Retractor 7% (19.7 cm)
1 ea.	11-102	Cushing Vein Retractor 9" (22.9 cm)
1 ea.	11-136	Ribbon Retractor 11/2" (3.8 cm) x 13" (33 cm)
2 ea.	11-224A	Deaver Retractor 3" (7.6 cm) x 12" (30.5 cm), w/handle
1 ea.	11-266	Kelly Retractor 2" (5.1 cm) x 11/2" (3.8 cm)
1 ea.	11-268	Kelly Retractor 3" (7.6 cm) x 2½" (6.4 cm)
1 ea.	11-314	Harrington Retractor 12" (30.5 cm), 21/2" (6.4 cm) wide
1 ea.	11-316	Harrington Retractor 12" (30.5 cm), 11/2" (3.8 cm) wide
2 ea.	11-548	Volkman Retractor 4 prong, sharp
1 ea.	18-1968	Freer Elevator, double end
1 ea.	19-574	Frazier Suction Tube 9 Fr. (3 mm)
1 ea.	24-550	DeBakey Thoracic Tissue Forceps 6" (15.2 cm)
1 ea.	24-554	DeBakey Thoracic Tissue Forceps 91/2" (24.1 cm)
2 ea.	24-702	Glover Bulldog Clamp 8 cm, str.
2 ea.	24-752	Glover Bulldog Clamp 8 cm, cvd.
1 ea.	24-1300	DeBakey-Semb Ligature Carrier 101/4" (26 cm)
1 ea.	24-1408	DeBakey Vascular Clamp 121/4" (31.1 cm), curved, 10 cm
1 ea.	24-1600	DeBakey Multipurpose Clamp 91/2" (24.1 cm), 60° angle, 4 cm
1 ea.	24-1604	DeBakey Multipurpose Clamp 91/2" (24.1 cm), 60° angle, 6.5 cm
1 ea.	24-2020	DeBakey Tangential Occlusion Clamp 101/4" (26 cm)

AORTIC-CORONARY GRAFT SET (add Basic Open Heart Set)

25-784

Quantity	Catalog No.	Description	
1 ea.	5-246	Potts-Smith Scissors 71/2" (19.1 cm), 25° angle	
1 ea.	5-248	Potts-Smith Scissors 71/2" (19.1 cm), 60° angle	
6 ea.	7-10	Halsted Mosquito Forceps 5" (12.7 cm), cvd., delicate	

Cont'd on next page

Rumel Thoracic Forceps 9" (22.9 cm), full curve



Surgical Instrument Sets **THORACIC AND CARDIAC SURGERY**



AORTIC-CORONARY GRAFT SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	8-6TC	Webster Needle Holder 5" (12.7 cm), smooth, Carb-N-Sert jaws
1 ea.	11-74	Senn Retractor, sharp
1 ea.	11-76	Senn Retractor, blunt
1 ea.	11-548	Volkman Retractor, 4 prong, blunt
1 ea.	11-549	Volkman Retractor, 6 prong, blunt
1 ea.	11-608	Weitlaner Retractor 51/2" (14 cm), sharp
1 ea.	19-576	Frazier-Ferguson Suction Tube, 11 Fr. (3.6 mm)
1 ea.	24-600	DeBakey Bulldog Clamp, 7.5 cm, straight
1 ea.	24-604	DeBakey Bulldog Clamp, 10.5 cm, straight
1 ea.	24-652	DeBakey Bulldog Clamp, 8.5 cm, curved
1 ea.	24-656	DeBakey Buildog Clamp, 12 cm, curved
2 ea.	24-800	DeBakey Ring Handle Bulldog Clamp 5" (12.7 cm), str.



Surgical Instrument Sets **OPHTHALMIC SURGERY**

BASIC EYE SET

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-9	Knife Handle No. 7
1 ea.	4-401	Miniature Blade Handle w/improved self-locking chuck 3¾" (9.5 cm), (3K type)
1 ea.	5-300	Iris Scissors 4" (10.2 cm), str.
1 ea.	5-302	Iris Scissors 4" (10.2 cm), cvd.
4 ea.	7-8	Halsted Mosquito Forceps 5" (12.7 cm), str. delicate
6 ea.	7-10	Halsted Mosquito Forceps 5" (12.7 cm), cvd. delicate
4 ea.	7-504	Backhaus Towel Forceps 31/2" (8.9 cm)
1 ea.	9-110	O'Brien Stitch Scissors 3¾" (9.5 cm)
1 ea.	18-2	Lancaster Eye Speculum
1 ea.	18-12	Castroviejo Eye Speculum, medium size
1 ea.	18-112	Desmarres Lid Retractor, size 2
1 ea.	18-452	Von Graefe Strabismus Hook 51/2" (14 cm), medium
1 ea.	18 - 564	Castroviejo Spatula 51/2" (14 cm)
1 ea.	18-650	Castroviejo Caliper 31/4" (8.3 cm)
1 ea.	18-780	Eye Dressing Forceps 4" (10.2 cm), str.
1 ea.	18-782	Eye Dressing Forceps 4" (10.2 cm), cvd.
1 ea.	18-786	Iris Forceps 4" (10.2 cm), 1 x 2 T, str.
1 ea.	18-787	Iris Forceps 4" (10.2 cm), 1 x 2 T, cvd.
1 ea.	18-850	Graefe Fixation Forceps 4%" (11.2 cm)
1 ea.	18-951	Castroviejo Suturing Forceps 4" (10.2 cm), 1 x 2 T
1 ea.	18-1000	Arruga Capsule Forceps 4" (10.2 cm)
1 ea.	18-1202	Desmarres Chalazion Forceps 31/2" (8.9 cm), medium
1 ea.	18-1464	Stevens Tenotomy Scissors 41/8" (10.5 cm), cvd. blunt
1 ea.	18-1560	Castroviejo Corneal Scissors 4½" (11.4 cm), right
1 ea.	18-1562	Castroviejo Corneal Scissors 41/2" (11.4 cm), left
1 ea.	18-1618	Micro Iris Scissors 4" (10.2 cm), str.
1 ea.	18-1620	Micro Iris Scissors 4" (10.2 cm), cvd.
1 ea.	18-1652	Eye Suture (Gradle) Scissors 3¾" (9.5 cm), slightly cvd.
1 ea.	18-1828TC	Castroviejo Needle Holder 5½" (14 cm), str. w/lock, Carb-N-Sert jaws
2 ea.	18-2200	Fox Eye Shields, aluminum



OPHTHALMIC SURGERY

ERY SET (add Basic Eye Se	t)
Catalog No.	Description
18-33	Maumenee-Park Eye Speculum, w/lock screws
18-172	Troutman-Barraquer Blade Breaker
18-452	Von Graefe Strabismus Hook 51/2" (14 cm), medium size
18-486	Wilder Lens Loop 5¼" (13.3 cm), smooth
18-613	Kirby Hook and Expressor 6" (15.2 cm)
18-827	Troutman-Barraquer Colibri Corneal Forceps 0.4 mm, 3" (7.6 cm)
18-904	O'Brien Fixation Forceps 4" (10.2 cm), 1 x 2 T
18-952	Castroviejo Suturing Forceps 4" (10.2 cm), 0.6 mm teeth
18-1049	Castroviejo Capsule Forceps 4" (10.2 cm), new curve
18-1487	Westcott Tenotomy Scissors 51/4" (13.3 cm), right
18-1524	DeWecker Iris Scissors 41/4" (10.8 cm), sharp points
18-1558	Castroviejo Corneal Scissors 4" (10.2 cm), right, 6 mm blades
18-1559	Castroviejo Corneal Scissors 4" (10.2 cm), left, 6 mm blades
18-1832TC	Castroviejo Needle Holder 5½" (14 cm), cvd. w/lock, Carb-N-Sert jaws
	Catalog No. 18-33 18-172 18-452 18-486 18-613 18-827 18-904 18-952 18-1049 18-1487 18-1558 18-1559

EYE MUSCLE AND ENUCLEATION SET (add Basic Eye Set)

LIL IN COUNTY THE LITE CAN AUGUS CON CONTROL OF CONTROL		
Quantity	Catalog No.	Description
1 ea.	5-312	Strabismus Scissors 4" (10.2 cm), str.
1 ea.	5-314	Strabismus Scissors 4" (10.2 cm), cvd.
1 ea.	7-256	Gemini-Mixter Forceps 7" (17.8 cm), delicate
2 ea.	7-280	Serrefine 11/2" (3.8 cm), str.
1 ea.	18-62	McKinney Eye Speculum
1 ea.	18-904	O'Brien Fixation Forceps 4" (10.2 cm), 1 x 2 T
1 ea.	18-1252	Prince Advancement Forceps 31/2" (9.5 cm)
1 ea.	18-1498	Enucleation Scissors 5" (12.7 cm), full curve

RETINAL SURGERY SET (add Basic Eye Set)

Quantity	Catalog No.	Description
2 ea.	7-280	Serrefine 11/2" (3.8 cm), str.
1 ea.	18-446	Jameson Strabismus Hook 5" (12.7 cm)
1 ea.	18-452	Von Graefe Strabismus Hook 5½" (14 cm), medium
1 ea.	18-910	Gifford Fixation Forceps 4" (10.2 cm), 1 x 2 T
1 ea.	18-966	McCullough Suturing Forceps 4" (10.2 cm), 1 x 2 T
1 ea.	18-1487	Westcott Tenotomy Scissors 51/4" (13.3 cm), right



OPHTHALMIC SURGERY

DACRYOCYSTORHINOSTOMY SET (add Basic Eye Set)

Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-180A-TC	Metzenbaum Scissors, 5½" (14 cm), delicate, Carb-N-Sert blades
1 ea.	9-124	Wire Cutting Scissors 49/4" (12.1 cm), angular
1 ea.	18-34	Wire Speculum
1 ea.	18-88	Goldstein Lachrimal Sac Retractor
1 ea.	18-412	Tyrell Hook 5" (12.7 cm), sharp
1 ea.	18-660	Stainless Ruler 6" (15.2 cm)
1 ea.	18-691	Castroviejo Lacrimal Dilator 51/2" (14 cm), double end
1 ea.	18-740	Bowman Probe size 4/0 and 3/0
1 ea.	18-742	Bowman Probe size 2/0 and 0
1 ea.	18-744	Bowman Probe size 1 and 2
1 ea.	18-746	Bowman Probe size 3 and 4
1 ea.	18-784	Eye Dressing Forceps 4" (10.2 cm), full cvd.
1 ea.	18-1960	Freer Chisel 61/2" (16.5 cm), cvd.
1 ea.	18-1962	West Chisel 61/2" (16.5 cm), str.
1 ea.	18-1968	Freer Elevator 7" (17.8 cm), double end
1 ea.	18-1990	Kerrison Rongeur, 3.5 mm bite, 3" (7.6 cm) shaft
1 ea.	19-192	Lathbury Applicator 51/2" (14 cm)
1 ea.	19-370	Lucae Bayonet Forceps 5½" (14 cm)
1 ea.	19-572	Frazier Ferguson Suction Tube 7 Fr. (2.3 mm)
1 ea.	19-580	Baron Suction Tube 3 Fr. (1 mm)
1 ea.	19-793	Miltex Mallet 71/2" (19.1 cm), stainless and nylon faces
1 ea.	19-1232	Lempert Rongeur 61/4" (15.9 cm), cvd.
1 ea.	20-24	Beckman Nasal Speculum 6" (15.2 cm)
1 ea.	20-36-50	Cottle Septum Speculum, 50 mm blades
1 ea.	20-570	Takahashi Nasal Forceps, 2.5 x 10 mm bite, 41/2" (11.4 cm) shaft

EAR, NOSE AND THROAT SURGERY

BASIC EAR SET		
Quantity	Catalog No.	Description
1 ea.	4-7	Knife Handle No. 3
1 ea.	5-122TC	Mayo Dissecting Scissors 51/2" (14 cm), cvd. Carb-N-Sert blades
1 ea.	5-180A-TC	Metzenbaum Scissors 51/2" (14 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-314	Strabismus Scissors 4" (10.2 cm), cvd.
1 ea.	6-118	Adson Dressing Forceps 4¾" (12.1 cm), serr.
1 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
6 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
2 ea.	7-36	Kelly Forceps 5½" (14 cm), str.
4 ea.	7-504	Backhaus Towel Clamp 31/2" (8.9 cm)
1 ea.	8-8A-TC	Halsey Needle Holder 5" (12.7 cm), serrated, Carb-N-Sert jaws
1 ea.	9-124	Wire Cutting Scissors 4¾" (12.1 cm), angular
1 ea.	11-2	Jansen Retractor 41/4" (10.8 cm), blunt
1 ea.	11-6	Allport Retractor 4" (10.2 cm), sharp
2 ea.	11-74	Senn Retractor 6%" (16.2 cm), sharp
1 ea.	11-160	Langenbeck Retractor %" (.95 cm) x 11/4" (3.2 cm)
1 ea.	11-164	Langenbeck Retractor 1/2" (1.3 cm) x 11/2" (3.8 cm)
2 ea.	16-6	Allis Tissue Forceps 6" (15.2 cm)
1 ea.	18-850	Graefe Tissue Forceps 4%" (11.2 cm)
1 ea.	18-1968	Freer Elevator 7" (17.8 cm), double end
1 ea.	19-50-1	Farrier Ear Speculum size 1
1 ea.	19-50-3	Farrier Ear Speculum size 3
1 ea.	19 - 50-5	Farrier Ear Speculum size 5
1 ea.	19-50-7	Farrier Ear Speculum size 7
1 ea.	19-260	Buck Ear Curette, size 00, str. sharp
1 ea.	19-264	Buck Ear Curette, size 1, str. sharp
1 ea.	19-371	Jansen Ear Forceps 61/4" (15.9 cm)
1 ea.	19-572	Frazier Ferguson Suction Tube 7 Fr. (2.3 mm)
1 ea.	19-574	Frazier Ferguson Suction Tube 9 Fr. (3 mm)
1 ea.	19-660	Wullstein Endaural Retractor
2 ea.	21-153	Joseph Single Hook 61/4" (15.9 cm)





EAR, NOSE AND THROAT SURGERY

MYRINGO [*]	TOMY	SET
----------------------	------	-----

Quantity	Catalog No.	Description
1 ea.	5-270	Plastic Surgery Scissors 4¾" (12.1 cm), str. sharp points
1 ea.	19-50-1	Farrier Ear Speculum, size 1
1 ea.	19-50-3	Farrier Ear Speculum, size 3
1 ea.	19-50-5	Farrier Ear Speculum, size 5
1 ea.	19-50-7	Farrier Ear Speculum, size 7
1 ea.	19-192	Lathbury Applicator 51/2" (14 cm)
1 ea.	19-270	Buck Ear Curette, size 00 angled, sharp
1 ea.	19-274	Buck Ear Curette, size 1 angled, sharp
1 ea.	19-371	Jansen Ear Forceps 61/4" (15.9 cm)
1 ea.	19-500	Sexton Ear Knife 7" (17.8 cm), str.
1 ea.	19-572	Frazier Ferguson Suction Tube 7 Fr. (2.3 mm)
1 ea.	19-582	Baron Suction Tube 5 Fr. (1.7 mm)

BASIC NASAL SET

Quantity	Catalog No.	Description
1 ea.	2-104SS	Yankauer Suction Tube
1 ea.	4-7	Knife Handle No. 3
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-122TC	Mayo Dissecting Scissors 51/2" (14 cm), cvd., Carb-N-Sert blades
1 ea.	5-180A-TC	Metzenbaum Scissors 5½" (14 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd., delicate, Carb-N-Sert blades
1 ea.	5-276	Plastic Surgery Scissors 4¾" (12.1 cm), cvd., S/S
2 ea.	6-120	Adson Tissue Forceps 4¾" (12.1 cm), 1 x 2 T
6 ea.	7-2	Halsted Mosquito Forceps 5" (12.7 cm), str.
4 ea.	7-504	Backhaus Towel Clamp 3¾" (8.9 cm)
1 ea.	8-8A-TC	Halsey Needle Holder 5" (12.7 cm), serr., Carb-N-Sert jaws
2 ea.	11-74	Senn Retractor, sharp
1 ea.	11-122	U.S. Army Retractor, Set of 2
2 ea.	16-4	Baby Allis Tissue Forceps 5" (12.7 cm)
1 ea.	18-850	Graefe Forceps 4%" (11.2 cm)
1 ea.	19-202	Farrell Applicator 61/2" (16.5 cm)
2 ea.	19-572	Frazier Suction Tube, 7 Fr. (2.3 mm)

Cont'd on next page



EAR, NOSE AND THROAT SURGERY

BASIC NASAL SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	20-3	Vienna Nasal Speculum, medium
1 ea.	20-30	Killian Septum Speculum, 3" (7.6 cm) blades
1 ea.	20-36-50	Cottle Septum Speculum, 50 mm delicate blades
1 ea.	20-120	Wilde Nasal Dressing Forceps 61/2" (16.5 cm)
1 ea.	21-153	Joseph Single Hook
1 ea.	21-154	Joseph Double Hook, 2 mm
1 ea.	21-158	Joseph Double Hook, 7 mm
1 ea.	21-188	Aufricht Nasal Retractor-Speculum 61/2" (16.5 cm)
1 ea.	21-630	Joseph Nasal Scissors 5¾" (14.6 cm), str.
1 ea.	21-632	Joseph Nasal Scissors 5¾" (14.6 cm), cvd.

MINOR NASAL SET

Quantity

MINOTI NACAL CE	••	
Quantity	Catalog No.	Description
1 ea.	10-220	Anthony Suction Tube
1 ea.	20-3	Vienna Nasal Speculum, medium
1 ea.	20-30	Killian Septum Speculum, 3" (7.6 cm) blades
1 ea.	20-36-50	Cottle Septum Speculum, 50 mm delicate blades
1 ea.	20-140	Noyes Alligator Forceps 5½" (14 cm), serr.
1 ea.	20-572	Takahashi Ethmoid Forceps 41/2" (11.4 cm), 2.5 mm bite

MAJOR NASAL SET (add Basic Nasal Set)

Catalog No.

Quantity	Catalog 140.	Description
1 ea.	5-180A-TC	Metzenbaum Scissors 51/2" (14 cm), cvd. B/B, Carb-N-Sert blades
1 ea.	5-182A-TC	Metzenbaum Scissors 7" (17.8 cm), cvd. B/B, Carb-N-Sert blades
1 ea.	5-300	Iris Scissors 4" (10.2 cm), str.
1 ea.	5-312	Strabismus Scissors 4" (10.2 cm), str.
1 ea.	5-314	Strabismus Scissors 4" (10.2 cm), cvd.
1 ea.	6-118	Adson Dressing Dorceps 4¾" (12.1 cm), serr.
1 ea.	6-197	Jansen (Gruenwald) Bayonet Dressing Forceps 61/4" (15.9 cm)
1 ea.	8-51TC	Crile-Wood Needle Holder 7" (17.8 cm), Carb-N-Sert jaws
1 ea.	18-110	Des Marres Lid Retractor
1 ea.	18-650	Castroviejo Caliper 31/4" (8.3 cm)
1 ea.	18-660	Stainless Steel Ruler
1 ea.	18-1968	Freer Elevator 7" (17.8 cm), double end
		Cont'd on next page

Description



EAR, NOSE AND THROAT SURGERY

MAJOR NASAL SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	19-172	Uebe Applicator 7" (17.8 cm), cross serrated
1 ea.	19-794	Miltex Mallet 71/2" (19.1 cm)
1 ea.	20-30	Killian Septum Speculum 6" (15.2 cm)
2 ea.	20-36-50	Cottle Septum Speculum 6" (15.2 cm)
1 ea.	20-124	Hartman Nasal Dressing Forceps 7" (17.8 cm)
1 ea.	20-218	Ballenger Swivel Knife, str., 3 mm
1 ea.	20-224	Ballenger Swivel Knife 73/4" (19.7 cm), bayonet
1 ea.	21-24	Converse Bistouri 6" (15.2 cm), cvd.
1 ea.	21-28	Fomon Knife 51/2" (14 cm), double edge, cvd.
1 ea.	21-32	Cottle Nasal Knife 51/2" (14 cm)
1 ea.	21-64	Joseph Elevator 7" (17.8 cm), cvd.
1 ea.	21-68	Fomon Periosteal Elevator 61/4" (15.9 cm)
1 ea.	21-89	Gillies Skin Hook 7" (17.8 cm)
1 ea.	21-144	Cottle Double Hook 51/2" (14 cm), sharp
1 ea.	21-148	Fomon Retractor 61/2" (16.5 cm), double ball end prongs
1 ea.	21-150	Cottle Retractor 51/2" (14 cm), 4 prong
1 ea.	21-153	Joseph Single Hook 61/4" (15.9 cm)
2 ea.	21-156	Joseph Double Hook 61/4" (15.9 cm), sharp, 5 mm
1 ea.	21-163	Converse Alar Retractor 4" (10.2 cm)
1 ea.	21-165	Cottle Knife Guide and Retractor 8" (20.3 cm)
1 ea.	21-176	Neivert Knife Guide and Retractor 61/2" (16.5 cm), double end
1 ea.	21-205-4	Sheehan Osteotome 61/4" (15.9 cm), 4 mm
1 ea.	21-205-8	Sheehan Osteotome 61/4" (15.9 cm), 8 mm
1 ea.	21-205-10	Sheehan Osteotome 61/4" (15.9 cm), 10 mm
1 ea.	21-205-12	Sheehan Osteotome 61/4" (15.9 cm), 12 mm
1 ea.	21-212	Cottle Septum Chisel 71/4" (18.4 cm), cvd.
1 ea.	21-215	Cinelli Osteotome 6½" (16.5 cm)
1 ea.	21-232	Neivert Osteotome 8" (20.3 cm)
1 ea.	21-280	Joseph Nasal Saw 71/2" (19.1 cm), str.
1 ea.	21-330	Joseph Rasp 61/4" (15.9 cm)
1 ea.	21-342	Maltz Rasp 7" (17.8 cm)
1 ea.	21-343	Maltz Rasp D/E 81/4" (21 cm)
1 ea.	21-348	Aufricht Glabella Rasp 81/4" (21 cm)

Cont'd on next page

EAR, NOSE AND THROAT SURGERY



MAJOR NASAL SET (Cont'd)

Quantity	Catalog No.	Description
1 ea.	21-440	Cottle Columella Forceps 41/2" (11.4 cm)
1 ea.	21-602	Fomon Lower Lateral Scissors 5" (12.7 cm), cvd.
1 ea.	21-604TC	Fomon Dorsal Scissors 51/2" (14 cm), angular, Carb-N-Sert blades
1 ea.	21-610	Cottle Dorsal Scissors 61/2" (16.5 cm)
1 ea.	21-660	Ash Septum Straightner 9" (22.9 cm)
2 ea.	22-662	Schnidt Hemostatic Forceps 7½" (19.1 cm)

TONSILLECTOMY AND ADENOIDECTOMY SET

Quantity	Catalog No.	Description
1 ea.	2-104SS	Yankauer Suction Tube
1 ea.	2-180	McIvor Mouth Gag w/3 blades
1 ea.	4-9	Knife Handle No. 7
1 ea.	5-182TC	Metzenbaum Scissors 7" (17.8 cm), cvd., Carb-N-Sert blades
1 ea.	7-190	Miculicz Tonsil Forceps 71/2" (19.1 cm)
4 ea.	7-506	Backhaus Towel Clamp 51/4" (13.3 cm)
1 ea.	7-604	Ballenger Sponge Forceps 7" (17.8 cm), str., serrated
1 ea.	8-50TC	Crile-Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws
1 ea.	8-52TC	Crile-Wood Needle Holder 8" (20.3 cm), Carb-N-Sert jaws
1 ea.	20-800	Barnhill Adenoid Curette 81/2" (21.6 cm), size 0
1 ea.	20-802	Barnhill Adenoid Curette 81/2" (21.6 cm), size 1
1 ea.	20-804	Barnhill Adenoid Curette 81/2" (21.6 cm), size 2
1 ea.	20-836	Love Uvula Retractor, 18 mm
1 ea.	20-850	LaForce Adenotome 91/2" (24.1 cm), extra small, size 0
1 ea.	20-854	LaForce Adenotome 91/2" (24.1 cm), medium, size 2
1 ea.	20-926	Meltzer Adenoid Punch, size 1, triangular jaws
1 ea.	22-4	Wieder Tongue Depressor, large
1 ea.	22-550	White Tonsil Seizing Forceps 71/2" (19.1 cm), cvd.
6 ea.	22-662	Schnidt Tonsil Forceps 71/2" (19.1 cm)
1 ea.	22-666	Sawtell Tonsil Forceps 71/2" (19.1 cm)
1 ea.	22-808	Fischer Tonsil Knife and Dissector 8½" (21.6 cm)
1 ea.	22-854	Hurd Tonsil Dissector and Pillar Retractor 9" (22.9 cm)
1 ea.	22-1010	Beck-Schenck Tonsil Snare with 1 straight and 3 vedder tips
1 ea.	22-1030	Sage Automatic Ratchet Tonsil Snare



SUTURE

E.K. 50	IUKE SEI
Quantity	y

Quantity	Catalog No.	
1 ea.	4-7	
1 ea.	5-290	
1 ea.	5-16	
1 ea.	6-118	
1 ea.	6-120	
2 ea.	7-2	
2 ea.	7-4	
1 ea.	7-36	
1 ea.	7-38	
2 ea.	7-504	
1 ea.	8-50TC	

Description

Knife Handle No. 3

Ragnell Dissecting Scissors 5" (12.7 cm) Operating Scissors 51/2" (14 cm), str. S/B Adson Dressing Forceps 43/4" (12.1 cm), serr. Adson Tissue Forceps 43/4" (12.1 cm), 1 x 2 T Halsted Mosquito Forceps 5" (12.7 cm), str. Halsted Mosquito Forceps 5" (12.7 cm), cvd. Kelly Forceps 51/2" (14 cm), str. Kelly Forceps 51/2" (14 cm), cvd.

Backhaus Towel Clamp 31/2" (8.9 cm)

Crile-Wood Needle Holder 6" (15.2 cm), Carb-N-Sert jaws

SMALL SUTURE SET (For Doctor's Office or E.R.)

	/
Quantity	Catalog No.
1 ea.	4-7
1 ea.	5-312
1 ea.	5-314
2 ea.	6-120
4 ea.	7-4
2 ea.	7-26
2 ea.	8-6TC

Description

Knife Handle No. 3

Strabismus Scissors 41/4" (10.8 cm), str. Strabismus Scissors 41/4" (10.8 cm), cvd. Adson Tissue Forceps 4%" (12.1 cm), 1 x 2 T Halstead Mosquito Forceps 5" (12.7 cm), cvd.

Hartman Mosquito Forceps 31/2" (8.9 cm), cvd.

Webster Needle Holder 5" (12.7 cm), smooth Carb-N-Sert jaws

Index Index



Abbey Needle Holders113, 341
Abdominal Scissors23
Adair Tissue Tenaculum208, 501
Adenoid Curettes304
Adenoid Punches306
Adenotomes305
Adson Aneurysm Needle143,425
Adson Bayonet
Forceps73, 89, 275, 292, 429
Adson Cannula425
Adson Dressing Forceps 68, 69, 88, 220,
322, 427, 537, 538
Adson Dura Hooks192, 431
Adson Ganglion Scissors423
Adson Hemostatic Forceps 103, 430
Adson Needle Holders116
Adson Periosteal Elevators434
Adson Retractors 187, 412, 438
Adson Rongeurs414
Adson Suction Tubes149, 426
Adson Tissue Forceps 68, 70, 88, 220,
323, 427, 537, 538
Aebli Corneal Scissors48, 261
Agricola Lacrimal Sac Retractor 153, 228
Alexander Chisels284, 449
Alexander Gouges
Alexander-Farabeuf Periosteotomes 388
Allen Intestinal Clamp104, 210
Allen-Schiotz Tonometers
Allen-Schiotz Tonometers230
Allen-Schiotz Tonometers230 Allis Tissue Forceps206, 400, 477
Allen-Schiotz Tonometers230 Allis Tissue Forceps206, 400, 477 Allis-Adair Tissue Forceps206, 476
Allen-Schiotz Tonometers
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273
Allen-Schiotz Tonometers
Allen-Schiotz Tonometers
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixalion Forceps 245, 541
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Alport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Armiotome 518
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Armniotome 518 Amputating Knives 465
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixalion Forceps 245, 541 Alvis Foreign Body Spud 233 Armiotome 518 Amputating Knives 465 Anastomosis Clamps 210, 378, 384, 385
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Armiotome 518 Amputating Krives 465 Anastomosis Clamps 210, 378, 384, 385 Anderson Osteotomes 315
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Alport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Amniotome 518 Amputating Knives 465 Anastomosis Clamps 210, 378, 384, 385 Anderson Osteotomes 315 Andrews Tongue Depressor 347
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Amniotome 518 Amputating Knives 465 Anastomosis Clamps 210, 378, 384, 385 Anderson Stongue Depressor 347 Andrews-Pynchon Suction 347
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis-Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Amniotome 518 Amputating Knives 465 Anastomosis Clamps 210, 378, 384, 385 Anderson Osteotomes 315 Andrews Tongue Depressor 347 Andrews-Pynchon Suction Tubes 6, 147, 351, 357
Allen-Schiotz Tonometers 230 Allis Tissue Forceps 206, 400, 477 Allis Adair Tissue Forceps 206, 476 Allison Lung Retractor 171, 394 Allport Mastoid Searchér 190, 273 Allport Retractors 151, 280 Alm Retractors 152, 281, 545 Alvin Fixation Forceps 245, 541 Alvis Foreign Body Spud 233 Amniotome 518 Amputating Knives 465 Anastomosis Clamps 210, 378, 384, 385 Anderson Stongue Depressor 347 Andrews-Pynchon Suction 347

Angiotribes1	07
Anoscopes4	71
Anthony Suction Tube 147, 149, 279, 2	99
Antrum Cannula298, 2	
Antrum Curettes2	97
Antrum Needles2	
Antrum Rasps2	
Antrum Rongeurs3	
Applicators, Ear2	73
Applicators, Laryngeal3	
Applicators, Nasal	
Approximator Set For Clips2	
Arruga Capsule Forceps2	
Arruga Lacrimal Trephines2	
Ash Septum Straightening Forceps'3	
Atraumatic Instruments	
Attic Cannula2	
Aufricht Nasal Retractors157, 3	
Aufricht Rasps3	
Auvard Vaginal Specula4	
Aver Chalazion Forceps2	Tr.
Ayer Gharazion Porceps2	06
В	
Babcock Suture Wire1	ac
Babcock Tissue Forceps207, 4	
Baby Crile Forceps92, 3	
Baby Ochsner Forceps93, 3	
Baby Pean Forceps94, 3	
Backhaus Towel Clamps1	
Bacon Rongeur414, 4	
Badgley Retractors186, 4	31
Badgley Rib Spreader3	
	91
Bailey Gigli Saw Guide413, 4	91 64
Bailey Gigli Saw Guide413, 4 Bailey Rib Contractors3	91 64 92
Bailey Gigli Saw Guide413, 4 Bailey Rib Contractors	91 64 92 92
Bailey Gigli Saw Guide	91 64 92 92 22
Bailey Gigli Saw Guide	91 64 92 92 22
Bailey Gigli Saw Guide	91 64 92 92 22 99
Bailey Gigli Saw Guide	91 64 92 92 22 99 79
Bailey Gigli Saw Guide	91 64 92 92 22 99 79
Bailey Gigli Saw Guide 413, 4 Bailey Rib Contractors 3 Bailey-Gibbon Rib Contractor 3 Bailey-Williamson Obstetrical Forceps. 5 Bainbridge Forceps 103, 197, 2 Bakes Dilators 1 Ballour Retractors 177-1 Ballantine Hysterectomy Forceps 96, 5 Ballenger Chisels 286, 3	91 64 92 92 22 99 79 13
Bailey Gigli Saw Guide 413. 4 Bailey Rib Contractors 3 Bailey Gibbon Rib Contractor 3 Bailey-Williamson Obstetrical Forceps. 5 Bainbridge Forceps 103, 197. 2 Bakes Dilators 1 Ballour Retractors 177-1 Ballantine Hysterectomy Forceps 96, 5 Ballenger Chisels 296, 3 Ballenger Sponge Forceps 3	91 64 92 92 22 09 79 13
Bailey Gigli Saw Guide	91 64 92 92 92 99 79 13 15
Bailey Gigli Saw Guide	91 64 92 92 92 99 79 15 15 56
Bailey Gigli Saw Guide	91 64 92 92 92 79 79 15 56 66
Bailey Gigli Saw Guide	91 64 92 92 92 79 79 15 56 66

Barkan Goniotomy Knife	232
Barnhill Adenoid Curettes	304
Baron Suction Tubes 149, 279,	
Barr Anal Retractor174,	469
Barr Crypt Hook	194
Barr Probes145, 146,	474
Barr-Shuford Rectal Speculum	473
Barraquer Blade Breaker	231
Barraquer Cilia Forceps80,	251
Barraquer Eye Specula	227
Barraquer Iris Scissors32,	259
Barraquer Needle	
Holders 121, 265, 346,	536
Barrett Uterine Tenaculum Forceps	.501
Baumgartner Needle	
Holders114, 124, 342,	535
Beacham Amniotome	
Beaupre Cilia Forceps80,	251
Beck-Schenck Tonsil Snares	360
Becker Cystotome	.232
Beckman Nasal Speculum	.289
Beckman Retractors187, 196,	
Beckman-Eaton Retractor173,	
Beckman-Weitlaner	
Retractors 184, 195, 282, 412,	438
Beer Cilia Forceps79,	
Beeson Cast Spreader	
Bellucci Micro Ear Scissors 49, 217,	288
Bennett Cilia Forceps80	25
Bennett Bone Elevator	.446
Bennett Tibia Retractors	. 182
Benson Pylorus Separator	.213
Bergh Cilia Forceps80,	25
Berke Ptosis Forceps252	325
Best Colon Clamps	.214
Bethune Rib Shears	.394
Beyer Rongeur286, 353, 415,	46
Bihrle Needle Holder Clamp128,	488
Bill Obstetrical Traction Handle	.519
Billeau Ear Loops	.274
Billroth Tumor Forceps	.508
Binocular Loupe	4
Biopsy Forceps, Rectal479	
Biopsy Forceps, Uterine 486, 509	
Biopsy Needles	
Biopsy Punches, Cervical510	
Biopsy Punches, Disposable	. 526
•	

Index E





Biro Dermal Naevus Scissors . 38, 332, 531
Bishop Mastoid Gouges284
Bishop-Harmon Dressing
Forceps76, 244, 322, 541
Bishop-Harmon Tissue
Forceps76, 221, 244, 322, 541
Bistouri, Tracheal369
Bladder Retractors486
Blade Breakers231
Blade Remover15
Biade-Safe Blade Remover15
Blades, Surgical15-17
Blair Cleft Palate Elevators351
Blake Dressing Forceps81
Blake Ear Forceps275
Blake Gall Stone Forceps205
Blakesley Ethmoid Forceps302
Blepharoplasty Scissors59, 335, 531
Blount Knee Retractor 183, 446
Blumenthal Rongeurs352, 459
Bodenhammer Rectal Speculum 469
Boettcher Scissors28, 359
Bonaccolto Utility Forceps79, 248
Bone Chisels449, 450
Bone Clamps441, 442
Bone Curettes318, 436, 456, 457
Bone Cutting
Forceps338, 395, 462, 463, 553
Bone Elevators446
Bone File454
Bone Gouges449, 451, 452
Bone Holding Forceps440, 441
Bone Hooks188, 447
Bone Mallets453
Bone Rasps 454
Bone Retractors446-448
Bone Rongeurs459-462
Bone Saws465
Bone Skid454
Bone Wire Guide143, 442
Bonn Iris Scissors31, 254, 330, 528
Bonn Suturing Forceps 75, 76, 220, 244
Bonney Tissue Forceps70
Bosworth Tongue Depressor347
Botvin iris Forceps
Boucheron Ear Specula269
Bower Sphenoid Cannula299

Bowman Lacrimal Probes		
Bowman Needle		.23
Boynton Needle Holder 114, 265	, 342	53
Boys-Allis Tissue Forceps		.20
Bozeman Sponge Forceps		
Braasch Biopsy Forceps		
Braces, Neurosurgical		
Bracken Iris Forceps		
Braun Episiotomy Scissors	23	514
Braun Tenaculum Forceps	···· ∠u,	EO(
Breast Tenaculum Forceps	•••••	. 500
Bridge Forceps		400
Brinkerhoff Rectal Specula	25,	400
Bronchial and Esophageal		.4/2
Instruments		
instruments	365,	366
Bronchoscopic Forceps		.365
Bronchus Clamps		
Brophy Scissors35		
Brown Applicators		
Brown Needle Holders		
Brown Tissue Forceps 71, 324	, 399,	429
Brown Tonsil Snare		361
Brown-Adson Tissue		
Forceps 68, 70, 71, 88, 324	, 428,	538
Bruening Ear Snare		277
Bruening Laryngeal Tip		364
Bruening Nasal Snare		294
Bruening Otoscope		
Bruening Septum Forceps		
Brun Bone Curettes		
Brun Mastoid Curettes		
Brunner Intestinal Forceps		
Bruns Cast Shears		
Brush For Rasps and Files		
Brushes, Instrument Cleaning	14	210
Brushes, Trachea Tubes		
Buck Curettes		
Buck Ear Probes		
Buck Neuro Hammers		
Buie Anal Retractors		
Buie Biopsy Forceps		
Bule Fistula Probe		
Bule Pile Clamp		
Buie Suction Tube		
Bulldog Clamp Applying Forceps		
Bulldog Clamps 108, 224, 267,		
Julidog Scissors	37,	337

Bumm Uterine Curettes506
Bunge Evisceration Spoons236
Bunnell Hand Drill445
Bunnell Tendon Strippers458
Bunt Forceps Holder8
Duff rolleps notices
Burford-Finochietto Rib Spreaders 390
Burke Cervical Biopsy Punch Forceps .510
Burs, Neurosurgical409, 410
Busch Scissors23, 523
Busch Umbilical Scissors523
Button End Knives307
Buxton Uterine Forceps505
С
C-N-S Pin & Wire Pullers444
C-N-S Wire Cutting Scissors59
CTS Sterilization Trays10
Calipers, E.k.G3
Calipers, Eye239
Calipers, Plastic Surgery317
Canfield Tonsil Knife358
Cannula Occlusion Clamp12
Cannula, Antrum298, 299
Cannula, Attic279
Cannula, Laryngeal364
Cannula, Sinus298
Cannula, Sphenoid299
Caplan Nasal Bone Scissors 38, 337
Capsule Forceps249, 250, 487
Capsulotomy Scissors262
Carb-N-Sert Dressing Forceps88, 89
Carb-N-Sert Needle
Holders 122-130, 222, 266
Carb-N-Sert Pin Cutter443, 444
Carb-N-Sert Pin Puller142, 444
Carb-N-Sert Scissors54-59
Carb-N-Sert Tissue Forceps88, 89
Cardiovascular Clamps 378-381, 384-387
Cardiovascular Scissors26, 405
Carmalt Splinter Forceps83
Cartilage Clamps458
Cartilage Knives457
Cartilage Scissors25, 458
Cases, Sterilizing230
Cast & Utility Shears53, 466, 467
Cast Instruments53, 466-468
Cast Spreaders 468

Index transfer and the second second



Casianares Facellit Scissors336
Castroviejo Blade Breaker231
Castroviejo Caliper239, 317
Castroviejo Capsule Forceps249
Castroviejo Colibri Corneal
Forceps77, 245
Castroviejo Corneal Scissors43, 47, 216, 260, 261
Scissors 43, 47, 216, 260, 261
Castroviejo Corneal Section
Scissors47, 260
Castroviejo Cyclodialysis Spatula237
Castroviejo Eye Specula225
Castroviejo Fixation Forceps 245, 541
Castroviejo Keratoplasty Scissors 261
Castroviejo Lacrimal Dilators240
Castroviejo Needle Holders 120, 128,
129, 264, 266, 345, 346, 536
Castroviejo Spatula237
Castroviejo Suturing Forceps 247, 541
Cataract Knives231
Catheter Guide485
Catheter Introducing Forceps5
Catheters, Eustachian274
Cathelers, Urethral484
Cautery 239
Cerebellar Retractors412
Cervix Holding Forceps499
Chalazion Curettes236
Chalazion Forceps252, 253, 325
Chandler Elevators455
Chandler Iris Forceps75, 243
Charriere Bone Saws465
Cheatle Sterilizer Forceps
Cheek Retractors162, 347, 348
Chelsea-Eaton Anal Specula473
Cheron Uterine Dressing Forceps 498
Child Intestinal Forceps
Childrens Hospital Spatula425
Chisels, Bone449, 450
Chisels, Lacrimal Sac268
Chisels, Mastoid284
Chisels, Plastic Surgery313-315
Chisels, Septum296, 313
Cicherelli Rongeur459
Cilia Forceps
Cinelli Osteotomes314
Clamps, Anastomosis 210, 384, 385
Clamps, Bone441, 442

Clamps, Bronchus404
Clamps, Bulldog108, 374-376, 406
Clamps, Cannula Occlusion12
Clamps, Cardiovascular . 378-381, 384-387
Clamps, Cartilage458
Clamps, Coarctation377, 404
Clamps, Colon214
Clamps, Cystic Duct Catheter204
Clamps, Intestinai210, 212-214
Clamps, Kidney488
Clamps, Multipurpose 381, 382, 385
Clamps, Muscle197
Clamps, Occlusion12, 383
Clamps, Patent Ductus377, 378
Clamps, Peripheral Vascular378
Clamps, Peritoneal515
Clamps, Pile477
Clamps, Pylorus213
Clamps, Thoracic404
Clamps, Towel109
Clamps, Umbilical
Clamps, Vascular377-381, 384-387
Clamps, Vascular, Pediatric386, 387
Clayman Intraocular Lens Forceps 250
Cleft Palate Elevators351
Clip Applying Forceps 136, 137, 424
Clip Removing Forceps136, 137
Clips, Vessel223, 224, 430
Clips, Wound135
Cloward Pituitary Rongeurs422
Coakley Antrum Curetles
Coakley Antrum Trocar298
Coarctation Clamps
Colibri Forceps
Coller Forceps
Collier Needle Holder112, 340
Collin Dynamometer
Collin Intestinal Forceps
Collin Retractors
Collin Tongue Seizing Forceps5, 350
Collin Uterine Elevating Forceps505
Collin Vaginal Specula491
Collin (Judd Demartel) Gall Bladder
Forceps204
Collin-Duval Tissue Forceps209, 401
Collyer Pelvimeter517
Colon Clamps214

Columbia Lip Retractor		
Colver Tonsil Retractor	163,	358
Colver Tonsil Seizing Forceps	******	355
Combo Neuro-Diagnostic Sets	*>+>+	1
Comedone Extractors	547,	548
Common Duct Dilators		199
Common Duct Scoops		200
Conjuctival Scissors	32,	257
Converse Alar Retractor		
Converse Bistouri	******	307
Converse Blade Retractors	157,	313
Converse Nasal Scissors		332
Converse Needle Holders	113,	341
Converse Skin Hooks190,	309,	542
Cook Eye Specula		227
Cook Rectal Retractor	174,	470
Cooley Needle Holders		
Cooley Pediatric Clamps 384,		
Cooley Tissue Forceps		
Cooley-Derra Clamp		
Cooley-Satinsky Clamp		
Cornea-Scieral Forceps		.238
Corneal Dissector		
Corneal Forceps243	-245,	249
Corneal Hooks		234
Corneal Knives		
Corneal Scissors43, 47, 48,	260,	261
Corneal Section Scissors . 46, 47,	258,	260
Corneal Spuds		
Cottle Bulldog Scissors		
Cottle Cartilage Crusher		.317
Cottle Chisel		
Cottle Dorsal Scissors		.337
Cottle Knife	******	.307
Cottle Knife Guide Retractor	. 156,	312
Cottle Osteotomes		
Cottle Rasp		
Cottle Retractor156		
Cottle Septum Elevator		
Cottle Septum Specula		.290
Cottle Skin Elevator	*****	.308
Cottle Skin Hooks . 190, 191, 309	, 311	, 543
Cottle Thumb Retractor	. 156	310
Cottle-Kazanjian Nasal Cutting		
Forceps		.338
Cotton Carrier, Rectal		
Coupland Suction Tube	148	, 299





Graidord Goardation Gramp
Crego Elevators458
Crile Forceps92, 327
Crile Needle Holders116, 343
Crile Retractor 158, 195
Crile-Wood Needle
Holders 115, 125, 126, 130, 343, 407
Crypt Hooks474
Cup Jaw Biopsy Forceps510
Curettes, Adenoid304
Curettes, Antrum297
Curettes, Bone318, 436, 456, 457
Curettes, Dermal524, 525
Curettes, Ear274, 284
Curettes, Endaural287
Curettes, Endocervical508
Curettes, Endometrial508
Curettes, Ethmoid
Curettes, Eye
Curettes, Mastoid284
Curettes, Nasal296, 297, 304
Curettes, Neurosurgical436
Curettes, Plastic Surgery318
Curettes, Uterine506, 507
Curettes, Uterine Biopsy508
Cusco Vaginal Specula490
Cushing "Little Joker" Elevator434
Cushing Decompression
Cushing Decompression Retractor161, 43
Cushing Dressing
Forceps 72, 73, 89, 292, 428, 429
Cushing Dura Hook192, 43
Cushing Perforator Drill409
Cushing Pituitary Rongeurs42
Cushing Spatula Spoons425
Cushing Tissue Forceps 72, 73, 428, 429
Cushing Vein Retractor16
Cuticle Nippers550, 55
Cuticle Scissors50
Cutter, Finger Ring546
Cyclodialysis Spoon and Spatula23
Cystic Duct Forceps101, 201, 202
Cystotomes23
System of management of the state of the sta
D
D'Errico Perforator Bur409
Dandy Forceps

Dandy Nerve Hooks	193,	431
Davidson Scapula Retractor	171,	394
Daviel Lens Spoon		
Davis Spatula		
Davis Foreign Body Spud		
Davis Mouth Gags		
Davis Retractors151,		
Davis Sterilizer Forceps		
Davis Tissue Forceps	.73.	399
Day Attic Cannula		279
Day Ear Hooks		
De Lee Stethoscope		
De Lee-Hillis Stethoscope		
De Lee Obstetrical Forceps		
De Lee Universal Retractor	165	493
Dean Knife Needle		
Dean Rongeurs		
Dean Scissors		
Dean Tonsil Knife		
Deaver Retractors		
Deaver Scissors		
DeBakey Bulldog Clamps		
DeBakey Multi Purpose Clamps		
DeBakey Needle Holders		
DeBakey Occlusion Clamps		
DeBakey Pediatric Clamps		
DeBakey Peripheral Vascular Clan		
DeBakey Penpheral Vascular Cian DeBakey Tissue	ips .	3/6
Forceps68, 82, 324,	074	
Parameter Clause 222	J/4,	535
DeBakey Vascular Clamps 378,	J8U,	381
DeBakey-Adson Tissue	074	rno
Forceps70, 324,	3/4,	538
DeBakey-Derra Clamps		
DeBakey-Potts Scissors		
DeBakey-Satinsky Clamps		
DeBakey-Semb Ligature Carrier		
Debridement Scissors		
Dejerine Percussion Hammer		
Denhardl Mouth Gags		
Dennis Anastomosis Clamp		
Depilatory Forceps		
Derf Needle Holder 112, 123,	340,	534
Dermal Curettes		
Dermal Instruments		
Dermal Punches	******	526
Deschamps Ligature Carriers		
Deschamps Needles		142

Desjardin Gall Stone Forceps	205
Desjardin Gall Stone Scoops	200
Desmarres Chalazion Forceps 253,	325
Desmarres Corneal Dissector	232
Desmarres Retractors 154, 229, 312,	546
DeVilbiss Cranial Rongeur	417
DeVilbiss Vaginal Specula	491
DeWecker Iris Scissors33,	259
Dewey Obstetrical Forceps	520
Diagnostic Instruments	.1-4
Diamond Pin Cutter	
Dieffenbach Serrefines108,	327
Diethrich Micro Bulldog	
Clamps108, 224, 267, 327,	406
Dilators, Common Duct	
Dilators, Lacrimal	240
Dilators, Tracheal	369
Dilators, Urethral484,	485
Dilators, Uterine496-	498
Dingman Bone Clamp	441
Directors	144
Disposable Dermal Punches	526
Disposable Scalpels	17
Dissecting Scissors 20-23, 28, 30, 42	50,
334, 335, 398,	405
Dissecting Tenaculum194,	
Dissectors, Tonsil	
Dittel Urethral Sounds	
Dix Foreign Body Needle and Gouge	
Dix Foreign Body Needle and Spud	
Dix Foreign Body Spud	
Dorsal Clamp/Needle Holder128,	
Dorsal Scissors336,	
Double Action Nail Nipper	549
Double Action Parallel Pliers 142,	
Douglas Cilia Forceps80,	
Downing Cartilage Knife	
Doyen Intestinal Forceps 105, 210,	457
	457 211
Doyen Rib Elevators-Raspatories	457 211 388
Doyen Rib Elevators-Raspatories	457 211 388 515
Doyen Rib Elevators-Raspatories	457 211 388 515 505
Doyen Rib Elevators-Raspatories	457 211 388 515 505 493
Doyen Rib Elevators-Raspatories	457 211 388 515 505 493 348
Doyen Rib Elevators-Raspatories	457 211 388 515 505 493 348 244,
Doyen Rib Elevators-Raspatories	457 211 388 515 505 493 348 244, 429,
Doyen Rib Elevators-Raspatories	457 211 388 515 505 493 348 244, 429, 541

Index Market



Drills, Hand445
Drills, Neurosurgical409, 410
Dubois Scissors25, 523, 552
Duckbill Ear Forceps277
Duckbill Rongeurs416, 417, 461, 462
Duke Trocar
Dunning Elevator and Curette296
Duplay Uterine Tenaculum Forceps500
Dupuys-Weiss Ligature Carriers 143, 358
Dupuys-Weiss Needles143
Dura Hooks431
Dura Separators431
Durai Scissors36, 335, 423
Duval Lung Grasping Forceps401
Dynamometers, Hand3
Dynamometers, nand
E
E-Z Off Knife Handles14
Ear Applicators273
Ear Curettes
Ear Forceps
Ear Hooks
Ear Instruments
Ear Knives278
Ear Loops274
Ear Probes145
Ear Scissors48, 49
Ear Snare Wire277
Ear Snares277
Ear Specula269, 270
Ear Spoons273
Ear Suction Tubes 148, 149, 279, 287
Ear Syringes275
Eastman Vaginal Retractors 165, 493
Echlin Duckbill Rongeurs416, 462
Ekg Calipers3
Elevators, Bone446
Elevators, Cleft Palate351
Elevators, Endaural287
Elevators, Nasal295, 296, 308
Elevators, Neurosurgical434
Elevators.
Periosteal .268, 308, 351, 434, 454, 455
Elevators, Rib
Elevators, Septurn295, 296, 308
Elevators, Skin
Elliott Obstetrical Forceps519
Emon Obstetitori i orocpo

Ellis Foreign Body Spud2	
Elsberg Cannula4	
Elschnig Fixation Forceps2	:46
Elschnig-O'Connor Fixation Forceps 2	46
Emmett Uterine Tenaculum	
Hooks194, 5	600
Endaural Curettes	
Endaural Elevators	:87
Endaural Retractors2	:83
Endaural Rongeurs	
Endocervical Curettes5	
Endocervical Specula5	
Endometrial Curettes5	
Endoscopic Forceps	
Engel Cast Saw4	
Enterotomy Scissors	
Entropium Forceps252, 3	
Enucleation Scissors46, 2	50
Enucleation Spoons	
Episiotomy Scissors23, 5	
Eppendorfer Utering Biopsy	
Forceps511, 5	:12
Erhardt Lid Forceps	
Esmarch Bandage Shears	
Esmarch Cast Shears4	
Esophageal Scissors	
Esophagoscopic Forceps	
Ethmoid Curettes	
Ethmoid Forceps	
Eustachian Catheters	
Evans Tissue Forceps	
Eve Tonsil Snares	
Evisceration Spoons	
Ewald Dressing Forceps	
Ewald Tissue Forceps	
Excavator Curettes	
Extra-Capsular Forceps249, 2	
Eye Calipers	
Eye Curettes	
Eye Dressing Forceps 64, 74, 241, 2	
Eye Hooks234, 2	
Eye Magnets	
Eye Needles232, 2	
Eye Scissors34, 256, 330-3	.0°
Eye Shield	
Eye Specula225-2	
Eye Sulure Scissors 41 263 3	

Facelift Scissors336, 337
Fansler Operating Speculum473
Fansler-Ives Rectal Speculum471
Farabeuf Retractors161
Farabeuf-Lambotte Bone Holding
Forceps441
Farr Wire Retractors
Farrell Applicators273, 291
Farrior Ear Specula270
Faulkner Ethmoid Curette297
Fehland Intestinal Clamp212
Feilchenfeld Splinter Forceps83
Female Catheters484
Fenestrated Handle Forceps 64-65,
374, 427, 428, 537, 538
Fenger Gall Duct Probe146, 199
Ferguson Angiotribes107, 477
Ferguson Gall Stone Scoops200
Ferguson Scissors23, 515
Ferguson Suction Tubes148
Ferguson Tenaculum Forceps477, 501
Ferris-Smith Pituitary Rongeurs423
Ferris-Smith Tissue Forceps 82, 429
Finger Goniometer442
Finger Nail Drills548
Finger Retractors158
Finger Ring Cutter546
Finochietto Ligature Carrier
Finochietto Needle Holders118
Finochietto Rib Spreaders389
Finochietto Thoracic Forceps 99, 404
Finochietto-DeBakey Rib Spreader 389
First Aid Shears52
Fisher Spoon And Needle237
Fisher Tonsil Knife-Dissector358
Fistula Scissors27
Fixation Forceps77, 238, 244-246
Fixation Hooks
Flash Autoclaving Pans9
Fletcher Sponge Forceps111, 499
Flexible Neck Retractors 154, 310, 544
Fluoride Coated Scissors52, 466
Foerster Eye Dressing Forceps 74, 243
Foerster Iris Forceps74, 243

Index ____



Fomon Dorsal Scissors		33
Fomon Elevator		30
Fomon Knife		
Fomon Rasps		31
Fomon Retractor	156,	31
Fomon Lateral Scissors		
Fomon Saber Back Scissors	. 36,	33
Forceps Holders		8,
Forceps, Angiotribe	107,	47
Forceps, Atraumatic82,	373	38
Forceps, Biopsy, Bladder		48
Forceps, Biopsy, Rectal	479-	48
Forceps, Bone Cutting 338,	395,	55
Forceps, Bone Holding	440,	44
Forceps, Breast Tenaculum		50
Forceps, Bronchoscopic		36
Forceps, Bulldog Clamp Applying		
Forceps, Capsule249,		
Forceps, Carb-N-Sert		
Forceps, Catheter Introducing		
Forceps, Cervical Biopsy	510-	-51
Forceps, Cervix Holding		49
Forceps, Chalazion	252.	25
Forceps, Cilia79, 80,	250.	25
Forceps, Clip Applying136,	137.	42
Forceps, Clip Removing	136.	13
Forceps, Colibri75, 77,	243.	24
Forceps, Columella		
Forceps, Cornea Scleral		23
Forceps, Corneal243-	245.	24
Forceps, Cystic Duct101,		
Forceps, Depilatory	.82.	54
Forceps, Dressing 64-82, 220,		
406, 427-429, 498, 499, 537,		
Forceps, Ear221,		
Forceps, Endoscopic		36
Forceps, Entroplum		
Forceps, Esophagoscopic		
Forceps, Ethrnoid		
Forceps, Extra-Capsular		
Forceps, Eye Dressing64,		
243, 244,		
Forceps. Fenestrated Handles		
68-70, 76, 82,		
374, 427, 428,		
Forceps, Fixation		
1 Oloopo, 1 Ixanori	,	

Forceps,	Foreign Body	7!
Forceps,	Gall Bladder	.202-20
Forceps,	Gall Duct	201, 202
Forceps,	Gall Stone	20
Forceps,	Goiter	196
Forceps,	Hemostatic 90-108,	402, 430
Forceps,	Hysterectomy	513, 514
	Intestinal	
	Intraocular lens	
Forceps,	Iris 66, 74, 75, 242-	244, 540
	Jeweler86, 87,	
	Kidney Stone	
	Laryngeal	
	Lid	
Forceps,	Lung Grasping	40
Forceps,	Membrane Puncturing	518
	Micro69, 70, 74	
	218-223,	427, 430
	Mosquito90-91,	
Forceps,	Nail Pulling	550
	Nasal Bone Cutting	
	Nasal Dressing	
	Nasal Polyp	
	Obstetrical	
	Placenta	
	Polyp Uterine	
	Ptosis	
	Recession	
	Septum	
	Sequestrum	
Forceps,	Skin Hook	194, 542
	Skin Staple Removing	
	Splinter	
	Splinter, Jewelers Type .	86, 87
	Sponge	
	j110, 111,	
	Staple Removing	
	Sterilizer	
	Strabismus	
	Strut	
	Suturing 78, 80, 88,	
	244, 247, 248,	
Forceps,	Tenaculum	477
Forceps,	Tendon Pulling	458
⊢orceps,	Thoracic399-	404, 408
Forceps,	Tissue 64-82, 88,	89, 220

Forceps, Tissue 321-324, 399, 400,	
477, 504, 537, 538, 540,	54
Forceps, Tongue Seizing5,	350
Forceps, Tonsil 102, 355	-357
Forceps, Tonsil	109
Forceps, Trachoma	259
Forceps, Traction	106
Forceps, Tube Occluding	
Forcess Tumor	I
Forceps, Tumor	. 000
Coronn Turbinate	າດຄ
Forceps, Turbinate	
Forceps, Undermining	
Forceps, Urethral	
Forceps, Uterine Artery	.488
Forceps, Uterine Biopsy 486, 509	512
Forceps, Uterine Dressing498,	
Forceps, Uterine Elevating	
Forceps, Uterine Tenaculum 500,	501
Forceps, Uterine Vulsellum 107, 502,	503
Forceps, Utility 78, 79, 243, 247,	248
Forceps, Vaginal	515
Ford Stethoscope	
Foreign Body Forceps	
Foreign Body Gouges	
Foreign Body Loops	
Foreign Body Needles	
Foreign Body Spuds233,	
Foss Intestinal Clamp	
Foss Resection Clamp	
Foss Retractors	
Fox Dermal Curettes	
Fox Eye Shield	
Fraenkel Laryngeal Forceps	
Francis Chalazion Forceps	
Frankel Head Band	
Franz Retractors	
Frazier Dura Hooks 190, 309,	
Frazier Dura Separator	
Frazier Needle	
Frazier Retractors186,	
Frazier Suction Tubes 148, 279, 299,	
Freer Chisels268, 296,	
Freer Gouges296,	313
Freer Periosteal Elevator 268, 351,	434
Freer Septum Elevator295,	308
Freer Septum Knife	
Froer Skin Hooks 101 300	

Index Index



Freer Submucous Knife		295
French Pattern Lacrimal Probes	145.	240
Friedman Rongeurs224,	287.	353
Fritsch Retractors		170
Fuchs Extra-Capsular Forceps		
Fullon Rongeurs		
. a.o		
G		
Gags, Mouth6, 7,		
Gall Bladder Forceps		
Gall Bladder Instruments	199	-205
Gall Duct Forceps	201,	202
Gall Stone Scoops		200
Gallaher Rasps		319
Galt Skull Trephines		411
Garland Hysterectomy Forceps	96,	513
Gauze Scissors		
Gaylor Uterine Biopsy Forceps		.509
Gellhorn Uterine Biopsy Forceps		
Gelpi Retractors173,		
Gemini Forceps		
Gemini-Mixter Forceps 105,	203.	402
Gerald Dressing Forceps	72.	429
Gerald Tissue Forceps		
Gerzog Mallet284,		
Giertz Rib Shears		
Gigli Saws		
Gigli-Strully Saws		
Gill Corneal Knife		
Gill Iris Forceps		
Gill-Hess Iris Forceps		
Gillies Skin Hooks190,		
Gillies-Sheehan Needle Holders		
Glover Bulldog Clamps		
Glover Coarctation Clamps		
Glover Patent Ductus Clamps		
Gluck Rib Shears		
Goelet Retractor		. 163
Goldman-Fox Scissors35,	334,	532
Goldstein Lacrimal Sac Retractor.	153,	228
Goniometer, Finger		
Gonlotomy Knives		
Goodell Uterine Dilators		.498
Gordon Uterine Vulsellum Forceps	s	.503
Gorney Scissors		
Gosset Retractors		
Gouges, Bone449,	451.	452

Gouges, Lacrimal Sac	
Gouges, Mastoid	.284
Gouges, Septum	.296
Gradle Cilia Forceps80,	
Gradle Scissors 41, 140, 263, 333,	531
Graefe Eye Dressing Forceps75,	243
Graefe Fixation Forceps76, 244,	324
Graefe Iris Forceps 75, 76, 243,	540
Graham Nerve Hook193,	
Grasping Forceps	
Graves Vaginal Specula	
Gray Cystic Duct Forceps	
Green Cystic Duct Forceps 101,	
Green Fixation Forceps77,	
Green Needle Holder 120, 264, 345,	
Green Retractor	
Green Spatula	546
Green Strabismus Hook189	
Greene Uterine Curettes	
Greg Comedone Extractor	
Grooved Directors	
Gross Hook And Spoon1990,	
Gruber Ear Specula	
Gruenwald Bayonet Forceps73,	
Gruenwald Ear Forceps73,	
Gruenwald Nasal Forceps	
Gum Scissors	
Guthrie Fixation Hook 188, 234, 309	
Guttman Vaginal Speculum	
Guyon Vessel Clamp	
Guyton-Noyes Fixation Forceps	
Guyton-Park Eye Specula	. 226
Gynecological And Obstetrical	
Instruments489	J-523
н	
Haight Rib Spreaders	
Hajek Antrum Rongeur	
Hajek-Ballenger Septum Elevator	
Halle Ethmoid Curette	
Halle Nasal Speculum	
Halsey Needle Holders . 112, 123, 340	, 534
Halsted Mosquito	
Forceps90, 91, 267, 326	
Hammer, Post Mortem	553
Handles, Knife1	3, 14
Hook Utoring Dilators	404

larken Rib Spreaders	391
larmon Fixation Forceps	
larms Suture Forceps	78, 247
larrington Retractors	170, 392
Harrington Scissors	.24, 398
lartman Ear Dressing Forceps	276
lartman Ear Specula	270
lartman Eustachian Catheters	
lartman Lee Mosquito Forceps	91
lartman Mosquito Forceps91, 2	
Hartman Nasal Forceos	292
lartman Rongeur285,	353, 459
lartman Tonsil Punch	361
dartman-Herzfeld Ear Forceps	277
lartman-Noyes Ear Forceps	277
lawks-Dennen Obstetrical Forcep	s522
Hays Hand Retractor	183, 447
lead Band, Mirror	3
lead Extractor, Obstetrical	522
Head Mirrors	3
Heaney Hysterectomy Forceps	.96, 513
leaney Needle Holders	
leaney Retractors	164, 494
Heaney Tissue Forceps	
leaney Uterine Biopsy Curette	508
Heaney-Ballantine Hysterectomy	
Forceps	.96, 513
leaney-Simon Retractors	165, 494
Heath Chalazion Forceps	252
Heath Curettes	
Heath Suture Scissors39,	
Hebra Curettes	
Hegar Uterine Dilators	
Hegenbarth Clip Applying Forceps	
Hemorrhoidal Ligator	
Hemostatic Forceps90-108,	
Hennig Cast Spreader	
Hercules Cast Shears	
Hercules Heavy Duty Scissors	
Herrick Kidney Pedicle Clamp	
Herzfeld Ear Forceps	277
Heyman Nasal Scissors	
Hibbs Chisels	450
Hibbs Gouges	
Hibbs Osteotomes	
Hibbs Retractors	182, 393
Hibbs Spinal Fusion Curettes	436

Index 🗀



Hill-Ferguson Hetractors 175, 473
Hirschman Anoscopes471
Hirschman Proctoscope473
Hoffman Ear Punch Forceps277
Hohmann Retractors183, 446
Hoke Osteotomes450
Hollh Corneoscleral Punch262
Hooks188-194, 309-311,
Hooks. Bone
Hooks, Corneal189
Hooks, Crypt
Hooks, Dura
Hooks, Ear
Hooks, Eye
Hooks, I.U.D. Extractor
Hooks, Nerve192, 193, 431, 432
Hooks, Skin188-192, 194,
309-311, 542, 543
Hooks, Strabismus189
Hooks, Sympathectomy193
Hooks, Tenotomy189
Hooks, Tracheal192, 367, 368
Hooks, Uterine Tenaculum 194, 500
Hosford Foreign Body Spuds233
Hosford Lacrimal Dilators240
House Strut Forceps221, 288
House Suction/Irrigating
Tubes148, 279, 426
Hudson Braces410
Hudson Burs410
Hudson Cerebellar Extension
Hudson Dressing Forceps69, 322
Hudson Tissue Forceps
Hunt Chalazion Forceps252
Hupp Trachea Hook192, 367
Hupp Trachea Hook192, 367
Hupp Trachea Retractor 160, 312, 367
Hurd Dissector-Retractor162, 358
Hysterectomy Forceps513, 514
Hysterectomy Retractors494
1
I.U.D. Extractor Hook194, 500
Inge Lamina Spreaders173, 439
Instrument Care Products11
Instrument Care System11
Instrument Care Video11

Instrument Cleaning Brushes11,	319
Instrument Holders	9
Instrument Sterilizing Cases	10
Instrument Sterilizing Trays	
Intestinal Instruments206	
Intraocular Lens Forceps	
Iowa Membrane Forceps	
lowa Trumpet Needle Guide	
Iridectomy Scissors	
Iris Forceps .66, 74, 75, 242-244, 321,	540
Iris Scissors 30-33, 254-256, 259,	วรณ
262, 329, 330, 527	
Iris Spatula	000
Israel Retractors172, 393,	
Ivan Laryngeal Applicator	ამხა
Jackson Bronchoscopic Forceps	966
Jackson Endoscopic Forceps	
Jackson Esophagoscopic Forceps	
Jackson Laryngeal Forceps	
Jackson Trachea Bistouri	
Jackson Trachea Retractor 160, 313,	367
Jackson Trachea Tenaculum	
Hook192,	
Jackson Trachea Tubes	
Jackson Vaginal Retractors 165,	493
Jacobs Uterine Vulsellum	
Forceps107,	503
Jacobson Hemostatic	
Forceps90, 103, 223, 402,	430
Jacobson Micro Mosquito	
Forceps90, 223, 267, 326,	
Jacobson Needle Holders 129, 222,	408
Jaeger Lid Plate	229
Jameson Caliper239,	317
Jameson Recession Forceps	253
Jameson Strabismus Hook189,	235
lansen Baynnet	
Forceps73, 81, 275, 292,	429
Jansen Ear Forceps	275
Jansen Retractors151, 280,	
Jansen-Middleton Septum Forceps	302
Jansen-Wagner Retractors 151, 280,	
Jarcho Uterine Tenaculum Forceps	
Jennings Mouth Gags7,	
Definings Mount Gags	349

Johns Hopkins Bulldog Clamps 108, 406
Johns-Hopkins Gall Bladder
Forceps101, 202
Jones Forceps93
Jones Towel Clamps109
Jorgenson Scissors25, 398
Joseph Button End Knives307
Joseph Hooks191, 311, 543
Joseph Nasal Knives307
Joseph Nasal Saws316, 317
Joseph Nasal Scissors38, 337, 531
Joseph Periosteal Elevator308
Joseph Rasp318
Judd-Allis Tissue Forceps207
Judd-Demartel Gall Bladder Forceps 204
Juers-Lempert Rongeur287
Julian Needle Holders127
Jung Mini-Chole Retractors 166, 204
К
Kahler Bronchoscopic Forceps365
Kahn Dissecting Scissors28
Kahn Uterine Tenaculum Forceps 500
Kalt Needle
Holders 120, 128, 264, 266, 345, 535
Kane Umbilical Clamp523
Kantrowitz Thoracic Forceps99, 403
Kazanjian Nasal Forceps338
Kelly Anoscope472
Kelly Hemostatic Forceps92
Kelly Placenta Forceps522
Kelly Procloscope472
Kelly Retractors169
Kelly Scissors27, 423, 478, 516
Kelly Sphincteroscopes472
Kelly Sponge Forceps
Kelly Tissue Forceps73, 399, 504
Kelly Uterine Vulsellum Forceps502
Keratoplasty Scissors47
Kern Bone Holding Forceps440
Kerrison Rongeurs268, 285, 303
Kevorkian-Younge Endocervical
Curettes508
Kevorkian-Younge Uterine Biopsy
Forceps511, 512
Key Elevators455
Keyes Punches526

Index Marketine



Kidney Clamps	
Kidney Stone Forceps	.487
Kielland Obstetrical Forceps	.520
Kielland-Luikart Obstetrical Forceps	.520
Killian Antrum Cannula	.299
Killian Septum Forceps	
Killian Septum Specula	.290
Kirby Cornea-Scieral Forceps	
Kirby Fixation Forceps	.238
Kirby Hook And Expressor	
Kirby Iris Forceps74	
Kirschner Wires	
Kleinert-Kutz Bone Cutting Forceps	
Kleinert-Kutz Hooks190, 309	
Kleinert-Kutz Rongeurs	
Kleinert-Kutz Scissors35, 335	
Knapp Cataract Knife	
Knapp Eye Speculum	
Knapp Iris Scissors 32, 256, 330	
Knapp Knife Needle	
Knapp Retractors155, 229, 310	
Knife Guides	
Knife Needles	
Knight Nasal Forceps	
14 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	
Knight Nasal Scissors49	294
Knight Septum Forceps	294
Knight Septum Forceps Knives, Abscess	294 300 527
Knives, Abscess	. 294 . 300 . 527 . 465
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Brain	. 294 . 300 . 527 . 465 . 552
Knight Septum Forceps	. 294 . 300 . 527 . 465 . 552 . 307
Knight Septum Forceps	. 294 . 300 . 527 . 465 . 552 . 307 . 457
Knight Septum Forceps	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Brain Knives, Button End Knives, Cartilage Knives, Cataract Knives, Ear	. 294 .300 .527 .465 .552 .307 .457 .231
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Bratin Knives, Button End Knives, Catrilage Knives, Cataract Knives, Ear Knives, Meniscus	. 294 .300 .527 .465 .552 .307 .457 .231 .278
Knight Septum Forceps	. 294 .300 .527 .465 .552 .307 .457 .231 .278 .457
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Brain Knives, Button End Knives, Cartilage Knives, Cataract Knives, Ear Knives, Hassal Knives, Meniscus Knives, Masal Massal M	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 278 . 457 . 307 . 466
Knight Septum Forceps	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 278 . 457 . 307 . 466
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Brain Knives, Button End Knives, Cartilage Knives, Cataract Knives, Ear Knives, Hassal Knives, Meniscus Knives, Masal Massal M	. 294 .300 .527 .465 .552 .307 .457 .278 .457 .307 .466 .295
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Bratin Knives, Button End Knives, Cartilage Knives, Cataract Knives, Cataract Knives, Meniscus Knives, Meniscus Knives, Plaster Knives, Plaster Knives, Septum	, 294 , 300 , 527 , 465 , 552 , 307 , 457 , 231 , 278 , 457 , 307 , 466 , 295 , 358
Knight Septum Forceps Knives, Abscess Knives, Ampulating Knives, Brain Knives, Button End Knives, Catrliage Knives, Cataract Knives, Cataract Knives, Maniscus Knives, Maniscus Knives, Plaster Knives, Plaster Knives, Septum Knives, Tonsil	, 294 , 300 , 527 , 465 , 552 , 307 , 457 , 231 , 278 , 457 , 307 , 466 , 295 , 358 ,51
Knight Septum Forceps	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231 . 278 . 457 . 307 . 466 295
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Button End Knives, Catrilage Knives, Catrilage Knives, Catracat Knives, Ear Knives, Button End Knives, Ear Knives, Button End Knives, Ear Knives, Button End Knives, Masal Less Knives, Nasal Less Knives, Plaster Knives, Septum	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231 . 278 . 457 . 466 . 295 . 358
Knight Septum Forceps Knives, Abscess Knives, Ampulating Knives, Brain Knives, Button End Knives, Catrilage Knives, Cataract Knives, Cataract Knives, Meniscus Knives, Meniscus Knives, Plaster Knives, Plaster Knives, Tonsil Knowles Bandage Scissors Kocher Hemostatic Forceps Kogan Endocervical Specula	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231 . 278 . 457 . 307 . 466 295
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Amputating Knives, Button End Knives, Catrilage Knives, Catrilage Knives, Catrilage Knives, Ear Knives, Meniscus Knives, Nasal 295 Knives, Plaster Knives, Septum Knives, Septum Knives, Tonsil Knowles Bandage Scissors Kocher Hemostalic Forceps 92 Kocher Intestinal Forceps Kogan Endocervical Specula Kollmann Dilators	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231 . 278 . 457 . 307 . 466 295
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Brain Knives, Button End Knives, Catrilage Knives, Catrilage Knives, Catrilage Knives, Catrilage Knives, Ear Knives, Meniscus Knives, Nasal 295 Knives, Plaster Knives, Septum Knives, Septum Knives, Tonsil Knowles Bandage Scissors Kocher Hemostatic Forceps 92 Kocher Intestinal Forceps Kogan Endocervical Specula Kollmann Dilators Krause Laryngeal Tips	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 278 . 457 . 307 . 466 . 295
Knight Septum Forceps Knives, Abscess Knives, Amputating Knives, Amputating Knives, Button End Knives, Catrilage Knives, Catrilage Knives, Catrilage Knives, Ear Knives, Meniscus Knives, Nasal 295 Knives, Plaster Knives, Septum Knives, Septum Knives, Tonsil Knowles Bandage Scissors Kocher Hemostalic Forceps 92 Kocher Intestinal Forceps Kogan Endocervical Specula Kollmann Dilators	. 294 . 300 . 527 . 465 . 552 . 307 . 457 . 231 . 278 . 457 . 307 . 466 . 295

L		
Laborde Trachea Dilator		369
Lacrimal Dilators		240
Lacrimal Probes		145
Lacrimal Sac Chisels		268
Lacrimal Sac Gouges		
Lacrimal Sac Retractors		
Lacrimal Trephines		
Laforce Adenotomes		305
LaGrange Scissors35, 257,		
Lahey Forceps92, 100, 196,	197.	201
Lahey Gall Duct Forceps	100.	201
Lahey Goiter Scissors		
Lahey Hemostatic Forceps		
Lahey Ligature Carriers		
Lahey Retractor		
Lahey Traction Forceps		
Lambert Chalazion Forceps		
Lambert-Kay Aortic Clamp		
Lambert-Lowman Bone Clamps		
Lambotte Bone Holding Forceps		
Lambotte Osteotomes		
Lamina Spreaders		
Laminectomy Retractors		
Lancaster Eye Speculum		
Lane Bone Holding Forceps		
Lange Eye Speculum		
Langenbeck Elevators		
Langenbeck Metacarpal Saw		
Langenbeck Retractors		
Larry Rectal Probe		
Laryngeal Applicator		
Laryngeal Extension Cannula		
Laryngeal Forceps		
Laryngeal Instruments		
Laryngeal MirrorsLaryngeal Mirrors		
Laryngeal Tips		
Laryngectomy Tubes		
Lateral Vaginal Retractor		
Lathbury Applicators		
Lebsche Sternum Knife		
Lee Bronchus Clamp		
Lee Hartman Mosquito Forceps		
Lefort Urethral Sounds		483
Leksell Rongeurs397,		
Lempert Curettes		287

Lempert Rongeurs	.286
Lens Loop	
Lerner Skin Hook Forceps 194,	542
Lester Fixation Forceps77,	
Lewin Bone Holding Forceps	.441
Lewis Laryngectomy Tubes	.372
Lewis Lens Loops	
Lewis Rasps	.319
Lewis Septum Forceps	
Lewis Tonsil Screw	
Lichtwicz Antrum Needle	
Lid Forceps252	325
Lid Plate	
Lid Retractors154, 229	546
Life Saving Tubes	.369
Ligature Carriers 142, 143	
Lilly Dissecting Scissors	
Lilly Tonsil Scissors	
Lindner Cyclodialysis Spatula And	
Spoon	.237
Lip Retractor	
Lister Bandage Scissors	
Liston Amputating Knives	
Liston Bone Cutting Forceps	
Littauer Cilia Forceps79	
Littauer Dressing Forceps	
Littauer Nasal Forceps	
Littauer Stitch	
Scissors39, 138, 139, 332, 333	533
Littauer-Liston Bone Cutting Forceps	
Little Joker Elevator	
Little Retractor	
Littler Suture Carrying	. 101
Scissors41, 140	333
Lobectomy Scissors	308
Locke Elevators	
Locklin Scissors	
Lockwood Tissue Forceps207	
Lockwood-Allis Tissue Forceps	
Long Hysterectomy Forceps97	
Long Stitch Scissors40, 139	
Loops, Ear	
Lorna Towel Clamps	
Lothrop Uvula Retractor	
Love Nerve Retractors 159, 313	
Love Uvula Retractors	
Love-Gruenwald Pituitary Rongeurs	

Index 🔤



Love-Kerrison Cervical Rongeurs	,419
Lovelace Hemostatic Forceps 93	197
Lovelace Gall Bladder Forceps 103	
Lovelace Lung Grasping Forceps	
Lowe-Breck Cartilage Knife	
Lower Gall Duct Forceps100	
Lower Lateral Scissors	
Lowman Bone Clamps	
Lowman-Hoglund Bone Clamps	.442
Luc Septum Forceps	
Lucae Dressing Forceps	
Lucae Ear Forceps	
Lucae Maliet284, 316	453
Luer Rongeurs414	459
Luer S-Shape Retractor	.158
Luikart Obstetrical Forceps	.520
Lukens Trachea Retractor 160, 312	367
Lung Grasping Forceps	.401
Lung Retractor171	394
-	
M	
M.G.H. Needle Holders	.116
M.G.H. Vulsellum Forceps	
Magili Catheter Forceps	
Magnifying Loupe	
Maier Ulerine Dressing Forceps	
Malik Cystic Duct Catheter Clamp	
Mallets, Bone	
Mallets, Mastoid284	
Mallets, Plastic Surgery	
Maltz Cartilage Knife	
Maltz Rasps318	, 319
Manhattan Eye And Ear Suturing	
Forceps	.248
Martin Cartilage Clamps	.458
Martin Cartifage Scissors25	458
Martin Laryngectomy Tubes	.372
Martin Pelvimeter	
Martin Tissue Forceps71, 324	504
Martini Bone Curette318	
Masson Needle Holders115	
Masson-Judd Retractors174	486
Mastin Muscle Clamp	
Mastoid Chisels	
Mastoid Curettes	
Mastoid Gouges	.284

Mastoid Mallets......284,285

Masiola Hetraciors	
Mastoid Rongeurs285,	286
Mastoid Searcher	
Mathieu Needie Holders . 14, 124, 342,	
Mathieu Retractors 159, 195, 312, 368,	545
Mathieu Urethral Forceps	
Matson Rib Stripper-Elevator	.388
Matson-Alexander Rib Elevator-	
Stripper	.388
Maumenee-Park Eye Speculum	.226
Maunder Oral Screw5,	350
Maxillo-Facial Instruments 143,	338
Mayo-Collins Retractor	.163
Mayo Common Duct Probes 146,	
Mayo Common Duct Scoops	
Mayo Dissecting Scissors20	
Mayo Gall Stone Scoop	
Mayo Operating Scissors	
Mayo Retractors	
Mayo Russian Tissue Forceps 73, 399,	504
Mayo Uterine Scissors	
Mayo Vessel Clamp	
Mayo-Adams Retractors	
Mayo-Harrington Scissors	.398
Mayo-Hegar Needle	
Holders115, 125, 130,	342
Mayo-Noble Dissecting Scissors	
Mayo-Robson Intestinal Forceps., 105,	
McClure Iris Scissors33,	
McCrea Infant Sounds	
McCullough Suturing Forceps 78,	
McCullough Utility Forceps78,	
McGannon Retractor 155, 235, 310,	
McGivney Grasping Forceps	475
McGivney Hemorrhoidal Ligator	
McGuire Corneal Scissors47.	
McIvor Mouth Gag7,	
McKenty Septum Elevator	
McKenzie Clip Forceps137,	
McKenzie Clip Rack135,	
McKenzie Enlarging Bur	
McKenzie Perforator Drill	
McKinney Eye Speculum	
McLean Obstetrical Forceps	
McLean-Luikart Obstetrical Forceps	
McLean-Tuckert-Luikart Obstetrical	
Forceps	520

McPherson Corneal Suturing Forceps .	244
McPherson Needle Holder 121, 265,	
McPherson Suluring Forceps .78, 220,	
McPherson Tissue Forceps	
McPherson-Vannas Iris	
Scissors	262
Mead Mallet285, 316,	
Mead Rongeurs352,	
Measuring Tape	
Mechanical Finger	
Meller Lacrimal Sac Retractor 153,	
Mellinger Eye Specula	
Meltzer Adenoid Punch	
Menghini Biopsy Needles	198
Meniscus Knives	157
Metzenbaum Needle Holders	116
Metzenbaum Scissors24, 334, 359,	
Meyerding Finger Retractors 160,	
Meyerding Retractors 168, 312, 368,	
Meyerding Self-Retaining Retractors	
Meyerding Skin Hook	
Meyerding Skin Retractors 160, 312,	
Meyhoefer Curettes	
Michel Clip Applying Forceps	
Michel Wound Clips	
Michele Trephines	
Micro Adson Dressing	410
Forceps	197
Micro Adson Tissue	461
Forceps	197
Micro Bulldog Clamps 108, 224,	
Micro Clip Approximator	
Micro Dissecting Scissors	
Micro Ear Forceps	
Micro Ear Scissors	
Micro Forceps69, 70, 74, 76, 78,	
427, 430, 537, 538,	5U,
427, 430, 337, 336, . Micro Mosauito	J4 I
Forceps90, 103, 267, 326,	420
Micro Rongeurs224,	40U 0E0
Micro Scissors33, 42-44, 49.	
Micro Surgery	00ء
Instruments 108, 215-224, 427,	420
Micro Surgery Needle Holders 119, 1	ZU, oer
Micro Vessel Clips21, 128, 129, 222, 3	
Micro-Friedman Rongeur224, 287, 3	
wicto-rhealtan nongeut 224, 287, a	ುವರ

Index Index



Miculicz Forceps95,	515
Miculicz Tonsil Hemostat	356
Millard Thimble Hook	
Retractor156, 191,	
Miller Rectal Scissors27,	
Millin Capsule Forceps	
Millin Urethral Forceps	
Miltex Blepharoplasty Scissors	
Miltex Cannulated Pin Cutter 443,	
Miltex Cast Spreader	468
Miltex Comedone Extractor	548
Miltex Cross Action	
Retractors 152, 281,	545
Miltex Depilatory Forceps82,	548
Miltex Dermai Curettes524,	525
Miltex Double Action Pin Cutter	
Miltex Hooks 191, 234, 311, 500,	543
Miltex I.U.D. Extractor Hook 194,	500
Miltex Mallet285, 316,	
Miltex Pin Puller142,	
Miltex Rhytidectomy Scissors	336
Miltex Stitch Scissors 40, 139, 333,	
Miltex Tissue Twist Hook 192, 431,	
Miltex Undermining Forceps	403
Miltex Uterine Tenaculum Hook 194,	500
Miltex-Kerrison Cervical Rongeurs	500 420
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum	500 420 507
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles	500 420 507 14
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal	500 420 507 14 .362
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 .362
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 .362 199
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 .362 199
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 .362 199 202, 403
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 .362 199 202, 403
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 .362 199 202, 403 .202
Miltex-Kerrison Cervical Rongeurs	500 420 507 14 362 362 199 202, 403 .202 .201
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mixer Dilating Probe 146, Mixer Forceps 98-100, 201, 205, 327, Mixter Gall Bladder Forceps Mixter Gall Duct Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps	500 420 507 14 362 362 199 202, 403 .202 .201 .205 20
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixter Dilating Probe 146, Mixter Forceps 98-100, 201, 205, 327, Mixter Gall Bladder Forceps Mixter Gall Stone Forceps Mixter Toracic Forceps Mixter Toracic Forceps	500 420 507 14 362 199 202, 403 .202 .201 205 20
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mixter Dilating Probe Mixter Forceps Mixter Forceps Mixter Gall Bladder Forceps Mixter Gall Stone Forceps Mixter Scissors Mixter Thoractc Forceps Mixter Mixte	500 420 507 14 362 199 202, 403 .202 .201 205 20
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miriature Blade Handles Mirror Handles, Laryngeal Mixter Dilating Probe Mixter Forceps Mixter Gall Bladder Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Mixter Mixter Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mott Mouth Gags 7, Morrison-Hurd Dissector-	500 420 507 14 362 362 199 202, 403 202 .201 .205 20 .403 349
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixter Dilating Probe 146, Mixter Forceps 98-100, 201, 205, 327, Mixter Gall Bladder Forceps Mixter Gall Duct Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Molt Mouth Gags 7, Morrison-Hurd Dissector- Retractor 163,	500 420 507 14 362 362 199 202, 403 .202 .201 .205 20 403 349
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miriature Blade Handles Mirror Handles, Laryngeal Mixter Dilating Probe Mixter Forceps Mixter Gall Bladder Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Mixter Mixter Mixter Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Mix	500 420 507 14 362 362 199 202, 403 .202 .201 .205 20 403 349
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirror Laryngeal Mixter Dilating Probe 146, Mixter Forceps 98-100, 201, 205, 327, Mixter Gall Bladder Forceps Mixter Gall Btone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Molt Mouth Gags 7, Morrison-Hurd Dissector Retractor 163, Mosher Ethmoid Curette Mosher Life Saving Tubes	500 420 507 14 362 362 199 202, 403 .202 .201 .205 20 20 20 349 358 297 369
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixter Dilating Probe 146, Mixter Forceps 98-100, 201, 205, 327, Mixter Gall Bladder Forceps Mixter Gall Blout Forceps Mixter Gall Stone Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mixter Horacic Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mixter Thoracic Forceps Motth Mouth Gags 7, Morrison-Hurd Dissector Retractor 163, Mosher Ethmoid Curette Mosher Life Saving Tubes Mosquito Cardiovascular Clamps	500 420 507 14 362 362 199 202, 403 .202 201 .205 20 403 349 358 .297 .369 .402
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixter Dilating Probe 146, Mixter Dilating Probe 205, 327, Mixter Gall Bladder Forceps Mixter Gall Duct Forceps Mixter Gall Duct Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mott Mouth Gags 7, Morrison-Hurd Dissector- Retractor 163, Mosher Ethmoid Curette Mosher Life Saving Tubes Mosquito Cardiovascular Clamps Mosquito Forceps Mosquito Forceps Mosquito Forceps Mosquito Forceps 90, 91, 103, 267,	500 420 507 14 362 362 199 202, 403 .202 .201 20 403 349 358 .297 .369 .402 326
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixer Dilating Probe 146, Mixer Forceps 98-100, 201, 205, 327, Mixer Gall Bladder Forceps Mixter Gall Duct Forceps Mixter Gall Stone Forceps	500 420 507 14 362 362 199 202, 403 .202 .201 20 403 349 358 .297 .402 326 349
Miltex-Kerrison Cervical Rongeurs Miltex-Kogan Endocervical Speculum Miniature Blade Handles Mirror Handles, Laryngeal Mirrors, Laryngeal Mixter Dilating Probe 146, Mixter Dilating Probe 205, 327, Mixter Gall Bladder Forceps Mixter Gall Duct Forceps Mixter Gall Duct Forceps Mixter Gall Stone Forceps Mixter Thoracic Forceps Mott Mouth Gags 7, Morrison-Hurd Dissector- Retractor 163, Mosher Ethmoid Curette Mosher Life Saving Tubes Mosquito Cardiovascular Clamps Mosquito Forceps Mosquito Forceps Mosquito Forceps Mosquito Forceps 90, 91, 103, 267,	500 420 507 14 362 362 199 202, 403 .202 .201 20 403 349 358 .297 .402 326 349

maidour Edonina Dilator mann		******	P
Multipurpose Clamps38	31.	382.	38
Murdock Eye Speculum			22
Murless Head Extractor			
Murphy Retractors			
Murphy-Lane Bone Skid			
Muscle Clamps			
Museholdt Nasal Forceps	••••		29
Museux Uterine Vulsellum			
Forceps		107,	50
N			
Nabatoff Vein Stripper			
Nail Nippers			
Nail Scissors			
Nail Splitters		50,	55
Nasal Applicators			29
Nasal Bone Scissors		38,	33
Nasal Curettes29			
Nasal Elevators29	95.	296.	30
Nasal Forceps 292, 293, 3	nn.	302	33
Nasal Instruments			
Nasal Knives			
Nasal Polyp Forceps			
Nasal Probes1			
Nasal Retractors			
Nasal Saws			
Nasal Scissors49, 29			
Nasal Snare Wire			
Nasal Snares			
Nasal Specula			
Nasal Suction Tubes1	40	440	20
Nasopharyngeal Specula			
Needle Holders112-13			
340-346, 407, 408, 4	87	534	-53
Needle Holders, Micro			
Surgery 119, 120, 121, 12			
265, 3			
Needle Holders, Urethral			
Needle Sterilizer Rack			
Needles, Aneurysm			
Needles, Antrum			
Needles, Biopsy			
Needles, Eye		232,	23
Needles, Surgical		. 131	-13
Neivert Knife Guide-Retractor .			
Neivert Needle Holders			

Neivert Osteotomes
Neivert-Anderson Osteolomes 315
Nelson Scissors
Nerve Hooks431, 432
Nerve Retractors 433
Nerve Separators431, 432
Neuro-Diagnostic Sets1
Neurological Hammers1
Neurological Pinwheel
Neurological Instruments409-439
New Orleans Needle Holders 116, 126
New's Suture Scissors39, 138
News Tracheal Hook192, 367
Noise Apparatus271
Northbent Stitch Scissors 39, 139, 333, 533
Novak Endometrial Biopsy Curettes 508
Noyes Alligator Ear Forceps 276, 293
Noyes Iridectomy Scissors47, 260
Noyes Iris Scissors32, 259, 530
Noyes Nasal Forceps293
Nugent Hook189, 234, 309
Nugent Utility Forceps78, 247
Nurse's Scissors52
0
O O'Brien Fixation Forcens 246
O'Brien Fixation Forceps246
O'Brien Fixation Forceps246 O'Brien Slitch
O'Brien Fixation Forceps246 O'Brien Slitch Scissors40, 139, 263, 333, 533
O'Brien Fixation Forceps 246 O'Brien Stitch Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234
O'Brien Fixation Forceps 246 O'Brien Stlich Scissors SCISSOR 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242
O'Brien Fixation Forceps 246 O'Brien Slitch 333,533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234
O'Brien Fixation Forceps 246 O'Brien Slitch 352 Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176
O'Brien Fixation Forceps 246 O'Brien Slitch Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal
O'Brien Fixation Forceps 246 O'Brien Stlich 305 Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Tis Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 1816267 Retractor 176, 494
O'Brien Fixation Forceps 246 O'Brien Slitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Obstetrical Forceps 518-522
O'Brien Fixation Forceps 246 O'Brien Slitch 333,533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Destetrical Forceps 518-522 Obwegeser Awls 143, 338
O'Brien Fixation Forceps 246 O'Brien Slitch 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Obstetrical Forceps 518-522 Obwegser Awls 143, 338 Occlusion Clamps 11
O'Brien Fixation Forceps 246 O'Brien Stlitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Retractor 176, 494 Obstetrical Forceps 518-522 Obwegaser Awls 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199
O'Brien Fixation Forceps 246 O'Brien Slitch 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Obstetrical Forceps 518-522 Obwegser Awls 143, 338 Occlusion Clamps 11
O'Brien Fixation Forceps 246 O'Brien Stlitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Retractor 176, 494 Obstetrical Forceps 518-522 Obwegaser Awls 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199
O'Brien Fixation Forceps 246 O'Brien Stlitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Tis Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Obstetrical Forceps 518-522 Obwegser Awls 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198
O'Brien Fixation Forceps 246 O'Brien Slitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Iris Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Destetrical Forceps 518-522 Obwegeser Awls 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198 Oldberg Pituliary Rongeur 422
O'Brien Fixation Forceps 246 O'Brien Stlitch 333, 333, 533 O'Connor Hook 189, 234 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Destetrical Forceps 518-522 Obseterical Forceps 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198 Oldberg Pituliary Rongeur 422 Ollier Retractor 172, 393, 448 Olsen-Heiger Needle Holders 113, 124,
O'Brien Fixation Forceps 246 O'Brien Stlich Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Tis Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Destetrical Forceps 518-522 Obstetrical Forceps 113, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198 Oldier Retractor 172, 393, 448 Olsen-Hegar Needle Holders 113, 124, 341, 535 341, 535
O'Brien Fixation Forceps 246 O'Brien Slitch Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Tis Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Obstetrical Forceps 518-522 Obwegser Awis 143, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198 Oldberg Pituitary Rongeur 422 Olier Retractor 172, 393, 448 Olsen-Hegar Needle Holders 113, 124, Operating Scissors 18-50, 54-63
O'Brien Fixation Forceps 246 O'Brien Stlich Scissors 40, 139, 263, 333, 533 O'Connor Hook 189, 234 O'Connor Tis Forceps 74, 242 O'Connor Tenotomy Hook 189, 234 O'Sullivan-O'Connor Retractors 176 O'Sullivan-O'Connor Vaginal 176, 494 Destetrical Forceps 518-522 Obstetrical Forceps 113, 338 Occlusion Clamps 11 Ochsner Gall Duct Probe 146, 199 Ochsner Trocars 198 Oldier Retractor 172, 393, 448 Olsen-Hegar Needle Holders 113, 124, 341, 535 341, 535





Oral Surgery Rongeurs	352-354
Orthopedic Instruments	440-465
Osleotomes, Orthopedic	449-452
Osteotomes, Plastic Surgery	314, 315
Otis Urethrotome	485
Otoscopes	
• • • • • • • • • • • • • • • • • • • •	
P	
Packer Mosquito Forceps	90
Parker Mott Retractors	164
Parker Retractors	161
Patent Ductus Clamps	
Paton Needle Holders 121, 264	1, 346, 536
Patterson Trocars	198
Payr Pylorus Clamps	214
Pean Forceps	
Pean Hysterectomy Forceps	514
Pederson Vaginal Specula	490
Pediatric Vascular Clamps	
Pedicle Clamps98	
Peet Splinter Forceps	85
Pelvimeters	517
Pencil-Type Micro Forceps	215
Pencil-Type Micro Scissors	42, 215
Penfield Dissectors	435
Pennington Forceps	102, 476
Pennington Rectal Speculum	471
Pennington Septum Elevator	296
Percussion Hammers	1
Percy Intestinal Forceps	208
Perineal Retractors	495
Periosteal Elevators 268, 308	, 351, 388,
Periosteal Elevators 389, 434	454, 455
Periosteotomes	388
Peripheral Vascular Clamps	378
Peritoneal Clamp	
Perkins Endaural Retractor	185, 283
Phaneuf Uterine Artery Forceps .	97, 514
Physicians Splinter Forceps	84
Pierce Antrum Trocar	298
Pierce Elevator	
Pierse Colibri Corneal Forceps	
Piffard Dermal Curettes	
Pillar Retractors	162, 358
Pin Cutters	
Pinard Stethoscope	
Piper Obstetrical Forceps	519

Plain Wire Eye Speculum226,	
Plaster Knives	.466
Plaster Saws	
Plaster Shears 51-53, 466,	
Plastic Surgery Calipers	.317
Plastic Surgery Curettes	
Plastic Surgery Instruments 307	
Plastic Surgery Mallets	
Plastic Surgery Rasps318,	
Plastic Surgery Scissors 29, 328-337,	339
Platypus Nail Forceps	
Pliers, Wire Cutting142,	
Politzer Ear Knives	
Polk Finger Goniometer	442
Polyp Forceps, Uterine	
Poole Suction Tubes	
Post Mortem Hammer	
Post Morlem Instruments552,	
Potts-Smith Dressing Forceps 71, 89,	
Potts-Smith Scissors26,	
Potts-Smith Tissue Forceps72,	
Pratt Crypt Hooks	
Pratt Ethmoid Curette	
Pratt Rectal Probe146,	
Pratt Rectal Scissors	
Pratt Rectal Speculum	469
Pratt T-Shape Forceps94, 208,	
Pratt Uterine Dilators	
Presbyterian Hosp. Occluding Clamp	
Prince Forceps	
Prince Trachoma Forceps	
Probes, Ear145,	
Probes, General Operating144	-146
Probes, Lacrimal	145
Probes, Nasal144, 145,	291
Probes, Rectal144, 146,	
Probes, Sinus	
Proctoscopes	4/3
Providence Hospital Forceps91, 196,	320
Ptosis Forceps252, Pudendal Needle Guide	
Punches, Adenoid	
Punches, Corneoscleral	
Punches, Dermal Punches, Dermal Disposable	
Punches, Dermai Disposable Punches, Tonsil	
Putti Bone Rasps	
Lam Date Uggby	404

Pylorus Clamps	
Pylorus Separator	213
Pynchon Suction Tubes 6, 147, 350,	357
O	
Quick Release Needle Holders	190
Quimby Scissors36,	
Quire Mechanical Finger	
Quire Mechanica: Friger	210
R	
Rabiner Neuro Hammer	
Rack For Surgical Needles8,	135
Rack, Clip Holding135,	
Ragnell Dissecting Scissors 30, 334,	532
Ragnell Retractor158, 312,	368
Rake Retractors154, 155, 310,	
Randall Endometrial Biopsy Curette	508
Randall Kidney Stone Forceps	
Raney Scalp Clip Forceps143,	
Raney Scalp Clips143,	
Raney-Crutchfield Skull Tong	411
Rankin Forceps	
Raspatories, Rib	
Rasps, Antrum	297
Rasps, Plastic Surgery318,	319
Recession Forceps	253
Rectal Cotton Carrier	
Rectal Instruments469-	
Rectal Probes144-146,	
Rectal Scissors	
Rectal Specula469-	
Reiner Cast Knife	
Reiner Rongeur	
Remover, Blade	
Renal Clamps	
Retractors, Anal174, 175,	
Retractors, Appendectomy	
Retractors, Bladder174,	
Retractors, Bone446-	
Retractors, Cerebellar	
Retractors, Decompression161,	
Retractors, Endaural	
Retractors, Gall Bladder	
Retractors, General Operating 151-	
446-	
Retractors, Hyslerectomy	
Retractors, Lacrimal Sac228,	229

Index



Retractors, Laminectomy 173,	437	-439
Retractors, Lid154,	229,	546
Retractors, Lung	171,	394
Retractors, Mastoid		280
Retractors, Nasal		304
Retractors, Nerve		432
Retractors, Nerve Root 159,	313,	433
Retractors, Perineal		495
Retractors, Pillar	162,	358
Retractors, Rectal 172, 174,	175,	470
Retractors, Scalp		
Retractors, Scapula171,		
Retractors, Spinal	171.	412
Retractors, Splanchnik		392
Retractors, Thoracic	.392	-394
Retractors, Tonsil		358
Retractors, Tracheal160,	367,	368
Retractors, Uvula		304
Retractors, Vaginal 176, 187, 493	495.	507
Retractors, Vein		432
Reynolds Dissecting Scissors 23,		
Rhoton Micro Forceps		
Rhytidectomy Scissors		
Rib Elevators		
Rib Shears394, 395,	552.	553
Rib Spreaders		
Rib Strippers		388
Ribbon Retractors		
Richardson Retractors		
Richardson-Eastman Retractors		169
Richter Dissecting Scissors		
Richter-Heath Clip Forceps		
Rigby Retractors	173,	495
Ring Tonque Refractor Blades		
Rochester Gall Stone Forceps	******	.205
Rochester Pean Forceps		
Rochester-Carmalt Forceps		
Rochester-Ochsner Forceps		
Roeder Towel Clamp		
Roger Septum Elevator		
Rollet Retractor 155,		
Rongeurs, Antrum		
Rongeurs, Bone		
Rongeurs, Duckbill 416, 417,	461,	462
Rongeurs, Endaural		
Rongeurs, Mastoid		
Rongeurs, Neurosurgical		

Rongeurs, Oral Surgery352-354
Rongeurs, Thoracic
Ronis Adenoid Punches306
Rosser Crypt Hook194, 474
Roux Retractors164
Rowland Nasal Hump Forceps
Rubin Septal Morselizer301
Ruedemann Lacrimal Dilator240
Ruler, Flexible239, 546
Rumel Thoracic Forceps
Ruskin Antrum Needles298
Ruskin Bone Cutting
Forceps338, 395, 463
Ruskin Rongeurs286, 397, 415, 460
Russian Tissue Forceps71, 399, 504
Ryder Needle Holders 126, 128, 407, 408
11ydel Needle Holdela 120, 120, 407, 400
S
Saalfield Cornedone Extractors547
Saber Back Scissors
Sachs Nerve Separator
Sage Tonsil Snare
Sarot Artery Forceps
Sarot Bronchus Clamps404 Sarot Needle Holders 117, 127, 343, 408
Satinsky Clamps
Satinsky Scissors27, 406
Satinsky-DeBakey Clamps
Satterlee Bone Saw465
Sauer Eye Speculum227
Sauer Suturing Forceps248
Sauerbruch Bone Rongeurs397
Sauerbruch Retractors161, 393
Saws, Bone465
Saws, Cast
Saws, Nasal316, 317
Saws, Wire413, 464
Sawtell Tonsil Hemostats
Sawyer Rectal Retractors172, 470
Sayre Elevators455
Scalp Retractors411, 412
Scalpel Blades 15-17
Scalpel Handles13-14
Scapula Retractor391
Schaaf Foreign Body Forceps79, 248
Schambaugh Adenotomes305
Schamberg Comedone Extractors 547

Scheie-Westcott Corneal Section
Scissors46, 258
Schink Retractor184
Schiotz Tonometers230
Schlesinger Cervical Rongeurs418
Schlesinger Pituitary Rongeurs422
Schnidt Tonsil Hemostats102, 356
Schroeder Uterine Vulsellum Forceps502
Schroeder-BraunTenaculum Forceps 500
Schubert Uterine Biopsy Forceps509
Schuhknecht Endaural Retractors 185, 283
Schwartz Vessel Clips224, 430
Schwartz Clip Applying Forceps 224, 430
Schweiger Extra-Capsular Forceps 249
Scissors, Abdominal23
Scissors, Bandage51, 52
Scissors, Blepharoplasty 59, 335, 531
Scissors, Brophy35, 331, 351
Scissors, Bulldog37, 337
Scissors, Capsulotomy48, 262
Scissors, Carb-N-Sert54-59
Scissors, Cardiovascular26, 405
Scissors, Cartilage
Scissors, Conjunctival32, 257
Scissors, Corneal 43, 47, 48, 260, 261
Scissors, Corneal
Section46, 47, 258, 260
Scissors, Cuticle50
Scissors, Debridement35, 532
Scissors, Decapitation523
Scissors, Dermal38, 528-533
Scissors, Dissecting 20-23, 28, 30, 42,
50, 334, 335, 398, 405, 532, 552
Scissors, Dorsal37, 336, 337
Scissors, Dural36, 335, 423
Scissors, Ear48, 49
Scissors, Enterotomy50, 552
Scissors, Enucleation46, 259
Scissors, Episiotomy23, 516
Scissors, Esophageal398
Scissors, Eye34, 41
Scissors, Facelift37, 336, 337
Scissors, Fistula27
Scissors, Gall Bladder205
Scissors, Gauze51, 52
Scissors, Gum331
Scissors, Iridectomy260





	Milio M	
Scissors	, Iris 259, 260, , Keratople	262, 329 sty

Scissors, Iris30-33, 44, 254-256,
259, 260, 262, 329, 330, 527-530
Scissors, Keratoplasty47
Scissors, Lobectomy398
Scissors, Lower Lateral36, 336
Scissors, Micro33, 44, 49, 288
Scissors, Nail50, 551
Scissors, Nasal38, 49, 294, 336, 337
Scissors, Nasal Bone38
Scissors, Neurosurgical423
Scissors, Operating 18-50, 54-63
Scissors, Plastic Surgery 29, 328-337, 339
Scissors, Rectal 478
Scissors, Rhylidectomy36, 336
Scissors, Saber Back36, 336
Scissors, Septal38
Scissors, Serratex19, 21
Scissors, Stitch39, 40, 138-140, 258,
263, 332, 333, 351, 533
Scissors, Strabismus34, 257, 529
Scissors, SuperCut60-63
Scissors, Suture39-41, 138-140, 258
Scissors, Suture Carrying 41, 140, 333
Scissors, Tenotomy35, 45, 46, 257, 258,
Scissors, Tonsil355, 359
Scissors, Tungsten Carbide54-59
Scissors, Turbinate49
Scissors, Úmbilical23, 523
Scissors, Undermining 38, 334, 531, 532
Scissors, Uterine515, 516
Scissors, Utility31, 32, 46, 52, 255,
257, 258, 332, 529
Scissors, Wire Cutting 40, 41, 141, 339
Scissors, Wound Debridement35
Scoops, Common Duct200
Scoops, Gall Stone200
Scoville Spatula425
Scoville Clip Forceps137, 424
Scoville Disc Curettes436
Scoville Nerve Retractors 159, 433
Screw, Tonsil358
Scudder Intestinal Forceps 105, 211
Sedillot Periosteal Elevator389, 454
Seiler Turbinate Scissors49, 294
Seldin Periosteal Elevator308
Seltzer Dermal Hook 194, 234, 311, 542

Semken Dressing Forceps	67 322
Semken Dressing Forceps	
Semken Tissue Forceps .67, 322	
Senn Retractors 158, 195, 312	
Septal Morselizer	
Septum Elevators295	
Septum Forceps	
Septum Gouges	
Septum Knives	
Septum Scissors	
Septum Specula	
Septum Straightners	
Serratex Operating Scissors	
Serratex-Mayo Scissors	
Serrefines108	
Sexton Ear Knives	
Shallcross Cystic Duct Forceps	
Shapleigh Ear Curettes	
Sheehan Osteotomes	
Sheehan-Gillies Needle Holders .	
Shepard Intraocular Lens Forceps	
Shortbent Stitch Scissors	39, 138
263	332, 53
Siegel Knife Handle	
Siegel Otoscope	
Silver Osteotomes	
Simpson Obstetrical Forceps	518
Simpson Uterine Sound	498
Simpson-Luikart Obstetrical Force	ps520
Sims Rectal Speculum	
Sims Scissors	
Sims Uterine Biopsy Curette	
Sims Uterine Curettes	
Sims Uterine Sound	
Sims Vaginal Specula	
Singley Tissue Forceps 73,	209 399
Sinus Cannula	
Sinus Probes	
Sistrunk Scissors	
Skeele Curettes	
Skene Uterine Tenaculum Forcep	
Skene Uterine Vulsellum Forceps	5 5U
Skillern Sinus Probe	145 604
Skin Elevators	
Skin Hook Forceps	IMA DAY
Skin Hooks188-192, 194,	309-311

Skin Staple Removing Forceps		
Skull Traction Tong		
Skull Trephines		
Smillie Meniscus Knives		
Smillie Meniscus Hook Retractor .		
Smith (Buie) Anal Retractors		
Smith Eye Speculum		
Smith-Fisher Iris Spatula		
Smith-Petersen Gouges		. 452
Smith-Petersen Osteolomes		
Smith-Petersen Rongeurs		
Smithwick Clip Forceps	137,	424
Smithwick Nerve Hook	193,	432
Smithwick Sympathectomy Hook .	193,	432
Snare Wire, Ear		
Snare Wire, Nasal		
Snares, Ear		
Snares, Nasal		
Snares, Tonsil		
Snellen Entropium Forceps		
Somer Uterine Elevating Forceps		
Sonnenschein Nasal Speculum		
Sounds, Urethral		
Spatula, Brain		
Specula, Ear		
Specula, Endocervical	209,	210
Specula, Eye		.007
Specula, Nasal289,		
Specula, Nasopharyngeal		
Specula, Rectal		
Specula, Septum		
Specula, Vaginal		
Spencer Chalazion Forceps		
Spencer Ear Probes		
Spencer Stitch Scissors39,	170, 120	263
openier olien obiaana	332	533
Sphenoid Cannula		
Sphincteroscopes		
Spinal Retractors		
Splanchnik Retractors		
Splinter Forceps		
Splinter Forceps With Magnifier		
Splinter Liberator		
Splinter Removal Kit		
Splitters, Nail		



Sponge Holding Forceps110, 111,
356, 499
Spoon And Spatula, Cyclodialysis
Spoons, Ear273
Spoons, Enucleation236
Spoons, Evisceration
Spratt Bone Curettes436, 456
Spratt Mastoid Curettes
Spratt Spinal Fusion Curettes436
Spray Lube11
Spring Wire Retractors152
Spurling Pituitary Rongeurs421
Spurling-Kerrison Cervical Rongeurs 419
St. Vincents Tube Occluding Forceps 11
Staple Removing Forceps136
Staude-Moore Tenaculum Forceps501
Steiner-Auvard Vaginal Speculum 492
Sterilizer Forceps8
Sterilizer Rack For Needles8
Sterilizing Aids8-10
Sterilizing Trays
Sternum Knife
Stelhoscope4
Stethoscopes, Obstetrical517
Stevens Iris Forceps
Stevens Needle Holder121, 265
Stevens Tenotomy Hook
Stevens Tenolomy Scissors45, 257,
Stevenson Lacrimal Sac Retractor 153, 228
Stewart Crypt Hooks193, 474
Stille Cranial Drills410, 445
Stille Patt. Vessel Clamp410, 445
Stille Cast Shears
Stille Rib Shears
Stille Tissue Forceps69
Stille-Horsley Bone Cutting
Forceps
Stille-Liston Bone Cutting
Forceps
Stille-Luer Rongeurs396, 417, 461
Stitch Scissors 39, 40, 138-140, 258,
Stone Intestinal Clamps213
Stone Tissue Forceps71
Stookey Rongeur414
Strabismus Forceps76, 244

Strabismus Hooks189
Strabismus Scissors 34, 257, 331, 529
Stratte Needle Holders 117, 127, 487
Struempel Ear Forceps277
Strully Scissors
Strut forceps
Struyken Turbinate Forceps302
Suction Tubes, Ear148, 149
Suction Tubes, General Operating 5. 6.
147-149, 350, 351, 426
Suction Tubes, Nasal 148, 149, 299
Suction Tubes, Rectal149, 481
SuperCut Scissors60-63
Surgi-Guide Template317
Surgical Blades15-17
Surgical Instrument Cleaner11
Surgical Instrument Stain Remover 11
Surgical Needle Rack8, 135
Surgical Needles131-134
Suture Carrying Scissors41, 140, 333
Suture Needles131-134
Suture Scissors 39-41, 138-140, 258, 263
Suture Sundries135-143
Suture Wire
Suturing Forceps
248, 251, 541
Sweet Scissors28, 398
Swiss Cilia And Suture Forceps 80, 138,
251, 325, 541
Swiss Jeweler Forceps218, 219
Swiss Pattern Osteotomes452
Syringes, Ear
т
T-C Carb-N-Sert Dressing
And Tissue Forceps88, 89
T-C Carb-N-Sert Needle
Holders 122-130, 222, 266
T-C Carb-N-Sert Scissors54-59
Takahashi Ethmoid Forceps302
Tape Measure4
Taylor Dural Scissors36, 335, 423
Taylor Percussion Hammers1
Taylor Spinal Retractor171, 412
Teale Uterine Vulsellum Forceps 502
Tenaculum Forceps, Breast501
Tenaculum Forceps, Uterine 500, 501

Tenaculum, Dissecting553
Tendon Pulling Forceps458
Tendon Strippers458
Tenotomy Hooks189
Tenotomy Scissors 35, 45, 46, 257, 258,
Theobald Sinus Probe145, 291
Thomas Uterine Curettes506
Thoms-Ailis Tissue Forceps .207, 476, 502
Thoms-Gaylor Uterine Biopsy Forceps .509
Thoracic And Cardiovascular
Instruments368-408
Thoracic Clamps404
Thoracic Forceps
Thoracic Tissue Forceps
Thorek Scissors
Thorpe Caliper
Thorpe Corneal And Splinter
Forceps87, 249
Thyroid-Goiter Instruments 195-197
Tieck-Halle Nasal Speculum289
Tischler Cervical Biopsy Punch
Forceps510, 512
Tissue Forceps 64-82, 88, 89, 321-324,
399, 400, 427-429, 476, 477.
Tissue Nippers
Tivnen Tonsil Seizing Forceps355
Tonque Depressors347
Tongue Retractor Blades6
Tongue Seizing Forceps
Tonometer Sterilizing Case230
Tonometers
Tonsil Dissectors
Tonsil Forceps
Tonsil Instruments
Tonsil Knives
Tonsil Punches
Tonsil Retractors
Tonsil Scissors
Tonsil Screw
Tonsil Snare Wire360, 361
Tonsil Snares
Tooke Corneal Knife232
Towel Clamps109 Townsend Cervical Biopsy Punch
Forceps510





roynbee Ear Specula269
Trachea Tube Brushes372
Trachea Tubes371, 372
Tracheal Bistouri369
Tracheal Dilators369
Tracheal Hooks367, 368
Tracheal Instruments367-372
Tracheal Retractors367, 368
Trachoma Forceps253
Trays, Sterilizing
Trephines, Lacrimal268
Trephines, Neurosurgical411, 413
Trocars198, 298
Tromner Percussion Hammer1
Trousseau Trachea Dilators
Troutman-Barraquer Blade Breaker231
Troutman-Barraquer Colibri
Forceps75, 243
Troutman-Barraquer Corneal
Forceps75, 243
Tube, Life Saving369
Tube Occluding Clamps12
Tubes, Laryngectomy372
Tubing Clamps12
Tuffier Retractor186, 437
Tuffier Rib Spreaders391
Tumor Forceps505
Tumor Screw505
Tungsten Carbide Dressing Forceps 88, 89
Tungsten Carbide Needle
Holders122-130, 222, 266
Tungsten Carbide Scissors54-59
Tungsten Carbide Tissue Forceps 88, 89
Tuning Forks
Turbinate Forceps302
Turbinate Scissors49
Turrell Biopsy Forceps480
Tuttle Tissue Forceps399
Twist Drills, Neurosurgical410, 445
Twist Hook, Tissue 192, 431, 542
Tyding Tonsil Snare361
Tyrell Hooks188, 234, 309, 542
IJ
U.S. Army Chisels
U.S. Army Gouges449
U.S. Army Osteotomes

U.S. Patt, Cast Knife	466
U.S.A. Gauze Scissors	51
Uebe Applicators273	
Umbilical Clamps	
Umbilical Scissors23	
Undermining Forceps	
Undermining Scissors38, 334, 531	
Universal Handle	
Universal Laryngeal Tips	
University of Minnesota Cheek	304
Retractor162	0.47
Unna Comedone Extractor	
Urethral Catheters	
Urethral Dilators484	
Urethral Forceps	
Urethral Sounds482	
Urethrotome	
Urological Instruments482	
Ulerine Arlery Forceps	
Uterine Biopsy Curettes	508
Uterine Biopsy Forceps486, 509	
Uterine Curettes506	, 507
Uterine Dilators496	
Uterine Dressing Forceps498	, 499
Uterine Elevating Forceps	505
Uterine Polyp Forceps	. 499
Uterine Scissors515	, 516
Uterine Sounds	. 498
Uterine Tenaculum Forceps500	, 501
Uterine Tenaculum Hooks	
Uterine Vulsellum Forceps 107, 502	. 503
Utility Forceps 78, 79, 243, 247	
Utility Scissors31, 32, 46, 52,	
258, 332	
Utility Shears51-53, 466	
Utility-Sterilizer Forceps	
Uvula Retractors	
orale Hobactory	
V	
Vaginal Forceps	.515
Vaginal Retractors493-495	
Vaginal Specula489	
Van Alyea Antrum Cannula 298	
Van Alyea Sinus Cannula	
Van Buren Sequestrum Forceps	
Van Buren Sounds	
Vannas Capsulotomy Scissors 48	
cap only coldoor on the	,

Vannas Scissors44, 217, 26	12
Vascular Clamps 377-381, 384-38	
Vascular Tissue Forceps	
Vaughan Abscess Knife	
Vein Strippers15	
Vessel Clips223, 224,43	
Video, Instrument Care1	11
Vienna Nasal Specula28	19
Virchow Brain Knives55	
Virchow Skull Breaker55	
Virtus Splinter Forceps8	
Vogel Adenoid Curette30	
Volkman Finger Retractors15	
Volkman Hooks181, 188, 44	
Volkman Retractors158, 181, 182	
447, 44	
Volkmann Bone Curettes318, 45	
Von Eicken Antrum Cannula29	
Von Graefe Cystotome23	
Von Graefe Strabismus Hooks 189, 23	
Von Graeffe Cataract Knife23 Vulsellum Forceps, Uterine 107, 502, 50	
vuisenum rorceps, Olerine 107, 502, 50	U
w	
Wachenfeld Clip Forceps13	6
Wachenfeld Clip Forceps13 Wadsworth-Todd Cautery23	
Wachenfeld Clip Forceps	9
Wadsworth-Todd Cautery23	1
Wadsworth-Todd Cautery23 Wagner Retractors151, 280, 41	1 5
Wadsworth-Todd Cautery23 Wagner Retractors151, 280, 41 Waldeau Fixation Forceps24	1915
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Watter Corneal Spuds 23 Walter Splinter Forceps 83, 8	19 1 5 18 12 15
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48	9 1 5 8 2 5 4
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 63, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48	19 15 18 12 15 14 13
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23	19 15 18 12 15 14 13
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle	9 1 5 8 2 5 4 3 3
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle Holders 118, 128, 40	9 1 5 8 2 5 4 3 3
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 63, 8 Walter Dilator-Catheters 48 Walther Dilator-Gatheters 48 Walther Forcign Body Gouge 23 Wangensteen Needle 118, 128, 40 Wangensteen Tissue 118, 128, 40	9 1 5 8 2 5 4 3 3 8
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walter Dilator-Catheters 48 Walther Dilator-Catheters 48 Walther Dilator-Gatheters 48 Walton Foreign Body Gouge 23 Wangensteen Needle Holders Holders 118, 128, 40 Wangensteen Tissue Forceps Forceps 73, 209, 39	9 1 5 8 2 5 4 3 3 8
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle Holders Holders 118, 128, 40 Wangensteen Tissue Forceps Forceps 73, 209, 39 Wartenberg Pinwheel 39	915825433
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle Holders 118, 128, 40 Wangensteen Tissue Forceps 73, 209, 39 Wartenberg Pinwheel Watson-Williams Potyp Forceps 30	915825433 8 930
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 63, 8 Walter Dilator-Catheters 48 Walther Dilator-Catheters 48 Walther Diethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle 118, 128, 40 Wangensteen Tissue Forceps 73, 209, 39 Wartenberg Pinwheel 30 Watson-Williams Polyp Forceps 30 Webster Needle Holders 112, 123, 130	19 1 15 18 12 15 14 13 13 18 19 3 10 0,
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 63, 8 Walter Dilator-Catheters 48 Walther Dilator-Catheters 48 Walter Foreign Body Gouge 23 Wangensteen Needle Holders Holders 118, 128, 40 Wangensteen Tissue Forceps Forceps 73, 209, 39 Watson-Williams Polyp Forceps 30 Webster Needle Holders 112, 123, 130 Webster Needle Holders 112, 123, 130 340, 53 340, 53	19 1 15 18 12 15 14 13 13 18 19 3 10 0, 14
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Dilator-Catheters 48 Walther Dilator-Sunds 48 Walton Foreign Body Gouge 23 Wangensteen Needle 118, 128, 40 Wangensteen Tissue Forceps Forceps 73, 209, 39 Wartenberg Pinwheel 30 Webster Needle Holders 112, 123, 130 340, 53 36 Wecker Iris Spatula 33	19 1 15 18 12 15 14 13 13 18 19 3 10 0,14 17
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Dilator-Catheters 48 Walther Urethral Sounds 48 Walton Foreign Body Gouge 23 Wangensteen Needle Holders 118, 128, 40 Wangensteen Tissue Forceps 73, 209, 39 Wartenberg Pinwheel Watson-Williams Polyp Forceps 30 Webster Needle Holders 112, 123, 130 Webster Needle Holders 112, 123, 130 Wecker Iris Spatula 23 Weingartner Ear Forceps 27	19 1 15 18 12 15 14 13 13 18 19 3 10 0 7 14 7 7
Wadsworth-Todd Cautery 23 Wagner Retractors 151, 280, 41 Waldeau Fixation Forceps 24 Walsham Septum Straightners 33 Walter Corneal Spuds 23 Walter Splinter Forceps 83, 8 Walther Dilator-Catheters 48 Walther Dilator-Catheters 48 Walther Dilator-Sunds 48 Walton Foreign Body Gouge 23 Wangensteen Needle 118, 128, 40 Wangensteen Tissue Forceps Forceps 73, 209, 39 Wartenberg Pinwheel 30 Webster Needle Holders 112, 123, 130 340, 53 36 Wecker Iris Spatula 33	19 1 1 5 18 12 15 14 13 13 18 19 3 10 0 14 7 7 18

Index Index




Wertheim-Cullen Pedicle
Clamp98, 487, 515
West Lacrimal Gouge268
West Lacrimal Sac Chisel268
Westcott Suture Scissors41, 140,
258, 333
Westcott Tenotomy Scissors 46, 258, 530
Westcott Utility Scissors46, 258
Wheeler Discission Knife
White Tonsil Seizing Forceps355
Wieder Tongue Depressors347
Wiener Antrum Rasps297
Wiener Eye Speculum227
Wiener Hooks
Wies Chalazion Forceps253
Wigmore Cast Saw
Wilde Dressing Forceps81, 291
Wilde Ear Forceps274
Wilde Ethmoid Forceps
Wilde Nasal Forceps291
Wilder Cyslotome
Wilder Lacrimal Dilator240
Wilder Lens Loops235
Willauer Scissors
Williams Eye Specula225
Williams Intestinal Forceps208
Williams Lacrimal Probes145, 240
Williams Splinter Forceps84
Williams Tissue Forceps400
Willinger Bone Curette318, 457
Wills Hospital Utility Forceps79, 248
Wilmer Conjunctival And Utility
Scissors32, 257, 529
Winter Placenta Forceps523
Wire Brush For Rasps And Files 12, 319
Wire Cutters444
Wire Cutting Pliers142, 443
Wire Cutting Scissors 40, 41, 141, 339
Wire Guide143, 442
Wire Retractors152
Wire Saws413
Wire Twisters
Wire, Suture
Wire, Tonsil Snare360, 361
Wires, Kirschner445
Wittner Uterine Biopsy Forceps511
Wolf-Bookler Cast Breaker 468

135
288
48, 288
185, 283
383
culum 305
5, 6, 147,
351, 357
479, 480
488
350
36, 415, 460
79, 250
232
240



CAT. NO	CAT. NO	CAT. NO
1-44	2-144	4-7
1-200	2-145	4-7R
1-202	2-146	4-8
1-208	2-147	4-8R ,,
1-2101	2-150	4-9
1-2111	2-154	4-10
1-212	2-160	4-11
1-2131	2-162	4-12
1-2141	2-164	4-16
1-220	2-180	4-17
1-222	2-182	4-18
1-226	3-2	4-20
1-228	3-4	4-30
1-240	3-10	4-31
1-242	3-14	4-100
1-243	3-20	4-110
1-2443	3-22	4-111
1-3000	3-24	4-112
1-30053	3-26	4-1128
1-30153	3-50	4-115
1-32004	3-60	4-115C
1-40004	3-70	4-120
1-40064	3-71	4-121
1-40104	3-72	4-122
2-20	3-74	4-123
2-50	3-400	4-125
2-60	3-401	4-310
2-62	3-402	4-311
2-705	3-404	4-312
2-725	3-405	4-312B
2-100	3-499	4-315
2-104	3-500	4-315C
2-104SS	3-501	4-320
2-106	3-502	4-321
2-106SS	3-550	4-322
2-109	3-552	4-323
2-110	3-554	4-325
2-110SS	3-570	4-400
2-120	3-572	4-401
2-122	3-574	4-410
2-130	3-576	4-411
2-132	3-620	4-412
2-132A	3-624	4-415
2-134	3-628	4-420
2-1344	3-640	4-421
2-134A	3-700	4-422
2-136A	3-720	4-423
2-138	3-740	5-2
2-138A	3-800	
2-140	3-810	5-6
2-141	3-1000	5-6TC
2-142	3-1001 11, 319	
2-143	3-2000	5-8







CAT. NO	CAT. NO	CAT. NO
5-271	5-SC-12060	6-64
5-272	5-SC-122	6-66
5-273	5-SC-124	6-72
5-274	5-SC-126	6-74
5-275	5-SC-128	6-76
5-276	5-SC-130	6-82
5-277	5-SC-180	6-84
5-278	5-SC-182	6-86
5-279	5-SC-182TC	6-90
5-280	5-SC-188	6-92
5-281	5-SC-283	6-94
5-282	5-SC-284	6-96
5-282A	5-SC-304	6-100
5-283	5-SC-304TC	6-104 67, 322, 428, 537
5-284	5-SC-306	6-106 67, 322, 428, 537
5-286	5-SC-306TC	6-107 67, 322, 428, 537
5-288	5-SC-32062	6-108 67, 322, 428, 537
5-290	5D-251 27, 423, 478, 516	6-109 67, 322, 428
5-291	5D-253 27, 423, 478, 516	6-109A 67, 322, 428
5-300 30, 254, 329, 527	5D-255	6-110
5-302 30, 254, 329, 527	5D-303 30, 254, 329, 527	6-112
5-304 30, 254, 329, 527	5D-322	6-114
5-304TC58	5D-324	6-118 69, 322, 427, 537
5-306	5D-326	6-118A 69, 322, 427, 537
5-306TC58	5D-328	6-118TC
5-308	5D-329	6-118XL 68, 69, 322, 427, 537
5-310	6-2	6-120 70, 323, 427, 537
5-312	6-4	6-120A 70, 323, 427, 537
5-312TC58	6-6	6-120XL 68, 70, 323, 427, 538
5-314	6-8	6-121 70, 323, 427, 538
5-314TC58	6-8XL	6-122
5-316	6-10	6-123
5-318	6-12	6-123TC88
5-320	6-14	6-123XL 68, 70, 323, 427, 538
5-320TC59	6-16	6-124
5-510	6-18	6-124TC 88
5-512	6-24	6-124XL 68, 70, 324, 428, 538
5-514	6-26	6-125
5-516	6-28	6-126
5-547	6-28XL	6-128
5-550	6-30	6-132
5-560	6-30TC	6-134
5-561	6-40	6-136
5-564	6-42	6-138
5-568	6-44	6-140
5-600	6-44XL	6-141
5-700	6-46	6-142
5-702	6-47	6-144
5-800	6-48	6-145
5-998	6-50	6-148
5-1000	6-52	6-150
5-SC-16	6-60	6-154
5-SC-46	6-62	6-154TC



CAT. NO	CAT, NO	CAT. NO
6-156	7-14	7-162
6-156TC	7-16	7-164
6-158	7-20	7-166
6-15BTC	7-24	7-167
6-159	7-24A	7-170
6-160	7-26	7-172
6-162	7-30	7-180
6-164	7-32	7-182
6-164A	7-36	7-190 95, 356, 515
6-165TC89	7-38	7-200
6-166	7-42	7-202
6-166TC89	7-44	7-206
6-168	7-46	7-210
6-180	7-48	7-210A
6-182	7-50	7-220
6-184	7-52	7-230
6-186	7-60	7-232 103, 197, 209
6-188	7-62	7-234
6-189	7-66	7-236
6-190	7-68	7-240
6-192	7-70	7-242
6-197	7-72	7-250
6-198	7-82	7-252
6-198TC	7-84	7-254
6-200	7-90	7-255
6-202	7-91	7-256
6-204	7-92	7-257 105, 203, 402
6-210	7-94	7-258 105, 203, 402
6-212	7-98	7-259
6-214	7-102	7-260
6-300	7-106	7-262
6-302	7-110	7-264
6-304	7-112	7-266
6-308	7-118	7-280 108, 267, 406
6-310	7-120	7-282
6-312	7-122	7-284
6-314	7-124	7-286
6-316	7-126	7-288
6-318	7-128	7-290
6-330	7-130	7-292
6-336	7-136	7-294 108, 267, 406
6-337	7-138	7-296
6-338	7-140	7-298
6-340	7-142	7-300
6-342	7-144	7-302
6-350	7-146	7-304
6-400	7-148	7-305
7-4	7-152	7-315
	7-152 95	7-500
7-8	7-156	7-501
7-11	7-157	7-504
7-12	7-160	7-506



CAT. NO	CAT. NO	CAT. NO
7-510	8-44TC125	8-136
7-514	8-46	8-706TC
7-515	8-46TC	8-744TC
7-600		
7-602	8-48	8-746TC
		8-748TC
7-604 110, 356	8-49TC 125,	8-750TC
7-606	8-50	8-751TC
7-610 110, 499	8-50A-TC	8-752TC
7-612 110, 499	8-50TC125	8-753TC
7-614	8-51	8D-36
7-616	8-51TC	9-1
7-620	8-52	9-2
7-622	8-52TC	9-4
7-624	8-53	9-6
7-630	8-53TC	9-8
7-632	8-55TC	9-10
8-2	8-57TC	9-32
8-4	8-58	9-40-18,
8-4TC123	8-59	9-40-20
8-6	8-62TC	9-40-22
8-6TC123	8-63TC	9-40-24
8-7	8-65TC	9-40-25
8-7A-TC	8-66TC	9-40-26
8-7TC123	8-67TC	9-40-28
8-8	8-69TC	9-40-30
8-8A-TC	8-70	9-40-32
8-8TC123	8-74	9-40-34
8-10	8-76	9-40-36
8-12	8-80	9-44-18
8-14	8-81TC126	9-44-20
8-14A	8-82TC127	9-44-22
8-14A-TC124	8-83TC127	9-44-24
8-14TC124	8-84	9-44-25
8-15	8-94	9-44-26
8-15TC124	8-96 117, 343	9-44-28
8-16	8-96TC	9-44-30
8-16TC124	8-100	9-44-32
8-17	8-100TC	9-44-34
8-17TC124	8-101TC	9-50
8-20	8-102	9-52
8-22	8-102TC	9-54
8-23	8-103TC	9-60
8-26	8-104TC	9-62
8-28	8-105TC	9-74
8-28TC	8-106TC	9-80
8-30	8-107TC	9-90
8-30TC	8-114	9-94
8-32	8-114TC	9-96
8-40	8-122	9-98
8-40TC124	8-130	9-100
8-42	8-130TC	9-101
8-42TC125	8-134	9-102
8-44	8-134TC	9-103
	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	2




DAT NO BACE	CAT. NO	CAT. NO
CAT. NO	10-74	11-82
9-104		
9-106	10-76	11-88
9-108	10-79	11-90 160, 312
9-110	10-80	11-92
9-120	10-82	11-94
9-122	10-84	11-96
9-124	10-220	11-100
9-124TC 59, 141, 339	10-290	11-102
9-128	10-291	11-110
9-140	10-310	11-112
9-142	10-312	11-114
9-146	11-2	11-116
9-148	11-4 151, 280, 411	11-118
9-154	11-6	11-120
9-156	11-8	11-122 163
9D-121	11-10 151, 281, 412	11-130
9D-123	11-12 152, 281, 545	11-132
9D-132	11-13 152, 281, 545	11-134
9D-133	11-14 152, 281, 545	11-136
9D-134	11-16	11-140 163
9D-135 41, 141, 339	11-17	11-144
10-2-NS	11-18	11-146 163
10-2-ST144	11-19	11-150
10-4-NS	11-20 154, 195	11-152
10-4-ST	11-22 154, 195	11-153
10-6-NS	11-23	11-154
10-6-ST 144	11-30 154, 229, 310, 368, 544	11-155
10-8-NS	11-32 154, 229, 310, 368, 544	11-160
10-8-ST144	11-34 154, 229, 310, 368, 544	11-162
10-10-NS	11-36 154, 310, 368, 544	11-164
10-10-ST144	11-38 154, 310, 368, 544	11-166
10-12-NS	11-40 154, 310, 368, 544	11-180
10-12-ST144	11-42	11-182
10-14-NS	11-43	11-184
10-14-ST144	11-44	11-185
10-22-SS144	11-45	11-192
10-22-ST	11-46	11-194
10-24-SS144	11-47 158	11-196
10-24-ST	11-48	11-200
10-26-SS	11-49	11-202
10-26-ST144	11-50 155, 310, 368, 544	11-204 166
10-28-SS	11-52 155, 310, 368, 544	11-206
10-28-ST144	11-54	11-208 167
10-30-SS144	11-56	11-210
10-30-ST144	11-58 155, 310, 368, 544	11-218
10-32-SS	11-60 155, 310, 368, 544	11-219
10-32-ST 144	11-73 158, 312, 368	11-220 167
10-34-SS	11-74	11-221 167
10-34-ST144	11-76	11-222
10-38	11-78	11-223 167
10-69	11-79 158	11-224
10-70	11-80	11-224A
10-72	11-81	11-225



CAT. NO	CAT. NO	CAT. NO
11-240	11-442 179	11-634-Bi
11-242	11-444 ,	11-634-SH
11-244	11-446	11-642-BL 184, 195, 412, 438
11-246	11-448	11-642-SH 184, 195, 412, 438
11-247	11-470	11-650
11-248	11-500 181, 188, 447	11-652
11-249	11-502	11-654
11-250	11-504	11-656
11-252	11-506	11-660
11-254	11-508	11-672
11-256	11-512	11-676 186, 391, 437
11-260	11-514 181, 447	11-684 187, 196, 438
11-260-A 169	11-516	11-686
11-260-B169	11-530	11-700
11-266	11-532	11-702
11-267	11-534	11D-27 154, 229, 310, 368, 544
11-268	11-536	12-80
11-269	11-538	12-90
11-270	11-540	12-92
11-272	11-542	12-170
11-274	11-544	13-8
11-276	11-546	13-9
11-290	11-547	13-10
11-292	11-548	13-11198
11-300	11-549 182	13-20
11-314	11-551 182	13-40
11-316 170, 392	11-552	13-42
11-318	11-553	13-44
11-342	11-554	13-46
11-346	11-560	13-48
11-349	11-562	13-50
11-350 172, 393, 448	11-564	13-60
11-352	11-566 182	13-62
11-360	11-568	13-64
11-362	11-570	13-66
11-365 173, 439, 495	11-572 182	13-68
11-366 173, 439, 495	11-574 182	13-70
11-368	11-580 182, 446	13-72
11-374 173, 495	11-582	13-74
11-390	11-590	13-120
11-392	11-600	13-122
11-400	11-602	13-124
11-402	11-608	13-126
11-404	11-610 184, 282	13-128
11-410	11-611 184	13-130
11-414	11-614	13-132
11-420	11-616	13-134
11-421	11-617	13-150
11-424	11-618	13-152
11-428	11-620	13-154
11-430	11-622	13-156
11-432	11-630-BL	14-2
11-440	11-630-SH	14-4-15146, 199








18-33	CAT. NO	CAT. NO	CAT. NO
18-56 227 18-372 223 18-564 237 18-46 227 18-376 233 18-564 237 18-47 227 18-389 233 18-565 237 18-40 227 18-389 233 18-565 237 18-50 227 18-384 233 18-565 237 18-50 227 18-384 233 18-570 237, 546 18-54 227 18-389 233 18-570 237, 546 18-54 227 18-389 233 18-576 237 18-63 227 18-399 233 18-576 227 18-64 227 18-399 233 18-576 227 18-66 228 18-399 233 18-576 227 18-66 228 18-400 234 18-60 237 18-66 228 18-400 234 18-60 74, 238 18-69 153, 228 18-402 234 18-62 74, 238 18-68 153, 228 18-402 234 18-62 74, 238 18-68 153, 228 18-402 234 18-62 74, 238 18-68 153, 228 18-404 234 18-62 74, 238 18-68 153, 228 18-404 234 18-62 74, 238 18-68 153, 228 18-404 234 18-62 74, 238 18-69 155, 229, 310, 544 18-41 188, 234, 309, 542 18-65 229, 317 18-100 155, 229, 310, 544 18-41 188, 234, 309, 542 18-65 229, 317 18-110 154, 229, 312, 546 18-418 189, 234, 309, 542 18-65 229, 317 18-111 154, 229, 312, 546 18-418 189, 234, 309, 542 18-66 229, 317 18-112 154, 229, 546 18-42 18-62, 309, 542 18-66 229, 317 18-114 154, 229, 546 18-42 18-52, 309, 542 18-66 229, 317 18-114 154, 229, 546 18-42 18-92, 309, 542 18-66 229, 317 18-114 154, 229, 546 18-42 18-92, 309, 542 18-66 229, 317 18-114 154, 229, 546 18-42 18-92, 309, 542 18-66 229, 317 18-115 154, 229, 312, 544 34-42 34, 309, 542 34-66 229, 317 18-116 154, 229, 314 34-42 34, 309, 542 34-66 229, 317 18-116 154, 229, 314 34-42 34, 309, 542 34-66 229, 317 18-116 154, 229, 546 34-42 34, 309, 542 34-66 229, 317 18-116 154, 229, 546 34-42 34, 309, 542 34-66 229, 34, 34, 34, 34, 34, 34, 34, 34, 34, 34		18-350	18-542
18-36	18-34	18-370	18-550
18-46	18-35	18-372	18-562
18-47	18-36	18-374	18-564
18-50	18-46	18-376	18-566
18-54	18-47	18-380	18-568
18-62	18-50	18-384	18-570
18-63	18-54	18-388	18-573A
18-64	18-62	18-390	18-574
18-66	18-63	18-398	18-576
18-86	18-64	18-398M	18-584
18-80	18-66	18-399	18-613
18-84		18-400	18-620
18-88	18-80	18-402	18-624
18-96	18-84	18-402A	18-628
18-100	18-88	18-402B	18-632
18-102	18-96	18-404	18-650
18-110	18-100	18-412	18-657
18-112			
18-114 154, 229, 546 18-420 189, 234, 309 18-680 240 18-116 154, 229, 546 18-424 189, 234, 309 18-680 240 18-120 229 18-428 189, 234, 309 18-691 240 18-132 230 18-439 155, 235, 310, 544 18-692 240 18-132A 230 18-440 189, 235 18-694 240 18-134 230 18-442 189, 235 18-696 240 18-134 230 18-446 189, 235 18-696 240 18-139 230 18-450 189, 235 18-700 240 18-140 231 18-452 189, 235 18-700 240 18-142 231 18-452 189, 235 18-700 240 18-144 231 18-452 189, 235 18-700 240 18-140 231 18-452 189, 235 18-710 145, 240 18-144 232 18-544		18-416	
18-116 154, 229, 546 18-424 189, 234 18-688 240 18-120 229 18-428 180, 234, 309 18-691 240 18-132 230 18-439 155, 235, 310, 544 18-692 240 18-132A 230 18-440 189, 235 18-694 240 18-134A 230 18-446 189, 235 18-698 240 18-134A 230 18-450 189, 235 18-698 240 18-139 230 18-450 189, 235 18-698 240 18-140 231 18-452 189, 235 18-700 240 18-142 231 18-452 189, 235 18-710 145, 240 18-144 231 18-452 189, 235 18-710 145, 240 18-144 231 18-452 189, 235 18-710 145, 240 18-148 231 18-480 235 18-718 145, 240 18-148 231 18-480 235	18-112	18-418 189, 234, 309	18-666
18-120 229			
18-132			18-688
18-132A 230 18-440 189, 235 18-694 240 18-134 230 18-442 189, 235 18-696 240 18-134A 230 18-446 189, 235 18-698 240 18-139 230 18-450 189, 235 18-700 240 18-140 231 18-452 189, 235 18-710 145, 240 18-142 231 18-454 189, 235 18-712 145, 240 18-144 231 18-472 236 18-714 145, 240 18-144 231 18-480 235 18-716 145, 240 18-148 231 18-480 235 18-716 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-720 145, 240 18-170 231 18-498 236 18-720 145, 240 18-171 231 18-499 236 18-726			
18-134 230 18-442 189, 235 18-696 240 18-134A 230 18-446 189, 235 18-698 240 18-139 230 18-450 189, 235 18-700 240 18-140 231 18-452 189, 235 18-710 145, 240 18-142 231 18-454 189, 235 18-712 145, 240 18-144 231 18-472 235 18-714 145, 240 18-144 231 18-480 235 18-714 145, 240 18-146 231 18-480 235 18-714 145, 240 18-148 231 18-480 235 18-714 145, 240 18-152 231 18-480 235 18-720 145, 240 18-154 231 18-480 235 18-720 145, 240 18-157 231 18-486 235 18-720 145, 240 18-164 231 18-498 236 18-722			18-692
18-134A 230 18-446 189, 235 18-698 240 18-139 230 18-450 189, 235 18-700 240 18-140 231 18-452 189, 235 18-710 145, 240 18-142 231 18-454 189, 235 18-712 145, 240 18-144 231 18-472 235 18-714 145, 240 18-146 231 18-480 235 18-716 145, 240 18-148 231 18-482 235 18-718 145, 240 18-148 231 18-482 235 18-718 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-722 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728			18-694
18-139 230 18-450 189, 235 18-700 240 18-140 231 18-452 189, 235 18-710 .145, 240 18-142 231 18-454 189, 235 18-712 .145, 240 18-144 231 18-472 235 18-714 .145, 240 18-146 231 18-480 235 18-716 .145, 240 18-148 231 18-482 235 18-718 .145, 240 18-152 231 18-484 235 18-720 .145, 240 18-154 231 18-486 235 18-722 .145, 240 18-170 231 18-498 236 18-722 .145, 240 18-171 231 18-499 236 18-724 .145, 240 18-172 231 18-499 236 18-726 .145, 240 18-292 232 18-500 236 18-728 .145, 240 18-292 232 18-502 236 18-7			
18-140 231 18-452 189, 235 18-710 145, 240 18-142 231 18-454 189, 235 18-712 145, 240 18-144 231 18-472 235 18-714 145, 240 18-146 231 18-480 235 18-716 145, 240 18-148 231 18-482 235 18-720 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-720 145, 240 18-170 231 18-498 236 18-722 145, 240 18-171 231 18-499 236 18-722 145, 240 18-171 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-728 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-504 236 18-732			
18-142 231 18-454 189,235 18-712 145,240 18-144 231 18-472 235 18-714 145,240 18-146 231 18-480 235 18-718 145,240 18-148 231 18-482 235 18-718 145,240 18-152 231 18-484 235 18-720 145,240 18-164 231 18-486 235 18-722 145,240 18-170 231 18-498 236 18-722 145,240 18-171 231 18-499 236 18-726 145,240 18-172 231 18-50 236 18-728 145,240 18-259 232 18-502 236 18-728 145,240 18-260 232 18-504 236 18-730 145,240 18-262 232 18-504 236 18-732 145,240 18-262 232 18-506 236 18-740 145,240 <td></td> <td></td> <td></td>			
18-144 231 18-472 235 18-714 145, 240 18-146 231 18-480 235 18-716 145, 240 18-148 231 18-482 235 18-718 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-726 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-502 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-262 232 18-510 236 18-740 145, 240 18-264 232 18-514 236 18-744 1			
18-146 231 18-480 235 18-716 145, 240 18-148 231 18-482 235 18-718 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-728 145, 240 18-260 232 18-502 236 18-730 145, 240 18-262 232 18-504 236 18-732 145, 240 18-262 232 18-504 236 18-740 145, 240 18-264 232 18-510 236 18-740 145, 240 18-270 232 18-512 236 18-744 1			
18-148 231 18-482 235 18-718 145, 240 18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-728 145, 240 18-260 232 18-504 236 18-730 145, 240 18-262 232 18-504 236 18-730 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-740 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-744 1			
18-152 231 18-484 235 18-720 145, 240 18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-730 145, 240 18-262 232 18-504 236 18-730 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748			
18-154 231 18-486 235 18-722 145, 240 18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-726 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-740 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-746 145, 240 18-304 232 18-520 236, 525 18-748 145, 240 18-308 232 18-524 236, 525 18-760 <td></td> <td></td> <td></td>			
18-170 231 18-498 236 18-724 145, 240 18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-504 236 18-732 145, 240 18-264 232 18-510 236 18-740 145, 240 18-270 232 18-512 236 18-742 145, 240 18-292 232 18-514 236 18-742 145, 240 18-292 232 18-512 236 18-744 145, 240 18-302 232 18-514 236 18-746 145, 240 18-304 232 18-522 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-778			
18-171 231 18-499 236 18-726 145, 240 18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-52 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525			
18-172 231 18-500 236 18-728 145, 240 18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-760 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236			
18-259 232 18-502 236 18-730 145, 240 18-260 232 18-504 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-748 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-331 232 18-536 236, 525 18-781 64, 241, 320, 539			
18-260 232 18-504 236 18-732 145, 240 18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 226, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-304 232 18-522 236, 525 18-779 64, 241, 320, 539 18-308 232 18-524 236, 525 18-789 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-330 232 18-532 236, 525 18-781 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-262 232 18-506 236 18-740 145, 240 18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-781 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-264 232 18-510 236 18-742 145, 240 18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-270 232 18-512 236 18-744 145, 240 18-292 232 18-514 236 18-746 145, 240 18-302 233 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-781 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-292 232 18-514 236 18-746 145, 240 18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-302 232 18-520 236, 525 18-748 145, 240 18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-304 232 18-522 236, 525 18-750 145, 240 18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-308 232 18-524 236, 525 18-779 64, 241, 320, 539 18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-314 232 18-530 236, 525 18-780 64, 241, 320, 539 18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-320 232 18-532 236, 525 18-781 64, 241, 320, 539 18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-330 232 18-534 236, 525 18-782 64, 241, 320, 539 18-331 232 18-536 236, 525 18-783 64, 241, 320, 539			
18-331			
10-040			
	10-0-1	10-040	10-704







CAT. NO	CAT. NO	CAT. NO
18-1554	18-1950-8	19-70
18-1558 43, 216, 260	18-1952-7	19-71
18-1559	18-1952-8	19-71-1
18-1560	18-1952-9	19-71-2
18-1562 47, 260	18-1952-10	19-71-3
18-1570 47, 261	18-1960	19-74
18-1572 47, 261	18-1962	19-75
18-1576 47, 261	18-1964	19-78
18-1578	18-1968 268, 295, 308, 351, 434	19-100
18-1594 47, 261	18-1990	19-102
18-1596 47, 261	18-1992 268, 285, 303	19-104
18-1600	18-1994	19-106
18-1602	18-1996	19-108
18-1604	18-2200	19-110
18-1618	18-SC-1472	19-112
18-1619	18-SC-1476	19-120 2, 272
18-1620	19-2	19-120A
18-1622	19-2-1	19-170
18-1623	19-2-2	19-172
18-1624	19-2-3	19-180
18-1625	19-2-4	19-182
18-1626	19-10	19-190
18-1630	19-10A	19-192
18-1631	19-12	19-200
18-1633	19-14	19-202
18-1634	19-16	19-210
18-1652 41, 140, 263, 333, 531	19-18	19-220
18-1653 41, 140, 263, 333, 531	19-20	19-222
18-1716	19-21	19-226
18-1718	19-22	19-228
18-1800 120, 264, 345, 535	19-23	19-240
18-1800TC	19-24	19-242
18-1810 120, 264, 345, 535	19-25	19-244
18-1816TC	19-30	19-246
18-1820 120, 264, 345, 536	19-32	19-248
18-1824 120, 264, 345, 536	19-34	19-260 274, 525
18-1828 120, 264, 345, 536	19-36	19-262
18-1828A-TC 129, 266	19-40	19-264
18-1828TC	19-48	19-266 274, 525
18-1831	19-48-A	19-268
18-1832 120, 264, 345, 536	19-48-B	19-270 274, 525
18-1832A-TC	19-48-C	19-272 274, 525
18-1832TC	19-48-D	19-274
18-1833	19-48-E	19-276 274, 525
18-1834	19-50-1	19-278
18-1837	19-50-2	19-280
18-1838	19-50-3	19-282
18-1839	19-50-4	19-284
18-1858	19-50-5	19-286
18-1934	19-50-7	19-288
18-1936	19-50-8	19-290
18-1950-7	19-50-9	19-294
	10 00 0	10-2042/4






CAT. NO. PAGE	CAT. NO	CAT. NO
19-296	19-580 149, 279, 426	19-1234
19-298	19-582 149, 279, 426	19-1330
19-300	19-584	19-1334
19-302	19-589	19-1336
19-314	19-590	19-1340
19-316	19-655	19-1350
19-318	19-660	19-1352
19-350	19-668	19-1354
19-356	19-669	19-1356
19-360	19-700 284, 436, 456	19-1820
19-370	19-702 284, 436, 456	19-2060
19-371	19-704 284, 436, 456	19-2081-B221, 288
19-380	19-706 284, 436, 456	19-2090-B221, 288
19-381	19-708	19-2100
19-382	19-710	19-2100-B221, 288
19-383	19-712 284, 436, 456	19-2110
19-384	19-714 284, 436, 456	19-2110-B221, 288
19-385	19-716	19-2120
19-390	19-718	19-2130
19-394	19-720	19-2138
19-396	19-750	19-2150-B 49, 217, 288
19-400	19-752	20-2
19-404	19-754	20-3
19-410	19-756 284, 449	20-4
19-414	19-758	20-6
19-429	19-760 284, 449	20-7
19-430	19-770	20-12
19-434	19-772	20-20
19-437	19-774	20-24
19-438	19-776	20-28
19-439	19-778	20-29
19-442	19-780 284, 449	20-30
19-464	19-784	20-32
19-465	19-785	20-36-35
19-468	19-790 284, 316, 453	20-36-50
19-500	19-792	20-36-75
19-502	19-793	20-36-90
19-504	19-794	20-37-45
19-510	19-795	20-37-75
19-512	19-796	20-100
19-520	19-797	20-102
19-524	19-804	20-120
19-562	19-808 285, 353, 414, 459	20-124
19-563	19-840	20-126
19-564	19-844	20-130
19-565	19-848	20-140
19-570	19-849	20-142
19-572	19-850	20-144
19-573	19-852	20-148
19-574	19-854	20-150
19-575	19-856	20-160
19-575	19-1230	20-164
19-576	19-1232	20-168
19-077 140, 279, 420	19-1202	20-10048, 234



CAT. NOPAGE	CAT. NO	CAT. NO
20-177	20-500	20-926
20-178	20-504	20-928
20-180	20-506	21-2
20-182	20-510	21-4
20-184	20-522	21-8
20-186-4	20-530	21-16
20-186-5	20-532	21-18
20-200	20-544	21-20
20-202	20-546	21-24
20-216	20-552	21-28
20-218	20-553	21-32
20-220	20-554	21-40
20-222	20-556	21-60
20-224	20-570	21-62
20-230	20-572	21-63
20-238	20-580	21-64
20-242	20-582	21-68
20-246	20-584	21-72
20-258	20-594	21-85
20-270	20-596	21-86
20-274	20-598	21-87
20-278	20-630	21-88
20-280	20-660	21-89
20-282	20-662	21-90
20-284	20-800	21-92 190, 309, 543
20-300	20-802	21-94 190, 309, 543
20-308	20-804	21-102 191, 309, 543
20-312	20-806	21-103
20-330	20-808	21-107
20-332	20-809	21-108
20-334	20-835	21-144 191, 311, 543
20-336	20-836	21-148
20-338	20-838	21-150
20-340	20-840	21-152
20-342	20-850	21-153 191, 311, 543
20-360	20-852	21-154 191, 311, 543
20-378	20-854	21-156
20-380	20-856	21-158
20-382	20-858-0	21-160 191, 311, 543
20-400	20-858-1	21-163
20-404	20-858-2	21-165
20-420	20-858-3	21-176
20-424	20-872	21-180 159, 313, 433
20-428	20-874	21-182
20-430	20-876	21-183
20-436	20-878-1305	21-184
20-438	20-878-2305	21-185
20-440	20-878-3305	21-188
20-442	20-890	21-190
20-450	20-910	21-205-2
20-456	20-912	21-205-4
20-460	20-922	21-205-6
20-466	20-924	21-205-8







CAT. NO PAGE	CAT. NO PAGE	CAT. NO
23-160	23-1210SS371	24-804
23-162	23-1212SS371	24-806
23-200	23-1214SS	24-808
23-230	23-1216SS	24-900
23-232	23-1218SS371	24-1000
23-282	23-1219SS371	24-1002
23-286	23-1243SS371	24-1100
23-320364	23-1244SS371	24-1102
23-322	23-1245SS371	24-1150
23-324	23-1246SS371	24-1152
23-326	23-1247SS371	24-1180
23-399	23-1248SS371	24-1190
23-400	23-1249SS371	24-1200
23-402	23-1250SS371	24-1250
23-410	23-1251SS371	24-1300
23-420	23-1350SS372	24-1400
23-600	23-1351SS	24-1404
23-622	23-1352SS	24-1408
23-640	23-1360ST372	24-1412380
23-642	23-1361ST372	24-1500381
23-644	23-1402369	24-1600
23-650	23-1420	24-1602
23-680	23-1421369	24-1604
23-690	23-1430	24-1606
23-720	23-1600372	24-1650
23-755	23-1602372	24-1654
23-756	24-550	24-1656
23-804	24-55282, 374	24-2000
23-1047	24-554	24-2010
23-1048	24-556	24-2020
23-1049	24-570	24-2030
23-1050	24-570XL	24-2040
23-1051	24-57282, 374	24-2050383
23-1060	24-574	24-2100384
23-1139SS371	24-58082, 374	24-2110
23-1140SS371	24-582	24-2120384
23-1142SS	24-600	24-2500
23-1144SS371	24-602	24-2700
23-1146SS	24-604	24-2710
23-1148SS	24-606	24-2720
23-1150SS	24-650	24-2960
23-1152SS	24-652	24-2961
23-1154\$\$	24-654	24-2962
23-1156SS	24-656	24-3000
		24-3010
23-1159SS	24-702	24-3020
		24-3200
23-1200SS	24-750	24-3210
23-1202SS	24-752	24-3220
23-120455	24-756	24-3230387 24-3240387
23-1208SS	24-736	
LO-120000	24-802	24-3250387 24-3260387
	۵۹-۵02	24-020038/



CAT. NO PAGE	CAT. NO	CAT. NO
24-3270387	25-782	26-62411
25-2	25-784	26-91171, 412
25-4	25-786	26-92171, 412
25-6	25-800	26-119
25-14	25-81299, 403	26-120
25-18	25-814	26-122
25-22	25-82199, 403	26-124
25-24	25-83499, 403	26-130
25-28	25-83599, 403	26-132
25-98	25-836	26-134
25-100	25-83799, 403	26-138
25-102389	25-838	26-140413
25-104389	25-840	26-141413
25-106	25-848	26-142413
25-130	25-880404	26-143413
25-132	25-882	26-190
25-134	25-1100	26-192
25-140	25-1104	26-198
25-142	25-1108	26-200414
25-144	25-1160	26-220414
25-150	25-1240B	26-222
25-152	25-1240C	26-255
25-154	25-1280	26-256416 26-258416
25-158	25-1420	26-261
25-170	26-16	26-262
25-170	26-18	26-264
25-200	26-20	26-265
25-202	26-24	26-266
25-204	26-28A	26-270417
25-272	26-28B	26-326
25-274	26-28C411	26-327418
25-302	26-28D	26-328
25-306	26-28E411	26-329418
25-350394	26-30410	26-330419
25-360	26-32	26-331419
25-362	26-33-1	26-332419
25-378	26-33-2	26-334419
25-384	26-33-3	26-335419
25-385	26-33-4	26-336419
25-390	26-35	26-337419
25-392	26-40410, 445	26-338419
25-396	26-42	26-339A419
25-440	26-44-1	26-339B
25-454	26-44-2	26-340419
25-456	26-44-3	26-370
25-458	26-44-4	26-388
25-480	26-46-1	26-389
25-510	26-46-2	26-390
25-511	26-46-4	26-392A
25-760	26-48	26-392B
25-780	26-60	26-393
20-700100	£0-00411	E0.000



CAT. NO	CAT. NO PAGE	CAT. NO
26-393A420	26-1042	27-7
26-393B	26-1043	27-11
26-394421	26-1054431	27-12441
26-395421	26-1080	27-13
26-396421	26-1084	27-14
26-397421	26-1086	27-16
26-398421	26-1090	27-20
26-399	26-1092	27-21
26-400	26-1150	27-22
26-402	26-1156	27-23
26-410	26-1170	27-24
26-412	26-1174	27-30
26-420	26-1326	27-32
26-430	26-1354	27-34
26-435	26-1360	27-44
26-436	26-1362	27-44
26-437	26-1403	27-47
26-440	26-1410	27-60
26-461	26-1412	27-110
26-462	26-1430	27-144
26-466	26-1432	27-158
26-470	26-1434	27-159TC
26-471	26-1436	27-160
26-472	26-1450435	27-161TC444
26-474423	26-1451	27-163TC
26-480422	26-1452	27-164TC
26-481422	26-1453435	27-165TC
26-482422	26-1454	27-175TC444
26-483	26-1620	27-176TC142, 444
26-484	26-1624	27-180445
26-485	26-1627	27-202445
26-590	26-1628	27-204
26-591	26-1629	27-206
26-660	26-1630	27-208
26-700	26-1632	27-210
26-710	26-1634	27-212
26-724	26-1636	27-213
26-730	26-1640	27-214
26-732	26-1650	27-216
26-760	26-1652	27-217
26-762	26-1654	27-218
26-770	26-1656	27-219
26-772	26-1658	27-223
26-774	26-1660	27-224
26-778	26-1973	27-225
26-780	26-2000	27-230
26-782	26-2001	27-232
26-95682, 429	26-2005	27-260
26-95882, 429	27-2	27-262
26-1040	27-4	27-264
26-1041	27-6	27-270449




CAT. NO	CAT. NO	CAT. NO
27-272449	27-496451	27-753
27-274449	27-498451	27-754455
27-280449	27-499451	27-820455
27-282449	27-500451	27-821455
27-331450	27-502	27-822
27-332	27-504	27-824455
27-333	27-506	27-826
27-334	27-508451	27-920
27-335	27-510451	27-924
27-336450	27-512	27-926
27-418450	27-514	27-928456
27-420450	27-516451	27-970457
27-422450	27-518451	27-972457
27-424450	27-520451	27-974457
27-426450	27-522451	27-980457
27-427450	27-530451	27-984457
27-428450	27-532451	27-990
27-430450	27-534451	27-992
27-432450	27-536451	27-1000
27-438450	27-538	27-1016458
27-438C	27-540451	27-1020458
27-440450	27-542451	27-1022
27-440C450	27-549	27-1024
27-442450	27-550	27-1026458
27-442C450	27-551	27-1028458
27-444	27-553	27-1295
27-444C	27-554	27-1297
27-446	27-556	27-1390462
27-446C	27-558	27-1392
27-447	27-599	27-1394
27-447C	27-600	27-1400
27-448	27-601	27-1402
27-448C	27-601A	27-1404
27-449	27-602A	27-1560
27-4490	27-603	27-1561
27-450	27-604	27-1564
27-450	27-605	27-1565-12
27-452C	27-609	27-1565-18
27-458	27-610	27-1600
27-45BC	27-660	27-1602
27-460	27-662	27-3002
27-460C	27-680	27-3004
27-462	27-700	27-3020
27-462C	27-702	27-3022
27-464C451	27-703	27-3059
27-466	27-744	27-3060
27-466C451	27-746	27-3062
27-468	27-747	27-3070
27-488451	27-748	27-3080
27-490	27-749	27-3090
27-492451	27-750455	27-3094
27-494451	27-752455	27-3100



CAT. NO	CAT. NO PAGE	CAT. NO PAGE
28-1	28-296	29-8-14
28-2	28-298	29-8-16
28-4	28-300	29-8-18
28-6	28-302	29-8-20
28-8174, 470	28-303	29-8-22
28-10	28-304	29-8-24
28-11	28-306	29-8-26
28-12		
		29-8-28
28-14	28-309	29-8-30
28-20	28-310	29-8-32
28-22	28-320	29-8-34
28-24	28-322	29-10-8
28-30	28-323	29-10-10
28-32	28-324	29-10-12
28-34	28-326480	29-10-14
28-36	28-328	29-10-16
28-38	28-329	29-10-18
28-40	28-330480	29-10-20
28-44	28-340	29-10-22
28-48	28-341	29-10-24
28-50	28-360	29-10-26
28-52	28-370	29-10-28
28-64	28-380	29-10-30
28-65	28-382	29-14-8
28-68	28-384	29-14-10
28-69175, 473	28-386	29-14-12
28-70	28-400	29-14-14
28-72	28-410	29-14-16
28-74	28-412	29-14-18
28-79	28-414	29-14-20
28-80	28-416	29-14-22
28-100	29-2	29-14-24
28-104	29-4	29-14-24
28-104	29-6-8	29-14-26
28-108	29-6-10	29-14-20
28-109	29-6-12	29-14-32
28-120	29-6-14	29-14-34
28-122	29-6-16	29-14-36
28-124	29-6-18	29-18-16
28-128	29-6-20	29-18-18
28-130	29-6-22	29-18-20
28-131	29-6-24	29-18-22
28-140	29-6-26	29-18-24
28-150475	29-6-28	29-18-26
28-151475	29-6-30	29-18-28
28-152475	29-6-32	29-18-30
28-153	29-6-34	29-18-32
28-196	29-6-36	29-18-34
28-200	29-6-38	29-18-36
28-240	29-6-40	29-32-8
28-244	29-8-8	29-32-10
28-280	29-8-10	29-32-12
28-294	29-8-12	29-32-14







CAT. NO PAGE	CAT. NO	CAT. NO
30-835	30-1270-3	30-2230517
30-840	30-1270-4	30-2235
30-845	30-1275-1	30-2250518
30-850	30-1275-4	30-2260518
30-855	30-1350507	30-2300
30-940	30-1351	30-2305518
30-950	30-1355	30-2315
30-951	30-1380508	30-2320519
30-952	30-1382	30-2330519
30-953	30-1383	30-2335519
30-954	30-1385	30-2345
30-955	30-1390508	30-2350
30-965	30-1391508	30-2355
30-970500	30-1392	30-2390520
30-975	30-1395508	30-2400520
30-985	30-1400	30-2410
30-988	30-1405	30-2415
30-990	30-1410509	30-2430
30-997	30-1420509	30-2435
30-998	30-1425	30-2445
30-1000502	30-1430	30-2450
30-1005	30-1435	30-2505
30-1010502	30-1440	30-2530
30-1020	30-1442510	30-2535523
30-1025	30-1443510	30-2560
30-1030	30-1445510	30-2575
30-1050-6	30-1480510	30-2580
30-1050-8	30-1482	30-2605
30-1050-10	30-1485	30-2632523
30-1077503	30-1486	33-2
30-1150	30-1487	33-3
30-1160	30-1490	33-4
30-1185	30-1491	33-5
30-1190	30-1492	33-5A
30-1205-0506	30-1493	33-5B524
30-1205-00	30-1494	33-6
30-1205-1506	30-1495	33-8
30-1205-2	30-1700	33-10
30-1205-3	30-1710	33-11
30-1205-5	30-1715	33-11A
30-1205-6	30-1716	33-12
30-1225-1	30-1718	
30-1225-2	30-1730	33-14
30-1225-3	30-1730	33-16-0
30-1225-4	30-1740	33-16-0
30-1225-5	30-1805	33-16-1
30-1225-6	30-1815	33-16-2
30-1255-2	30-2190	33-16-3
30-1255-4	30-2195	33-20
30-1255-6	30-2200	33-22
30-1270-1	30-2205	33-24
30-1270-2	30-2210	33-25



CAT. NO	CAT. NO	CAT. NO
33-26	40-58/0	MS117-2
33-28	40-58/1	MS117-3
33-30	40-58/1-2525	MS117-4
33-30-1	40-58/1-3	MS117-5
33-30-2	40-58/2	MS117-6
33-30-3	40-58/3	MS117-7
33-30-4	40-94	MS117-8
33-30-5	40-96	MS126-2
33-30-6	40-200	MS126-4
33-30-7	40-205	MS126-6
33-31	40-207	MS126-8
33-31A	40-210-CH549	MS126-10
33-32	40-210-SS	MS126-12
33-33	40-212-CH549	MS126-14
33-34	40-212-SS	MS126-16
33-35	40-215-CH549	MS126-18
33-36	40-215-SS	MS126-20
33-37	40-217	MS126-22
33-40	40-219549	MS140-0
33-121	40-220	MS140-1
33-140546	40-224550	MS140-2
33-142546	40-225-CH550	MS140-3
33-200547	40-225-SS	MS140-4
33-201547	40-225A550	MS140-5
33-202	40-226-CH550	MS140-6
33-203547	40-226-SS	MS140-7
33-204547	40-226A	MS140-8
33-206	40-227-CH550	MS140-9
33-210547	40-227-SS	MS140-10132
33-214548	40-228	MS140-11
33-215	40-230	MS140-12
33-216	40-234	MS140-13
33-230548	40-245-CH	MS140-14
33-232	40-245-SS	MS140-16
33-240	40-250-CH	MS140-17
33-242	40-250-SS	MS140-1B
33-280	40-251-CH551	MS140-19
33-500	40-251-SS	MS140-20
34-50	40-252	MS140A-1
34-52	40-260	MS140A-2
34-80	40-425	MS140A-3
34-9250, 552	40-435	MS140A-4
34-98	40-445	MS140A-5
34-210	40-1235	MS140A-6
34-220553	40-2000	MS140A-7
34-250	63-10	MS140A-8
34-256553	63-15	MS140A-9
40-54524	63-20	MS140A-10
40-57/0	65-25	MS140A-11
40-57/1	73-62308	MS140A-12
40-57/2	DELS23	MS140A-13
40-57/3	MS117-1133	MS140A-14



CAT. NO	CAT. NO PAGE	CAT. NO
MS140A-15	MS155-3	MS192A-7
MS140A-16	MS155-4	MS192A-8
MS140A-17132	MS155-5	MS192A-9
MS140A-17	MS156-1	MS192A-9
MS140A-19	MS156-2	MS192A-11
		MS192A-12
MS140S-0	MS156-4	MS192A-13
MS140S-1	MS156-5	MS192A-14
MS140S-2	MS172-1	MS192A-15
MS140S-3	MS172-2	MS192A-16
MS140S-5	MS172-4	MS192A-18
MS140S-7	MS175-1	
MS140S-8	MS175-3	MS192A-20
MS1405-8	MS175-4	MS192S-0
MS140S-10132	MS176-1	MS192S-2
MS140S-11	MS176-2	MS192S-3
MS140S-12	MS176-3	MS192S-4
MS140S-13	MS176-4	MS192S-5
MS140S-14	MS177	MS192S-6
M\$140\$-15	MS179-1132	MS192S-7
MS140S-16	MS179-2	MS192S-8
		MS192S-9
MS140S-18	MS179-4	MS192S-10
MS140S-20	MS192-0	MS192S-11
MS141-0	MS192-1	MS192S-12
MS141-1		
MS141-2	MS192-3	MS192S-14
MS141-3	MS192-4	MS1925-15
MS141-3	MS192-6	MS1925-16
MS141-5	MS192-7	MS192S-18
MS141-6	MS192-8	MS192S-19
MS141-7	MS192-9	MS1925-19
MS141-8	MS192-10	MS431-4
MS141-9	MS192-10	MS431-3
MS141-10	MS192-12	MS431-2%
MS141-10	MS192-12	MS431-2½
MS141-12	MS192-14	MS431-2
MS150-1	MS192-14	MS431-134
MS150-2	MS192-16	MS431-1½
MS150-3	MS192-16	MS450-2/0
MS150-4	MS192-17	MS450-2/0
MS150-5	MS192-19	MS450-1
MS150-6	MS192-20	MS450-2
MS150-7	MS192A-0	MS450-6
MS153-1¾	MS192A-1	MS450-7
M\$153-174	MS192A-2	MS450-7
MS153-1	MS192A-3	MS450-9
MS153-4	MS192A-4	MS1274-2
MS155-1	MS192A-5	MS1274-2
MS155-2	MS192A-6	MS1274-4
MOTOUTE	WIG102A-U	134 July 1274-0





CAT. NO	CAT. NO	CAT. NO
MS1281-1134	MS1291-1	MS1295-8
MS1281-2	MS1291-2	MSRC-1413
MŠ1281-3	MS1291-3	MSRC-1813
MS1281-4	MS1291-4	V96-329
MS1281-5	MS1291-5	V96-410
MS1281-6	MS1291-6	V97-584
MS1281-8	MS1295-1	V97-585
MS1282-1	MS1295-2	V97-5881
MS1282-3	MS1295-3	V97-5901
MS1282-4	MS1295-4	V97-5911
MS1282-5	MS1295-5	V930-85650
MS1282-6	MS1295-6	V930-99649
MS1989-8 194	MS1905-7 194	



.

.



589 Davies Drive
York, PA 17402
toll-free phone 866 854-8300
toll-free fax 866 854-8400
phone 717 840-9335
fax 717 840-9347
email surgical@miltex.com
www.miltex.com

Please call Customer Service today to receive other Miltex catalogs.



+HB34CAI SIMSUNGUN