

CaviWipes AF & CaviCide AF

(Alcohol Free)



CaviWipes™ AF and CaviCide™ AF are here. No other company offers quat-based alcohol free wipes and sprays. The CaviWipes AF and CaviCide AF family of products provide convenient packaging and ease of use in various sizes. They are recommended for use on hard, non-porous surfaces. When used as directed, they will effectively clean and disinfect surfaces and can help reduce the risk of cross contamination.

Features:

- Alcohol & Fragrance Free
- Can be used on materials that are sensitive to alcohol
- Clinically relevant efficacy claims including TB, Staph, Pseudomonas & MRSA
- Conveniently packaged and ready to use

You've trusted Metrex for over 30 years, now put your trust in CaviWipes AF and CaviCide AF.

For more information on CaviWipes AF and CaviCide AF, visit metrex.com or call 800-841-1428.




CaviWipes™ AF & CaviCide™ AF (Alcohol Free)

Selected Pathogens:	CaviWipes AF	CaviCide AF
Mycobacterium bovis BCG (Tuberculosis)	5 minutes	3 minutes
Staphylococcus aureus	5 minutes	3 minutes
Pseudomonas aeruginosa	5 minutes	3 minutes
Salmonella (choleraesuis) enterica	5 minutes	3 minutes
Staphylococcus aureus (MRSA)	5 minutes	3 minutes
Vancomycin Resistant Enterococcus faecalis (VRE)	5 minutes	3 minutes
Escherichia coli (E. coli)	5 minutes	3 minutes
Influenza A (H1N1) Virus	5 minutes	3 minutes
Hepatitis B Virus (HBV)	5 minutes	X
Hepatitis C Virus (HCV)	5 minutes	X
Human Immunodeficiency Virus Type 1 (HIV-1)	2 minutes	1 minute
Streptococcus pyogenes	X	3 minutes
Enterococcus faecium	X	3 minutes
Paramyxovirus (Mumps)	X	3 minutes
Rhinovirus type 39	X	3 minutes
Rotavirus	×	3 minutes
Norovirus	X	30 seconds

Ordering Information

Part No.	Product Description	Packaging
13-8008	CaviCide AF - 8 ounce spray	12 bottles/case
13-8024	CaviCide AF - 24 ounce spray	12 bottles/case
13-8000	CaviCide AF - 1 gallon	4 gallons/case
13-8025	CaviCide AF - 2.5 gallon	2 bottles/case
13-8100	CaviWipes AF Canister (6"x12") 80 wipes/canister	12 canisters/case
13-1175	CaviWipes Canister Wall Bracket	12 brackets/case

Now with CaviWipes AF and CaviCide AF, you get the same formulation in an easy to use wipe or liquid format designed for use in each of the seven critical touch points in the healthcare facility: from surgery to housekeeping.



For more information on CaviWipes AF and CaviCide AF, visit metrex.com or call 800-841-1428.