# the sleek, powerful efficient precise light



# **OUTPATIENT LED**EXAMINATION LIGHT

# The next generation of examination lighting has arrived.

Outpatient LED is the perfect low maintenance, top-performing light for examination applications in the hospital, physician's office and veterinary field.

## Designed for patient comfort and ease of use

- LED sources provide virtually heat-free light for patient comfort
- High quality light for accurate color rendition
- Designed to reduce shadows on the examination area
- Modularity for ease of installation and assembly

### **Enhanced performance**

- Increased illuminance to 90,000 lux (8361 fc) at 24" (61 cm)
- CRI (Color Rendering Index) > 95
- High R9 > 87, to accurately render skin tones
- Anti-microbial paint to inhibit the growth of microbes
- · Lower profile to reduce space impact


### Lower cost of ownership

- LED illumination offers 10x (or 20,000 hours average) the operational life of traditional halogen solutions
- Provides a 75% energy savings over current halogen option
- · Smooth finish aids routine cleaning
- · Reduced maintenance




# OUTPATIENT LED EXAMINATION LIGHT


# the next generation of exam lighting



**Floorstand** 



**Wall Mount** 




Single Ceiling Mount

### Summary of technical data 90,000 lux (8361 fc) at 24" (61 cm) Illuminance: Color Temperature: CRI (Color Rendering Index): 95 Light Source: Light Emitting Diode (LED) Rated Life of LED: 20,000 hours average LED life Power: 42 Watts Power Supply: Integral electronic driver Floorstand, Wall, Ceiling Mounting Options: Maximum Arm Reach: 37.5" (95.25 cm) Certifications/Approvals: AAMI ES60601-1, IEC 60601-1-2, EN/IEC 60601-2-41 and CAN/CSA C22.2 NO. 60601-1 certified 5-year warranty **Product Warranty:**

### Ordering Information

120V Model	Description
OPLEDFL	Floorstand
OPLEDSC	Single ceiling mount
OPLEDDC	Double ceiling mount
OPLEDW	Wall mount



**Double Ceiling Mount** 

© 2016 Koninklijke Philips N.V. All rights reserved. Philips reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Philips Burton 11500 Melrose Avenue Franklin Park, IL 60131 Burtonmedical.com

