Material Natural rubber from latex in compliance with ASTM D1076-06 containing dispersed bismuth oxide.

Packaging 1 Sterile Pair/Package - 5 Pairs/Case

Product Numbers

SIZE (identified inside glove)	ASPN (Ansell Standard Product Number)
6 1/2	20873065
7	20873070
7 1/2	20873075
8	20873080
8 1/2	20873085
9	20873090

Dimensional Specifications (minimum length 280mm)

SIZE	Average Length mm (inches)	Average Palm Width mm (inches)
6 1/2	280 (11.0)	83 (3.3)
7	280 (11.0)	89 (3.5)
7 1/2	285 (11.2)	95 (3.7)
8	285 (11.2)	102 (4.0)
8 1/2	290 (11.4)	108 (4.3)
9	290 (11.4)	114 (4.5)

Thickness (single)

Data from Sampling	Average Thickness
Fingertip, mm (mils)	0.35 (14)
Cuff, mm (mils)	0.30 (12)

Physical Properties (ASTM Requirements)

	Before Aging	After Aging
Freedom from Holes (inspection level I)	1 AQL*	N/A
Tensile, MPa (p.s.i.) minimum	12 (1,740)	12 (1,740)
Elongation (%) minimum	500	500
Force at Break (N)	21	18

Ansell manufacturing specifications meet or exceed all the above.

COIOI	Miaki
Accelerators	Zinc diethyldithiocarbamate (ZDEC) Zinc dibutyldithiocarbamate (ZDBC)
Donning Agent	None (powder-free)
Residual	< 2 mg/glove

Certification Manufacturing facility registered to ISO 9001 and ISO 13485

Contains 50 micrograms or less of total water-extractable protein per gram per glove

Sterilization Gamma radiation; Sterility Assurance Level (SAL) 1/1,000,000

References ISO 2859-1 Sampling Procedures and Tables for Inspection by Attributes

Country of Origin Malaysia

Storage Cool environment, shielded from ozone and UV light

*Acceptable Quality Level

Powder Content



For additional information and samples, call Ansell at 800-952-9916 or visit us at www.ansellhealthcare.com

Ansell Healthcare Products LLC I ansellhealthcare.com
200 Schulz Drive I Red Bank, NJ 07701 U.S.A. I Phone: (800) 952-9916 (U.S.A. only)
Ansell is a registered trademark owned by Ansell Limited or one of its affiliates. ©2011 Ansell Healthcare Products LLC. All Rights Reserved.

