Powered by

VALHALLA SCIENTIFIC









ADVANTAGES

15 Second Measurement Time

On screen display delivers rapid access to weight, body fat %, body water %, lean mass, BMR & BMI with a full illustrated printout in 60 seconds.

Large Platform

The spacious platform allows patients to find a comfortable and secure stance.

Few Pretest Requirements

Analysis can be performed at any time of day and under almost any circumstances, resulting in the greatest flexibility of any body composition analyzer.

Safe to Use

This scale is FDA registered.

High Capacity

The scale features an industry leading 1000 lb capacity, that provides the flexibility and accuracy to weigh any patient.

Patient Report Templates*

The G6 Series scales offer users 3 graphic printout formats to use with their patients: Adult/Bariatric, Pediatric and Fitness.

* User must connect scale to a printer to provide personalized graphic printout. Printer is optional and sold separately, not available through Health o meter® Professional scales.





Model	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
BCS-G6	46" x 19" x 10" / 10	37 lb	854896003702

Specifications

• Capacity: 1000 lb / 453.5 Kg • Resolution: 0.1 lb / 0.1 kg

• Platform Size: 14" (w) x 16" (d) x 1.875" (h)

• Display: 6-window display (.5" digits) includes weight, % body fat, % body water, % muscle mass, BMI, REE

• Power Source: 110-220V power adapter (included)

• Printer Connectivity: USB or WiFi (optional printer, sold separately, not available through Health o meter® Professional scales)

• 1 Year Limited Warranty

Bio-Impedance Specifications

• Electrode Method: Tetrapolar System

• Frequency: 50kHz • Test Current: 500µA

Conversion Rate: 750msec



Model	Description	UPC Code
BCS-G6-adult Body Comp Scale	The whole-body, adult only scale is ideal for practices that exclusively see patients 18 and over	854896003702
BCS-G6-limbs Body Comp Scale	The adult segmented version provides patients with body composition information about their torso as well as specific information about the composition of their arms and legs	854896003719
BCS-G6-DUO Body Comp Scale	The whole-body adult/pediatric version is a good choice for providers who see adults and children	854896003733
BCS-G61-upper Body Comp Scale	The adult upper-body-only configuration provides a partial body composition analysis for patients who do not want to remove their shoes and socks	854896003726
IPO-ADU	Patient Report Template - Adult (5-page report, 25 reports/package)	854896003740
IPO-PED	Patient Report Template - Pediatric (7-page report, 25 reports/package)	854896003757
IPO-ATH	Patient Report Template - Athletic (3-page report, 50 reports/package)	854896003764
BCS-WIRELESS	Body Comp Scale Wireless Adapter	854896003771

Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330



