



MODERATE TO HEAVY



FRAGRANCE FREE



HEAVILY SOILED

## Three Ways to Protect

THE FIRST AND MOST IMPORTANT STEP IN THE STERILIZATION PROCESS

Metrex now offers 3 different enzymatic detergents to meet your various needs in cleaning reusable medical instruments and devices. Product offering ranges from moderate to heavily soiled instruments and for the respiratory-sensitive.

All three detergents clean away blood, tissue, mucous and other protein rich fluids, while leaving your instruments clean and shiny. Neutral pH formulas provide advanced corrosion protection. They are compatible with stainless steel, soft metals, hard plastics, rubber and glass.



**NEW**

EmPower® Fragrance Free offers the same great enzymatic cleaning formula without the fragrance for customers who are respiratory-sensitive or prefer a non-scented detergent.

## EmPower® & EmPower® Fragrance-Free Dual-Enzymatic Detergents

Formulated for moderate to heavily soiled instruments. Removes protein rich debris in both manual and automated reprocessing systems. Contains powerful cleaning agents that will not get cloudy at the use temperature.

Order Information	Description	Packaging
10-4102	EmPower 2-ounce	48 bottles/case
10-4100	EmPower 1-gallon	4 bottles/case
10-4150	EmPower 5-gallon	1 drum
10-4115	EmPower 15-gallon	1 drum
10-4130	EmPower 30-gallon	1 drum
10-4400	EmPower Fragrance Free 1-gallon	4 bottles/case



## EmPower Foam Enzymatic Spray & EmPower® Foam Touch-n-Spray

The AAMI Guideline on "Safe Handling and Biological Decontamination of Reusable Medical Devices in Healthcare Facilities and in Non-Clinical Settings" recommend using products like EmPower Foam on instruments immediately after use to keep them moist in transport containers.\* EmPower Foam Touch-n-Spray battery-operated sprayer head provides a continuous spray with just one touch.

Order Information	Description	Packaging
10-4224	EmPower Foam - 24oz. spray	12 bottles/case
46-1510	Touch-n-Spray Trigger	10 sprayers/case

\*AAMI Guideline Transportation, section 6.2



## MetriZyme® Dual-Enzymatic Detergent

Offers powerful cleaning, yet gentle on instruments. MetriZyme is designed for superior performance in the presence of high organic debris while offering the broadest cleaning action on a variety of protein soils. Surfactants provide additional cleaning on carbohydrates, lipids and protein soils. MetriZyme includes an enzyme stabilizer, a buffer to maintain pH, a preservative, an antifoaming agent, and solution dye to differentiate it from water and other solutions.

Order Information	Description	Packaging
10-4005	MetriZyme - 1 quart	4 bottles/case
10-4010	MetriZyme - 1/2 gallon	4 bottles/case
10-4000	MetriZyme - 1 gallon	2 bottles/case
10-4050	MetriZyme - 5 gallon	1 drum
10-4015	MetriZyme - 15 gallon	1 drum
10-4030	MetriZyme - 30 gallon	1 drum
10-4055	MetriZyme - 55 gallon	1 drum



## MetriSponge®

MetriSponge is a contoured sponge saturated with extra strength MetriZyme. According to the 2008 SGNA guidelines, manual cleaning of endoscopes immediately after patient use is the first and most important step in removing the bioburden prior to automated or manual disinfection.<sup>1</sup> Easily forms around endoscopes, simplifying effective precleaning of surgical instruments and endoscopes.

Order Information	Description	Packaging
10-4025	MetriSponge	4 boxes/case

<sup>1</sup>Kelsey, Loralee. Et al. Standards of Infection Control in Reprocessing of Flexible Gastrointestinal Endoscopes. Society of Gastroenterology Nurses and Associates, Inc. 2008



## About Metrex Research

Metrex offers free static cling wall charts that provide step-by-step instructions on how to clean your instruments properly. To order your free wall chart, call 800.841.1428 and reference literature number 77-1001.

**Metrex**  
Protecting People