



WE CARE LIKE FAMILY™




NEW

See side of box for full instructions for use and safety information

A Complete Line of Syringes and IV Administration Sets

- *Safety Insulin Syringes*
- *Safety Syringes with Needle*
- *Dialysis Safety Syringes with Needle*
- *Non-Safety Syringes with Needle*
- *Non-Safety Insulin Syringes*
- *Syringes Without Needle*
- *Hypodermic Needles*
- *Pen Needles*
- *Blood Collection Tube Holders*



Find us on:



dynarex® is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative.

Or call **888-396-2739** to find your nearest distributor.

10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

dynarex COMPLETE LINE OF SYRINGES

• Insulin Syringes • Safety Syringes with Needle • Dialysis Safety Syringes with Needle • Syringes Without Needle • Non-Safety Insulin Syringes • Hypodermic Needles • Non-Safety Syringes with Needle • Blood Collection Tube Holders • Pen Needles • IV Administration Sets • IV Extension Sets

Safety Insulin Syringes With Needle



Reorder No.	Description	Quantity
6901	Safety Insulin Syringe With Needle, 0.5cc, 29G x 1/2"	5/100/Cs
6902	Safety Insulin Syringe With Needle, 0.5cc, 30G x 5/16"	5/100/Cs
6903	Safety Insulin Syringe With Needle, 1cc, 27G x 1/2"	5/100/Cs
6904	Safety Insulin Syringe With Needle, 1cc, 29G x 1/2"	5/100/Cs
6905	Safety Insulin Syringe With Needle, 1cc, 30G x 1/2"	5/100/Cs

Safety Syringes With Needle (Luer Lock or Fixed Needle)



Reorder No.	Description	Quantity
6915	With Needle, Luer Lock, 10cc, 20G x 1 1/2"	6/100/Cs
6916	With Needle, Luer Lock, 10cc, 21G x 1 1/2"	6/100/Cs
6917	With Needle, Luer Lock, 10cc, 22G x 1 1/2"	6/100/Cs
6918	With Needle, Luer Lock, 5cc, 20G x 1 1/2"	4/100/Cs
6919	With Needle, Luer Lock, 5cc, 21G x 1 1/2"	4/100/Cs
6920	With Needle, Luer Lock, 5cc, 22G x 1 1/2"	4/100/Cs
6921	With Needle, Luer Lock, 3cc, 20G x 1 1/2"	8/100/Cs
6922	With Needle, Luer Lock, 3cc, 20G x 1"	8/100/Cs
6923	With Needle, Luer Lock, 3cc, 21G x 1 1/2"	8/100/Cs
6924	With Needle, Luer Lock, 3cc, 21G x 1"	8/100/Cs
6925	With Needle, Luer Lock, 3cc, 22G x 1 1/2"	8/100/Cs
6926	With Needle, Luer Lock, 3cc, 22G x 1"	8/100/Cs
6927	With Needle, Luer Lock, 3cc, 23G x 1"	8/100/Cs
6928	With Needle, Luer Lock, 3cc, 25G x 5/8"	8/100/Cs
6929	With Needle, Fixed Needle, 1cc, 25G x 5/8"	8/100/Cs
6930	With Needle, Fixed Needle, 1cc, 27G x 1/2"	8/100/Cs

Dialysis Safety Syringes with Needle, Low Dead Space



Reorder No.	Description	Quantity
6931	Fixed Needle, 1cc, 22G x 1 1/2"	12/100/Cs
6932	Luer Lock, 3cc, 22G x 1 1/2"	12/100/Cs

Blood Collection Tube Holders



Reorder No.	Description	Quantity
6933	Without Needle, Non-Sterile	4/250/Cs
6934	Luer Lock with Needle	200/Box Sold by inner box
6935	Luer Slip with Needle	200/Box Sold by inner box

Syringes Without Needle



Reorder No.	Description	Quantity
6987	Syringes Without Needle, Luer Slip, 1mL/cc	10/100/Cs
6988	Syringes Without Needle, Luer Lock, 3mL/cc	10/100/Cs
6989	Syringes Without Needle, Luer Lock, 5mL/cc	10/100/Cs
6990	Syringes Without Needle, Luer Lock, 10mL/cc	10/100/Cs
6991	Syringes Without Needle, Luer Lock, 20mL/cc	10/50/Cs
6992	Syringes Without Needle, Luer Lock, 30mL/cc	10/50/Cs
6993	Syringes Without Needle, Luer Lock, 60mL/cc	10/25/Cs

Non-Safety Insulin Syringes



Reorder No.	Description	Quantity
6940	Non-Safety Insulin Syringe With Needle, 0.5cc, 27G x 1/2", Individually Packaged	5/100/Cs
6941	Non-Safety Insulin Syringe With Needle, 0.5cc, 28G x 1/2", Individually Packaged	5/100/Cs
6942	Non-Safety Insulin Syringe With Needle, 0.5cc, 29G x 1/2", Individually Packaged	5/100/Cs
6943	Non-Safety Insulin Syringe With Needle, 0.5cc, 30G x 5/16", Individually Packaged	5/100/Cs
6944	Non-Safety Insulin Syringe With Needle, 0.5cc, 31G x 5/16", Individually Packaged	5/100/Cs
6945	Non-Safety Insulin Syringe With Needle, 1cc, 27G x 1/2", Individually Packaged	5/100/Cs
6946	Non-Safety Insulin Syringe With Needle, 1cc, 28G x 1/2", Individually Packaged	5/100/Cs
6947	Non-Safety Insulin Syringe With Needle, 1cc, 29G x 1/2", Individually Packaged	5/100/Cs
6948	Non-Safety Insulin Syringe With Needle, 1cc, 30G x 1/2", Individually Packaged	5/100/Cs
6949	Non-Safety Insulin Syringe With Needle, 1cc, 31G x 5/16", Individually Packaged	5/100/Cs
6950	Non-Safety Insulin Syringe With Needle, 0.5cc, 27G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6951	Non-Safety Insulin Syringe With Needle, 0.5cc, 28G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6952	Non-Safety Insulin Syringe With Needle, 0.5cc, 29G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6953	Non-Safety Insulin Syringe With Needle, 0.5cc, 30G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs
6954	Non-Safety Insulin Syringe With Needle, 0.5cc, 31G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs
6955	Non-Safety Insulin Syringe With Needle, 1cc, 27G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6956	Non-Safety Insulin Syringe With Needle, 1cc, 28G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6957	Non-Safety Insulin Syringe With Needle, 1cc, 29G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6958	Non-Safety Insulin Syringe With Needle, 1cc, 30G x 1/2", 10/Bag, 10 Bags/Box	5/100/Cs
6959	Non-Safety Insulin Syringe With Needle, 1cc, 31G x 5/16", 10/Bag, 10 Bags/Box	5/100/Cs

Hypodermic Needles



Reorder No.	Description	Quantity
6960	Hypodermic Needles, 18G x 1"	10/100/Cs
6961	Hypodermic Needles, 18G x 1 1/2"	10/100/Cs
6962	Hypodermic Needles, 19G x 1"	10/100/Cs
6963	Hypodermic Needles, 19G x 1 1/2"	10/100/Cs
6964	Hypodermic Needles, 20G x 1"	10/100/Cs
6965	Hypodermic Needles, 20G x 1 1/2"	10/100/Cs
6966	Hypodermic Needles, 21G x 1"	10/100/Cs
6967	Hypodermic Needles, 21G x 1 1/2"	10/100/Cs
6968	Hypodermic Needles, 22G x 1"	10/100/Cs
6969	Hypodermic Needles, 22G x 1 1/2"	10/100/Cs
6970	Hypodermic Needles, 23G x 1"	10/100/Cs
6971	Hypodermic Needles, 23G x 1 1/2"	10/100/Cs
6972	Hypodermic Needles, 25G x 5/8"	10/100/Cs
6973	Hypodermic Needles, 25G x 1 1/2"	10/100/Cs
6974	Hypodermic Needles, 26G x 1/2"	10/100/Cs
6975	Hypodermic Needles, 27G x 1 1/2"	10/100/Cs
6976	Hypodermic Needles, 27G x 1/2"	10/100/Cs

Non-Safety Syringes With Needle

Reorder No.	Description	Quantity
6994	With Needle, Luer Slip, 1cc, 25G x 5/8"	8/100/Cs
6995	With Needle, Luer Lock, 3cc, 20G x 1"	8/100/Cs
6996	With Needle, Luer Lock, 3cc, 20G x 1 1/2"	8/100/Cs
6997	With Needle, Luer Lock, 3cc, 21G x 1"	8/100/Cs
6998	With Needle, Luer Lock, 3cc, 21G x 1 1/2"	8/100/Cs
6999	With Needle, Luer Lock, 3cc, 22G x 1"	8/100/Cs
7000	With Needle, Luer Lock, 3cc, 22G x 1 1/2"	8/100/Cs
7001	With Needle, Luer Lock, 3cc, 23G x 1"	8/100/Cs
7002	With Needle, Luer Lock, 3cc, 25G x 5/8"	8/100/Cs
7003	With Needle, Luer Lock, 3cc, 25G x 1"	8/100/Cs
7004	With Needle, Luer Lock, 5cc, 20G x 1"	4/100/Cs
7005	With Needle, Luer Lock, 5cc, 20G x 1 1/2"	4/100/Cs
7006	With Needle, Luer Lock, 5cc, 21G x 1"	4/100/Cs
7007	With Needle, Luer Lock, 5cc, 21G x 1 1/2"	4/100/Cs
7008	With Needle, Luer Lock, 5cc, 22G x 1"	4/100/Cs
7009	With Needle, Luer Lock, 5cc, 22G x 1 1/2"	4/100/Cs
7010	With Needle, Luer Lock, 10cc, 20G x 1"	4/100/Cs
7011	With Needle, Luer Lock, 10cc, 20G x 1 1/2"	4/100/Cs
7012	With Needle, Luer Lock, 10cc, 21G x 1"	4/100/Cs
7013	With Needle, Luer Lock, 10cc, 21G x 1 1/2"	4/100/Cs
7014	With Needle, Luer Lock, 10cc, 22G x 1"	4/100/Cs
7015	With Needle, Luer Lock, 10cc, 22G x 1 1/2"	4/100/Cs

Pen Needles



Reorder No.	Description	Quantity
7026	Pen Needles, 30G x 8mm	12/100/Cs
7027	Pen Needles, 31G x 5mm	12/100/Cs
7028	Pen Needles, 32G x 4mm	12/100/Cs

IV Administration Sets (sold by inner box)



Reorder No.	Description	Quantity
7030	10 Drop/mL, 83", 1 inj. site, 1 Luer lock site	50/Box
7031	15 Drop/mL, 83", 1 inj. site, 1 Luer lock site	50/Box
7032	15 Drop/mL, 78", 1 inj. site	50/Box
7033	15 Drop/mL, 106", 2 injection sites	50/Box
7034	15 Drop/mL, 73", 1 injection site	50/Box
7035	15 Drop/mL, 84", 2 Luer lock needle-free sites	50/Box
7036	15 Drop/mL, 85", no injection sites	50/Box
7037	15 Drop/mL, 92", 2 injection sites	50/Box
7048	20 Drop/mL, 100", 1 injection site with micron filter	50/Box
7049	20 Drop/mL, 92", 1 injection site with micron filter and regulator	50/Box
7050	20 Drop/mL, 92", 1 injection site with micron filter	50/Box
7051	60 Drop/mL, 83", 1 injection site, 1 Luer lock site	50/Box
7052	60 Drop/mL, 78", 2 injection sites	50/Box
7053	60 Drop/mL, 67", no injection sites	50/Box

IV Extension Sets (sold by inner box)



Reorder No.	Description	Quantity
7060	6", Needle-Free Luer Lock Connector	100/Box
7061	7", Needle-Free Luer Lock Connector	100/Box
7062	8", Needle-Free Luer Lock Connector	100/Box
7063	15", Needle-Free Luer Lock Connector	100/Box

See our website for all our medical products! www.dynarex.com