

Advantages

Durable Construction

The scale has an easy-to-clean stainless steel solid construction and durable, nonskid mat and can withstand frequent use

600 lb Capacity Floor Scale

The high capacity of this floor scale provides more versatile use

AC Adapter Included

Scale is powered by a 120V AC adaptor (included) or 6 C-Cell batteries (not included) which allows for the convenience of dual power sources

Motion Detection

Scale has the ability to weigh accurately despite movement of the subject on the scale



See reverse side for additional specifications



Model	Master Pack	Shipping Dims. (L x W x H) / Cube	Weight (LB)	UPC Code
2842KL	1	47.5" x 25.625" x 5.25" / 3.70	61.75 lb	892076002480



Specifications

- Capacity: 600 lb / 270 kg
- Resolution: 0.2 lb / 0.1 kg
- 7/8" LCD display
- Functions: LB / KG Conversion, Zero Out / Tare, Motion Detection
- Platform Size: 22¼" (w) x 42" (d) x 2½" (h) / 565 mm (w) x 1067 mm (d) x 64 mm (h)
- 120V adapter included, 6 C-cell batteries (not included)
- 2 Year Limited Warranty



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615 Fax: 1-800-654-7330